

## Improved HPLC Separation of Steroids

### Overview

The chromatograms below demonstrate the complete separation of a mixture of 12 steroids using an Advanced Chromatography Technologies ACE C18-AR column. Under these conditions, ACE C18 and ACE Phenyl columns both show two incompletely resolved critical pairs of peaks. The ACE C18-AR column combines the benefits of the hydrophobic interactions associated with a standard C18 phase with the additional aromatic ( $\pi$ - $\pi$ ) interactions commonly exhibited by Phenyl bonded phases.

### HPLC Conditions\*

Column: ACE 3 C18-AR (3 $\mu$ m, 150 x 4.6mm)

Mobile Phase: A. Water B. Acetonitrile

Gradient: Time (mins)	%B
0	25
24	46
26	46
27	25

Flow Rate: 1.0ml/min

Column Temperature: 20°C


Detection: UV, 214nm

Sample Concentration: 25 $\mu$ g/ml methanol

\* Conditions optimised using Drylab 2000 Plus method development software

### Separation of 12 Steroid Standards

App No. 1160



- |                          |                                   |
|--------------------------|-----------------------------------|
| 1. Estriol               | 7. Cortisone-21-acetate           |
| 2. Prednisolone          | 8. 17 $\alpha$ -Estradiol         |
| 3. Hydrocortisone        | 9. 19-Norethindrone               |
| 4. Cortisone             | 10. 17 $\alpha$ -Ethinylestradiol |
| 5. Corticosterone        | 11. 21-Hydroxyprogesterone        |
| 6. 17 $\beta$ -Estradiol | 12. Estrone                       |