




HPLC Chromatogram Index

COSMOSIL HILIC



Selection guide of mobile phase

COMOSIL HILIC column generates retention and separation by hydrophilic interaction (mainly hydrogen bond) and anion-exchange. Refer to following recommendations to select an appropriate mobile phase condition.

(1) The effect of organic solvent type and content

- In general, acetonitrile/water is used as mobile phase.
- Retention increases as water content in the mobile phase decreased. (Fig.1)
- Use acetonitrile content in the mobile phase within the range of 0-95% (in general 50-95%).
- Methanol/water generates shorter retention than acetonitrile/water. (Fig.2)
- Use only HPLC grade solvents.




Fig.1 The effect of acetonitrile content on retention

Column; COSMOSIL HILIC

Mobile phase; Acetonitrile/ 10mmol/l $\text{CH}_3\text{COONH}_4$




Fig.2 Difference of acetonitrile and methanol on retention

Column; COSMOSIL HILIC

Mobile phase; Organic solvent/ 10mmol/l $\text{CH}_3\text{COONH}_4$ = 90/10

(2) The effect of buffer pH

- Keep pH of the mobile phase within the range of 2-7.5.
- The buffer around neutrality generates larger retention. (Fig.3)




Fig.3 The effect of buffer pH on retention

Column; COSMOSIL HILIC

Mobile phase; Acetonitrile / 10mmol/l buffer = 90/10

(3) The effect of salt type and concentration

- When analyze ionic compounds, add salts or buffers in the mobile phase.
- When mobile phase has high acetonitrile content, note dissolubility of the salt. The dissolubility of phosphate buffers used widely in reversed phase chromatography is low in acetonitrile. Therefore use of phosphate buffers is not recommended. Keep the concentration of acetonitrile under 70% if use a phosphate buffer. Check that the salt does not precipitate when mixed with acetonitril before use.
- Ammonium acetate or formic acid ammonium buffers are recommended because they are soluble in high acetonitrile content.

- Use the buffer concentration within the range of 5 - 100mmol/l. Moreover, Check that the salt does not precipitate after mixing buffer and acetonitrile.
- High salt concentration inhibits ion exchange and decreases retention. (Fig.4)
- Run mobile phase through membrane filter (less than 0.45µm in pore size) before use.




Fig.4 The effect of salt concentration on retention

Column: COSMOSIL HILIC

Mobile phase; Acetonitrile / 10mmol/l CH₃COONH₄ = 50/50

(4) Selection of mobile phase

Following are the recommended mobile phases for different compound types.

Neutral compounds	→ Acetonitrile / Water = 90/10
Basic compounds	→ Acetonitrile / 10mmol/l CH ₃ COONH ₄ = 90/10
Amphoteric compounds	→ Acetonitrile / 10mmol/l CH ₃ COONH ₄ = 70/30
Acidic compounds	→ Acetonitrile / 10mmol/l CH ₃ COONH ₄ = 50/50
	↓ not eluted
	Acetonitrile / 10mmol/l Phosphate buffer (pH7.0)= 50/50

(5) Improvement of peak shape

Try following if peak shape is tailing. The peak shape might improve.

- Add 5mmol/l EDTA to mobile phase.
- Change to citrate buffer. (i. e. 10 mmol/l citrate buffer pH7.0)




Fig.5 Improvement of peak shape

Column: COSMOSIL HILIC(4.6mmI.D.-250mm)

Sample: Tryptophan(1ng), Flow rate: 1.0ml/min

Detection: UV254nm, Temperature: 30°C

(6) Others


- Use scrupulously degassed mobile phase. Air bubbles generate detection noise and accelerate column deterioration.
- We recommend keeping the chromatography conditions constant, since frequent changes of mobile phase shorten column life.

Index

Sample name	Cas No.	Page	Sample name	Cas No.	Page
A Acetamide	60-35-5	1	L Homoserine	672-15-1	14
Acetazolamide	59-66-5	1	Hydantoic Acid	462-60-2	14
Acetrizoic Acid	85-36-9	1	Hydantoin	461-72-3	14
Acrylic Acid	79-10-7	1	Hydroxylamine O sulfonic Acid	2950-43-8	14
L- Alanine	56-41-7	1	cis-4-Hydroxy-D-proline	2584-71-6	14
Alanine	107-95-9	1	L Hydroxyproline	51-35-4	15
Allantoic Acid	99-16-1	2	N Hydroxysuccinimide	6066-82-6	15
Allantoin	97-59-6	2	I Indigo carmine	860-22-0	15
p Aminobenzamidine	3858-83-1	2	D-Isoascorbic Acid	89-65-6	15
p Aminobenzoic Acid	150-13-0	2	Isocinchomeronic Acid	100-26-5	15
4 Amino n butyric Acid (GABA)	56-12-2	2	L-Isoleucine	73-32-5	15
6 Aminohexanoic Acid	60-32-2	2	Isonicotinic Acid	55-22-1	16
5 Aminolevulinic Acid	106-60-5	3	Isonicotinohydrazide	54-85-3	16
2-Aminopyridine	504-29-0	3	Isopropyl -D-1-thiogalactopyranoside (IPTG)	367-93-1	16
3-Aminopyridine	462-08-8	3	K Kojic Acid	501-30-4	16
5-Amino-1H-tetrazole	4418-61-5	3	L L-Leucine	61-90-5	16
3 Amino 1H 1,2,4 triazole	61-82-5	3	N-D-Leucyl-L-tyrosine	3303-29-5	16
5-Aminouracil	932-52-5	3	L Lysine	56-87-1	17
Amphotericin B	1397-89-3	4	M Maleic Acid	110-16-7	17
Angiotensin II(Human)	484-42-4	4	L (-) Malic Acid	97-67-6	17
Angiotensin II(Human)	4474-91-3	4	Malonic Acid	141-82-2	17
L Arginine	74-79-3	4	Mecobalamin	13422-55-4	17
L(+) Ascorbic Acid	50-81-7	4	Metanilic Acid	121-47-1	17
L-Asparagine	70-47-3	4	L Methionine	63-68-3	18
L Aspartic Acid	56-84-8	5	6-Methyl-2-thiouracil	56-04-2	18
6 Azauracil	461-89-2	5	N Methylglucamine	6284-40-8	18
Aztreonam	78110-38-0	5	N Methylhydroxylamine	593-77-1	18
B Benzamidine	618-39-3	5	Mucic Acid	526-99-8	18
Benzenesulfonic Acid	98-11-3	5	Murexide	3051-09-0	18
Benzoic Acid	65-85-0	5	N Nicotinamide	98-92-0	19
Bromoacetic Acid	79-08-3	6	Nicotinic Acid	59-67-6	19
C Cacotheline	561-20-6	6	L Noradrenaline	51-41-2	19
Camostat	59721-28-7	6	DL-Norleucine	616-06-8	19
L Carnitine	541-15-1	6	DL Norvaline	760-78-1	19
Ceftriaxone	73384-59-5	6	O L Ornithine	70-26-8	19
Chloroacetic Acid	79-11-8	6	Orotic Acid	65-86-1	20
Citrazinic Acid	99-11-6	7	Oxalic Acid	144-62-7	20
Creatine	57-00-1	7	Oxamic Acid	471-47-6	20
Creatinine	60-27-5	7	Oxytocin	50-56-6	20
Cyanoacetic Acid	372-09-8	7	P D Pantothenic Acid	79-83-4	20
Cyanuric Acid	108-80-5	7	L-(-)-Phenylalanine	63-91-2	20
L Cysteine	52-90-4	7	p Phenylenediamine	106-50-3	21
L (-) Cystine	56-89-3	8	L (+) Phenylglycine	2935-35-5	21
Cytidine	65-46-3	8	Phosphocreatine	67-07-2	21
Cytosine	71-30-7	8	O Phospho L serine	407-41-0	21
D 3,4-Diaminobenzoic Acid	619-05-6	8	Picolinic acid	98-98-6	21
3,5-Diaminobenzoic Acid	535-87-5	8	Pivalic Acid	75-98-9	21
2,4-Diaminophenol	95-86-3	8	Procaterol	72332-33-3	22
DL 2,6 Diaminopimelic Acid	583-93-7	9	L-Proline	147-85-3	22
DL 2,3 Diaminopropionic Acid	54897-59-5	9	Propionic Acid	79-09-4	22
Diatrizoic Acid	117-96-4	9	Pyruvic Acid	127-17-3	22
Dipicolinic acid	499-83-2	9	R Ribose 5 phosphate	4300-28-1	22
Dithiouracil	2001-93-6	9	S D Saccharic Acid	87-73-0	22
L DOPA	59-92-7	9	Sarcosine	107-97-1	23
Dopamine	51-61-6	10	Sebacic Acid	111-20-6	23
meso Erythritol	149-32-6	10	L-Serine	56-45-1	23
F Famotidin	76824-35-6	10	Sinigrin	3952-98-5	23
Folic Acid	59-30-3	10	Succinic Acid	110-15-6	23
Folinic Acid	58-05-9	10	Sulbactam	68373-14-8	23
Formamide	75-12-7	10	Sulfanilic acid	121-57-3	24
D Fructose 6 phosphate	643-13-0	11	T L (+) Tartaric Acid	87-69-4	24
Fuchsine, Acid (Rubin S)	3244-88-0	11	Taurine	107-35-7	24
Fumaric Acid	110-17-8	11	L-Theanine	3081-61-6	24
G Gluconic Acid	526-95-4	11	2 Thiobarbituric Acid	504-17-6	24
D Glucose 1 phosphate	59-56-3	11	2 Thiouracil	141-90-2	24
D Glucose 6 phosphate	56-73-5	11	L-Threonine	72-19-5	25
D-Glucuronic Acid	6556-12-3	12	Todralazine	14679-73-3	25
L Glutamic Acid	56-86-0	12	Trichloroacetic Acid	76-03-9	25
L Glutamine	56-85-9	12	Tris(hydroxymethyl)aminomethane	77-86-1	25
Glutaric Acid	110-94-1	12	L-Tryptophan	73-22-3	25
DL Glyceric Acid	600-19-1	12	L-Tyrosine	60-18-4	25
Glycinamide	598-41-4	12	U Uracil	66-22-8	26
Glycine	56-40-6	13	Urea	57-13-6	26
Glycolic Acid	79-14-1	13	Uridine	58-96-8	26
Glycylglycine	556-50-3	13	V L-Valine	72-18-4	26
Guanidoacetic Acid	352-97-6	13			
H 1,2,6-Hexanetriol	106-69-4	13			
L-Histidine	71-00-1	13			
L Homocystine	626-72-2	14			


COSMOSIL Chromatogram Index

Sample: Acetamide
 CAS No.: [60-35-5]
 Molecular formula: C₂H₅NO
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ H₂O=95/5
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV220 nm
 Attenuation: 0.128 aufs
 Sample conc.: 10.0mg/ml
 Injection volume: 0.5µl
 Retention time: 4.75min
 Capacity factor: 0.57




COSMOSIL Chromatogram Index

Sample: Acetazolamide
 CAS No.: [59-66-5]
 Molecular formula: C₄H₆N₄O₃S₂
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =90/10
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV254 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.2mg/ml
 Injection volume: 0.5µl
 Retention time: 5.99min
 Capacity factor: 1.05




COSMOSIL Chromatogram Index

Sample: Acetrizoic Acid
 CAS No.: [85-36-9]
 Molecular formula: C₇H₆I₂NO₃
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV254 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.8mg/ml
 Injection volume: 1.0µl
 Retention time: 16.39min
 Capacity factor: 4.76




COSMOSIL Chromatogram Index

Sample: Acrylic Acid
 CAS No.: [79-10-7]
 Molecular formula: C₃H₄O₂
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV220 nm
 Attenuation: 0.128 aufs
 Sample conc.: 1.0mg/ml
 Injection volume: 1.0µl
 Retention time: 15.05min
 Capacity factor: 4.28




COSMOSIL Chromatogram Index

Sample: L-α-Alanine
 CAS No.: [56-41-7]
 Molecular formula: C₃H₇NO₂
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV 210nm
 Attenuation: 0.128 aufs
 Sample conc.: 5.0mg/ml
 Injection volume: 2.0µl
 Retention time: 6.67min
 Capacity factor: 153




COSMOSIL Chromatogram Index

Sample: β-Alanine
 CAS No.: [107-95-9]
 Molecular formula: C₃H₇NO₂
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV 210nm
 Attenuation: 0.128 aufs
 Sample conc.: 10.0mg/ml
 Injection volume: 0.5µl
 Retention time: 7.38min
 Capacity factor: 1.81




COSMOSIL Chromatogram Index

Sample: Allantoic Acid
 CAS No.: [99-16-1]
 Molecular formula: $C_2H_4N_4O_4$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV220 nm
 Attenuation: 0.128 aufs
 Sample conc.: 5.0mg/ml
 Injection volume: 0.5µl
 Retention time: 4.45min
 Capacity factor: 0.69




COSMOSIL Chromatogram Index

Sample: Allantoin
 CAS No.: [97-59-6]
 Molecular formula: $C_4H_6N_4O_3$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =90/10
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV220 nm
 Attenuation: 0.128 aufs
 Sample conc.: 1.0mg/ml
 Injection volume: 1.0µl
 Retention time: 8.75min
 Capacity factor: 2.02




COSMOSIL Chromatogram Index

Sample: p-Aminobenzamidine
 CAS No.: [3858-83-1]
 Molecular formula: $C_7H_9N_3$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH7.0)=70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV254 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.5mg/ml
 Injection volume: 1.0µl
 Retention time: 6.07min
 Capacity factor: 1.31




COSMOSIL Chromatogram Index

Sample: p-Aminobenzoic Acid
 CAS No.: [150-13-0]
 Molecular formula: $C_7H_7NO_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV254 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.4mg/ml
 Injection volume: 1.0µl
 Retention time: 16.97min
 Capacity factor: 4.91




COSMOSIL Chromatogram Index

Sample: 4-Amino-n-butyric Acid
 CAS No.: [56-12-2]
 Molecular formula: $C_4H_9NO_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 10.0mg/ml
 Injection volume: 1.0µl
 Retention time: 7.67min
 Capacity factor: 1.92





COSMOSIL Chromatogram Index

Sample: 6-Aminohexanoic Acid
 CAS No.: [60-32-2]
 Molecular formula: $C_6H_{13}NO_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV 210nm
 Attenuation: 0.128 aufs
 Sample conc.: 10.0mg/ml
 Injection volume: 1.0µl
 Retention time: 7.98min
 Capacity factor: 2.03



COSMOSIL Chromatogram Index


Sample: 5-Aminolevulinic Acid
 CAS No.: [5451-09-2]
 Molecular formula: $C_5H_9NO_3$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV 210nm
 Attenuation: 0.128 auFS
 Sample conc.: 5.0mg/ml
 Injection volume: 1.0µl
 Retention time: 12.80min
 Capacity factor: 3.87



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index


Sample: 2-Aminopyridine
 CAS No.: [504-29-0]
 Molecular formula: $C_5H_6N_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =95/5
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV254 nm
 Attenuation: 0.128 auFS
 Sample conc.: 1.0mg/ml
 Injection volume: 0.5µl
 Retention time: 4.25min
 Capacity factor: 0.39



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index


Sample: 3-Aminopyridine
 CAS No.: [462-08-8]
 Molecular formula: $C_5H_6N_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =90/10
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV254 nm
 Attenuation: 0.128 auFS
 Sample conc.: 0.1mg/ml
 Injection volume: 1.0µl
 Retention time: 4.05min
 Capacity factor: 0.51



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index


Sample: 5-Amino-1H-tetrazole
 CAS No.: [4418-61-5]
 Molecular formula: CH_3N_5
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV220 nm
 Attenuation: 0.128 auFS
 Sample conc.: 0.5mg/ml
 Injection volume: 1.0µl
 Retention time: 15.76min
 Capacity factor: 4.49



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index


Sample: 3-Amino-1H-1,2,4-triazole
 CAS No.: [61-82-5]
 Molecular formula: $C_2H_4N_4$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =90/10
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV220 nm
 Attenuation: 0.128 auFS
 Sample conc.: 0.2mg/ml
 Injection volume: 1.0µl
 Retention time: 7.01min
 Capacity factor: 1.42



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index



Sample: 5-Aminouracil
 CAS No.: [932-52-5]
 Molecular formula: $C_4H_5N_3O_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =90/10
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV260 nm
 Attenuation: 0.128 auFS
 Sample conc.: 0.5mg/ml
 Injection volume: 0.5µl
 Retention time: 7.01min
 Capacity factor: 1.42



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index



Sample: Amphotericin B
 CAS No.: [1397-89-3]
 Molecular formula: $C_{47}H_{73}NO_{17}$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.25mg/ml
 Injection volume: 0.5µl
 Retention time: 5.34min
 Capacity factor: 0.99



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index



Sample: Angiotensin I(Human)
 CAS No.: [484-42-4]
 Molecular formula: $C_{62}H_{89}N_{17}O_{14}$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV220 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.4mg/ml
 Injection volume: 0.5µl
 Retention time: 5.28min
 Capacity factor: 0.84



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index



Sample: Angiotensin II(Human)
 CAS No.: [4474-91-3]
 Molecular formula: $C_{59}H_{71}N_{13}O_{12}$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV220 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.4mg/ml
 Injection volume: 0.5µl
 Retention time: 6.56min
 Capacity factor: 1.29



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index



Sample: L-Arginine
 CAS No.: [74-79-3]
 Molecular formula: $C_6H_{14}N_4O_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH7.0)=60/40
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 10.0mg/ml
 Injection volume: 1.0µl
 Retention time: 7.97min
 Capacity factor: 1.95



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index



Sample: L(+)-Ascorbic Acid
 CAS No.: [50-81-7]
 Molecular formula: $C_6H_8O_6$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV245nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.2mg/ml
 Injection volume: 3.0µl
 Retention time: 17.80min
 Capacity factor: 5.31



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index



Sample: L-Asparagine
 CAS No.: [70-47-3]
 Molecular formula: $C_4H_8N_2O_3$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH7.0)=60/40
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 5.0mg/ml
 Injection volume: 1.0µl
 Retention time: 4.88min
 Capacity factor: 0.80



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index


Sample: L-Aspartic Acid
 CAS No.: [56-84-8]
 Molecular formula: $C_4H_7NO_4$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV 210nm
 Attenuation: 0.128 aufs
 Sample conc.: 10.0mg/ml
 Injection volume: 2.0µl
 Retention time: 19.79min
 Capacity factor: 6.01



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index



Sample: 6-Azauracil
 CAS No.: [461-89-2]
 Molecular formula: $C_4H_4N_2O_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =95/5
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV260 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.3mg/ml
 Injection volume: 0.5µl
 Retention time: 9.65min
 Capacity factor: 2.19



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index

Sample: Aztreonam
 CAS No.: [78110-38-0]
 Molecular formula: $C_{13}H_{17}N_5O_8S_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 20mmol/l Phosphate buffer(pH7)=50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV280 nm
 Attenuation: 0.128 aufs
 Sample conc.: 2.5mg/ml
 Injection volume: 1.0µl
 Retention time: 17.57min
 Capacity factor: 5.18



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index


Sample: Benzamidine
 CAS No.: [618-39-3]
 Molecular formula: $C_7H_9N_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =90/10
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV254 nm
 Attenuation: 0.128 aufs
 Sample conc.: 1.5mg/ml
 Injection volume: 0.5µl
 Retention time: 7.16min
 Capacity factor: 1.46



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index


Sample: Benzenesulfonic Acid
 CAS No.: [98-11-3]
 Molecular formula: $C_6H_6O_3S$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV254 nm
 Attenuation: 0.128 aufs
 Sample conc.: 5.0mg/ml
 Injection volume: 1.0µl
 Retention time: 11.54min
 Capacity factor: 3.05



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index


Sample: Benzoic Acid
 CAS No.: [65-85-0]
 Molecular formula: $C_7H_6O_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV254 nm
 Attenuation: 0.128 aufs
 Sample conc.: 5.0mg/ml
 Injection volume: 0.5µl
 Retention time: 15.19min
 Capacity factor: 4.29



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index



Sample: Bromoacetic Acid
 CAS No.: [79-08-3]
 Molecular formula: $C_2H_3BrO_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 1.0mg/ml
 Injection volume: 1.0µl
 Retention time: 15.04min
 Capacity factor: 4.31



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index



Sample: Cacotheleine
 CAS No.: [561-20-6]
 Molecular formula: $C_{21}H_{21}N_3O_7$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV220 nm
 Attenuation: 0.128 aufs
 Sample conc.: 1.0mg/ml
 Injection volume: 1.0µl
 Retention time: 12.19min
 Capacity factor: 3.23



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index



Sample: Camostat
 CAS No.: [59721-28-7]
 Molecular formula: $C_{29}H_{22}N_4O_5$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =90/10
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV265 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.5mg/ml
 Injection volume: 0.5µl
 Retention time: 7.16min
 Capacity factor: 1.47



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index



Sample: L-Carnitine
 CAS No.: [541-15-1]
 Molecular formula: $C_7H_{15}NO_3$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: ELSD
 Attenuation: Gain=6, Atten=8
 Sample conc.: 2.0mg/ml
 Injection volume: 1.5µl
 Retention time: 6.96min
 Capacity factor: 1.78



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index


Sample: Ceftriaxone
 CAS No.: [73384-59-5]
 Molecular formula: $C_{18}H_{18}N_6O_7S_3$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV254 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.5mg/ml
 Injection volume: 1.0µl
 Retention time: 11.36min
 Capacity factor: 3.05



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index


Sample: Chloroacetic Acid
 CAS No.: [79-11-8]
 Molecular formula: $C_2H_2ClO_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 10.0mg/ml
 Injection volume: 1.0µl
 Retention time: 14.69min
 Capacity factor: 4.15



NACALAI TESQUE, INC


COSMOSIL Chromatogram Index

Sample: Citrazinic Acid
 CAS No.: [99-11-6]
 Molecular formula: $C_7H_5NO_4$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV220 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.5mg/ml
 Injection volume: 0.5µl
 Retention time: 14.16min
 Capacity factor: 3.98




COSMOSIL Chromatogram Index

Sample: Creatine
 CAS No.: [57-00-1]
 Molecular formula: $C_4H_9N_3O_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV220 nm
 Attenuation: 0.128 aufs
 Sample conc.: 1.0mg/ml
 Injection volume: 1.0µl
 Retention time: 6.35min
 Capacity factor: 1.40




COSMOSIL Chromatogram Index

Sample: Creatinine
 CAS No.: [60-27-5]
 Molecular formula: $C_4H_7N_3O$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =90/10
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV220 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.5mg/ml
 Injection volume: 0.5µl
 Retention time: 8.93min
 Capacity factor: 2.08




COSMOSIL Chromatogram Index

Sample: Cyanoacetic Acid
 CAS No.: [372-09-8]
 Molecular formula: $C_3H_3NO_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: ELSD
 Attenuation: Gain=6, Atten=8
 Sample conc.: 10.0mg/ml
 Injection volume: 0.5µl
 Retention time: 11.78min
 Capacity factor: 3.56




COSMOSIL Chromatogram Index

Sample: Cyanuric Acid
 CAS No.: [108-80-5]
 Molecular formula: $C_3H_3N_3O_3$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV220 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.5mg/ml
 Injection volume: 1.0µl
 Retention time: 7.61min
 Capacity factor: 1.68




COSMOSIL Chromatogram Index

Sample: L-Cysteine
 CAS No.: [52-90-4]
 Molecular formula: $C_3H_7NO_2S$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH7.0)=60/40
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 2.0mg/ml
 Injection volume: 1.0µl
 Retention time: 4.69min
 Capacity factor: 0.73




COSMOSIL Chromatogram Index

Sample: L-(-)-Cystine
 CAS No.: [56-89-3]
 Molecular formula: $C_4H_8N_2O_4S_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH7.0)=60/40
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 5.0mg/ml
 Injection volume: 0.5µl
 Retention time: 6.42min
 Capacity factor: 1.38




COSMOSIL Chromatogram Index

Sample: Cytidine
 CAS No.: [65-46-3]
 Molecular formula: $C_9H_{13}N_3O_5$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =80/20
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV260 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.5mg/ml
 Injection volume: 0.5µl
 Retention time: 7.22min
 Capacity factor: 1.58




COSMOSIL Chromatogram Index

Sample: Cytosine
 CAS No.: [71-30-7]
 Molecular formula: $C_4H_5N_3O$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =90/10
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV260 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.5mg/ml
 Injection volume: 0.5µl
 Retention time: 11.22min
 Capacity factor: 2.87




COSMOSIL Chromatogram Index

Sample: 3,4-Diaminobenzoic Acid
 CAS No.: [619-05-6]
 Molecular formula: $C_7H_7N_2O_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV254 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.10mg/ml
 Injection volume: 4.0µl
 Retention time: 16.13min
 Capacity factor: 4.62




COSMOSIL Chromatogram Index

Sample: 3,5-Diaminobenzoic Acid
 CAS No.: [535-87-5]
 Molecular formula: $C_7H_7N_2O_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV254 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.1mg/ml
 Injection volume: 4.0µl
 Retention time: 16.54min
 Capacity factor: 4.76




COSMOSIL Chromatogram Index

Sample: 2,4-Diaminophenol
 CAS No.: [95-86-3]
 Molecular formula: $C_6H_7N_2O$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 20mmol/l Ammonium acetate =90/10
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV254 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.5mg/ml
 Injection volume: 0.5µl
 Retention time: 4.40min
 Capacity factor: 0.51



COSMOSIL Chromatogram Index


Sample: DL-2,6-Diaminopimelic Acid
 CAS No.: [583-93-7]
 Molecular formula: $C_7H_{14}N_2O_4$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH7.0)=60/40
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 auFS
 Sample conc.: 10.0mg/ml
 Injection volume: 1.5µl
 Retention time: 6.93min
 Capacity factor: 1.56



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index


Sample: DL-2,3-Diaminopropionic Acid
 CAS No.: [54897-59-5]
 Molecular formula: $C_3H_6N_2O_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: ELSD
 Attenuation: Gain=6, Atten=8
 Sample conc.: 5.0mg/ml
 Injection volume: 2.0µl
 Retention time: 14.38min
 Capacity factor: 4.52



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index


Sample: Diatrizoic Acid
 CAS No.: [117-96-4]
 Molecular formula: $C_{11}H_9I_3N_2O_4$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV254 nm
 Attenuation: 0.128 auFS
 Sample conc.: 0.8mg/ml
 Injection volume: 1.0µl
 Retention time: 14.98min
 Capacity factor: 4.26



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index


Sample: Dipicolinic acid
 CAS No.: [499-83-2]
 Molecular formula: $C_7H_7NO_4$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH7.0)=50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV254 nm
 Attenuation: 0.128 auFS
 Sample conc.: 0.5mg/ml
 Injection volume: 1.0µl
 Retention time: 6.37min
 Capacity factor: 1.23



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index


Sample: Dithiouracil
 CAS No.: [2001-93-6]
 Molecular formula: $C_4H_4N_2S_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =95/5
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV260 nm
 Attenuation: 0.128 auFS
 Sample conc.: 0.2mg/ml
 Injection volume: 1.5µl
 Retention time: 15.60min
 Capacity factor: 4.15



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index


Sample: L-DOPA
 CAS No.: [59-92-7]
 Molecular formula: $C_9H_{11}NO_4$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH7.0)=70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV254 nm
 Attenuation: 0.128 auFS
 Sample conc.: 3.0mg/ml
 Injection volume: 3.0µl
 Retention time: 7.12min
 Capacity factor: 1.72



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index


Sample: Dopamine
 CAS No.: [51-61-6]
 Molecular formula: $C_8H_{11}NO_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH7.0)=70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV254 nm
 Attenuation: 0.128 auFS
 Sample conc.: 1.0mg/ml
 Injection volume: 4.0µl
 Retention time: 7.73min
 Capacity factor: 1.96



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index


Sample: meso-Erythritol [meso-Erythrite]
 CAS No.: [149-32-6]
 Molecular formula: $C_4H_{10}O_4$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ H_2O =90/10
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: ELSD
 Attenuation: Gain=6, Atten=8
 Sample conc.: 1.0mg/ml
 Injection volume: 1.0µl
 Retention time: 5.78min
 Capacity factor: 1.18



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index


Sample: Famotidin
 CAS No.: [76824-35-6]
 Molecular formula: $C_{16}H_{15}N_7O_2S_3$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =90/10
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV254 nm
 Attenuation: 0.128 auFS
 Sample conc.: 0.25mg/ml
 Injection volume: 2.0µl
 Retention time: 5.99min
 Capacity factor: 1.06



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index


Sample: Folic Acid
 CAS No.: [59-30-3]
 Molecular formula: $C_{19}H_{19}N_7O_6$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 20mmol/l Phosphate buffer(pH7)=50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV254 nm
 Attenuation: 0.128 auFS
 Sample conc.: 0.25mg/ml
 Injection volume: 2.0µl
 Retention time: 14.09min
 Capacity factor: 3.95



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index


Sample: Folinic Acid
 CAS No.: [58-05-9]
 Molecular formula: $C_{20}H_{23}N_7O_7$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV254 nm
 Attenuation: 0.128 auFS
 Sample conc.: 0.25mg/ml
 Injection volume: 2.0µl
 Retention time: 13.36min
 Capacity factor: 3.68



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index



Sample: Formamide
 CAS No.: [75-12-7]
 Molecular formula: CH_3NO
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ H_2O =95/5
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV220 nm
 Attenuation: 0.128 auFS
 Sample conc.: 10.0mg/ml
 Injection volume: 0.5µl
 Retention time: 4.58min
 Capacity factor: 0.52



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index



Sample: D-Fructose-6-phosphate
 CAS No.: [643-13-0]
 Molecular formula: $C_6H_{13}O_9P$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: RI
 Attenuation: 4×10^3 RIU/FS
 Sample conc.: 10.0mg/ml
 Injection volume: 5.0µl
 Retention time: 12.16min
 Capacity factor: 3.64



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index



Sample: Fuch sine, Acid
 CAS No.: [3244-88-0]
 Molecular formula: $C_{20}H_{17}N_3Na_2O_9S_3$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV254 nm
 Attenuation: 0.128 aufs
 Sample conc.: 5.0mg/ml
 Injection volume: 1.5µl
 Retention time: 13.82min
 Capacity factor: 3.85



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index



Sample: Fumaric Acid
 CAS No.: [110-17-8]
 Molecular formula: $C_4H_4O_4$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.2mg/ml
 Injection volume: 0.5µl
 Retention time: 10.63min
 Capacity factor: 2.75



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index



Sample: Gluconic Acid
 CAS No.: [526-95-4]
 Molecular formula: $C_6H_{12}O_7$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 10.0mg/ml
 Injection volume: 1.0µl
 Retention time: 5.15min
 Capacity factor: 0.81



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index



Sample: α-D-Glucose-1-phosphate
 CAS No.: [59-56-3]
 Molecular formula: $C_6H_{13}O_9P$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: RI
 Attenuation: 4×10^3 RIU/FS
 Sample conc.: 10.0mg/ml
 Injection volume: 5.0µl
 Retention time: 12.26min
 Capacity factor: 3.68



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index


Sample: D-Glucose-6-phosphate
 CAS No.: [56-73-5]
 Molecular formula: $C_6H_{13}O_9P$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: RI
 Attenuation: 4×10^3 RIU/FS
 Sample conc.: 10.0mg/ml
 Injection volume: 5.0µl
 Retention time: 12.95min
 Capacity factor: 3.94



NACALAI TESQUE, INC


COSMOSIL Chromatogram Index

Sample: D-Glucuronic Acid
 CAS No.: [6556-12-3]
 Molecular formula: $C_6H_{10}O_7$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 10.0mg/ml
 Injection volume: 1.0µl
 Retention time: 5.45min
 Capacity factor: 0.92




COSMOSIL Chromatogram Index

Sample: L-Glutamic Acid
 CAS No.: [56-86-0]
 Molecular formula: $C_5H_9NO_4$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV 210nm
 Attenuation: 0.128 aufs
 Sample conc.: 10.0mg/ml
 Injection volume: 2.0µl
 Retention time: 19.38min
 Capacity factor: 5.87




COSMOSIL Chromatogram Index

Sample: L-Glutamine
 CAS No.: [56-85-9]
 Molecular formula: $C_5H_{10}N_2O_3$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 10.0mg/ml
 Injection volume: 0.5µl
 Retention time: 7.50min
 Capacity factor: 1.85




COSMOSIL Chromatogram Index

Sample: Glutaric Acid
 CAS No.: [110-94-1]
 Molecular formula: $C_5H_8O_4$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 10.0mg/ml
 Injection volume: 1.0µl
 Retention time: 10.45min
 Capacity factor: 2.68




COSMOSIL Chromatogram Index

Sample: DL-Glyceric Acid
 CAS No.: [600-19-1]
 Molecular formula: $C_3H_4O_4$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 6.0mg/ml
 Injection volume: 5.0µl
 Retention time: 15.68min
 Capacity factor: 4.50




COSMOSIL Chromatogram Index

Sample: Glycinamide
 CAS No.: [598-41-4]
 Molecular formula: $C_2H_5N_2O$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: ELSD
 Attenuation: Gain=6, Atten=8
 Sample conc.: 1.0mg/ml
 Injection volume: 3.0µl
 Retention time: 13.35min
 Capacity factor: 3.64




COSMOSIL Chromatogram Index

Sample: Glycine
 CAS No.: [56-40-6]
 Molecular formula: C₂H₅NO₂
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 5.0mg/ml
 Injection volume: 2.0µl
 Retention time: 7.29min
 Capacity factor: 1.77




COSMOSIL Chromatogram Index

Sample: Glycolic Acid
 CAS No.: [79-14-1]
 Molecular formula: C₂H₃O₃
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 10.0mg/ml
 Injection volume: 2.0µl
 Retention time: 15.28min
 Capacity factor: 4.39




COSMOSIL Chromatogram Index

Sample: Glycylglycine
 CAS No.: [556-50-3]
 Molecular formula: C₄H₈N₂O₃
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 1.0mg/ml
 Injection volume: 0.5µl
 Retention time: 6.40min
 Capacity factor: 1.27




COSMOSIL Chromatogram Index

Sample: Guanidoacetic Acid
 CAS No.: [352-97-6]
 Molecular formula: C₃H₅N₃O₂
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.5mg/ml
 Injection volume: 1.0µl
 Retention time: 6.61min
 Capacity factor: 1.51




COSMOSIL Chromatogram Index

Sample: 1,2,6-Hexanetriol
 CAS No.: [106-69-4]
 Molecular formula: C₆H₁₄O₃
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ H₂O=90/10
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: ELSD
 Attenuation: Gain=6, Atten=8
 Sample conc.: 1.0mg/ml
 Injection volume: 2.0µl
 Retention time: 5.19min
 Capacity factor: 0.80





COSMOSIL Chromatogram Index

Sample: L-Histidine
 CAS No.: [71-00-1]
 Molecular formula: C₆H₉N₃O₂
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH7.0)=70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV220 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.50mg/ml
 Injection volume: 1.0µl
 Retention time: 8.38min
 Capacity factor: 2.19



COSMOSIL Chromatogram Index



Sample: L-Homocystine
 CAS No.: [626-72-2]
 Molecular formula: $C_4H_{10}N_2O_4S_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l
 Ammonium acetate=70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV 210nm
 Attenuation: 0.128 aufs
 Sample conc.: 2.0mg/ml
 Injection volume: 1.0µl
 Retention time: 13.41min
 Capacity factor: 4.10



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index



Sample: L-Homoserine
 CAS No.: [672-15-1]
 Molecular formula: $C_4H_9NO_3$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l
 Ammonium acetate=70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV 210nm
 Attenuation: 0.128 aufs
 Sample conc.: 10.0mg/ml
 Injection volume: 1.0µl
 Retention time: 7.03min
 Capacity factor: 1.67



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index


Sample: Hydantoinic Acid
 CAS No.: [462-60-2]
 Molecular formula: $C_3H_6N_2O_3$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l
 Ammonium acetate =50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 5.0mg/ml
 Injection volume: 1.0µl
 Retention time: 16.33min
 Capacity factor: 4.72



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index


Sample: Hydantoin
 CAS No.: [461-72-3]
 Molecular formula: $C_3H_4N_2O_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l
 Ammonium acetate =95/5
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV220 nm
 Attenuation: 0.128 aufs
 Sample conc.: 1.0mg/ml
 Injection volume: 0.5µl
 Retention time: 5.01min
 Capacity factor: 0.66



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index


Sample: Hydroxylamine-O-sulfonic Acid
 CAS No.: [2950-43-8]
 Molecular formula: H_3NO_3S
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l
 Ammonium acetate =70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: ELSD
 Attenuation: Gain=6, Atten=8
 Sample conc.: 2.0mg/ml
 Injection volume: 3.0µl
 Retention time: 15.60min
 Capacity factor: 5.24



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index


Sample: cis-4-Hydroxy-D-proline
 CAS No.: [2584-71-6]
 Molecular formula: $C_5H_9NO_3$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l
 Ammonium acetate=70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 10.0mg/ml
 Injection volume: 1.0µl
 Retention time: 6.96min
 Capacity factor: 1.65



NACALAI TESQUE, INC


COSMOSIL Chromatogram Index

Sample: L-Hydroxyproline
 CAS No.: [51-35-4]
 Molecular formula: $C_5H_9NO_3$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV 210nm
 Attenuation: 0.128 aufs
 Sample conc.: 10.0mg/ml
 Injection volume: 1.0µl
 Retention time: 6.49min
 Capacity factor: 1.47




COSMOSIL Chromatogram Index

Sample: N-Hydroxysuccinimide
 CAS No.: [6066-82-6]
 Molecular formula: $C_4H_5NO_3$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV220 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.1mg/ml
 Injection volume: 1.5µl
 Retention time: 11.29min
 Capacity factor: 3.22




COSMOSIL Chromatogram Index

Sample: Indigo carmine
 CAS No.: [860-22-0]
 Molecular formula: $C_{16}H_8N_2Na_2O_6S_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV254 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.2mg/ml
 Injection volume: 1.0µl
 Retention time: 7.82min
 Capacity factor: 1.79




COSMOSIL Chromatogram Index

Sample: D-Isoascorbic Acid
 CAS No.: [89-65-6]
 Molecular formula: $C_6H_8O_6$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV 245nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.2mg/ml
 Injection volume: 3.0µl
 Retention time: 17.26min
 Capacity factor: 5.11




COSMOSIL Chromatogram Index

Sample: Isocinchomeronic Acid
 CAS No.: [100-26-5]
 Molecular formula: $C_7H_7NO_4$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH7.0)=50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV254 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.5mg/ml
 Injection volume: 0.5µl
 Retention time: 5.91min
 Capacity factor: 1.07




COSMOSIL Chromatogram Index

Sample: L-Isoleucine
 CAS No.: [73-32-5]
 Molecular formula: $C_6H_{13}NO_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 10.0mg/ml
 Injection volume: 1.0µl
 Retention time: 5.12min
 Capacity factor: 0.95



COSMOSIL Chromatogram Index


Sample: Isonicotinic Acid
 CAS No.: [55-22-1]
 Molecular formula: $C_6H_5NO_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV254 nm
 Attenuation: 0.128 aufs
 Sample conc.: 1.0mg/ml
 Injection volume: 0.5µl
 Retention time: 16.45min
 Capacity factor: 4.78



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index



Sample: Isonicotinohydrazide
 CAS No.: [54-85-3]
 Molecular formula: $C_6H_5N_3O$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =90/10
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV265 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.50mg/ml
 Injection volume: 0.5µl
 Retention time: 5.37min
 Capacity factor: 0.85



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index



Sample: Isopropyl β-D-1-thiogalactopyranoside
 CAS No.: [367-93-1]
 Molecular formula: $C_{18}H_{27}O_6S$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ H₂O=90/10
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: ELSD
 Attenuation: Gain=6, Atten=8
 Sample conc.: 0.1mg/ml
 Injection volume: 0.5µl
 Retention time: 6.23min
 Capacity factor: 1.15



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index



Sample: Kojic Acid
 CAS No.: [501-30-4]
 Molecular formula: $C_6H_6O_4$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH7.0)=70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV245 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.1mg/ml
 Injection volume: 1.0µl
 Retention time: 3.83min
 Capacity factor: 0.46



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index



Sample: L-Leucine
 CAS No.: [61-90-5]
 Molecular formula: $C_6H_{13}NO_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 10.0mg/ml
 Injection volume: 1.0µl
 Retention time: 4.91min
 Capacity factor: 0.87



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index


Sample: D-Leucyl-L-tyrosine
 CAS No.: [3303-29-5]
 Molecular formula: $C_{13}H_{23}N_2O_4$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV 254nm
 Attenuation: 0.128 aufs
 Sample conc.: 10.0mg/ml
 Injection volume: 1.0µl
 Retention time: 7.79min
 Capacity factor: 1.96



NACALAI TESQUE, INC


COSMOSIL Chromatogram Index

Sample: L-Lysine
 CAS No.: [56-87-1]
 Molecular formula: $C_6H_{14}N_2O_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l
 Citrate buffer(pH7.0)=60/40
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 10.0mg/ml
 Injection volume: 2.0µl
 Retention time: 8.46min
 Capacity factor: 2.13




COSMOSIL Chromatogram Index

Sample: Maleic Acid
 CAS No.: [110-16-7]
 Molecular formula: $C_4H_4O_4$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l
 Phosphate buffer(pH7)=50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.1mg/ml
 Injection volume: 0.5µl
 Retention time: 7.45min
 Capacity factor: 1.62




COSMOSIL Chromatogram Index

Sample: L-(-)-Malic Acid
 CAS No.: [97-67-6]
 Molecular formula: $C_4H_6O_5$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l
 Phosphate buffer(pH7)=50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 10.0mg/ml
 Injection volume: 0.5µl
 Retention time: 10.55min
 Capacity factor: 2.71




COSMOSIL Chromatogram Index

Sample: Malonic Acid
 CAS No.: [141-82-2]
 Molecular formula: $C_3H_4O_4$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l
 Phosphate buffer(pH7)=50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 10.0mg/ml
 Injection volume: 0.5µl
 Retention time: 10.78min
 Capacity factor: 2.81




COSMOSIL Chromatogram Index

Sample: Mecobalamin
 CAS No.: [13422-55-4]
 Molecular formula: $C_{63}H_{91}CoN_{13}O_{14}P$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l
 Ammonium acetate =70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV266 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.5mg/ml
 Injection volume: 1.0µl
 Retention time: 6.22min
 Capacity factor: 1.35




COSMOSIL Chromatogram Index

Sample: Metanilic Acid
 CAS No.: [121-47-1]
 Molecular formula: $C_6H_7NO_3S$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l
 Ammonium acetate =50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV254 nm
 Attenuation: 0.128 aufs
 Sample conc.: 1.0mg/ml
 Injection volume: 1.0µl
 Retention time: 13.68min
 Capacity factor: 3.80




COSMOSIL Chromatogram Index

Sample: L-Methionine
 CAS No.: [63-68-3]
 Molecular formula: $C_5H_{11}NO_2S$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH7.0)=60/40
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 1.0mg/ml
 Injection volume: 0.5µl
 Retention time: 4.15min
 Capacity factor: 0.54





COSMOSIL Chromatogram Index

Sample: 6-Methyl-2-thiouracil
 CAS No.: [56-04-2]
 Molecular formula: $C_5H_6N_2OS$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =95/5
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV260 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.1mg/ml
 Injection volume: 0.5µl
 Retention time: 5.58min
 Capacity factor: 0.84




COSMOSIL Chromatogram Index

Sample: N-Methylglucamine
 CAS No.: [6284-40-8]
 Molecular formula: $C_7H_{17}NO_5$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH7.0)=70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: RI
 Attenuation: 4×10^3 RIU/FS
 Sample conc.: 10.0mg/ml
 Injection volume: 2.0µl
 Retention time: 15.52min
 Capacity factor: 4.22





COSMOSIL Chromatogram Index

Sample: N-Methylhydroxylamine
 CAS No.: [593-77-1]
 Molecular formula: CH_3NO_2
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: ELSD
 Attenuation: Gain=6, Atten=8
 Sample conc.: 1.0mg/ml
 Injection volume: 2.0µl
 Retention time: 13.45min
 Capacity factor: 4.21





COSMOSIL Chromatogram Index

Sample: Mucic Acid
 CAS No.: [526-99-8]
 Molecular formula: $C_8H_{10}O_8$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 10.0mg/ml
 Injection volume: 1.0µl
 Retention time: 10.27min
 Capacity factor: 2.62




COSMOSIL Chromatogram Index

Sample: Murexide
 CAS No.: [3051-09-0]
 Molecular formula: $C_8H_8N_4O_6$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV220 nm
 Attenuation: 0.128 aufs
 Sample conc.: 1.0mg/ml
 Injection volume: 0.5µl
 Retention time: 13.47min
 Capacity factor: 3.69




COSMOSIL Chromatogram Index

Sample: Nicotinamide
 CAS No.: [98-92-0]
 Molecular formula: $C_6H_7N_2O$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =95/5
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV254 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.2mg/ml
 Injection volume: 1.0µl
 Retention time: 5.40min
 Capacity factor: 0.77




COSMOSIL Chromatogram Index

Sample: Nicotinic Acid
 CAS No.: [59-67-6]
 Molecular formula: $C_6H_5NO_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV254 nm
 Attenuation: 0.128 aufs
 Sample conc.: 1.0mg/ml
 Injection volume: 0.5µl
 Retention time: 16.67min
 Capacity factor: 4.87




COSMOSIL Chromatogram Index

Sample: L-Noradrenaline
 CAS No.: [51-41-2]
 Molecular formula: $C_8H_{11}NO_3$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH7.0)=60/40
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV254 nm
 Attenuation: 0.128 aufs
 Sample conc.: 5.0mg/ml
 Injection volume: 1.0µl
 Retention time: 5.47min
 Capacity factor: 1.07




COSMOSIL Chromatogram Index

Sample: DL-Norleucine
 CAS No.: [616-06-8]
 Molecular formula: $C_6H_{13}NO_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 5.0mg/ml
 Injection volume: 1.0µl
 Retention time: 4.89min
 Capacity factor: 0.86




COSMOSIL Chromatogram Index

Sample: DL-Norvaline
 CAS No.: [760-78-1]
 Molecular formula: $C_5H_{11}NO_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV 210nm
 Attenuation: 0.128 aufs
 Sample conc.: 10.0mg/ml
 Injection volume: 0.5µl
 Retention time: 5.43min
 Capacity factor: 1.07





COSMOSIL Chromatogram Index

Sample: L-Ornithine
 CAS No.: [70-26-8]
 Molecular formula: $C_5H_{12}N_2O_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH7.0)=60/40
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 10.0mg/ml
 Injection volume: 2.0µl
 Retention time: 8.39min
 Capacity factor: 2.10




COSMOSIL Chromatogram Index

Sample: Orotic Acid
 CAS No.: [65-86-1]
 Molecular formula: $C_4H_4N_2O_4$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV220 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.5mg/ml
 Injection volume: 1.0µl
 Retention time: 15.24min
 Capacity factor: 4.36




COSMOSIL Chromatogram Index

Sample: Oxalic Acid
 CAS No.: [144-62-7]
 Molecular formula: $C_2H_2O_4$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 5.0mg/ml
 Injection volume: 0.5µl
 Retention time: 12.08min
 Capacity factor: 3.27





COSMOSIL Chromatogram Index

Sample: Oxamic Acid
 CAS No.: [471-47-6]
 Molecular formula: $C_2H_3NO_3$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.1mg/ml
 Injection volume: 1.0µl
 Retention time: 4.83min
 Capacity factor: 0.71





COSMOSIL Chromatogram Index

Sample: Oxytocin
 CAS No.: [50-56-6]
 Molecular formula: $C_{43}H_{66}N_{12}O_{12}S_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV220 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.4mg/ml
 Injection volume: 0.5µl
 Retention time: 3.71min
 Capacity factor: 0.39





COSMOSIL Chromatogram Index

Sample: D-Pantothenic Acid
 CAS No.: [79-83-4]
 Molecular formula: $C_9H_{17}NO_5$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV220 nm
 Attenuation: 0.128 aufs
 Sample conc.: 10.0mg/ml
 Injection volume: 1.0µl
 Retention time: 13.21min
 Capacity factor: 3.60




COSMOSIL Chromatogram Index

Sample: L-(-)-Phenylalanine
 CAS No.: [63-91-2]
 Molecular formula: $C_9H_9NO_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV 254nm
 Attenuation: 0.128 aufs
 Sample conc.: 10.0mg/ml
 Injection volume: 0.5µl
 Retention time: 5.10min
 Capacity factor: 0.94




COSMOSIL Chromatogram Index

Sample: p-Phenylenediamine
 CAS No.: [106-50-3]
 Molecular formula: C₆H₈N₂
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =95/5
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV254 nm
 Attenuation: 0.128 auFS
 Sample conc.: 0.1mg/ml
 Injection volume: 0.5µl
 Retention time: 4.15min
 Capacity factor: 0.36




COSMOSIL Chromatogram Index

Sample: L-(+)-α-Phenylglycine
 CAS No.: [2935-35-5]
 Molecular formula: C₉H₉NO₂
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV254 nm
 Attenuation: 0.128 auFS
 Sample conc.: 5.0mg/ml
 Injection volume: 1.0µl
 Retention time: 5.96min
 Capacity factor: 1.27




COSMOSIL Chromatogram Index

Sample: Phosphocreatine
 CAS No.: [67-07-2]
 Molecular formula: C₄H₁₀N₃O₅P
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV220 nm
 Attenuation: 0.128 auFS
 Sample conc.: 1.0mg/ml
 Injection volume: 1.0µl
 Retention time: 11.42min
 Capacity factor: 3.00




COSMOSIL Chromatogram Index

Sample: O-Phospho-L-serine
 CAS No.: [407-41-0]
 Molecular formula: C₃H₉NO₆P
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 20mmol/l Phosphate buffer(pH7)=50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 auFS
 Sample conc.: 10.0mg/ml
 Injection volume: 3.0µl
 Retention time: 9.19min
 Capacity factor: 2.24




COSMOSIL Chromatogram Index

Sample: Picolinic acid
 CAS No.: [98-98-6]
 Molecular formula: C₆H₅NO₂
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH7.0)=70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV254 nm
 Attenuation: 0.128 auFS
 Sample conc.: 0.5mg/ml
 Injection volume: 0.5µl
 Retention time: 5.03min
 Capacity factor: 0.92





COSMOSIL Chromatogram Index

Sample: Pivalic Acid
 CAS No.: [75-98-9]
 Molecular formula: C₅H₁₀O₂
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 auFS
 Sample conc.: 10.0mg/ml
 Injection volume: 1.0µl
 Retention time: 12.14min
 Capacity factor: 3.28




COSMOSIL Chromatogram Index

Sample: Procaterol
 CAS No.: [72332-33-3]
 Molecular formula: $C_{16}H_{22}N_2O_3$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 20mmol/l Ammonium acetate =85/15
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV254 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.5mg/ml
 Injection volume: 0.5µl
 Retention time: 6.17min
 Capacity factor: 1.25




COSMOSIL Chromatogram Index

Sample: L-Proline
 CAS No.: [147-85-3]
 Molecular formula: $C_5H_9NO_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 10.0mg/ml
 Injection volume: 1.0µl
 Retention time: 5.83min
 Capacity factor: 1.22




COSMOSIL Chromatogram Index

Sample: Propionic Acid
 CAS No.: [79-09-4]
 Molecular formula: $C_3H_6O_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 10.0mg/ml
 Injection volume: 2.0µl
 Retention time: 14.85min
 Capacity factor: 4.24





COSMOSIL Chromatogram Index

Sample: Pyruvic Acid
 CAS No.: [127-17-3]
 Molecular formula: $C_3H_4O_3$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 20mmol/l Phosphate buffer(pH7)=50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV220 nm
 Attenuation: 0.128 aufs
 Sample conc.: 1.0mg/ml
 Injection volume: 2.0µl
 Retention time: 3.97min
 Capacity factor: 0.39





COSMOSIL Chromatogram Index

Sample: Ribose-5-phosphate
 CAS No.: [4300-28-1]
 Molecular formula: $C_5H_{11}O_8P$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 20mmol/l Phosphate buffer(pH7)=50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: RI
 Attenuation: 4×10^3 RIU/FS
 Sample conc.: 10.0mg/ml
 Injection volume: 5.0µl
 Retention time: 8.02min
 Capacity factor: 2.06





COSMOSIL Chromatogram Index

Sample: D-Saccharic Acid
 CAS No.: [87-73-0]
 Molecular formula: $C_6H_{10}O_8$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 10.0mg/ml
 Injection volume: 2.0µl
 Retention time: 10.48min
 Capacity factor: 2.69



COSMOSIL Chromatogram Index



Sample: Sarcosine
 CAS No.: [107-97-1]
 Molecular formula: $C_2H_5NO_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV 210nm
 Attenuation: 0.128 aufs
 Sample conc.: 10.0mg/ml
 Injection volume: 1.0µl
 Retention time: 6.30min
 Capacity factor: 1.40



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index



Sample: Sebacic Acid
 CAS No.: [111-20-6]
 Molecular formula: $C_{10}H_{18}O_4$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 5.0mg/ml
 Injection volume: 1.5µl
 Retention time: 9.43min
 Capacity factor: 2.28



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index



Sample: L-Serine
 CAS No.: [56-45-1]
 Molecular formula: $C_3H_7NO_3$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 10.0mg/ml
 Injection volume: 2.0µl
 Retention time: 7.92min
 Capacity factor: 2.01



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index



Sample: Sinigrin
 CAS No.: [3952-98-5]
 Molecular formula: $C_{10}H_{16}KNO_9S_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 1.0mg/ml
 Injection volume: 1.0µl
 Retention time: 12.57min
 Capacity factor: 3.38



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index



Sample: Succinic Acid
 CAS No.: [110-15-6]
 Molecular formula: $C_4H_6O_4$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 10.0mg/ml
 Injection volume: 0.5µl
 Retention time: 10.64min
 Capacity factor: 2.74



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index



Sample: Sulbactam
 CAS No.: [68373-14-8]
 Molecular formula: $C_8H_{11}NO_5S$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV220 nm
 Attenuation: 0.128 aufs
 Sample conc.: 5.0mg/ml
 Injection volume: 0.5µl
 Retention time: 10.86min
 Capacity factor: 2.81



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index



Sample: Sulfanilic acid
 CAS No.: [121-57-3]
 Molecular formula: $C_6H_7NO_3S$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV254 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.10mg/ml
 Injection volume: 1.0µl
 Retention time: 13.87min
 Capacity factor: 3.87



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index



Sample: L-(+)- Tartaric Acid
 CAS No.: [87-69-4]
 Molecular formula: $C_4H_6O_6$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 10.0mg/ml
 Injection volume: 1.5µl
 Retention time: 10.52min
 Capacity factor: 2.70



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index



Sample: Taurine
 CAS No.: [107-35-7]
 Molecular formula: $C_2H_7NO_3S$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: ELSD
 Attenuation: Gain=6, Atten=8
 Sample conc.: 1.0mg/ml
 Injection volume: 1.0µl
 Retention time: 6.40min
 Capacity factor: 1.25



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index


Sample: L-Theanine
 CAS No.: [3081-61-6]
 Molecular formula: $C_7H_{14}N_2O_3$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV220 nm
 Attenuation: 0.128 aufs
 Sample conc.: 5.0mg/ml
 Injection volume: 0.5µl
 Retention time: 5.89min
 Capacity factor: 1.21



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index


Sample: 2-Thiobarbituric Acid
 CAS No.: [504-17-6]
 Molecular formula: $C_4H_4N_2O_2S$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7)=50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV254 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.1mg/ml
 Injection volume: 0.5µl
 Retention time: 5.18min
 Capacity factor: 0.82



NACALAI TESQUE, INC

COSMOSIL Chromatogram Index



Sample: 2-Thiouracil
 CAS No.: [141-90-2]
 Molecular formula: $C_4H_4N_2OS$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =95/5
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV260 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.1mg/ml
 Injection volume: 0.5µl
 Retention time: 6.38min
 Capacity factor: 1.11



NACALAI TESQUE, INC



COSMOSIL Chromatogram Index

Sample: L-Threonine
 CAS No.: [72-19-5]
 Molecular formula: $C_4H_9NO_3$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 10.0mg/ml
 Injection volume: 2.0µl
 Retention time: 7.19min
 Capacity factor: 1.73




COSMOSIL Chromatogram Index

Sample: Todralazine
 CAS No.: [14679-73-3]
 Molecular formula: $C_{11}H_{12}N_4O_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =95/5
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV240 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.5mg/ml
 Injection volume: 0.5µl
 Retention time: 4.56min
 Capacity factor: 0.51





COSMOSIL Chromatogram Index

Sample: Trichloroacetic Acid
 CAS No.: [76-03-9]
 Molecular formula: $C_2HCl_3O_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =50/50
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV210 nm
 Attenuation: 0.128 aufs
 Sample conc.: 1.0mg/ml
 Injection volume: 1.0µl
 Retention time: 11.83min
 Capacity factor: 3.17





COSMOSIL Chromatogram Index

Sample: Tris(hydroxymethyl)aminomethane
 CAS No.: [77-86-1]
 Molecular formula: $C_4H_{11}NO_3$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =80/20
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: ELSD
 Attenuation: Gain=6, Atten=8
 Sample conc.: 2.0mg/ml
 Injection volume: 1.0µl
 Retention time: 6.47min
 Capacity factor: 1.48





COSMOSIL Chromatogram Index

Sample: L-Tryptophan
 CAS No.: [73-22-3]
 Molecular formula: $C_{11}H_{12}N_2O_2$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Citrate buffer(pH7.0) =70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV254 nm
 Attenuation: 0.128 aufs
 Sample conc.: 0.5mg/ml
 Injection volume: 0.5µl
 Retention time: 5.14min
 Capacity factor: 0.95




COSMOSIL Chromatogram Index

Sample: L-Tyrosine
 CAS No.: [60-18-4]
 Molecular formula: $C_9H_9NO_3$
 Column: HILIC
 Column size: 4.6mm I.D.-250mm
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30
 Flow rate: 1.0 ml/min
 Temperature: 30°C
 Detection: UV 254nm
 Attenuation: 0.128 aufs
 Sample conc.: 5.0mg/ml
 Injection volume: 1.0µl
 Retention time: 5.92min
 Capacity factor: 1.25




COSMOSIL Chromatogram Index

Sample: Uracil
CAS No.: [66-22-8]
Molecular formula: $C_4H_4N_2O_2$
Column: HILIC
Column size: 4.6mm I.D.-250mm
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =90/10
Flow rate: 1.0 ml/min
Temperature: 30°C
Detection: UV260 nm
Attenuation: 0.128 aufs
Sample conc.: 0.1mg/ml
Injection volume: 0.5µl
Retention time: 5.33min
Capacity factor: 0.84




COSMOSIL Chromatogram Index

Sample: Urea
CAS No.: [57-13-6]
Molecular formula: CH_4N_2O
Column: HILIC
Column size: 4.6mm I.D.-250mm
Mobile phase: Acetonitrile/ H_2O =90/10
Flow rate: 1.0 ml/min
Temperature: 30°C
Detection: UV210 nm
Attenuation: 0.128 aufs
Sample conc.: 10.0mg/ml
Injection volume: 2.0µl
Retention time: 6.12min
Capacity factor: 1.15




COSMOSIL Chromatogram Index

Sample: Uridine
CAS No.: [58-96-8]
Molecular formula: $C_9H_{12}N_2O_6$
Column: HILIC
Column size: 4.6mm I.D.-250mm
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate =90/10
Flow rate: 1.0 ml/min
Temperature: 30°C
Detection: UV260 nm
Attenuation: 0.128 aufs
Sample conc.: 0.1mg/ml
Injection volume: 1.0µl
Retention time: 8.30min
Capacity factor: 1.86



COSMOSIL Chromatogram Index

Sample: L-Valine
CAS No.: [72-18-4]
Molecular formula: $C_6H_{11}NO_2$
Column: HILIC
Column size: 4.6mm I.D.-250mm
Mobile phase: Acetonitrile/ 10mmol/l Ammonium acetate=70/30
Flow rate: 1.0 ml/min
Temperature: 30°C
Detection: UV210 nm
Attenuation: 0.128 aufs
Sample conc.: 10.0mg/ml
Injection volume: 1.0µl
Retention time: 5.63min
Capacity factor: 1.14



ナカライテスク株式会社

〒604-0855 京都市中京区二条通烏丸西入東玉屋町498

Web site

<http://www.nacalai.co.jp>

試薬はここに

価格・納期のご照会

0120-489-552

製品に関する技術的なご照会

E-mail: info-tech@nacalai.co.jp

TEL:075-211-2703 FAX:075-211-2673

■販売取扱店