

## ChromTech Application Bulletin No.7 (2001)

### Strategies for bioanalytical work on the **CHIRAL-AGP** column

The **CHIRAL-AGP** column has been used for the determination of many different kinds of enantiomers present at low concentration in biological material, i.e. plasma and urine. The updated reference list for **CHIRAL-AGP** contains more than 200 publications. Around 50% of these publications discuss bioanalytical methods.

In analysis of drugs from biological material there are at least two main sources for interferences in the chromatogram:

1. endogenous compounds
2. metabolites

Although the **CHIRAL-AGP** column shows high enantioselectivity for a wide variety of compounds, the selectivity for structure-related compounds, i.e. parent compound-metabolite, may be lower.

There are several methods that can be used to overcome these interference problems.

### Different ways to separate enantiomers in complex samples (plasma, urine)

The examples given for each way are taken from the **CHIRAL-AGP** reference list. The latest version of the reference list can be found in the new "Chiral Application Handbook" and on [www.chromtech.se](http://www.chromtech.se).

#### 1. **CHIRAL-AGP** column alone

Retention and enantioselectivity are regulated with the mobile phase composition, i.e. pH, organic modifier nature and concentration, buffer type and concentration.

MS detection is preferable due to the high detection selectivity, however also UV and fluorescence detectors are used.

Examples: Bupivacaine, Mepivacaine, Ibuprofen, Naproxen, Metoprolol, Verapamil, Alfuzosin, Oxamniquine, Methadone, Ketoprofen, Labetalol, Vamicamide, Mefloquine, Ifosfamide, Ketorolac, Ketamine

#### 2. **CHIRAL-AGP** column + preparation of enantiomeric derivatives

Can affect both chiral selectivity and selectivity between endogenous compounds and the drug. When preparing an enantiomeric derivative, the derivatization reagent used is nonchiral, as for example acetylation of amines by acetic acid anhydride or an acid chloride. Acids can be transformed to enantiomeric ester derivatives.

Examples: Atenolol, Verapamil

#### 3. Non-chiral column coupled in series with the **CHIRAL-AGP** column

Used to increase the resolution between endogenous compounds and the drug and the metabolites

Examples: Disopyramide, Methadone,

#### 4. Non-chiral reversed-phase column in off-line combination with the **CHIRAL-AGP** column

Collect the eluate from the reversed-phase column and inject it onto the **CHIRAL-AGP** column after concentration. Suitable if the mobile phase on the reversed phase column is not suitable for obtaining selectivity on the **CHIRAL-AGP** column.

Examples: Verapamil, Hydroxychloroquine, Terfenadine, Mosapride, Albendazole sulphoxide

#### 5. Non-chiral normal phase column in off-line combination with the **CHIRAL-AGP** column

Collect the eluate, containing the drug and the metabolite separately, from the normal phase column, evaporate the solvent, dissolve in the mobile phase used for the chiral chromatography and inject on the chiral column. The method is used for very complicated samples

Examples: Disopyramide, Chloroquine, Trimipramine, Mianserin

#### 6. Column Switching

Cut a band from a non-chiral reversed-phase column and switch it automatically over to the **CHIRAL-AGP** column using a valve. The method is used for very complicated samples.

Examples: Metoprolol, Verapamil, Bupivacaine, Warfarin, Terbutaline, Amlodipine