



Method Optimisation

Moving to Fortis SpeedCore
particles

Transfer of Methods

- Reduction of current HPLC methods
 - 1 – ‘Legacy’ method transferred to new technology
 - 2 – ‘New’ method back-scaled to prep or production scale
 - 3 – Method transfer between differing systems (Equivalence studies)


Our Goal




Trimethoprim

➤ Original Method:


- 5µm Fortis C18 250x4.6mm
- A: Phosphate pH3
- B: MeOH
- 70:30 → 100%B in 20min
- 1ml/min
- 280nm



250x4.6mm Original Method



Adjust k' for increased R_s



Efficiency HPLC → UHPLC

Column Length	Efficiency of 5µm	Efficiency of 3µm	Efficiency of 2.6µm Core-shell
250	22000		
150	12700	16800	26460
100	8300	10700	21000
50	4000	6000	11200
30		3200	7000
20			3000

↑
Increase Speed
Save Solvent
↓
8 fold

Moving from a 250mm 5µm column to a 150mm 3µm will provide the same efficiency

Therefore equivalent Separating Power

Calculations

- 3 Critical Equations
 - Alter Flow Rate – In line with i.d. change
 - Alter Gradient time
 - Calculate dwell volume

Calculations

Alter Flow Rate

$$F_2 = F_1 \times \frac{(dc_2^2)}{(dc_1^2)}$$

Alter Gradient time

$$t_{g2} = t_{g1} \times \frac{V_{m2} \times F_1}{V_{m1} \times F_2}$$

F = flow rate (ml/min)
 dc = column diameter (mm)
1 = Original column parameter
2 = New column parameter
 T_g = gradient time (min)
 V_m = interstitial column of column

Calculations

Interstitial volume

$$V_m = \pi \times r^2 \times L \times W$$

V_m = interstitial column of column

r = column radius (mm)

L = column length (mm)

W = column % interstitial porosity

Approximation

W for fully porous material = 68% = 0.68

W for core-shell material = 55% = 0.55

1st Change

250x4.6mm 5μ to 150x4.6mm 3μ

Alter Flow Rate

*N/a as both 4.6mm
i.d. columns*


Alter Gradient time

$$t_{g2} = t_{g1} \times \frac{V_{m2} \times F_1}{V_{m1} \times F_2}$$

$$t_{g2} = \frac{20 \times 1695 \times 1.0}{2825 \times 1.0}$$

$$t_{g2} = 12min$$

150x4.6mm First Method Change



Efficiency HPLC → UHPLC

Column Length	Efficiency of 5µm	Efficiency of 3µm	Efficiency of 2.6µm Core-shell
250	22000		
150	12700	16800	26460
100	8300	10700	21000
50	4000	6000	11200
30		3200	7000
20			3000

↑ Increase Speed
↓ Save Solvent
8 fold

Moving from a 150mm 3µm column to a 100mm 2.6µm will provide the same efficiency

Therefore equivalent Separating Power

2nd Change

150x4.6mm to 100x4.6mm SpeedCore

Alter Flow Rate

*N/a as both 4.6mm
i.d. columns*


Alter Gradient time

$$t_{g2} = t_{g1} \times \frac{V_{m2} \times F_1}{V_{m1} \times F_2}$$

$$t_{g2} = \frac{20 \times 914 \times 1.0}{2825 \times 1.0}$$

$$t_{g2} = 6.5min$$

100x4.6mm SpeedCore Change



Calculator



Method Transfer Calculator

Current HPLC Method

Adjust Column Length

Existing Column length	250	mm
Existing Particle Size	5	µm
Existing Column Diameter	4.6	mm

Adjust injection Volume

Existing Injection Volume	20	µl
---------------------------	----	----

Adjust Flow Rate

Existing flow rate	1.00	ml/min
--------------------	------	--------

Adjust Gradient Program

Existing Gradient Time	20	min
------------------------	----	-----

Backpressure

Existing Backpressure	105	bar
-----------------------	-----	-----

Solvent Consumption

Existing Solvent Used	20	ml
-----------------------	----	----

SpeedCore Core-Shell Meth

New Column Length	100	mm
New Particle Size	2.6	µm
New Column diameter	4.6	mm

New injection volume	8.00	µl
----------------------	------	----

New flow rate	1.00	ml/min
---------------	------	--------

New Gradient Time	6.5	min
If higher flow rate is required then please enter here	1.50	ml/min
New Gradient Time will be	4.3	min



New Backpressure @ 1.0ml/mir	155	bar
New Backpressure @ 1.2ml/mir	186	bar
New Backpressure @ 1.5ml/mir	233	bar

New Solvent Usage	12	ml
-------------------	----	----

Saving in Time (%) 67.65 %

Saving in Solvent (%) 38 %

Final optimised method



Allowable Changes

Type of Change	USP Adjustment	EP Adjustment
Column Length	± 70%	± 70%
Column i.d.	± 25%	± 25%
Particle size	- 50%	- 50%
Column Temperature	± 10%	± 10%
Flow rate	± 50%	± 50%
Mobile phase pH	± 0.2 units	± 0.2 units
Concentration of buffer	± 10%	± 10%
Ratio of components in mobile phase	± 30% of the minor component, but a change in any component cannot exceed ± 10% absolute	± 30% of the minor component, or 2% absolute of that component, whichever is greater
Injection Volume	Reduced as far as consistent with accepted precision and detection limits	Increased by as much as twice the volume specified, providing no adverse effects

Summary

- Moving to a Core-shell particle is simple as long as 3 equations are adhered to:
 - Change Flow rate
 - Alter gradient time
 - Calculate dwell volume

Conclusion

- Moving to Speedcore particles provides gains in:
 - Time
 - Sensitivity
 - Resolution
 - Less solvent usage