

Improve Your USP LC Gradient Separation Methods Following the New <621> Guidelines

Stephanie A. Schuster, Ph.D.

Advanced Materials Technology, Inc.



LCGC





February 28, 2023

- Benefits of converting to methods that use superficially porous particle (SPP) columns
- Review guidelines of USP Chapter <621> for changing isocratic and gradient methods
- Case study examples
 - Isocratic
 - Gradient
- Summary

The background of the slide is a dark blue gradient with a bokeh effect of light blue and white particles, creating a starry or particle-like appearance. The text is centered in the lower half of the image.

Benefits of Converting to Methods that Use SPP Columns

How SPP Benefits Separations



Effect of Particle Size and Type

Columns: 4.6 x 50 mm
 5 μm FPP C18
 3.5 μm FPP C18
 1.8 μm FPP C18
 2.7 μm HALO C18

Solute: naphthalene
 Mobile phase: 60% ACN/40% water
 Temperature: 24 °C

van Deemter Equation

H = height equivalent to theoretical plate

A = eddy diffusion term (particle size and how well bed was packed) **30 - 40% smaller**


B = longitudinal diffusion term **25 - 30% smaller**

C = resistance to mass transfer term (kinetics of the analyte b/w mobile phase and stationary phase)

μ = mobile phase linear velocity (L/t₀)


$$H = A + \frac{B}{\mu} + C\mu$$

Increased Speed & Reduced Mobile Phase Consumption



PEAK IDENTITIES:

1. homovanillic acid
2. caffeic acid
3. syringic acid
4. vanillic acid
5. chlorogenic acid
6. sinapic acid
7. ferulic acid
8. *p*-coumaric acid
9. *trans*-cinnamic acid
10. resveratrol




SPP method

- 4x faster
- 12x less solvent

Increased Efficiency Demonstrated Using Fat Soluble Vitamins


Sharper peaks and increased resolution with the HALO[®] C30 column!



PEAK IDENTITIES:

1. Retinyl acetate (A)
2. Delta tocopherol (E)
3. Ergocalciferol (D2)
4. Cholecalciferol (D3)
5. Alpha tocopherol (E)
6. DL-alpha-tocopherol acetate (E)
7. 2,3-*trans*-phyllloquinone (K)


Isocratic: 100% Methanol
Wavelength: 280nm
Injection: 2 μL
Temperature: 30 °C
Flow Rate: 1.5 mL/min
Columns: 4.6 x 150 mm



**USP Chapter <621> Guidelines
for Changing
Isocratic and Gradient Methods**

When Did Changes Go Into Effect

- For United States Pharmacopeia (USP), changes went into effect December 1, 2022



The screenshot shows the USP website's navigation bar with links for Global Presence, Select Country, Careers, Store, Help, FAQs, USP en Español, Contact Us, and Login. Below the navigation bar is the USP logo and a search bar. The main content area displays a notice titled "<621> Chromatography" with the following details:

- Type of Posting:** Notice of Adoption of Harmonized Standard
- Posting Date:** 19-Nov-2021
- Official Date:** 01-Dec-2022
- Expert Committee:** General Chapters–Chemical Analysis (GCCA)
- Coordinating Pharmacopeia:** European Pharmacopeia

The text below the details states: "The new harmonized standard for <621> Chromatography has been approved by the Pharmacopeial Discussion Group (PDG) as described in its PDG Stage 4 sign-off cover sheet. The chapter has been formally approved by the General Chapters–Chemical Analysis Expert Committee in accordance with the Rules and Procedures of the Council of Experts."

The changes from the current official version of the chapter are:

- For British Pharmacopoeia (BP), European Pharmacopoeia (EP), and Japanese Pharmacopoeia (JP), changes went into effect January 2023

Allowable Changes to USP Methods



Parameter	Isocratic Methods	Gradient Methods
Stationary Phase	Must keep same L category Change from totally porous particle (TPP) to superficially porous particle (SPP) is allowed	Same as isocratic
Column Dimensions	Particle size and/or length of the column may be modified, provided that ratio of the column length (L) to the particle size (d_p) remains constant or in the range between -25% to +50% of the prescribed L/d_p ratio*	Same as isocratic^
Internal Diameter	Change allowed	Change allowed
Flow Rate	$F_2 = F_1 \times [(dc_2^2 \times dp_1) \div (dc_1^2 \times dp_2)]$	Same as isocratic


*For changes from TPP to SPP, other combinations of L and d_p can be used, provided that the plate number (N) is within -25% to +50%, relative to the prescribed column. These changes are acceptable, provided that system suitability criteria are fulfilled, and selectivity and elution order of the specified impurities to be controlled are demonstrated to be equivalent

^For changes from TPP to SPP, other combinations of L and d_p can be used provided that the ratio $(t_R/W_h)^2$ is within -25% to +50%, relative to the prescribed column for all the peaks used to determine the system suitability parameters. These changes are acceptable provided system suitability criteria are fulfilled, and selectivity and elution order of the specified impurities to be controlled are demonstrated to be equivalent.


Allowable Changes to USP Methods - continued

Parameter	Isocratic Methods	Gradient Methods
Detector Wavelength	No change allowed	No change allowed
Mobile Phase Composition	Amount of the minor components of the mobile phase can be adjusted by $\pm 30\%$ relative. However, the change in any component cannot exceed $\pm 10\%$ absolute	Gradient time adjusted by equation $t_{G2} = t_{G1} \times (F_1/F_2) [(L_2 \times d_{c2}^2)/(L_1 \times d_{c1}^2)]$
pH	± 0.2 pH units, unless otherwise prescribed	Same as isocratic
Buffer	$\pm 10\%$	Same as isocratic
Injection Volume	$V_{inj2} = V_{inj1} \times (L_2 d_{c2}^2)/(L_1 d_{c1}^2)$	Same as isocratic
Column Temperature	$\pm 10^\circ\text{C}$	$\pm 5^\circ\text{C}$

- https://ispso.unige.ch/labs/fanal/hplc_calculator:en



- <https://www.acdlabs.com/resources/freeware/translator/index.php>



Translate Liquid Chromatography Methods Quickly and Easily

ACD/Labs have teamed up with experienced chromatographers Patrik Petersson and Mel Euerby to provide an application for translating liquid chromatography methods.

We also offer a full suite of commercially available chromatography software.
[Compare the functionality of the freeware and commercial products.](#)

INSTRUMENTAL AND COLUMN PARAMETERS				TYPICAL DWELL VOLUME ESTIMATES	
	Original	New Scaled	New Selected	Alliance	1.40 mL
Flow (mL/min)	<input type="text" value="0.20"/>		<input type="text" value="0.20"/>	Acquity UPLC "classic" with 50 µL mixer	0.10 mL
				Acquity UPLC "classic" with 385 µL mixer	0.44 mL

Words of Caution for Moving to Smaller Particle Sizes and Smaller Column I.D.s




- Smaller particle sizes and smaller column I.D.s are more susceptible to extra-column band broadening

Words of Caution for Moving to Smaller Particle Sizes and Smaller Column I.D.s

- Smaller particle sizes and smaller column I.D.s are more susceptible to extra-column band broadening


$$\sigma_{V,tot}^2 = \sigma_{V,pre-col}^2 + \sigma_{V,col}^2 + \sigma_{V,post-col}^2$$



Words of Caution for Moving to Smaller Particle Sizes and Smaller Column I.D.s

- Smaller particle sizes and smaller column I.D.s are more susceptible to extra-column band broadening

$$\sigma_{V,tot}^2 = \sigma_{V,pre-col}^2 + \sigma_{V,col}^2 + \sigma_{V,post-col}^2$$



Ways to Reduce Extra-Column Band Broadening HALO®



- Use the shortest and smallest I.D. connecting tubing between instrument components



- Use a smaller volume detector flow cell ($\leq 2.5 \mu\text{L}$)
- Increase the data acquisition rate ($\geq 20 \text{ Hz}$)
- Decrease the injection volume (more critical for isocratic than gradient separations)

MarvelXACT and MarvelX are registered trademarks of IDEX.

Extra-Column Band Broadening Impact on Resolution and Efficiency



50% average increase in plates is observed by reducing the excess volume in the system!




Isocratic Separations with HALO 90 Å C18, 2 μ m, 3.0 x 150 mm



Isocratic Case Studies

What Changes are Possible?



$L/dp = 150/.005 = 30,000$
For -25-50%, L/dp can be
22,500-45,000

$L/dp = 150/0.005 = 30,000 =$
criteria met

37% higher plates



$L/dp = 100/.0027 = 37,037 =$
criteria met

57% higher plates



$L/dp = 50/0.002 = 25,000 =$
criteria met


3x times faster
41% higher plates



Increased Speed without Loss of Resolution



FPP C18
4.6 x 250 mm, 5 μ m
1 mL/min



HPLC in Standard Configuration

- ECV \sim 35 μ L
- Standard flow cell, 14 μ L
- 0.5 sec. response time
- Standard length and ID tubing (0.007" ID x 750 mm)


$$L/dp = 250/0.005 = 50,000$$

For -25-50%, L/dp can be 37,500-77,500

$$N = 25,900$$

For -25-50%, N can be 19,425-38,850

HALO C18
4.6 x 100 mm, 2.7 μ m
2 mL/min



HPLC in Ultra-Low ECV Configuration

- ECV \sim 10 μ L
- Semi-micro flow cell, 5 μ L
- 0.5 sec. response time
- Reduced length and ID tubing (0.005" ID x 460 mm)

$$L/dp = 100/0.0027 = 37,037$$

= criteria not met, but is met for plates

$$N = 24,500$$


6 times faster run time!



USP Monograph for Estradiol – Modified



HALO 90 Å C18, 2.7 μm
1.5 x 100 mm



	Original Method	Modified Method
Column	3.9 x 300 mm	1.5 x 100 mm
Flow Rate (mL/min)	1	0.2
Time (min)	8	1.5
Volume (mL)	8	0.3

$$L/dp = 300/0.01 = 30,000$$

For -25-50%, L/dp can be 22,500-45,000

$$L/dp = 100/0.0027 = 37,037$$

96% solvent savings!



Gradient Case Studies

While Making Changes, Take Notice!

- Adjustment of conditions with gradient elution (HPLC) ...is more critical than with isocratic (HPLC) ...elution, since it may shift some peaks to a different step of the gradient ..., potentially causing partial or complete coelution of adjacent peaks or peak inversion, and, thus leading to the incorrect assignment of peaks and to the masking of peaks or a shift such that elution occurs beyond the prescribed elution time.




- For some parameters the adjustments are explicitly defined in the monograph to ensure the system suitability.

USP Monograph for Cobamamide – Original




$L/dp = 100/0.0026 = 38,462$

For -25-50%, L/dp can be 28,846-57,693




$L/dp = 100/0.0027 = 37,037$ ✓



- Cobamamide (labeled peak) is also known as adenosylcobalamin and is one of the biologically active forms of Vitamin B12
- Suitability requirements for cobamamide
 - Column efficiency: NLT 22,000 plates
 - Tailing factor: NMT 2.0


	Monograph Column	HALO C18
Efficiency	Pass	Pass
Tailing Factor	Pass	Pass

USP Monograph for Cobamamide – Original System Suitability



Peak Identities

1. hydroxocobalamin chloride
2. cyanocobalamin
3. cobamamide
4. methylcobalamin



Elution order is identical
on both columns


Modifying the Original Method

$$t_{G2} = t_{G1} \times (F_1/F_2) [(L_2 \times d_{c2}^2)/(L_1 \times d_{c1}^2)]$$

Flow rate 1 (F1)	1.2
Flow rate 2 (F2)	0.4
Length 1 (L1)	100
Length 2 (L2)	75
Column Diameter 1 (dc1)	4.6
Column Diameter 2 (dc2)	2.1

tG1	tG2	%B
0	0	15
2	0.94	15
2.5	1.17	20
8	3.75	23
10	4.69	25
22	10.32	30
24	11.25	35
30	14.07	35
31	14.54	70
40	18.76	70
40.5	18.99	15
45	21.10	15



- Flow rate increased to 0.4 mL/min instead of 0.25 mL/min (same linear velocity as 1.2 mL/min)
- Equilibration time increased by 4 minutes for 20 column volumes for a total run time of 25.5 min



- Injection volume reduced to 0.8 μ L from 5 μ L by this equation:

$$V_{inj2} = V_{inj1} \times (L_2 d_{c2}^2)/(L_1 d_{c1}^2)$$

USP Monograph for Cobamamide – Modified




Method was moved to a shorter SPP column with smaller i.d. and smaller particle size

Total run time is 1.8 times faster!


54 mL/run compared to 10 mL/run which is 5.4 times less solvent used!

USP Monograph for Cobamamide – Modified System Suitability



Peak Identities

1. hydroxocobalamin chloride
2. cyanocobalamin
3. cobamamide
4. methylcobalamin




Elution order is identical
on both columns


Particle Size Distribution Comparison



2.6 μm SPP




HALO[®] 2.7 μm




- No significant difference in particle size
- Slightly better standard deviation for HALO[®]

USP Monograph for Itraconazole – Original




- Itraconazole is an antifungal medication
- **Suitability Requirements**
 - Tailing factor: NMT 2.0
 - RSD: NMT 2.0%




	Monograph Column	HALO C18
Tailing Factor	Pass	Pass
RSD	Pass	Pass

USP Monograph for Itraconazole – Modified




Method was moved to a shorter SPP column with smaller i.d. and smaller particle size




Total run time is > 3 times faster with the HALO[®] column!

USP Monograph for Itraconazole: Mobile Phase Reduction

For 10 injections:



compared to



10X reduction
in mobile
phase used!

- HALO[®] columns are made for ruggedness with high efficiency and high speed separations.
- It is important to optimize the instruments in order to gain the most benefit from HALO[®] column technology.
- Following new USP guidelines for gradient method modernization enables both FPP and SPP methods to be improved for speed and mobile phase savings using HALO[®] Fused-Core[®] column technology.

Acknowledgements



- Advanced Materials Technology, Inc.
 - Peter Pellegrinelli
 - Conner McHale
 - Dr. Taylor Harmon
 - Ben Libert
 - William Miles
 - Dr. Chuping Luo
 - Dan Messick
 - Melissa Ash
 - Stephanie Rosenberg

support@advanced-materials-tech.com



Photo by [Alexander Grey](#) on [Unsplash](#)



advancedmaterialstechnology



halocolumns.com



3521 Silverside Road, Suite 1-K
Quillen Building
Wilmington, DE 19810



(302) 992-8060

