


MCI[®] GEL

Mitsubishi chemical's packed columns and packing material for HPLC



ISO 9001 Certified
Working to enhance quality management



excellent performance

spherical and sharp particle size distribution

persistence and highest quality

offers packing materials and packed columns,
under strict quality control

wide range of product line

MCI® GEL has been designed based on technology of
the world famous Diaion® and Sepabeads®,
specialized in polymeric packing medias including
from analytical to preparative use,
for ion exchange, reversed phase mode

abundant accumulation of technology and experience

for more than 20 years, MCI® GEL has been used for
HPLC applications

CONTENTS

1	Column selection guide	3 ~ 4
2	Ion exchange columns	5 ~ 20
	Column list	5 ~ 6
	[Applications] Amino acids 《CK10U》	7 ~ 8
	Sugars • Sugar alcohols • Organic acids 《CK08E series》	9 ~ 12
	Examples of peak retention time 《CK08E series》	13
	Oligo saccharides 《CK04S • CK02A series》	14 ~ 16
	Sugars • Organic acids 《CA08F》	17 ~ 18
	Nucleic acids, etc. • Sugars • Human urine 《CDR10》	19 ~ 20
3	Packed columns for ion chromatography	21 ~ 26
	Column list	22
	[Applications] Cations 《SCK01》	22 ~ 23
	Anions 《SCA04》	24 ~ 26
4	Bioseparation columns	27 ~ 38
	Column list	27
	Size exclusion chromatography columns	
	《CQP series》	
	Column list	28
	[Applications] Calibration curve 《CQP series》	29
	Proteins • Water soluble polymers 《CQP series》	29
	Ion exchange chromatography columns	
	《ProtEx series》	
	Column list	30
	[Applications] Proteins	31 ~ 34
	《CQA • CQK series》	
	Column list	35
	[Applications] Proteins	35 ~ 36
	Hydrophobic interaction chromatography columns	
	《CQH series》	
	Column list	37
	[Applications] Proteins	38
5	Reversed phase columns	39 ~ 51
	Polymeric packing material columns	
	Column list	40
	Column durability	41
	[Applications] 《CHP10M • CHP5C • CHP2MG • CHP2MGM》	42 ~ 51
6	Chiral separation columns	52 ~ 57
	Separation mechanism of CRS series	52
	[Applications] 《CRS10W • CRS15W》	
	DL- α -Amino acids, DL- α -Hydroxy carboxylic acids	53 ~ 56
	Separation conditions for various amino acids	57
7	Chromatography media for preparative use	58 ~ 64
8	MCI[®] GEL columns	65 ~ 66
9	MCI[®] GEL chromatography media	67 ~ 71
10	Compounds index	72 ~ 79

1

MCI® GEL

Column selection guide

Column selection guide

1

Ion exchange columns

2

Packed columns for ion chromatography

3

Bioseparation columns

4

Reversed phase columns

5

Chiral separation columns

6

Chromatography media for preparative use

7

MCI® GEL columns

8

MCI® GEL chromatography media

9

Compounds Index


10

Nature of sample	Separation mode	MCI® GEL column	pH range	Applications	pages		
Sample	Water Soluble	M.W. >2,000	Size Exclusion	CQP10 CQP30	2 ~ 12	Proteins, Biopolymers Water soluble polymers	28 ~ 29
			Ion Exchange	ProtEx Series CQA Series CQK Series	2 ~ 12	Proteins, Biopolymers Proteins, Biopolymers	30 ~ 34 35 ~ 36
		Hydrophobic	CQH Series	2 ~ 12	Proteins, Biopolymers	37 ~ 38	
		Reversed Phase	CHP2MG CHP2MGM CHP5C CHP10M	2 ~ 12 1 ~ 13	Proteins, Biopolymers Proteins, Biopolymers	47 ~ 50	
		Size Exclusion	CK02AS CK02A	6 ~ 7	Oligosaccharides	14 ~ 16	
			CK04SS CK04S	6 ~ 7	Oligosaccharides	14 ~ 16	
	M.W. <2,000	Ion Exchange	CQP06	2 ~ 12	Peptides	28 ~ 29	
			CK10U CDR10 CA08F	1 ~ 14 1 ~ 13	Amino acids Organic acids Saccharides	7 ~ 8 17 ~ 20	
		Reversed Phase	CDR10 SCA04 SCK01	1 ~ 13 3 ~ 7 1.5 ~ 12	Nucleotides Anions Cations	19 24 ~ 26 22 ~ 23	
			CHP2MG CHP2MGM CHP5C CHP10M	2 ~ 12 1 ~ 13	Organic Compounds Organic Compounds	40 ~ 50 40 ~ 50	
			CK08EH	1 ~ 7	Organic acids	11 ~ 12	
		Ligand Exchange	CK08E Series CRS10W CRS15W	5 ~ 7	Saccharides Optical isomers (α -amino acids · α -hydroxy carboxylic acids)	9 ~ 10 50 ~ 55	
Organic Solvent Soluble	M.W. <2,000	Reversed Phase	CHP5C CHP10M CHP207M CHP2MG CHP2MGM CHPOD1M	1 ~ 13 2 ~ 12	Organic Compounds Organic Compounds	40 ~ 51 40 ~ 51	

Column selection guide

Variation of MCI® GEL products line

Particle size [μm]	Analytical			Preparative		
	5	10	30	50	150	
Ion exchange	CK ProtEx	CK CA CDR10	CK CA CQA_S CQK_S	CK CA PrepEx CQA_P CQK_P		CK CA
Ion chromatography	SCA	SCK				
Size exclusion		CQP	CQP_P			
Hydrophobic		CQH3_S	CQH3_P			
Reversed phase	CHP10M CHP2MGM	CHP5C CHP2MG	CHP55A CHP20A	CHP55Y CHP20Y CHP2MGY	CHP20P	CHP20P
Ligand exchange	CRS_W					



Column selection guide

1

Ion exchange columns

2

Packed columns for ion chromatography

3

Bioprocession columns

4

Reversed phase columns

5

Chiral separation columns

6

Chromatography media for preparative use

7

MCI® GEL columns

8

MCI® GEL chromatography media

9

Compounds Index

10

Polystyrenic cation exchange resins

- CK series -

Polystyrenic anion exchange resins

- CA series -

Mitsubishi Chemical Ion Exchange Resins

MCI® GEL specializes in polymer based packing materials. Specifically, polystyrene polymer based ion exchange resins are derived from over 45 years of manufacturing experience of Diaion® product line. MCI® GEL ion exchange resins for HPLC have been developed with the same attention to performance and quality. For several decades, Mitsubishi Chemical has been providing MCI® GEL ion exchange columns are offered in a variety of chemistries, particle sizes and counter ions to support a broad range of applications.

Features


- Variety of products** gel type, porous type, DVB%, particle size, particle size distribution
analytical use, preparative use
- Persistence of high quality, excellent separation performance**
- Accumulation of abundant knowledge and experience of applications**

Ion exchange resins are generally used for analysis of amino acids, sugars, organic acids and amines, etc. MCI® GEL custom pre-packed columns are specifically designed for each application using the most appropriate packing material among our product line and using the most suitable column dimensions.


Typical application for each column is shown in this catalog. These data will suggest an appropriate column.

Chemical structure of ion exchange resin

Strongly acidic cation exchange resin



Strongly basic anion exchange resin



MCI® GEL Columns for HPLC

	Description					Typical usage					
	Product name	Column dimensions I.D × L (mm)	Packing material			Amino acid	Mono, small-saccharide	Oligo-saccharide	Carboxylic acid	Amine	Physiological fluid
			Cross linkage (%)	Counter ion	Particle size (μm)						
MCI® GEL Cation exchange columns	MCI® GEL CK10U	6 × 120	10	Na ⁺	5	○				○	
	MCI® GEL CK08S	8 × 500	8	Na ⁺	11		○				
	MCI® GEL CK08E	8 × 300	8	Na ⁺	9		○				
	MCI® GEL CK08EC	8 × 300	8	Ca ²⁺	9		○				
	MCI® GEL CK08ES	8 × 300	8	Ag ⁺	9		○	○			
	MCI® GEL CK08EH	8 × 300	8	H ⁺	9		○		○	○	
	MCI® GEL CK04S	10 × 200	4	Na ⁺	11			○			
	MCI® GEL CK04SS	10 × 200	4	Ag ⁺	11			○			
	MCI® GEL CK02A	20 × 250	2	Na ⁺	20			○			
	MCI® GEL CK02AS	20 × 250	2	Ag ⁺	20			○			
MCI® GEL Anion exchange columns	MCI® GEL CA08F	4.6 × 250	8	Cl ⁻	7		○		○		
	MCI® GEL CDR10	4.6 × 250	porous	AcO ⁻	7		○		○		○

Description of a gel type ion exchange column

MCI® GEL CK08EC

for HPLC use

Cation = K
Anion = A

Counter ion

(no letter = Na⁺, C = Ca²⁺)
(S = Ag⁺, H = H⁺)

Particle size(mode)

(A = 20μm, S = 11μm)
(E = 9μm, F = 7μm,
U = 5μm)

Note ; Pre-column and guard column


- Please consider using a guard column concerning purity of injection sample. Guard columns, are listed in the end of this catalog, should be selected in accordance with a main column.
- As for analysis of amino acids by MCI® GEL CK10U, MCI® GEL AFR2-PC is recommended as a pre-column. The AFR2-PC column is very effective to stabilize base line because it can trap ammonium ion in eluent. A peak caused of the ammonium ion may disturb base line stability.



CK10U 6 x 120

Separation of amino acids

Fig. 2-1 Protein hydrolyzates amino acids



As for analysis of amino acids by a cation exchange column such as MCI® GEL CK10U, MCI® GEL AFR2-PC is recommended as a pre-column. The AFR2-PC column is very effective to stabilize base line because ammonium in eluent is trapped in this column. The ammonium ion may disturb base line stability. The AFR2-PC should be installed between an outlet of HPLC pump and an inlet of sample injector. A gradient elution, commonly used for amino acid analysis, is influenced by HPLC instrument. So to obtain a satisfactory chromatogram, gradient conditions should be optimized in accordance with the HPLC equipment.

Separation of amino acids

Fig. 2-2 Valine, β -Alanine






Fig. 2-3 γ -Aminobutyric acid



CK08E series

Cation exchange columns applications ; sugars, carboxylic acids, (poly)alcohols, etc.



CK08EC 8 × 300

CK08EH 8 × 300

Column list

MCI® GEL Column	Counter ion	Application areas
MCI® GEL CK08S MCI® GEL CK08E	Na ⁺	General sugar separation columns
MCI® GEL CK08EC	Ca ²⁺	The most general sugar separation column Highly recommended for fructose and glucose
MCI® GEL CK08ES	Ag ⁺	Gel permeation chromatographic effect
MCI® GEL CK08EH	H ⁺	Organic acids with H ₃ PO ₄ eluent; sugars with distilled water eluent

Application data of CK08EC

Fig. 2-4 Sugars

Fig. 2-5 Apple juice

Conditions
 Column : MCI® GEL CK08EC
 8mm I.D. × 300mm
 Eluent : H₂O
 Flow rate : 0.6ml/min
 Column temp. : 75
 Detection : RI






Fig. 2-6 Lactulose syrup



Application data of CK08EC

Fig. 2-7 Sports drink A




Fig. 2-8 Sports drink B




Fig. 2-9 Honey




Fig. 2-10 Jam





Fig. 2-11 Sugars/Alcohols



Conditions
 Column : MCI® GEL CK08EC
 8mm I.D. x 300mm
 Eluent : H₂O
 Flow rate : 0.6mL/min
 Column temp. : 75
 Detection : RI

Application data of CK08EH

Fig. 2-12 Carboxylic acids




Fig. 2-13 Amino sugars






Fig. 2-14 Alcohols



Conditions
 Column : MCI® GEL CK08EH, 8mm I.D. x 300mm
 Eluent : 1% H₃PO₄(Fig.2-12,2-13), H₂O(Fig. 2-14)
 Flow rate : 0.6ml/min
 Column temp. : 45 (Fig. 2-12), ambient(Fig. 2-13), 60 (Fig. 2-14)
 Detection : 210nm(Fig. 2-12), RI(Fig. 2-13, 2-14)

Fig. 2-15 Sugar alcohols/Alcohols

Conditions
 Column : MCI® GEL CK08EH
 8mm I.D. x 300mm
 Eluent : H₂O
 Flow rate : 0.6ml/min
 Column temp. : 45
 Detection : RI



- 1 Column selection guide
- 2 Ion exchange columns
- 3 Packed columns for ion chromatography
- 4 Bioseparation columns
- 5 Reversed phase columns
- 6 Chiral separation columns
- 7 Chromatography media for preparative use
- 8 MCI® GEL columns
- 9 MCI® GEL chromatography media
- 10 Compounds Index

Application data of CK08EH

Fig. 2-16 Poly alcohols

Conditions
 Column : MCI® GEL CK08EH
 8mm I.D. x 300mm
 Eluent : 1% H₃PO₄
 Flow rate : 0.6ml/min
 Column temp. : 25
 Detection : RI




Fig. 2-17 Chloroacetic acids

Conditions
 Column : MCI® GEL CK08EH
 8mm I.D. x 300mm
 Eluent : 1% H₃PO₄
 Flow rate : 0.6ml/min
 Column temp. : 45
 Detection : 210nm





Fig. 2-18 Carboxylic acids

Conditions
 Column : MCI® GEL CK08EH
 8mm I.D. x 300mm
 Eluent : 2% H₃PO₄
 Flow rate : 0.6ml/min
 Column temp. : ambient
 Detection : 210nm



Elution time for Sugars and Sugar alcohols on various columns[min]

CK08EC Ca ²⁺	CK08E Na ⁺	CK08ES Ag ⁺
Stachyose 9	Stachyose 8	* Melezitose 12
Melezitose 10	Melezitose 9	* Stachyose 13
Raffinose 11	Raffinose 10	* Raffinose 14
Gentiobiose 12	Gentiobiose 11	* Sucrose 15
Cellobiose 13	Cellobiose 12	Trehalose 16
Trehalose 14	Trehalose 13	Cellobiose 17
Isomaltose 15	Sucrose 14	Gentiobiose 18
Sucrose 16	Isomaltose 15	Maltose 19
Maltose 17	Melibiose 16	Isomaltose 20
Melibiose 18	Maltose 17	Maltulose 21
Lactose 19	Maltulose 18	Lactose 22
Maltulose 20	Lactose 19	Melibiose 23
Lactulose 21	Lactulose 20	Lactulose 24
Glucose 22	Glucose 21	Adonitol 25
Xylose 23	Mannitol 22	Adonitol 26
Galactose 24	Rhamnose 23	Digitoxose 27
Mannose 25	Rhamnose 24	Rhamnose 28
Rhamnose 26	Adonitol 25	Glucose 29
Fructose 27	Sorbitol 26	Xylose 30
Fucose 28	Digitoxose 27	Xylitol 31
Inositol 29	Mannose 28	Erythritol 32
Arabinose 30	Xylose 29	Mannitol 33
Digitoxose 31	Galactose 30	Fructose 34
Adonitol 32	Fructose 31	Dulcitol 35
Erythritol 33	Inositol 32	Galactose 36
Mannitol 34	Xylitol 33	Sorbitol 37
Salicin 35	Fucose 34	Mannose 38
Dulcitol 36	Dulcitol 35	Arabinose 39
Xylitol 37	Arabinose 36	Fucose 40
Sorbitol 38	Erythritol 37	Ribose 41
Ribose 39	Ribose 38	Ribose 42
	Salicin 39	Inositol 43
		Salicin 52

Column temp : CK08EC...75 , CK08E...45 , CK08ES...75
 Column size : 8mm I.D. x 300mm
 Eluent : H₂O
 Flow rate : 0.6ml/min
 Sample : 1% aq. solution
 Injection vol. : 20µl

* ; These sugars, containing Fructose component, may partially be decomposed by CK08ES and CK08EH.

Column selection guide
 1
 Ion exchange columns
 2
 Packed columns for ion chromatography
 3
 Bioseparation columns
 4
 Reversed phase columns
 5
 Chiral separation columns
 6
 Chromatography media for preparative use
 7
 MC® GEL columns
 8
 MC® GEL chromatography media
 9
 Compounds Index
 10

CK04S, CK04SS,
CK02A, CK02AS

Cation exchange columns
applications ; oligosaccharides

The separation mechanism is based on gel filtration chromatography and elution is achieved via simple distilled water. A larger molecule elutes ahead.



CK02A 20 x 250



CK04S 10 x 200



CK04SS 10 x 200

Separation ability of each column

MCI® GEL column	Counter ion	Separation ability (degrees of polymerization)
MCI® GEL CK04S	Na ⁺	8 ~ 9
MCI® GEL CK04SS	Ag ⁺	12 ~ 13
MCI® GEL CK02A	Na ⁺	15 ~ 16
MCI® GEL CK02AS	Ag ⁺	19 ~ 20

Calibration curve of malto-oligosaccharides

Fig. 2-19





Fig. 2-20



Comparison data of malto-oligosaccharides

Fig. 2-21 MCI® GEL CK04S
10mm I.D. x 200mm




Fig. 2-22 MCI® GEL CK04SS
10mm I.D. x 200mm




Fig. 2-23 MCI® GEL CK02A
20mm I.D. x 250mm





Fig. 2-24 MCI® GEL CK02AS
20mm I.D. x 250mm



Conditions
 Eluent : H₂O
 Flow rate : 0.4 ml/min (Fig. 2-21, 2-22, 2-25, 2-26)
 1.0 ml/min (Fig. 2-23, 2-24, 2-27)
 Column temp. : 85
 Detection : RI

* On Fig. 2-21 to 2-27, the numbers indicate degree of polymerization.

Column selection guide
1

Ion exchange columns
2

Packed columns for ion chromatography
3

Bioseparation columns
4

Reversed phase columns
5

Chiral separation columns
6

Chromatography media for preparative use
7

MCI® GEL columns
8

MCI® GEL chromatography media
9

Compounds Index
10

Comparison data of authentic malto-oligosaccharides samples

Fig. 2-25 MCI® GEL CK04S
10mm I.D. x 200mm




Fig. 2-26 MCI® GEL CK04SS
10mm I.D. x 200mm





Fig. 2-27 MCI® GEL CK02AS
20mm I.D. x 250mm



Application data of CK04S

Fig. 2-28 Honey




Fig. 2-29 Jam






Fig. 2-30 Fructo-oligosaccharides



Conditions
 Column : MCI® GEL CK04S
 10mm I.D. x 200mm
 Eluent : H₂O
 Flow rate : 0.4ml/min (Fig. 2-28, 2-29) 0.3ml/min (Fig. 2-30)
 Column temp. : 85 (Fig. 2-28, 2-29) 45 (Fig. 2-30)
 Detection : RI

Application data of CA08F


Fig. 2-31 Sugars



Conditions

Column : MCI® GEL CA08F
4.6mm I.D. x 250mm
Eluent : A 0.15M Borate buffer pH7.5
B 0.5M Borate buffer pH9.5
C 0.6M Borate buffer pH9.5
D 0.7M Borate buffer pH8.5
4-step-gradient
Flow rate : 0.5 ml/min
Column temp.: 65
Post column reaction
Reagent : 2.5% Boric acid,
5% Monoethanolamine pH7.9
Flow rate : 1.0 ml/min
Temp. : 150
Detection : Fluorescence Ex 360nm, Em 440nm

Fig. 2-32 Carboxylic acids




Conditions

Column : MCI® GEL CA08F
4.6mm I.D. x 250mm
Eluent : 0.6M Na₂SO₄ pH3.0
Flow rate : 0.5 ml/min
Column temp.: 60
Detection : 210nm


Application data of CA08F

Fig. 2-33 Carboxylic acids



Conditions
Column : MCI® GEL CA08F
 : 4.6mm I.D × 250mm
Eluent : 0.6M Na₂SO₄ pH2.0
Flow rate : 0.5ml/min
Column temp. : 60
Detection : 210nm

Fig. 2-34 Carboxylic acids




Conditions
Column : MCI® GEL CA08F
 : 4.6mm I.D × 250mm
Eluent : 0.1M NaH₂PO₄ pH3.3
Flow rate : 0.4ml/min
Column temp. : 55
Detection : 210nm

Packing material of MCI® GEL CDR10 column is based on a highly porous polystyrene functionalized with a quaternary ammonium anion exchange resin. Since a high porous type ion exchange resin is rigid, CDR10 allows usage of aggressive gradient elution, for example water to 6M of acetate buffer gradient. MCI® GEL CDR10 is highly recommended for rapid analysis of physiological fluids like urine and blood.

Application data of CDR10


Fig. 2-35 Nucleic acids and related substances



Conditions

Column : MCI® GEL CDR10
4.6mm I.D. x 250mm
Eluent : A H₂O
B 6M Acetate buffer pH4.4
A B 30min linear gradient
Flow rate : 0.5ml/min
Column temp.: 60
Detection : 254nm

Fig. 2-36 Mono-nucleotides




Conditions

Column : MCI® GEL CDR10
4.6mm I.D. x 250mm
Eluent : 1M Acetate buffer pH3.3
Flow rate : 1.2ml/min
Column temp.: 60
Detection : 254nm


Application data of CDR10

Fig. 2-37 Sugars



Conditions
 Column : MCI® GEL CDR10
 4.6mm I.D. x 250mm
 Eluent : A 0.15M Borate buffer pH7.5
 B 0.6M Borate buffer pH9.5
 A B60min linear gradient
 Flow rate : 0.5ml/min
 Column temp.: 65
 Post column reaction
 Reagent : 2.5% Boric acid, 5% Monoethanolamine pH7.9
 Flow rate : 0.5ml/min
 Temp. : 150
 Detection : Fluorescence Ex 360nm, Em 440nm

Fig. 2-38 human urine



Conditions
 Column : MCI® GEL CDR10
 4.6mm I.D. x 250mm
 Eluent : A 0.006M Acetate buffer pH4.4
 B 6M Acetate buffer pH4.4
 A B 60min. linear gradient
 Flow rate : 1.0ml/min
 Column temp.: 60
 Detection : 254nm

Packed columns for ion chromatography

For cation analysis : MCI® GEL SCK01, 6mmI.D × 50mm

For anion analysis : MCI® GEL SCA04, 4.6mmI.D × 150mm

The MCI® GEL ion chromatography columns are based on surface functionalized cation and anion exchange resins designed for non-suppressed ion chromatography applications. The non-suppressed ion chromatography is an analysis technique of cations and anions with combination of a packed column of low capacity ion exchange resin and low concentration of electrolyte solution as an eluent. The advantage of the ion chromatography is that several ions can be analyzed by only one injection with free of complicated sample pre-treatment.

For cation analysis MCI® GEL SCK01

Packing material of MCI® GEL SCK01 is crosslinked polystyrene functionalized with sulfonic acid. This column is characterized by excellent resolution and rapid analysis for monovalent and divalent cations. Standard monovalent cations like Li⁺, Na⁺, NH₄⁺, K⁺, Rb⁺, Cs⁺ and simple amines such as mono-, di- and trimethylamine can be resolved using a nitric acid solution as eluent. Divalent cations, such as alkaline earth metals and transition metal elements, can be efficiently resolved using tartaric acid and complexing reagent such as ethylene diamine to selectively elute the metals from the column.

Note:

When using the MCI® GEL SCK01 column for monovalent cations, it is recommended that a pre-column, MCI® GEL SCK-PC, be used to trap heavy metals which might otherwise poison the SCK01 column resulting in a rapid loss of capacity and chromatographic performance.

For anion analysis MCI® GEL SCA04

Packing material of MCI® GEL SCA04 is based on a hydrophilic vinyl polymer matrix functionalized with quaternary ammonium group and particle size of 5 μm. A solution of potassium hydrogen phthalate and a vanilic acid (VA)/N-methyldiethanolamine (MDEA) solution both can be used as a mobile phase. The unique VA/MDEA eluent, is developed for the SCA04 column, which allows users to determine 7 standard anions in 14 minutes without system peak.

Note:

A pre-column, MCI® GEL SCA-PC is recommended for prevention of contamination to the SCA04 column when the VA/MDEA eluent is used. The SCA-PC is effectively prolong SCA04 column life. The SCA-PC should be installed between an outlet of HPLC pump and an sample injector.

Column list

Cation analysis	MCI® GEL SCK01	6mm I.D × 50mm	stainless steel column
Cation analysis	MCI® GEL SCK01	4.6mm I.D × 150mm	stainless steel column
Pre-column for cation analysis	MCI® GEL SCK-PC	6mm I.D × 50mm	stainless steel column
Anion analysis	MCI® GEL SCA04	4.6mm I.D × 150mm	PEEK column
Pre-column for anion analysis	MCI® GEL SCA-PC	8mm I.D × 10mm	stainless steel column

Application data of SCK01

Fig. 3-1 Monovalent cations




Fig. 3-2 Amines




Fig. 3-3 Monovalent cations in rain




Fig. 3-4 Monovalent cations in tap water




Fig. 3-5 Sports drink





Fig. 3-6 Nucleoside



Conditions
 Column : MCI® GEL SCK01 6mm I.D. × 50mm
 Eluent : 5mM HNO₃
 Flow rate : 1.0ml/min
 Column temp. : 40
 Detection : Conductivity(Fig. 3-1, 3-2, 3-3, 3-4, 3-5) 254nm(Fig. 3-6)

Application data of SCK01

Fig. 3-7 Alkaline earth metals




Fig. 3-8 Transition metals




Fig. 3-9 Divalent cations




Fig. 3-10 Sports drink A





Fig. 3-11 Sports drink B



Conditions

- Column : MCI® GEL SCK01 6mm I.D. × 50mm
(On Fig. 3-8, two columns are connected in series)
- Eluent : 2mM Tartaric acid, 1.5mM Ethylenediamine (Fig. 3-7, 3-9, 3-10, 3-11)
1.5mM Tartaric acid, 0.8mM Ethylenediamine (Fig. 3-8)
- Flow rate : 1.0mL/min
- Column temp. : 40
- Detection : Conductivity

Column selection guide
1
Ion exchange columns
2
Packed columns for ion chromatography
3
Bioseparation columns
4
Reversed phase columns
5
Chiral separation columns
6
Chromatography media for preparative use
7
MCI® GEL columns
8
MCI® GEL chromatography media
9
Compounds Index
10



SCA04 4.6 × 150 PEEK

Application data of SCA04

Fig. 3-12 Standard anions eluent ; VA/MDEA





Fig. 3-13 Standard anions eluent ; Potassium hydrogenphthalate



Application data of SCA04

Fig. 3-14 Rain

A; Beginning of rain fall

B; After 4 hours

C; After 38 hours





Fig. 3-15 River water



Conditions (Fig. 3-14, 3-15)

Column : MCI® GEL SCA04
4.6mm I.D. x 150mm

Eluent : 3mM VA 2.8mM MDEA pH6.2

Flow rate : 1.2ml/min

Column temp. : 40

Detection : Conductivity

- 1 Column selection guide
- 2 Ion exchange columns
- 3 Packed columns for ion chromatography
- 4 Bioseparation columns
- 5 Reversed phase columns
- 6 Chiral separation columns
- 7 Chromatography media for preparative use
- 8 MCI® GEL columns
- 9 MCI® GEL chromatography media
- 10 Compounds Index

Application data of SCA04

Fig. 3-16 Sulfur compounds





Fig. 3-17 Instant coffee



4

MCI® GEL

Bioseparation columns

For size exclusion chromatography

MCI® GEL CQP series

For ion exchange chromatography

MCI® GEL ProtEx series,

MCI® GEL CQA/CQK series

For hydrophobic interaction chromatography

MCI® GEL CQH series

MCI® GEL bioseparation columns

MCI® GEL bioseparation columns are based on a hydrophilic, wide pore and rigid polymer designed for analytical chromatography of proteins, peptides, enzymes and other biomolecules.

MCI® GEL CQP series is for size exclusion chromatography.

For ion exchange chromatography, MCI® GEL ProtEx series and MCI® GEL CQA/CQK series are used. MCI® GEL ProtEx series columns are unique and brilliant packed columns provide excellent separation of proteins, good protein selectivity and high protein recovery. Specifically, proteins of small structural differences (isoforms) can be effectively separated and small amount of proteins (less than several tens µg) can be quantitatively recovered without nonspecific adsorption. From that point of view, the ProtEx columns can be applied in the field of purification of small amount of protein to obtain sample for structural determination and quality control of proteinaceous pharmaceuticals.

MCI® GEL CQH series is for hydrophobic interaction chromatography.

Column name	Separation mode	Type
MCI® GEL CQP06	Size exclusion	Exclusion limit MW ~ 10 ³
MCI® GEL CQP10	"	" ~ 10 ⁴
MCI® GEL CQP30	"	" ~ 10 ⁶
MCI® GEL ProtEx-DEAE	Anion exchange	DEAE
MCI® GEL CQA31S	"	DEAE
MCI® GEL CQA35S	"	QA
MCI® GEL ProtEx-SP	Cation exchange	SP
MCI® GEL CQK30S	"	SP
MCI® GEL CQK31S	"	CM
MCI® GEL CQH3ES	Hydrophobic interaction	Ether
MCI® GEL CQH3BS	"	Butyl
MCI® GEL CQH3PS	"	Phenyl

Size exclusion chromatography

Size exclusion chromatography is a liquid chromatographic technique which separates solute molecules according to their size in solution. The column is packed with porous particles and separation takes place as a result of the differential solute distribution outside and within the pores of the packing material. Solute molecules which are larger than the pores of the packing material will be excluded and therefore will elute first and have a lower retention time than the smaller one. The CQP series columns based on a hydrophilic polymer are designed for analysis of water soluble polymers such as oligosaccharides and PEG, etc.

Column list

CQP Series

MCI® GEL column	Column dimensions	Packing materials		Theoretical plates number (TP / column)	Exclusion limit (PEG)
		Particle size[μm]	Pore size[nm]		
MCI® GEL CQP06	7.5mm I.D. x 600mm	10	12	10000	$\sim 1 \times 10^3$
MCI® GEL CQP10	7.5mm I.D. x 600mm	10	20	6000	$\sim 1 \times 10^4$
MCI® GEL CQP30	7.5mm I.D. x 600mm	10	60	6000	$\sim 1 \times 10^6$

Guard columns

MCI® GEL column	Column dimensions
MCI® GEL CQP06G	4.0mm I.D. x 50mm
MCI® GEL CQP10G	4.0mm I.D. x 50mm
MCI® GEL CQP30G	4.0mm I.D. x 50mm

Packing materials available
Package size; 10 g , 25 g , 50 g

Application data of CQP series

Fig. 4-1 Calibration curve

Conditions
 Column : MCI® GEL CQP06
 CQP10
 CQP30
 7.5mm I.D. x 600mm
 Eluent : H₂O
 Flow rate : 1.0ml/min
 Column temp.: ambient
 Detection : RI
 Sample : PEG 100µl inj.




Fig. 4-2 Separation of PEG mixture

Conditions
 Column : MCI® GEL CQP30 7.5mm I.D. x 600mm
 Eluent : H₂O
 Flow rate : 1.0ml/min
 Column temp.: 25
 Detection : RI
 Sample : 1. PEG 145,000
 2. 40,000
 3. 6,000




Fig. 4-3 Separation of protein mixture

Conditions
 Column : MCI® GEL CQP30 7.5mm I.D. x 600mm
 Eluent : 14mM Tris-HCl/O₄ buffer
 Flow rate : 1.0ml/min
 Column temp.: ambient
 Detection : 280nm
 Sample : 1. Ferritin(MW440,000)
 2. Ovalbumin(MW43,000)
 3. Myoglobin(MW17,500)
 4. Cytochrome α (MW12,400)





Fig. 4-4 Separation of gluconic acid and glucose

Conditions
 Column : MCI® GEL CQP06 7.5mm I.D. x 600mm
 Eluent : H₂O
 Flow rate : 0.8ml/min
 Column temp.: ambient
 Detection : RI
 Sample : 1. 5% Gluconic acid
 2. 5% Glucose




ProtEx series

Ion exchange chromatography columns



ProtEx-DEAE 4.6 x 50



ProtEx-SP 4.6 x 50



ProtEx-DEAE 7.5 x 100



ProtEx-SP 7.5 x 100

Separation mechanism and Characteristic of ProtEx columns

MCI® GEL ProtEx series packed columns are for ion exchange chromatography mode which separates sample proteins mainly via ionic interaction between packing material and sample molecules.

The packing materials for ProtEx series columns are based on 5 μm, mono disperse, porous type, methacrylate polymer, are specifically designed for separation of proteins.


On a conventional protein separation column, non-specific adsorption of sample proteins is sometimes occurs resulting in loss of valuable sample. But on the ProtEx columns, non-specific adsorption is eliminated because the surface of the packing material is surrounded by hydrophilic layer is chemically bonded to base material and ion exchange functional group are effectively increased.

Two types of ion exchange columns, weakly basic diethylaminoethyl (DEAE) type and strongly acidic sulfopropyl (SP) type are available.


Column list

Column name	Column dimensions	Column format	Packing material		pH range
			Particle size	Functional group	
MCI® GEL ProtEx-DEAE	4.6mm I.D. x 50mm	PEEK	5μm	diethylaminoethyl	2 ~ 12
	7.5mm I.D. x 100mm	PEEK	5μm	diethylaminoethyl	2 ~ 12
MCI® GEL ProtEx-SP	4.6mm I.D. x 50mm	PEEK	5μm	sulfopropyl	1 ~ 13
	7.5mm I.D. x 100mm	PEEK	5μm	sulfopropyl	1 ~ 13

Packing material of ProtEx-DEAE



Scanning electron micrograph



Surface of ProtEx-DEAE

Application data of ProtEx series

Fig. 4-5 Separation of proteins mixture

Conditions
 Column : MCI® GEL ProtEx-DEAE 4.6mm I.D. x 50mm
 Eluent : A 20mM Tris-HCl pH8.0
 B A + 0.5M NaCl
 A B 30min linear gradient
 Flow rate : 0.5ml/min
 Column temp. : ambient
 Detection : 280nm
 Sample : 1. Myoglobin 25µg
 2. Conalbumin 25µg
 3. Trypsin Inhibitor 50µg




Fig. 4-6 Separation of hemoglobin (Hb) isoforms

Conditions
 Column : MCI® GEL ProtEx-DEAE 4.6mm I.D. x 50mm
 Eluent : A 20mM Tris-HCl pH8.0
 B A + 0.5M NaCl
 A 10% B 30min linear gradient
 Flow rate : 0.5ml/min
 Column temp. : ambient
 Detection : 280nm
 Sample : 1. Hb A₂ 100µg
 2. Hb S 100µg
 3. Hb A₀ 100µg




Fig. 4-7 Separation of human growth hormone (hGH)

Conditions
 Column : MCI® GEL ProtEx-DEAE 4.6mm I.D. x 50mm
 Eluent : A 20mM Tris-HCl pH8.0
 B A + 0.5M NaCl
 5% B 70% B 30min linear gradient
 Flow rate : 0.5ml/min
 Column temp. : ambient
 Detection : 280nm
 Sample : recombinant hGH 10µg





Fig. 4-8 Protein recovery

Conditions
 Column : MCI® GEL ProtEx-DEAE 4.6mm I.D. x 50mm
 Eluent : A 20mM Tris-HCl pH8.15
 B A + 0.5M NaCl
 A 50% B 30min linear gradient
 Flow rate : 0.5ml/min
 Column temp. : ambient
 Detection : 280nm
 Sample : recombinant epidermal growth factor (EGF)



Application data of ProtEx series

Fig. 4-9 Separation of interleukin 2 (IL-2) coexisting large amount of bovine serum albumin (BSA) as a stabilizer

Conditions
 Column : MCI® GEL ProtEx-DEAE 4.6mm I.D. x 50mm
 Eluent : A 20mM Trimethylenediamine-HCl pH9.75
 B A + 0.5M NaCl
 Flow rate : 0.5ml/min
 Column temp. : ambient
 Detection : 280nm
 Sample : recombinant IL-2 1.5µg
 1. BSA (stabilizer) 400µg
 2. IL-2 (Met-ox)
 3. IL-2




Fig. 4-10 Separation of protein mixture

Conditions
 Column : MCI® GEL ProtEx-SP 4.6mm I.D. x 50mm
 Eluent : A 20mM Phosphate buffer pH6.0
 B A + 0.5M NaCl
 A B 20min linear gradient
 Flow rate : 0.5ml/min
 Column temp. : ambient
 Detection : 280nm
 Sample : 1. Ribonuclease A 10µg
 2. α-Chymotrypsinogen A 5µg
 3. Cytochrome C 5µg





Fig. 4-11 Separation of glycohemoglobin (GHb)

Conditions
 Column : MCI® GEL ProtEx-SP 4.6mm I.D. x 50mm
 Eluent : A 20mM Bis-Tris HCl buffer pH6.0
 B A + 0.5M NaCl
 7% B 40% B 20min linear gradient
 Flow rate : 0.5ml/min
 Column temp. : ambient
 Detection : 415nm
 Sample : GHb
 1. GHb A1c
 2. Hb A0



Application data of ProtEx series

Fig. 4-12 Separation of mouse brain sap




Fig. 4-13 Separation of RNA







Fig. 4-14 Separation of IgG2b, (mouse)



Application data of ProtEx series

Fig. 4-15 Separation of collagenase




Fig. 4-16 Separation of IgG1 MOPC21(mouse)

Conditions

- Column : MCI® GEL ProtEx-DEAE 4.6mm I.D. x 50mm
- Eluent : A 10mM HEPES buffer pH8.0
B A + 0.5M NaCl
A100% B100% in 30min. linear gradient
- Flow rate : 0.5ml/min
- Column temp. : 25
- Detection : 280nm
- Sample : IgG1 MOPC21 (mouse) 10µg





Fig. 4-17 Separation of pancreatin

Conditions

- Column : MCI® GEL ProtEx-DEAE 4.6mm I.D. x 50mm
- Eluent : A 20mM HEPES buffer pH8.0
B A + 1.0M NaCl
A100% B40% in 30min. linear gradient
- Flow rate : 0.5ml/min
- Column temp. : 25
- Detection : 280nm
- Sample : Pancreatin 200µg



CQA series,
CQK series

Ion exchange chromatography columns

CQA and CQK series packed columns are for ion exchange chromatography mode which separates sample proteins mainly via ionic interaction between packing material and sample molecules.

Four types of ion exchange columns, strongly basic quaternary ammonium (QA), weakly basic diethylaminoethyl (DEAE), strongly acidic sulfopropyl (SP) and weakly acidic carboxymethyl (CM) are available.

Column list

CQA, CQK Series

Column name	Column dimensions	Packing material		pH range
		Particle size [μm]	Functional group	
MCI® GEL CQA31S	7.5mm I.D. × 75mm	10	DEAE	2 ~ 12
MCI® GEL CQA35S	7.5mm I.D. × 75mm	10	QA	2 ~ 12
MCI® GEL CQK30S	7.5mm I.D. × 75mm	10	SP	1 ~ 13
MCI® GEL CQK31S	7.5mm I.D. × 75mm	10	CM	4 ~ 13

Packing materials available

Package size; 10 g , 25 g , 50 g

Application data of CQA and CQK series

Fig. 4-18 Separation of protein mixture

Conditions
 Column : MCI® GEL CQA31S 7.5mm I.D. × 75mm
 MCI® GEL CQA35S 7.5mm I.D. × 75mm
 Eluent : A 14mM Tris-HCl buffer pH8.2
 B A + 0.5M NaCl
 A B 30min linear gradient
 Flow rate : 1.0ml/min
 Column temp.: ambient
 Detection : 280nm
 Sample : 1. Myoglobin 60μg
 2. Ovalbumin 200μg
 3. Trypsin Inhibitor 200μg







Fig. 4-19 Separation of protein mixture

Conditions
 Column : MCI® GEL CQA31S 7.5mm I.D. × 75mm
 MCI® GEL CQA35S 7.5mm I.D. × 75mm
 Eluent : A 14mM Tris-HCl buffer pH8.2
 B A + 0.5M NaCl
 A B 30min linear gradient
 Flow rate : 1.0ml/min
 Column temp.: ambient
 Detection : 280nm
 Sample : 1. Myoglobin 120μg
 2. Transferrin 160μg
 3. β-Lactoglobulin 400μg





Application data of CQA and CQK series

Fig. 4-20 Separation of protein mixture




Conditions
 Column : MCI® GEL CQK30S 7.5mm I.D. x 75mm
 MCI® GEL CQK31S 7.5mm I.D. x 75mm
 Eluent : A 20mM Phosphate buffer pH6.8
 B A + 0.5M NaCl
 A B 30min linear gradient
 Flow rate : 1.0ml/min
 Column temp.: ambient
 Detection : 280nm
 Sample : 1. Myoglobin 50µg
 2. Ribonuclease A 200µg
 3. Cytochrome C 60µg
 4. Lysozyme 80µg

Fig. 4-21 Separation of protein mixture




Conditions
 Column : MCI® GEL CQK30S 7.5mm I.D. x 75mm
 MCI® GEL CQK31S 7.5mm I.D. x 75mm
 Eluent : A 20mM Phosphate buffer pH6.8
 B A + 0.5M NaCl
 A B 30min linear gradient
 Flow rate : 1.0ml/min
 Column temp.: ambient
 Detection : 280nm
 Sample : 1. Trypsinogen 80µg
 2. Ribonuclease A 90µg
 3. α-Chymotrypsinogen A 60µg

Fig. 4-22 Separation of lipoxidase



Conditions
 Column : MCI® GEL CQA31S 7.5mm I.D. x 75mm
 Eluent : A 14mM Tris-HCl buffer pH8.2
 B A + 0.5M NaCl
 A B 30min linear gradient
 Flow rate : 1.0ml/min
 Column temp.: ambient
 Detection : 280nm
 Sample : Lipoxidase



CQH series

Hydrophobic interaction columns

MCI® GEL CQH series packed columns are for hydrophobic interaction chromatography mode. Functional groups of the packing materials are phenyl, butyl and ether.

The relative hydrophobicity of the CQH series columns decrease in the following order. CQH3PS>CQH3BS>CQH3ES

Column list

CQH Series

Column name	Column dimensions	Particle size (µm)	Functional group
MCI® GEL CQH3ES	7.5mm I.D. x 75mm	10	Ether
MCI® GEL CQH3BS	7.5mm I.D. x 75mm	10	Butyl
MCI® GEL CQH3PS	7.5mm I.D. x 75mm	10	Phenyl

Packing materials available

Package size; 10 g , 25 g , 50 g

Application data of CQH series

Fig. 4-23 Separation of human serum

Conditions
 Column : MCI® GEL CQH3ES 7.5mm I.D. x 75mm
 MCI® GEL CQH3PS 7.5mm I.D. x 75mm
 Eluent : A B + 1.7M(NH₄)₂SO₄
 B 0.1M Phosphate buffer pH6.8
 A B 60min linear gradient
 Flow rate : 1 ml/min
 Column temp.: ambient
 Detection : 280nm
 Sample : Human serum




Fig. 4-24 Separation of colibacillus extract

Conditions
 Column : MCI® GEL CQH3ES 7.5mm I.D. x 75mm
 Eluent : A B + 1.7M(NH₄)₂SO₄
 B 0.1M Phosphate buffer pH6.8
 A B 30min linear gradient
 Flow rate : 1.0ml/min
 Column temp.: ambient
 Detection : 280nm
 Sample : Colibacillus extract




Fig. 4-25 Separation of colibacillus extract

Conditions
 Column : MCI® GEL CQH3PS 7.5mm I.D. x 75mm
 Eluent : A B + 1.7M(NH₄)₂SO₄
 B 0.1M Phosphate buffer pH6.8
 A B 30min linear gradient
 Flow rate : 1.0ml/min
 Column temp.: ambient
 Detection : 280nm
 Sample : Colibacillus extract





Fig. 4-26 Separation of mixture of peptides

Conditions
 Column : MCI® GEL CQH3PS 7.5mm I.D. x 75mm
 Eluent : A B + 1.7M(NH₄)₂SO₄
 B 0.1M Phosphate buffer pH6.8
 A B 30min linear gradient
 Flow rate : 1.0ml/min
 Column temp.: ambient
 Detection : 220nm
 Sample : 1. Met-Leu-Tyr
 2. Leu-Enkephalin
 3. Bacitracin



MCI® GEL CHP series

MCI® GEL Polymeric reversed phase columns

A partition chromatography, an adsorption chromatography, an ion exchange chromatography and a size exclusion chromatography are typical separation mechanisms of high performance liquid chromatography. The partition chromatography is most commonly used, separates solute samples in accordance with the difference of partition of the samples between a stationary phase and a mobile phase, can be applied to broad range of applications of organic compounds such as pharmaceuticals, agricultural chemicals and those intermediate substances. There are two separation mechanisms in the partition chromatography, one is a normal phase and the other is a reversed phase are discriminated by comparison of polarity of stationary phase and mobile phase. On the normal phase chromatography, a polarity of the stationary phase is stronger than that of the mobile phase. As for the reversed phase mode, the relationship of the polarities of the two phases reverses. The RP chromatography is the most popular separation mode is said that RP occupies 60-70 % of HPLC applications.

For packing materials of the RP chromatography, silica based and synthetic polymer based materials are widely used.

MCI® GEL specializes in polymer based packing materials. The use of polymeric based RP columns has become more widespread thanks to unique selectivity of the polymer matrix, no specific adsorption common with silica based packings and can be operated with a wide pH range, basic eluents and acidic eluents due to the chemical stability of the inert polymeric materials. The MCI® GEL reversed phase columns are based on a polystyrenic and polymethacrylate porous polymers are normally applied to the separation of aromatic and aliphatic based compounds in the isocratic and gradient elution modes. The applications include pharmaceuticals, steroids, small peptides, amphoteric molecules such as sulfonamides and cephalosporin antibiotics, plus basic drugs, simple amines, antihistamines and carbamate pesticides.

The MCI® GEL reversed phase packings are based on the same chemistries offered in the Diaion® and Sepabeads® synthetic adsorbents resins. These polymer chemistries, like Diaion® HP series and Sepabeads® SP series are widely used and documented in the biopharmaceutical industry for fermentation extraction, the food industry and industrial reversed phase separations. The MCI® GEL reversed phase packings are available as packed columns for analytical applications and as bulk packings for analytical, preparative and production chromatography applications.

MCI® GEL CHP20P bulk polymer resin is widely used as a reversed phase resin for preparative extractive separation in biopharmaceutical applications.

Columns packed with polymers, MCI® GEL CHP series, are suitable for reversed phase chromatography and there are four kinds of columns of various hydrophobicity; Porous polystyrene, Modified Porous polystyrene, Acrylates and Octadecyl-alkylated aliphatic Porous polymers. Thus proper kind of columns can be selected in accordance with the properties of the target compounds.

Modified polystyrene packing	: MCI® GEL CHP207M
Polystyrene packing	: MCI® GEL CHP10M, CHP5C
Octadecyl-alkylated aliphatic packing	: MCI® GEL CHPOD1M
Acrylates packing	: MCI® GEL CHP2MGM, CHP2MG

The hydrophobicities of the columns are in the following order:

MCI® GEL CHP207M > CHP10M=CHP5C > CHPOD1M ODS columns CHP2MGM

Polymer columns for HPLC, with superior chemical resistance, can be applied with various mobile phases of broad pH range, acidic through alkaline. They have the following advantages due to their high hydrophobicity:

- 1) In the reversed phase distribution chromatography to separate acidic or alkaline compounds, the eluents suppressing the ionic properties of such compounds are generally used. Polymer columns can be applied for the unsuitable compounds to ODS columns.
- 2) Some of high hydrophilic compounds, e.g. amino acids, can be separated with strong hydrophobic CHP207M column.
- 3) Polymer columns can be washed with acidic and/or basic solutions when deteriorated by contamination.

< Reversed phase distribution chromatography >

Acrylates, CHP2MGM and CHP2MG, can be applied not only for reversed phase distribution chromatography but also for normal phase one.



CHP10M 4.6 × 150

CHP5C 4.6 × 250

CHP2MGM 4.6 × 150

Column list

CHP series

Column Name	Chemical Structures	particle size	HPLC Column	Preparative Column	pH range
MCI® GEL CHP10M	Styrene Divinylbenzene	4μm	4.6mm* 150mmL	NEW 20mm* 150mmL	Whole range
MCI® GEL CHP5C	Styrene Divinylbenzene	10μm	4.6mm* 250mmL	—	Whole range
MCI® GEL CHP2MGM	Methacrylates	4μm	4.6mm* 150mmL	NEW 20mm* 150mmL	2 ~ 12
MCI® GEL CHP2MG	Methacrylates	10μm	4.6mm* 250mmL	—	2 ~ 12
MCI® GEL CHP207M	Brominated Styrene Divinylbenzene	4μm	NEW 4.6mm* 150mmL	(20mm* 200mmL)	Whole range
MCI® GEL CHPOD1M	C18-alkylated aliphatics	4μm	NEW 4.6mm* 150mmL	NEW 20mm* 200mmL	2 ~ 12


MCI® GEL CHP5C and CHP2MG are both available in bulk.

Durability of polymeric column

The polymeric RP columns are chemically stable. Specifically, the columns have resistance to an alkaline eluent. The following graphs demonstrate stability of the polymeric columns. After feeding a solution of pH 12 into the MCI® GEL CHP10M, there is no change of column performance.


Fig. 5-1 Column durability at pH12 comparison between CHP10M and an ODS column

Conditions
Column : MCI® GEL CHP10M 4.6mm I.D × 150mm
Eluent : 20mM Na₂HPO₄ pH12/CH₃CN=60/40
Flow rate : 0.4 mL/min
Column temp. : 25
Detection : 254nm
Sample : 1000ppm Dimethyl phthalate 5μ L




Application data of CHP series

Fig. 5-2 Separation of catecholamines



Conditions
Column : MCI® GEL CHP10M
4.6mm I.D. x 150mm
Eluent : 50mM Na-phosphate pH2.0,
1.5% Hexanesulfonic acid /
CH₃CN = 80/20
Flow rate : 0.25ml/min
Column temp. : ambient
Detection : 280nm
Sample : 1. Epinephrine
2. Dopamine
3. 5-Hydroxy tryptophan
4. Serotonin
5. Tryptophan


Fig. 5-3 Separation of phthalic acid esters



Conditions
Column : MCI® GEL CHP10M
4.6mm I.D. x 150mm
Eluent : H₂O/CH₃CN = 50/50
Flow rate : 0.75ml/min
Column temp. : 60
Detection : 254nm
Sample : 1. Dimethyl phthalate
2. Diethyl phthalate
3. Dipropyl phthalate


Application data of CHP series

Fig. 5-4 Purine alkaloids



Conditions
Column : MCI® GEL CHP10M
 : 4.6mm I.D. x 150mm
Eluent : H₂O/CH₃CN=10/90
Flow rate : 0.4ml/min
Column temp. : 25
Detection : 275nm
Sample : 1.Theophylline
 : 2.Theobromine
 : 3.Caffeine

Fig. 5-5 Cinchona alkaloids



Conditions
Column : MCI® GEL CHP10M
 : 4.6mm I.D. x 150mm
Eluent : 0.1M NaH₂PO₄ pH2.0
 : CH₃CN=88/12
Flow rate : 0.3ml/min
Column temp. : 25
Detection : 275nm
Sample : 1.Cinchonine
 : 2.Quinine

Column selection guide

1

Ion exchange columns

2

Packed columns for ion chromatography

3

Bioseparation columns

4

Reversed phase columns

5

Chiral separation columns

6

Chromatography media for preparative use

7

MCI® GEL columns

8

MCI® GEL chromatography media


9

Compounds Index

10

Application data of CHP series


Fig. 5-6 Uric acid and related compounds




Conditions
 Column : MCI® GEL CHP10M
 4.6mm I.D. × 150mm
 Eluent : 10mM TBA/CH₃CN=75/25
 Flow rate : 0.4ml/min
 Column temp. : 25
 Detection : 284nm
 Sample : 1.Hypoxanthine
 2.Xanthine
 3.Theophylline
 4.Uric acid
 5.Orotic acid

Comparison with an ODS column

Fig. 5-7 Bile acids



Conditions
 Column : MCI® GEL CHP10M
 4.6mm I.D. × 150mm
 Eluent : 50mM NaH₂PO₄ pH2.0
 /CH₃CN = 40/60
 Flow rate : 0.3ml/min
 Column temp. : 25
 Detection : 210nm
 Sample : 1.Cholic acid
 2.Ursodeoxycholic acid
 3.Deoxycholic acid



Conditions
 Column : ODS column (5micron)
 4.6mm I.D. × 150mm
 Eluent : 50mM NaH₂PO₄ pH6.0
 /CH₃CN = 40/60
 Flow rate : 0.3ml/min
 Column temp. : 25
 Detection : 210nm
 Sample : 1.Cholic acid
 2.Ursodeoxycholic acid
 3.Deoxycholic acid

Column
 selection
 guide

1

Ion
 exchange
 columns

2

Packed
 columns
 for
 ion
 chromatography

3

Bioseparation
 columns

4

Reversed
 phase
 columns

5

Chiral
 separation
 columns

6

Chromatography
 media
 for
 preparative
 use

7

MCI® GEL
 columns

8

MCI® GEL
 chromatography
 media


9

Compounds
 index

10


Application data of CHP series

Fig. 5-8 Glycyrrhizae radix



Conditions
Column : MCI® GEL CHP10M
4.6mm I.D. × 150mm
Eluent : 2.06% acetic acid/CH₃CN
=63/37
Flow rate : 0.5ml/min
Column temp. : 45
Detection : 254nm
Sample : extract of
glycyrrhizae radix

Fig. 5-9 Adrenal cortex hormones



Conditions
Column : MCI® GEL CHP10M
4.6mm I.D. × 150mm
Eluent : H₂O/CH₃CN
: 60/40
Flow rate : 0.5ml/min
Column temp. : 45
Detection : 280nm
Sample : 1. Hydrocortisone
2. Corticosterone
3. 11-Deoxycortisol

Column selection guide

1

Ion exchange columns

2

Packed columns for ion chromatography

3

Bioseparation columns

4

Reversed phase columns

5

Chiral separation columns

6

Chromatography media for preparative use

7

MCI® GEL columns

8

MCI® GEL chromatography media


9

Compounds Index

10

Application data of CHP series


Fig. 5-10 Ingredients of medicine




Conditions
 Column : MCI® GEL CHP2MGM
 4.6mm I.D. × 150mm
 Eluent : 50mM phosphoric acid(pH2.0)/CH₃OH
 = 60/40
 Flow rate : 0.5mL/min
 Column temp. : 45
 Detection : 254nm
 Sample : 1.4-Dimethylaminoantipyrine
 2.Antipyrine
 3.Caffeine
 4.Aspirin
 5.Phenacetin

Comparison with an ODS column

Fig. 5-11 Sulfa drugs




Conditions
 Column : MCI® GEL CHP2MGM
 4.6mm I.D. × 150mm
 Eluent : 20mM phosphate pH6.8/CH₃CN
 = 82/18
 Flow rate : 0.5mL/min
 Column temp. : 45
 Detection : 254nm
 Sample : 1.Succinylsulfathiazole
 2.Sulfanilamide
 3.Sulfathiazole
 4.Sulfamerazine



Conditions
 Column : ODS column
 4.6mm I.D. × 150mm
 Eluent : 20mM phosphate pH6.8/CH₃CN
 = 90/10
 Flow rate : 0.5mL/min
 Column temp. : 45
 Detection : 254nm
 Sample : 1.Succinylsulfathiazole
 2.Sulfanilamide
 3.Sulfathiazole
 4.Sulfamerazine


Application data of CHP series

Fig. 5-12 Peptides



Conditions
Column : MCI® GEL CHP2MGM
4.6mm I.D. × 150mm
Eluent : 0.1%TFA/CH₃CN
=70/30
Flow rate : 0.5ml/min
Column temp. : 25
Detection : 220nm
Sample : 1.Gly-Tyr
2.Met Enkephalin
3.Leu Enkephalin
4.Angiotensin II

Fig. 5-13 Proteins



Conditions
Column : MCI® GEL CHP2MGM
4.6mm I.D. × 150mm
Eluent : A;0.05%TFA/CH₃CN
=80/20
B;0.05%TFA/CH₃CN
=20/80
A B 30min.linear
Flow rate : 0.5ml/min
Column temp. : 25
Detection : 280nm
Sample : 1.Ribonuclease A
2.Cytochrome c
3. -Chymotrypsinogen A

Column selection guide

1

Ion exchange columns

2

Packed columns for ion chromatography

3

Bioseparation columns

4

Reversed phase columns

5

Chiral separation columns

6

Chromatography media for preparative use

7

MCI® GEL columns

8

MCI® GEL chromatography media

9

Compounds Index

10

Application data of CHP series

Fig. 5-14 Procainamide, Procaine






Fig. 5-15 Water-soluble vitamins




Application data of CHP series

Fig. 5-16 Tetracycline



Conditions
Column : MCI® GEL CHP5C
Eluent : 4.6mm I.D. x 250mm
: potassium phosphate buffer (pH9.0)
(the Japanese pharmacopoeia)
Flow rate : 0.5ml/min
Column temp.: 60
Detection : 254nm
Sample : tetracycline HCl heat-treatment 90
for 60min.

Fig. 5-17 Anticonvulsants



Conditions
Column : MCI® GEL CHP5C
4.6mm I.D. x 250mm
Eluent : A ; 0.01M Na₂HPO₄
pH7.0/CH₃CN = 70/20
B ; 100% CH₃CN
A B 25min linear gradient
Flow rate : 0.3ml/min
Column temp.: 25
Detection : 260nm
Sample : 1.Barbital
2.Phenobarbital
3.Carbamazepine
4.Diazepam

Column selection guide

1

Ion exchange columns

2

Packed columns for ion chromatography

3

Bioseparation columns

4

Reversed phase columns

5

Chiral separation columns

6

Chromatography media for preparative use

7

MCI® GEL columns

8

MCI® GEL chromatography media

9

Compounds Index

10

Application data of CHP series

Fig. 5-18 Oxine copper, Thiuram

Conditions
 Column : MCI® GEL CHP2MG
 Eluent : 4.6mm I.D. x 250mm
 : 25mM KH₂PO₄(pH3.0)/CH₃CN
 = 45/55
 Flow rate : 0.8ml/min
 Column temp. : 40
 Detection : 240nm
 Sample : 1. Oxine copper 3ppm
 2. Thiuram 3ppm




Fig. 5-19 Aromatic carboxylic acids

Conditions
 Column : MCI® GEL CHP2MG
 Eluent : 4.6mm I.D. x 250mm
 : 0.1M NaH₂PO₄(pH5.8)/CH₃CN
 = 90/10
 Flow rate : 0.5ml/min
 Column temp. : 25
 Detection : 260nm
 Sample : 1. Terephthalic acid 0.7%
 2. 4 -Carboxybenzaldehyde 100ppm
 3. Benzoic acid 100ppm
 4. p-Toluic acid 100ppm




Fig. 5-20 Proteins

Conditions
 Column : MCI® GEL CHP2MG
 4.6mm I.D. x 250mm
 Eluent : A 0.05% TFA/CH₃CN = 80/20
 B 0.05% TFA/CH₃CN = 30/70
 A B 45min linear gradient
 Flow rate : 0.5ml/min
 Column temp. : 25
 Detection : 280nm
 Sample : 1. Ribonuclease A
 2. Cytochrome C
 3. Transferrin
 4. α-Chymotrypsinogen A
 5. β-Lactoglobulin





Fig. 5-21 Rice Bran Oil

Conditions
 Column : MCI® GEL CHP2MG(10mm)
 4.6A I.D. x 150A
 Eluent : Hexane-EtOH = 98/2 (vol.)
 Flow rate : 0.5ml/min
 Detection : 295nm
 Sample : Rice Bran Oil, 50g/ml
 Injection : 10µL



Column selection guide

1

Ion exchange columns

2

Packed columns for ion chromatography

3

Bioseparation columns

4

Reversed phase columns

5

Chiral separation columns

6

Chromatography media for preparative use

7

MCI® GEL columns

8

MCI® GEL chromatography media

9

Compounds index

10

Comparison between CHP207M and CHP10M

(Polyphenon 60)

Fig. 5-22 Modified Styrene Divinylbenzene
CHP207M

Conditions
 Column : MCI®GEL CHP207M
 4.6A I.D. x 150mm
 Eluent : CH₃OH/10mM-Acetic acid = 60/40
 Flow rate : 0.46ml/min
 Column temp.: 0.46ml/min
 Detection : 60
 Sample : 280nm
 : Polyphenon 60(10mg/ml) each 10mL





Fig. 5-23 Styrene Divinylbenzene CHP10M

Conditions
 Column : MCI®GEL CHP10M
 4.6mm I.D. x 150mm
 Eluent : CH₃OH/10mM-Acetic acid = 60/40
 Flow rate : 0.46ml/min
 Column temp.: 60
 Detection : 280nm
 Sample : Polyphenon 60(10mg/ml) each 10μL



Comparison with an ODS column

(TritonX-100)

Fig. 5-24 CHPOD1M

Conditions
 Column : MCI®GEL CHP-OD1M
 4.6mm I.D. x 150mm
 Eluent : 50vol%CH₃CN
 Flow rate : 0.50ml/min
 Column temp.: 40
 Detection : 254nm
 Sample : Triton X-100
 (polyoxyethylen Octylphey Ether)
 1% each 10μL





Fig. 5-25 ODS-1HU(ODS)

Conditions
 Column : MCI®GEL ODS-1HU
 4.6mm I.D. x 250mm
 Eluent : 50vol%CH₃CN
 Flow rate : 1.00ml/min
 Column temp.: 40
 Detection : 254nm
 Sample : Triton X-100
 (polyoxyethylen Octylphey Ether)
 1% each 10μL



MCI® GEL CRS10W (DLAA) 4.6mm I.D × 50mm

MCI® GEL CRS15W (LDAA) 4.6mm I.D × 50mm




CRS10W 4.6 × 50



CRS15W 4.6 × 50

Separation mechanism and Characteristic of MCI® GEL CRS series

Separation mechanism




Separation mechanism

MCI® GEL CRS10W and its companion product MCI® GEL CRS15W (an optical isomer of CRS10W) are based on a $3\mu\text{m}$ with 10nm mean pore diameter of silica gel coated with N,N-Dioctyl -L(or D)-alanine which is a novel optical activated ligand. The chiral resolution mechanism is a combination of ligand exchange and hydrophobic interaction. A copper sulfate aqueous solution is used as an eluent. Elution samples are directly detected at wave length of 254 nm because complex compound, composed of sample molecule and copper in the eluent, are object of detection. With the CRS10W, D-isomers generally elute in front of L-isomers while L-isomers elute ahead of D-isomers on the CRS15W. The hydrophobic interaction mechanism allows hydrophilic samples to elute faster than hydrophobic molecules. Long alkyl chain or aromatic compounds will elute late or require an organic solvent (CH_3CN or CH_3OH , max. of $15\text{v}/\text{v}\%$) to prevent adsorption onto the stationary phase.

Application of CRS10W

Fig. 6-1 DL-Isoleucine



Strong points

1. The CRS series columns separate over 20 D,L- α -Amino acids by only single column. The columns separate not only α -Amino acids but also α -Hydroxy carboxylic acids and derivative amino acids such as Acetylated amino acids.
2. The columns provide excellent resolution operated at room temperature.
3. The columns show high durability.

Application data of CRS10W

For all chromatograms, column temperature are room temperature and wave length are 254nm.

All eluents are CuSO₄ aqueous solution except for Fig. 6-9 and Fig. 6-10.

Fig. 6-2 Separation of amino acids mixture




Fig. 6-3 Separation of amino acids mixture




Fig. 6-4 Separation of DL-Ser.




Fig. 6-5 Separation of DL-aspartic acid




Fig. 6-6 Separation of DL-glutamic acid




Fig. 6-7 Separation of DL-histidine




Fig. 6-8 Separation of DL-lysine




Fig. 6-9 Separation of DL-phenylalanine




Fig. 6-10 Separation of DL-tryptophan




Fig. 6-11 Separation of DL-lactic acid





Fig. 6-12 Durability test

The sample was continuously injected 800 times for approximately 500 hrs. Changes of retention times and separation ability are not observed.



Application data of CRS10W

Fig. 6-13 Separation of DL- α -Phenylglycine

Conditions
 Column : MCI® GEL CRS10W 4.6mm I.D. x 50mm
 Eluent : 2mM CuSO₄/CH₃OH = 85/15
 Flow rate : 1.0mL/min
 Column temp.: 25
 Detection : 254nm
 Sample : 1. D- α -Phenylglycine
 2. L- α -Phenylglycine




Fig. 6-14 Separation of methionine and acetylmethionine

Conditions
 Column : MCI® GEL CRS10W 4.6mm I.D. x 50mm
 Eluent : 2mM CuSO₄/CH₃CN = 90/10
 Flow rate : 1.0mL/min
 Column temp.: 25
 Detection : 254nm
 Sample : 1. D-Met
 2. L-Met
 3. Acetyl-D-Met
 4. Acetyl-L-Met





Fig. 6-15 Separation of diaminopimelic acid

Conditions
 Column : MCI® GEL CRS10W 4.6mm I.D. x 50mm
 Eluent : 2mM CuSO₄
 Flow rate : 0.5mL/min
 Column temp.: 25
 Detection : 254nm
 Sample : 1. D,D-2,6-Diaminopimelic acid
 2. meso-2,6-Diaminopimelic acid
 3. L,L-2,6-Diaminopimelic acid



Column selection guide

1

Ion exchange columns

2

Packed columns for ion chromatography

3

Bioseparation columns

4

Reversed phase columns

5

Chiral separation columns

6

Chromatography media for preparative use

7

MCI® GEL columns

8

MCI® GEL chromatography media


9

Compounds Index

10


Application data of CRS10W

Fig. 6-16 Separation of 2-hydroxy carboxylic acids



Conditions
 Column : MCI® GEL CRS10W 4.6mm I.D. x 50mm
 Eluent : 2mM CuSO₄/CH₃CN = 90/10
 Flow rate : 1.0ml/min
 Column temp. : ambient
 Detection : 254nm
 Sample : 1,2 = D,L-Lactic acid
 3,4 = D,L-2-Hydroxy-n-butyric acid
 5,6 = D,L-α-Hydroxy-n-valeric acid
 7,8 = D,L-α-Hydroxy isocaproic acid

Fig. 6-17 Separation of 2-hydroxy carboxylic acids



Conditions
 Column : MCI® GEL CRS10W 4.6mm I.D. x 50mm
 Eluent : 2mM CuSO₄/CH₃CN = 90/10
 Flow rate : 1.0ml/min
 Column temp. : ambient
 Detection : 254nm
 Sample : 1,2 = D,L-p-Hydroxymandelic acid
 3,4 = D,L-m-Hydroxymandelic acid
 5,6 = D,L-Mandelic acid

Comparison data of CRS10W and CRS15W

Fig. 6-18 Separation of DL-alanine





Fig. 6-19 Analysis of a trace of L-lactic acid in 50 ppm D-lactic acid

The CRS15W is recommended for analysis of a trace of L-isomer in a principal D-isomer when the CRS10W does not provide an adequate chromatogram.



Examples of chromatographic conditions and data

	Amino acids	CuSO ₄ aq. soln. (mM)	Flow rate (ml/min)	Retention time; L-isomers (min)	Separation factor (α)	Separation rate (Rs)
1	Orn•HCl	0.1	0.2	6.8	1.26	< 1
2	Lys•HCl	0.1	0.2	7.7	1.45	< 1
3	Ala	0.1	0.5	11.0	1.39	1.4
4	His•HCl	0.1	0.5	10.5	1.63	1.7
5	Ser	0.1	0.5	10.1	1.25	1.0
6	Thr	0.1	0.5	11.3	1.29	1.3
7	Cit	0.5	0.5	10.4	1.75	2.3
8	Hyp	1.0	0.2	23.8	1.23	1.1
9	Pro	1.0	1.0	7.3	2.13	4.5
10	Val	1.0	1.0	8.9	2.04	5.0
11	Nval	1.0	1.0	11.5	2.07	4.7
12	Asp	2.0	0.5	13.2	1.18	0.8
13	Glu	2.0	1.0	16.2	1.54	2.3
14	Ileu(DL)	2.0	0.5	30.4	2.14	6.5
15	Ileu(allo)	2.0	0.5	21.9	1.97	6.0
16	Leu	2.0	1.0	14.6	1.97	4.6
17	Nleu	2.0	1.0	24.1	2.16	6.5
18	Met	2.0	1.0	10.3	1.64	2.6
19	Tyr	2.0	1.0	22.5	1.85	5.3
20	Eth	2.0	1.0	26.4	1.69	5.0
21	Phe	2.0	1.0	37.8	1.84	6.3

1. Column temperatures; ambient Detection; 254nm
2. These are example data and do not guarantee the column specifications.
3. Improved resolution or appropriate chromatogram can be obtained by further investigating chromatographic conditions.
4. For each amino acid in the table, D-isomer elutes ahead of L-isomer except for Hydroxyproline.

Notes

1. It will take hours for equilibrium between ligand of stationary phase and copper ion of eluent. Two to three hours of conditioning the column with the eluent is advised before sample injection or after changing concentration of CuSO₄ of eluent.
2. For acidic amino acids, higher CuSO₄ concentration of eluent provides better resolution.
3. For weakly retained hydrophilic amino acids, low flow rate (0.2-0.5 mL/min) yields better resolution.
4. Peak area may decrease with continuous injection of samples, when the concentration of amino acids in sample solution is much higher than that of CuSO₄ in the eluent.
5. Please be careful not to flow both water soluble organic solvents (CH₃CN, CH₃OH, etc) and non water soluble organic solvents (n-hexane, chloroform, etc) into the column. The column will be fatally damaged and will never separate optical isomers. Please be particularly careful if HPLC equipment is used together with RP mode and NP mode.
6. Please do not use acid or alkali solutions to adjust pH of eluent. And also do not use buffer solutions. These solutions may cause forming precipitation, hence cause of blockage of the column.
7. For strongly retained hydrophobic amino acids, addition of CH₃CN or CH₃OH in the eluent enables faster elution. The concentration of these organic solvents should be below 15 v/v%.
8. DOPA and other non-polar amino acids will be strongly adsorbed on the packing material and will cause contamination of the column.
9. Regeneration of contaminated column is difficult.

Chromatography media for preparative uses

For ion exchange chromatography
MCI® GEL PrepEx series

For hydrophobic interaction chromatography
MCI® GEL CQH_P series

For reversed phase chromatography
MCI® GEL CHP series

MCI® GEL chromatography media for preparative use for application of various chromatography modes are available.

MCI® GEL offers analytical use and preparative use products both are the same in chromatographic characteristics. Experimental results such as separating conditions of HPLC columns can be applied to preparative chromatography because both packing media have almost the same chromatographic characteristics. Our proposal of the developing process will accelerate your R&D and reduce the time.

The chromatography media for preparative use are based on rigid polymers and their particle size distributions are narrow compare to conventional packing media. These characteristics provide some chromatographic merit. Rigid polymers are allowed to use at high flow rates and narrow particle size distribution can provide high separation chromatography. Productivity will be extremely improved by using MCI® GEL chromatography media for preparative use.

MCI® GEL chromatography media for preparative use

Separation mode	Name	Base polymer	Functional group	Particle size μm	Main application	Equivalent analytical packing media
Ion exchange	PrepEx-DEAE PrepEx-SP	HMA	DEAE SP	30	Biomolecules	ProtEx-DEAE ProtEx-SP
Hydrophobic interaction	CQH3BP CQH3PP	HMA	Butyl Phenyl	30	Biomolecules	CQH3BS CQH3PS
Reversed phase	CHP55A CHP55Y CHP20A CHP20Y CHP2MGY	ST/DVB ST/DVB ST/DVB ST/DVB MA	— — — — —	18 30 18 30 30	Organic compounds	CHP5C,CHP10M CHP5C,CHP10M CHP5C,CHP10M CHP5C,CHP10M CHP2MG,CHP2MGM

HMA ; polyhydroxymethacrylate ST/DVB ; styrene/divinylbenzene MA ; polymethacrylate

Separation mechanism and characteristics of PrepEx series

MCI® GEL PrepEx series is ion exchange chromatography media for preparative chromatography whose average particle size is 30 μ m.


Surface characteristic of MCI® GEL PrepEx series packing media is the same as analytical use MCI® GEL PrepEx series. Two types of ion exchanger, weak anion exchanger PrepEx-DEAE and strong cation exchange PrepEx-SP are available.

Because ProtEx series and PrepEx series have the same chromatographic characteristics, experimental results of separating conditions of ProtEx series can be applied to PrepEx series.

High flow rates (over 300cm/hr.) are able to be applied to the chromatography media because they are based on rigid polymer. Productivity will be extremely improved.


Application data of PrepEx series

Fig. 7-1 Proteins



Conditions
 Column : ProtEx-SP, 8.0mm I.D. x 100mL
 PrepEx-SP, 8.0mm I.D. x 100mL
 Eluent : A : 20mM Phosphate buffer (pH7.0)
 B : A+0.5M NaCl
 Gradient : A B, 60min, linear
 Flow rate : 0.5mL/min
 Column temp. : 25
 Detection : 280nm
 Sample : 1. Ribnuclease A 150 μ g
 2. Cytochrome C 60 μ g
 3. Lysozyme 60 μ g

Fig. 7-2 Cow milk acid whey



Conditions
 Column : PrepEx-SP, 8.0mm I.D. x 100mL
 Eluent : A : 20mM Phosphate buffer (pH7.11)
 B : A+1.0M NaCl (pH7.03)
 : 0-50min, Flow rate 1.0mL/min
 50-110min, Eluent A, Flow rate 1.0mL/min
 110-160min, linear gradient 0-50% B
 Flow rate 0.5mL/min
 160-190min, linear gradient 50-100% B
 Flow rate 0.5mL/min
 190-210min, Eluent B, Flow rate 0.5mL/min
 Column temp. : 25
 Detection : 280nm
 Sample : Cow milk acid whey (pH4.7) 50mL

Application data of PrepEx series

Fig. 7-3 Proteins






Fig. 7-4 Egg white



Chromatograms of analytical column



Separation mechanism and characteristics of CQH_P series

MCI® GEL CQH_P series is preparative chromatography media for separating biomolecules in the basis of differences of their hydrophobic properties. Average particle size is 30µm.


Two types of hydrophobic ligands, butyl type and phenyl type, these are differences in hydrophobic properties, are available. An appropriate packing media can be selected in accordance with properties of injection samples.

Because chromatographic characteristics of CQH_S series and CQH_P series are the same, experimental results of separating conditions of CQH_S series can be applied to CQH_P series.


Application data of CQH_P series

Fig. 7-5 Proteins

10µm CQH 3BS



30µm CQH 3BP




Conditions


Column	: CQH3BS, 7.5mm I.D. x 75mmL	
	: CQH3BP, 7.5mm I.D. x 75mmL	
Eluent	: A : B+1.7M (NH ₄) ₂ SO ₄	
	: B : 0.1M Phosphate buffer (pH6.8)	
Gradient	: A B, 30min, linear	
Flow rate	: 1.0mL/min	
Column temp.	: 25	
Detection	: 280nm,	
Sample	: 1. Ribonuclease A	112µg
	: 2. Transferrin	154µg
	: 3. α-Chymotrypsinogen A	60µg

Fig. 7-6 Proteins

10µm CQH 3PS



30µm CQH 3PP



Conditions

Column	: CQH3PS, 7.5mm I.D. x 75mmL	
	: CQH3PP, 7.5mm I.D. x 75mmL	
Eluent	: A : B+1.7M (NH ₄) ₂ SO ₄	
	: B : 0.1M Phosphate buffer (pH6.8)	
Gradient	: A B, 30min, linear	
Flow rate	: 1.0mL/min	
Column temp.	: 25	
Detection	: 280nm,	
Sample	: 1. Ribonuclease A	112µg
	: 2. Transferrin	154µg
	: 3. α-Chymotrypsinogen A	60µg

Separation mechanism and characteristics of CHP series

MCI® GEL CHP series is reversed phase chromatography media of porous type polymers.

Because polymeric media are chemically stable, wide pH range, from acidic to alkaline eluents are able to be applied to MCI® GEL CHP series.

MCI® GEL CHP55 series and CHP20 series are both ST/DVB polymers, but they differences in porosity. Pore size of CHP20 series is fairly larger than that of CHP55 series. Appropriate packing media can be selected in accordance with molecular size of injection samples.

Application of CHP series

Fig. 7-7 Phthalic acid esters




Fig. 7-8 Penicillin antibiotics





Fig. 7-9 Extract of green tea leaves



Application of CHP series

Fig.7-10 Senna pulv. extract

Conditions

	Chromatogram A	Chromatogram B	Chromatogram C
Column	MCI® GEL CHP5C 4.6mm I.D. × 250mm	MCI® GEL CHP20A 10.0mm I.D. × 250mm	MCI® GEL CHP20Y 10.0mm I.D. × 250mm
Eluent	: *1	: *1	: *1
Flow rate	: 0.5ml/min.	: 2.4ml/min.	: 2.4ml/min.
Detection	: 270nm	: 270nm	: 270 nm
Sample	: *2 10 µL	: *2 80 µL	: *2 80 µL

*1: CH₃OH + 1% Acetic acid = 60 + 40 (vol.)

*2: Extract of senna pulv.





Fig. 7-11 Elution profile of senna pulv. extract separated on MCI GEL CHP20Y

Conditions

Column	MCI® GEL CHP20Y 32mm I.D. × 490mm
Eluent	: CH ₃ OH + 1% Acetic acid = 60 + 40 (vol.)
Flow rate	: 7.88ml/min.
Detection	: 270 nm
Sample	: Extract of senna pulv., partially purified by Diaion HP20 injected amount : 39.4 ml



Application of CHP series

Fig. 7-12 Elution profile of gardenia fructus extract separated on MCI GEL CHP20Y





Fig. 7-13 Rice brain oil



8

MCI® GEL

MCI® GEL columns

main column			guard/pre-column		
Code No.	Name	Column dimensions (mm)	Code No.	Name	Column dimensions (mm)
Cation exchange resin 001901	For amino acids CK10U	6 × 120	003321	AFR2-PC	6 × 50
Cation exchange resin 000101	For oligo saccharides CK02A	20 × 250	000111	CK02AG	8 × 10
000102	CK02AS	20 × 250	000112	CK02ASG	8 × 10
000301	CK04S	10 × 200	001711 000311	CK10SG CK04SG	6 × 50 8 × 10
000302	CK04SS	10 × 200	001711 000312	CK10SG CK04SG	6 × 50 8 × 10
Cation exchange resin 000901	For sugars CK08S	8 × 500	000911	CK08SG	6 × 50
001001	CK08E	8 × 300	001011	CK08EG	6 × 50
001002	CK08EC	8 × 300	001012	CK08ECG	6 × 50
001003	CK08ES	8 × 300	001013	CK08ESG	6 × 50
Cation exchange resin 001005	For carboxylic acids CK08EH	8 × 300	001015	CK08EHG	6 × 50
Anion exchange resin 011101	For carboxylic acids · sugars CA08F	4.6 × 250	011111	CA08FG	4 × 10
011901	CDR10	4.6 × 250	011911	CDR10G	4 × 10
Ion chromatography 003401	For cations SCK01	6 × 50	001721	SCK-PC	6 × 50
003404	SCK01	4.6 × 150			
Ion chromatography 013302	For anions SCA04	4.6 × 150	013312	SCA-04G	4.6 × 30
			013021	SCA-PC	8 × 10
Bioseparation 021301	For size exclusion Hydrophilic polymer CQP06	7.5 × 600	021311	CQP06G	4 × 50
021401	CQP10	7.5 × 600	021411	CQP10G	4 × 50
021501	CQP30	7.5 × 600	021511	CQP30G	4 × 50

main column			guard/pre-column		
Code No.	Name	Column dimensions (mm)	Code No.	Name	Column dimensions (mm)
Bioseparation For anion exchange chromatography 014603 ProtEx-DEAE 4.6 x 50 014604 ProtEx-DEAE 7.5 x 100 012601 CQA31S 7.5 x 75 013001 CQA35S 7.5 x 75					
Bioseparation For cation exchange chromatography 003703 ProtEx-SP 4.6 x 50 003704 ProtEx-SP 7.5 x 100 003601 CQK30S 7.5 x 75 003801 CQK31S 7.5 x 75					
Bioseparation For hydrophobic interaction chromatography 021601 CQH3BS 7.5 x 75 021701 CQH3ES 7.5 x 75 021801 CQH3PS 7.5 x 75					
Reversed phase columns 040101 CHP10M 4.6 x 150 020401 CHP5C 4.6 x 250 040201 CHP2MG 4.6 x 150 020201 CHP2MG 4.6 x 250					
Chiral separation columns 021901 CRS10W 4.6 x 50 022001 CRS15W 4.6 x 50					



Characteristics

1 . Excellent performance

Sphere packing and sharp particle size distribution provide high performance.

2 . Persistence and highest quality

Produced with Mitsubishi Chemical 's excellent technology, experience and under strict quality control.

3 . Wide range of product line

MCI® GEL bulk packings include ion exchange resins (cation and anion), non-functionalized polymer used for reversed phase chromatography and other varieties of products. Also MCI® GEL offers mean particle size of 4 µm to approximately 300 µm packing materials, this means that MCI® GEL products are applied to analysis use and preparative use.

4 . Abundant experience

Mitsubishi Chemical has been supplying packing materials for more than 40 years.

MCI® GEL bulk packings

For application of each packing material, please refer application data of the same packed column.

For Ion chromatography [SCK series, SCA series]

Code No.	Name	Packing size (g)	Base material	Functional group	Counter ion	Mean particle size (µm)	Cross linkage (%)	Ion exchange capacity (µeq / g)	Remarks
1-034-01	SCK01	5	ST/DVB	RSO ₃ ⁻	H ⁺	11	-	25	Cation analysis
1-034-02	SCK01	10							
1-133-01	SCA04	5	HMA	QA	Cl ⁻	5	-	30	Anion analysis
1-133-02	SCA04	10							

Abbreviations; ST/DVB = Styrene-divinylbenzene copolymer HMA = Polyhydroxymethacrylate QA = Quaternary ammonium

Cation exchange resins [CK series, AFR series]

Code No.	Name	Packing size (g)	Base material	Functional group	Counter ion	Mean particle size (μm)	Cross linkage (%)	Ion exchange capacity (meq / ml)	Remarks
1-001-01	CK02A	10	ST/DVB	RSO ₃ ⁻	Na ⁺	20	2	> 0.5	Oligosaccharides
1-003-01	CK04S	10	ST/DVB	RSO ₃ ⁻	Na ⁺	11	4	> 0.8	Oligosaccharides
1-003-02	CK04S	25							
1-003-03	CK04S	50							
1-004-01	CK06S	10	ST/DVB	RSO ₃ ⁻	Na ⁺	11	6	> 1.5	Oligosaccharides
1-004-02	CK06S	25							
1-004-03	CK06S	50							
1-009-01	CK08S	10	ST/DVB	RSO ₃ ⁻	Na ⁺	11	8	> 1.9	Sugars, Carboxylic acids
1-009-02	CK08S	25							
1-009-03	CK08S	50							
1-010-01	CK08E	10	ST/DVB	RSO ₃ ⁻	Na ⁺	9	8	> 1.9	Sugars, Carboxylic acids
1-010-02	CK08E	25							
1-010-03	CK08E	50							
1-013-01	CK08Y	50	ST/DVB	RSO ₃ ⁻	Na ⁺	25	8	> 1.9	Sugars, Carboxylic acids
1-013-02	CK08Y	300							
1-014-01	CK08P	100ml	ST/DVB	RSO ₃ ⁻	H ⁺	75/150	8	> 1.9	Sugars, Carboxylic acids
1-017-01	CK10S	10	ST/DVB	RSO ₃ ⁻	Na ⁺	11	10	> 2.0	Carboxylic acids, Amino acids
1-017-02	CK10S	25							
1-017-03	CK10S	50							
1-018-01	CK10F	5	ST/DVB	RSO ₃ ⁻	Na ⁺	7	10	> 2.0	Amino acids
1-018-02	CK10F	10							
1-019-01	CK10U	3	ST/DVB	RSO ₃ ⁻	Na ⁺	5	10	> 2.0	Amino acids
1-019-03	CK10U	5						> 2.0	
1-019-04	CK10U	10							
1-020-05	CK10M	5	ST/DVB	RSO ₃ ⁻	Na ⁺	4	10	> 2.0	Amino acids
1-020-06	CK10M	3							
1-021-01	CK10Y	50	ST/DVB	RSO ₃ ⁻	Na ⁺	25	10	> 1.9	Amino acids
1-033-01	AFR2	5	ST/DVB	RSO ₃ ⁻	H ⁺	25	-	> 1.9	Ammonia trap

Abbreviation; ST/DVB = Styrene-divinylbenzene copolymer

Anion exchange resins[CA series, CDR series]

Code No.	Name	Packing size (g)	Base material	Functional group	Counter ion	Mean particle size (μm)	Cross linkage (%)	Ion exchange capacity (meq / mL)	Remarks
1-104-01	CA06S	10	ST/DVB	QA	Cl ⁻	11	6	> 1.2	Sugars, Carboxylic acids
1-104-02	CA06S	25							
1-104-03	CA06S	50							
1-109-01	CA08S	10	ST/DVB	QA	Cl ⁻	11	8	> 1.2	Sugars, Carboxylic acids
1-109-02	CA08S	25							
1-109-03	CA08S	50							
1-111-01	CA08F	5	ST/DVB	QA	Cl ⁻	7	8	> 1.2	Sugars, Carboxylic acids
1-111-02	CA08F	10							
1-112-01	CA08Y	50	ST/DVB	QA	Cl ⁻	25	8	> 1.2	Sugars, Carboxylic acids
1-113-01	CA08P	100mL	ST/DVB	QA	Cl ⁻	75 / 150	8	> 1.3	Sugars, Carboxylic acids
1-116-01	CA10S	10	ST/DVB	QA	Cl ⁻	11	10	> 1.2	Sugars, Carboxylic acids
1-116-02	CA10S	25							
1-116-03	CA10S	50							
1-119-01	CDR10	7	ST/DVB	QA	Cl ⁻	7	-	> 0.3	Nucleic acids, Sugars
1-119-02	CDR10	14							

Abbreviations; ST/DVB = styrene-divinyl benzene copolymer QA; Quaternary ammonium

Ion exchange resins for bioseparation [CQK series, CQA series]

Code No.	Name	Packing size (g)	Base material	Functional group	Counter ion	Mean particle size (μm)	Pore size (nm)	pH range	Remarks
1-036-01	CQK30S	10	HMA	SP	Na ⁺	10	60	1 ~ 13	Proteins
1-036-02	CQK30S	25							
1-036-03	CQK30S	50							
1-037-01	CQK30P	100mL							
1-038-01	CQK31S	10	HMA	CM	Na ⁺	10	60	> 4	Proteins
1-038-02	CQK31S	25							
1-038-03	CQK31S	50							
1-039-01	CQK31P	100mL							
1-126-01	CQA31S	10	HMA	DEAE	Cl ⁻	10	60	< 11	Proteins
1-126-02	CQA31S	25							
1-126-03	CQA31S	50							
1-127-01	CQA31P	100mL							
1-130-01	CQA35S	10	HMA	QA	Cl ⁻	10	60	2 ~ 12	Proteins
1-130-02	CQA35S	25							
1-130-03	CQA35S	50							
1-131-01	CQA35P	100mL							
1-501-01	PrepEx-DEAE	25mL	HMA	DEAE	Cl ⁻	30	100	< 11	Proteins
1-501-02	PrepEx-DEAE	100mL							
1-501-03	PrepEx-DEAE	1000mL							
1-502-01	PrepEx-SP	25mL	HMA	SP	Na ⁺	30	100	1 ~ 13	Proteins
1-502-02	PrepEx-SP	100mL							
1-502-03	PrepEx-SP	1000mL							

Abbreviations; HMA = Polyhydroxymethacrylate SP = Sulfopropyl CM = Carboxymethyl DEAE = Diethylaminoethyl
QA = Quaternary ammonium

Hydrophobic interaction chromatography resins for bioseparation [CQH series]

Code No.	Name	Packing size (g)	Base material	Functional group	Counter ion	Mean particle size (µm)	Pore size (nm)	Ion exchange capacity (meq / ml)	Remarks
1-216-01	CQH3BS	10	HMA	Butyl	-	10	60	-	Proteins
1-216-02	CQH3BS	25							
1-216-03	CQH3BS	50							
1-217-01	CQH3ES	10	HMA	Ether	-	10	60	-	Proteins
1-217-02	CQH3ES	25							
1-217-03	CQH3ES	50							
1-218-01	CQH3PS	10	HMA	Phenyl	-	10	60	-	Proteins
1-218-02	CQH3PS	25							
1-218-03	CQH3PS	50							

Abbreviation; HMA = Polyhydroxymethacrylate

Size exclusion chromatography resins for bioseparation [CQP series]

Code No.	Name	Packing size (g)	Base material	Functional group	Counter ion	Mean particle size (µm)	Pore size (nm)	Exclusion limit	Remarks
1-213-01	CQP06	10	HMA	-	-	10	12	1×10^3	Water soluble polymer
1-213-02	CQP06	25							
1-213-03	CQP06	50							
1-214-01	CQP10	10	HMA	-	-	10	20	1×10^4	Water soluble polymer
1-214-02	CQP10	25							
1-214-03	CQP10	50							
1-215-01	CQP30	10	HMA	-	-	10	60	1×10^6	Water soluble polymer
1-215-02	CQP30	25							
1-215-03	CQP30	50							
1-222-01	CQP30P	100ml	HMA	-	-	30	60	1×10^6	

Abbreviation; HMA = Polyhydroxymethacrylate

For reversed phase chromatography [CHP series]

Code No.	Name	Packing size (g)	Base material	Functional group	Counter ion	Mean particle size (μm)	Pore size (nm)	pH range	Remarks
1-202-01	CHP2MG	10	MA	-	-	10	25	2 ~ 12	Reversed phase chromatography
1-202-02	CHP2MG	1000ml							
1-204-01	CHP5C	10	ST/DVB	-	-	10	25	1 ~ 13	Reversed phase chromatography
1-201-01	CHP20P	100	ST/DVB	-	-	37 / 75	45	1 ~ 13	Reversed phase chromatography
1-302-01	CHP20P	100ml	ST/DVB	-	-	75 / 150	45	1 ~ 13	Reversed phase chromatography
1-302-02	CHP20P	500ml							
1-302-03	CHP20P	1000ml							
1-304-01	CHP55A	25ml	ST/DVB	-	-	18	25	1 ~ 13	Reversed phase chromatography
1-304-02	CHP55A	100ml							
1-304-03	CHP55A	1000ml							
1-303-01	CHP55Y	25ml	ST/DVB	-	-	30	25	1 ~ 13	Reversed phase chromatography
1-303-02	CHP55Y	100ml							
1-303-03	CHP55Y	1000ml							
1-307-01	CHP20A	25ml	ST/DVB	-	-	18	45	1 ~ 13	Reversed phase chromatography
1-307-02	CHP20A	100ml							
1-307-03	CHP20A	1000ml							
1-305-01	CHP20Y	25ml	ST/DVB	-	-	30	45	1 ~ 13	Reversed phase chromatography
1-305-02	CHP20Y	100ml							
1-305-03	CHP20Y	1000ml							
1-306-01	CHP2MGY	25ml	MA	-	-	30	25	2 ~ 12	Reversed phase chromatography
1-306-02	CHP2MGY	100ml							
1-306-03	CHP2MGY	1000ml							

Abbreviations; MA = Polymethacrylate ST/DVB = Styrene-divinylbenzene copolymer

Synthetic adsorbents for enrichment organic compounds in environmental water

Code No.	Name	Packing size (ml)	Base material	Functional group	Counter ion	Mean particle size (μm)	Specific surface area (m ² / g)	Ion exchange capacity (meq / ml)	Remarks
1-219-01	CSP800	50	ST/DVB	-	-	75 / 150	600	-	Synthetic adsorbents for non-ionic substances
1-132-01	CHPA25	20	ST/DVB	QA	Cl ⁻	150 / 300	20	> 2.0	Synthetic adsorbents for anionic substances

Abbreviations; ST/DVB = Styrene-divinylbenzene copolymer QA = Quaternary ammonium

MCI® GEL CSP800 and MCI® GEL CHPA25 are used for enrichment traces of organic compounds in environmental water with high concentration ratio and high recovery, are recommended for sample preparation for mutagenicity study and GC-MS analysis. The CSP800 is for non-ionic substances such as trichloroethylene. The CHPA25 is for anionic substances such as humin. It is advised combined use these adsorbents.

Cation exchange resin

Name	Packing size	Functional group	Mean particle size (μm)	Cross linkage (%)	Ion exchange capacity (meq / ml)	Remarks
CSP800	50g	RSO ₃ ⁻	25	8	> 1.9	Sugars, Carboxylic acids
CK08P	100ml	RSO ₃ ⁻	75/ 150	8	> 1.9	Sugars, Carboxylic acids

Anion exchange resin

Name	Packing size	Functional group	Mean particle size (μm)	Cross linkage (%)	Ion exchange capacity (meq / ml)	Remarks
CA08Y	50g	QA	25	8	> 1.2	Sugars, Carboxylic acids
CA08P	100ml	QA	75/ 150	8	> 1.3	Sugars, Carboxylic acids

Chelating resin for ion exchange chromatography

Name	Packing size	Functional group	Mean particle size (μm)	Cross linkage (%)	Ion exchange capacity (meq / ml)	Remarks
CHL10P	100g	Iminodiacetic acid	75/ 150	-	> 1.5	Metal
CHL 20P	100g	Polyamine	75/ 150	-	> 1.8	Metal
CLB10P	100g	Glucamine	75/ 150	-	> 1.0	Boron

Reversed phase packing

Name	Packing size	Functional group	Mean particle size (μm)	Pore size()	pH range	Remarks
CHP5C	10g	-	10	250	Whole range	Reversed phase
CHP2MG	10g	-	10	250	2 ~ 12	Reversed phase
CHP20P	100g	-	37/ 75	450	Whole range	Reversed phase
CHP20P	100ml	-	75/ 150	450	Whole range	Reversed phase
CSP207P	100ml	-	75/ 150	250	Whole range	Reversed phase
CHP2MGP	100ml	-	75/ 150	500	2 ~ 12	Reversed phase
CHPOD1Y	100ml	-	30	250	Whole range	Reversed phase
CHP50P			75/ 150	<14000		Reversed phase
CHP75P			75/ 150	14000		Reversed phase

Protein exclusion limit MW

We have an assortment of MCI® GEL cation ion exchange resins, CK series, anion exchange resins, CA series, and chelating resins, CHL series, as solid phase adsorbents for the pretreatment in analyzing rare earth metals. We can also provide solid phase adsorbents with various micro-pore sizes and hydrophobicity, i.e. chemical structures, for R&D of new pharmaceuticals.

CHP50P and CHP75P with the controlled micro-pore size, in particular, have a distinctive advantage not to adsorb high molecular weight proteins but to adsorb only low molecular weight organic compounds.

	Compound	Classification	MCI® GEL Column	Figure	Page
1	Acetic acid	Carboxylic acid	CK08EH	2-12	11
2	Acetic acid	Carboxylic acid	CK08EH	2-17	12
3	Acetic acid	Carboxylic acid	CK08EH	2-18	12
4	Acetic acid	Carboxylic acid	CA08F	2-32	17
5	Acetic acid	Carboxylic acid	CA08F	2-33	18
6	Acetic acid	Carboxylic acid	CA08F	2-34	18
7	Acetylgalactosamine	Amino sugar	CK08EH	2-13	11
8	Acetylglucosamine	Amino sugar	CK08EH	2-13	11
9	Acetyl-D-Met.	Acetyl-D-amino acid	CRS10W	6-14	54
10	Acetyl-L-Met.	Acetyl-L-amino acid	CRS10W	6-14	54
11	5'-ADP	Nucleotide	CDR10	2-35	19
12	2'-AMP	Nucleotide	CDR10	2-35	19
13	2'-AMP	Nucleotide	CDR10	2-36	19
14	3'-AMP	Nucleotide	CDR10	2-36	19
15	5'-AMP	Nucleotide	CDR10	2-35	19
16	5'-AMP	Nucleotide	CDR10	2-36	19
17	5'-ATP	Nucleotide	CDR10	2-35	19
18	Adenine	Nucleic base	CDR10	2-35	19
19	Adenosine	Nucleoside	SCK01	3-6	22
20	Adonitol	Sugar alcohol	CK08EC	2-4	9
21	Alanine	Amino acid	CK10U	2-1	7
22	-Alanine	Amino acid	CK10U	2-2	8
23	D-Alanine	D-Amino acid	CRS10W/CRS15W	6-18	56
24	L-Alanine	L-Amino acid	CRS10W/CRS15W	6-18	56
25	-Aminobutyric acid	Amino acid	CK10U	2-3	8
26	6-Aminopenicillanic acid	Penicillin antibiotic	CHP55A	7-8	62
27	Ammonia	Ammonia	SCK01	3-2	22
28	Ammonium ion	Cation	SCK01	3-1	22
29	Ammonium ion	Cation	SCK01	3-2	22
30	Ammonium ion	Ammonia	SCK01	3-3	22
31	Angiotensin II	Peptide	CHP2MGM	5-12	47
32	Antipyrine	Ingredients of medicine	CHP2MGM	5-10	46
33	Arginine	Amino acid	CK10U	2-1	7
34	Aspirin	Ingredients of medicine	CHP2MGM	5-10	46
35	Aspartic acid	Amino acid	CK10U	2-1	7
36	D-Aspartic acid	D-Amino acid	CRS10W	6-5	53
37	L-Aspartic acid	L-Amino acid	CRS10W	6-5	53
38	Bacitracin	Peptide	CQH3PS	4-26	38
39	Barbital	Anticonvulsant	CHP5C	5-17	49
40	Barium ion	Cation	SCK01	3-7	23
41	Benzoic acid	Aromatic carboxylic acids	CHP2MG	5-19	50
42	Bovine Serum Albumin	Protein	ProtEx-DEAE	4-9	32
43	Bromide	Anion	SCA04	3-12	24
44	Bromide	Anion	SCA04	3-13	24
45	n-Butyl alcohol	Alcohol	CK08EH	2-14	11
46	sec-Butyl alcohol	Alcohol	CK08EH	2-14	11
47	Cadmium ion	Cation	SCK01	3-8	23
48	Caffeine	Purine alkaloid	CHP10M	5-4	43
49	Caffeine	Ingredients of medicine	CHP2MGM	5-10	46
50	Caffeine	Purine alkaloid	CHP55A	7-9	62
51	Calcium ion	Cation	SCK01	3-7	23
52	Calcium ion	Cation	SCK01	3-8	23
53	Calcium ion	Cation	SCK01	3-9	23
54	Calcium ion	Cation	SCK01	3-10	23
55	Calcium ion	Cation	SCK01	3-11	23
56	Carbamazepine	Anticonvulsant	CHP5C	5-17	49
57	Carbonate	Anion	SCA04	3-15	25
58	4-Carboxybenzaldehyde	Aromatic carboxylic acids	CHP2MG	5-19	50
59	Catechin	Catechin	CHP55A	7-9	61
60	Cellobiose	Disaccharide	CA08F	2-31	17

	Compound	Classification	MCI® GEL Column	Figure	Page
61	Cellobiose	Disaccharide	CDR10	2-37	20
62	Cesium ion	Cation	SCK01	3-1	22
63	Chloride	Anion	SCA04	3-12	24
64	Chloride	Anion	SCA04	3-13	24
65	Chloride	Anion	SCA04	3-14	25
66	Chloride	Anion	SCA04	3-15	25
67	Chloride	Anion	SCA04	3-17	26
68	Chloroacetic acid	Carboxylic acid	CK08EH	2-17	12
69	Chloroacetic acid	Carboxylic acid	CK08EH	2-18	12
70	Cholic acid	Bile acid	CHP10M	5-7	44
71	-Chymotrypsinogen A	Protein	ProtEx-SP	4-10	32
72	-Chymotrypsinogen A	Protein	CQK31S	4-21	36
73	-Chymotrypsinogen A	Protein	CQK30S	4-21	36
74	-Chymotrypsinogen A	Protein	CHP2MGM	5-13	47
75	-Chymotrypsinogen A	Protein	CHP2MG	5-20	50
76	-Chymotrypsinogen A	Protein	CQH3BP	7-5	61
77	-Chymotrypsinogen A	Protein	CQH3PP	7-6	61
78	2'-CMP	Nucleotide	CDR10	2-36	19
79	3'-CMP	Nucleotide	CDR10	2-36	19
80	5'-CMP	Nucleotide	CDR10	2-35	19
81	5'-CMP	Nucleotide	CDR10	2-36	19
82	Cinchonine	Cinchona alkaloid	CHP10M	5-5	43
83	Citric acid	Carboxylic acid	CK08EH	2-12	11
84	Citric acid	Carboxylic acid	CA08F	2-32	17
85	Citric acid	Carboxylic acid	CA08F	2-33	18
86	Cobalt ion	Cation	SCK01	3-8	23
87	Colibacillus extract	Protein	CQH3ES	4-24	38
88	Colibacillus extract	Protein	CQH3PS	4-25	38
89	Collagenase	Protein	ProtEx-DEAE	4-15	34
90	Conalbumin	Protein	ProtEx-DEAE	4-5	31
91	Conalbumin	Protein	PrepEx-DEAE	7-3	60
92	Cow milk acid whey	Cow milk acid whey	PrepEx-SP	7-2	59
93	Crocin	Crude drug	CHP20Y	7-12	64
94	Corticosterone	Adrenal cortex hormone	CHP10M	5-9	45
95	3',5'-Cyclic AMP	Nucleotide	CDR10	2-35	19
96	Cystine	Amino acid	CK10U	2-1	7
97	Cytidine	Nucleoside	SCK01	3-6	22
98	Cytochrome C	Protein	CQP30	4-3	29
99	Cytochrome C	Protein	ProtEx-SP	4-10	32
100	Cytochrome C	Protein	CQK30S	4-20	36
101	Cytochrome C	Protein	CQK31S	4-20	36
102	Cytochrome C	Protein	CHP2MGM	5-13	47
103	Cytochrome C	Protein	CHP2MG	5-20	50
104	Cytochrome C	Protein	PrepEx-SP	7-1	59
105	Cytosine	Nucleic base	CDR10	2-35	19
106	Deoxycholic acid	Bile acid	CHP10M	5-7	44
107	11-Deoxycortisol	Adrenal cortex hormone	CHP10M	5-9	45
108	Deoxyribose	Deoxysugar	CA08F	2-31	17
109	Deoxyribose	Deoxysugar	CDR10	2-37	20
110	D,D-2,6-Diaminopimelic acid	D,D-Diamino carboxylic acid	CRS10W	6-15	54
111	L,L-2,6-Diaminopimelic acid	L,L-Diamino carboxylic acid	CRS10W	6-15	54
112	meso-2,6-Diaminopimelic acid	meso-Diamino carboxylic acid	CRS10W	6-15	54
113	Diazepam	Anticonvulsant	CHP5C	5-17	49
114	Dibutyl phthalate	Phthalic acid ester	CHP10M	5-3	42
115	Dibutyl phthalate	Phthalic acid ester	CHP55A	7-7	62
116	Dichloroacetic acid	Carboxylic acid	CK08EH	2-17	12
117	Dichloroacetic acid	Carboxylic acid	CK08EH	2-18	12
118	Diethylene glycol	Polyalcohol	CK08EH	2-16	12
119	Dimethylamine	Amine	SCK01	3-2	22
120	4-Dimethylaminoantipyrine	Ingredients of medicine	CHP2MGM	5-10	46
121	Diethyl phthalate	Phthalic acid ester	CHP10M	5-3	42
122	Dimethyl phthalate	Phthalic acid ester	CHP55A	7-7	62
123	Dimethyl phthalate	Phthalic acid ester	CHP10M	5-3	42
124	Dimethyl phthalate	Phthalic acid ester	CHP55A	7-7	62
125	Dopamine	Catecholamine	CHP10M	5-2	42

Column selection guide

1

Ion exchange columns

2

Packed columns for ion chromatography

3

Bioseparation columns

4

Reversed-phase columns

5

Chiral separation columns

6

Chromatography media for preparative use

7

MCI® GEL columns

8

MCI® GEL chromatography media

9

Compounds Index

10

	Compound	Classification	MCI®GEL Column	Figure	Page
126	Egg white	Egg white	PrepEx-DEAE	7-4	60
127	Epicatechin	Catechin	CHP55A	7-9	62
128	Epicatechin gallate	Catechin	CHP55A	7-9	62
129	Epigallocatechin	Catechin	CHP55A	7-9	62
130	Epigallocatechin gallate	Catechin	CHP55A	7-9	62
131	Epinephrine	Catecholamine	CHP10M	5-2	42
132	4-Epitetracycline	Antibiotics	CHP5C	5-16	49
133	Erythritol	Sugar alcohol	CK08EC	2-4	9
134	Erythritol	Sugar alcohol	CK08EC	2-11	10
135	D-Ethionine	D-Amino acid	CRS10W	6-2	53
136	L-Ethionine	L-Amino acid	CRS10W	6-2	53
137	Ethyl alcohol	Alcohol	CK08EC	2-11	10
138	Ethyl alcohol	Alcohol	CK08EH	2-14	11
139	Ethyl alcohol	Alcohol	CK08EH	2-15	11
140	Ethylene glycol	Polyalcohol	CK08EH	2-15	11
141	Ethylene glycol	Polyalcohol	CK08EH	2-16	12
142	Extract of green tea leaves	Catechins	CHP55A	7-9	62
143	Ferritin	Protein	CQP30	4-3	29
144	Fluoride	Anion	SCA04	3-12	24
145	Fluoride	Anion	SCA04	3-13	24
146	Formic acid	Carboxylic acid	CK08EH	2-12	11
147	Formic acid	Carboxylic acid	CA08F	2-32	17
148	Fructose	Sugar	CK08EC	2-4	9
149	Fructose	Sugar	CK08EC	2-5	9
150	Fructose	Sugar	CK08EC	2-7	10
151	Fructose	Sugar	CK08EC	2-8	10
152	Fructose	Sugar	CK08EC	2-9	10
153	Fructose	Sugar	CK08EC	2-10	10
154	Fructose	Sugar	CK08EC	2-11	10
155	Fructose	Sugar	CK04S	2-28	16
156	Fructose	Sugar	CK04S	2-29	16
157	Fructose	Sugar	CK04S	2-30	16
158	Fructose	Sugar	CA08F	2-31	17
159	Fructose	Sugar	CDR10	2-37	20
160	Fructo-oligosaccharide	Fructo-oligosaccharide	CK04S	2-30	16
161	Fucose	Sugar	CA08F	2-31	17
162	Galactose	Sugar	CK08EC	2-6	9
163	Galactose	Sugar	CA08F	2-31	17
164	Galactose	Sugar	CDR10	2-37	20
165	Gallocatechin	Catechin	CHP55A	7-9	62
166	Gentiobiose	Disaccharide	CK08EC	2-4	9
167	Gluconic acid	Carboxylic acid	CA08F	2-33	18
168	Gluconic acid	Carboxylic acid	CA08F	2-34	18
169	Gluconic acid	Carboxylic acid	CQP06	4-4	29
170	Glucose	Sugar	CK08EC	2-4	9
171	Glucose	Sugar	CK08EC	2-5	9
172	Glucose	Sugar	CK08EC	2-7	10
173	Glucose	Sugar	CK08EC	2-8	10
174	Glucose	Sugar	CK08EC	2-9	10
175	Glucose	Sugar	CK08EC	2-10	10
176	Glucose	Sugar	CK08EC	2-11	10
177	Glucose	Sugar	CK08EH	2-13	11
178	Glucose	Sugar	CK04S	2-28	16
179	Glucose	Sugar	CK04S	2-29	16
180	Glucose	Sugar	CK04S	2-30	16
181	Glucose	Sugar	CA08F	2-31	17
182	Glucose	Sugar	CDR10	2-37	20
183	Glucose	Sugar	CQP06	4-4	29
184	Glutamic acid	Amino acid	CK10U	2-1	7
185	D-Glutamic acid	D-Amino acid	CRS10W	6-6	53
186	L-Glutamic acid	L-Amino acid	CRS10W	6-6	53
187	Glycerol	Polyalcohol	CK08EC	2-11	10
188	Glycerol	Polyalcohol	CK08EH	2-15	11
189	Glycine	Amino acid	CK10U	2-1	7
190	Glycohemoglobin	Protein	ProtEx-SP	4-11	32

Column selection guide
1

Ion exchange columns
2

Packed columns for ion chromatography
3

Bioseparation columns
4

Reversed phase columns
5

Chiral separation columns
6

Chromatography media for preparative use
7

MCI®GEL columns
8

MCI®GEL chromatography media
9

Compounds Index
10

	Compound	Classification	MCI® GEL Column	Figure	Page
191	Glycolic acid	Carboxylic acid	CK08EH	2-12	11
192	Glycolic acid	Carboxylic acid	CK08EH	2-18	12
193	Glycyrrhizic acid	Chinese medicinal drug	CHP10M	5-8	45
194	Gly-Tyr	Peptide	CHP2MGM	5-12	47
195	3'-GMP	Nucleotide	CDR10	2-36	19
196	5'-GMP	Nucleotide	CDR10	2-36	19
197	5'-GTP	Nucleotide	CDR10	2-35	19
198	Guanosine	Nucleoside	SCK01	3-6	22
199	Hemoglobin A0	Protein	ProtEx-DEAE	4-6	31
200	Hemoglobin A2	Protein	ProtEx-DEAE	4-6	31
201	Hemoglobin S	Protein	ProtEx-DEAE	4-6	31
202	Histidine	Amino acid	CK10U	2-1	7
203	D-Histidine	D-Amino acid	CRS10W	6-7	53
204	L-Histidine	L-Amino acid	CRS10W	6-7	53
205	Human growth hormone	Hormone	ProtEx-DEAE	4-7	31
206	Human serum	Serum	CQH3ES	4-23	38
207	Human serum	Serum	CQH3PS	4-23	38
208	Hydrocortisone	Adrenal cortex hormone	CHP10M	5-9	45
209	5-Hydroxytryptophan	Amino acid	CHP10M	5-2	42
210	D-2-Hydroxy-n-butyric acid	D- -Hydroxycarboxylic acid	CRS10W	6-16	55
211	L-2-Hydroxy-n-butyric acid	L- -Hydroxycarboxylic acid	CRS10W	6-16	55
212	D- -Hydroxy isocaproic acid	D- -Hydroxycarboxylic acid	CRS10W	6-16	55
213	L- -Hydroxy isocaproic acid	L- -Hydroxycarboxylic acid	CRS10W	6-16	55
214	D- -Hydroxy-n-valeric acid	D- -Hydroxycarboxylic acid	CRS10W	6-16	55
215	L- -Hydroxy-n-valeric acid	L- -Hydroxycarboxylic acid	CRS10W	6-16	55
216	D-m-Hydroxymandelic acid	D- -Hydroxycarboxylic acid	CRS10W	6-17	55
217	L-m-Hydroxymandelic acid	L- -Hydroxycarboxylic acid	CRS10W	6-17	55
218	D-p-Hydroxymandelic acid	D- -Hydroxycarboxylic acid	CRS10W	6-17	55
219	L-p-Hydroxymandelic acid	L- -Hydroxycarboxylic acid	CRS10W	6-17	55
220	Hypoxanthine	Uric acid related compound	CHP10M	5-6	44
221	IgG1 MOPC21(mouse)	monoclonal antibody	ProtEx-DEAE	4-16	34
222	IgG2b, (mouse)	monoclonal antibody	ProtEx-DEAE	4-14	33
223	5'-IMP	Nucleotide	CDR10	2-36	19
224	Interleukin 2	Protein	ProtEx-DEAE	4-9	32
225	Isoleucine	Amino acid	CK10U	2-1	7
226	allo-D-Isoleucine	D-Amino acid	CRS10W	6-1	52
227	allo-L-Isoleucine	L-Amino acid	CRS10W	6-1	52
228	D-Isoleucine	D-Amino acid	CRS10W	6-1	52
229	L-Isoleucine	L-Amino acid	CRS10W	6-1	52
230	Isopropyl alcohol	Alcohol	CK08EH	2-14	11
231	Isopropyl alcohol	Alcohol	CK08EH	2-15	11
232	Lactic acid	Carboxylic acid	CK08EH	2-12	11
233	Lactic acid	Carboxylic acid	CA08F	2-32	17
234	Lactic acid	Carboxylic acid	CA08F	2-34	18
235	-Lactoglobulin	Protein	CQA31S	4-19	35
236	-Lactoglobulin	Protein	CQA35S	4-19	35
237	-Lactoglobulin	Protein	CHP2MG	5-20	50
238	Lactose	Disaccharide	CK08EC	2-4	9
239	Lactose	Disaccharide	CK08EC	2-6	9
240	Lactose	Disaccharide	CA08F	2-31	17
241	Lactose	Disaccharide	CDR10	2-37	20
242	Lactulose	Disaccharide	CK08EC	2-6	9
243	D-Lactic acid	D- -Hydroxycarboxylic acid	CRS10W	6-11	53
244	L-Lactic acid	L- -Hydroxycarboxylic acid	CRS10W	6-11	53
245	D-Lactic acid	D- -Hydroxycarboxylic acid	CRS10W	6-16	55
246	L-Lactic acid	L- -Hydroxycarboxylic acid	CRS10W	6-16	55
247	D-Lactic acid	D- -Hydroxycarboxylic acid	CRS10W/CRS15W	6-19	56
248	L-Lactic acid	L- -Hydroxycarboxylic acid	CRS10W/CRS15W	6-19	56
249	Leucine	Amino acid	CK10U	2-1	7
250	D-Leucine	D-Amino acid	CRS10W	6-2	53
251	L-Leucine	L-Amino acid	CRS10W	6-2	53
252	Leu Enkephalin	Peptide	CQH3PS	4-26	38
253	Leu Enkephalin	Peptide	CHP2MGM	5-12	47
254	Lipoxidase	Enzyme	CQA31S	4-22	36
255	Lithium ion	Cation	SCK01	3-1	22

Column selection guide

1

Ion exchange columns

2

Packed columns for ion chromatography

3

Bioseparation columns

4

Reversed phase columns

5

Chiral separation columns

6

Chromatography media for preparative use

7

MCI® GEL columns

8

MCI® GEL chromatography media

9

Compounds Index

10

	Compound	Classification	MCI® GEL Column	Figure	Page
256	Lysine	Amino acid	CK10U	2-1	7
257	D-Lysine	D-Amino acid	CRS10W	6-8	53
258	L-Lysine	L-Amino acid	CRS10W	6-8	53
259	Lysozyme	Protein	PrepEx-SP	7-1	59
260	Lysozyme	Protein	CQK30S	4-20	36
261	Lysozyme	Protein	CQK31S	4-20	36
262	Magnesium ion	Cation	SCK01	3-7	23
263	Magnesium ion	Cation	SCK01	3-9	23
264	Magnesium ion	Cation	SCK01	3-10	23
265	Magnesium ion	Cation	SCK01	3-11	23
266	Malic acid	Carboxylic acid	CK08EH	2-12	11
267	Malic acid	Carboxylic acid	CA08F	2-32	17
268	Malonic acid	Carboxylic acid	CK08EH	2-12	11
269	Malonic acid	Carboxylic acid	CA08F	2-32	17
270	Maltose	Disaccharide	CA08F	2-31	17
271	Maltose	Disaccharide	CDR10	2-37	20
272	D-Mandelic acid	D- -Hydroxycarboxylic acid	CRS10W	6-17	55
273	L-Mandelic acid	L- -Hydroxycarboxylic acid	CRS10W	6-17	55
274	Manganese ion	Cation	SCK01	3-8	23
275	Mannitol	Sugar alcohol	CK08EC	2-4	9
276	Mannitol	Sugar alcohol	CK08EC	2-11	10
277	Mannose	Sugar	CK08EC	2-4	9
278	Mannose	Sugar	CA08F	2-31	17
279	Mannose	Sugar	CDR10	2-37	20
280	D-Methionine	D-Amino acid	CRS10W	6-3	53
281	L-Methionine	L-Amino acid	CRS10W	6-3	53
282	D-Methionine	D-Amino acid	CRS10W	6-14	54
283	L-Methionine	L-Amino acid	CRS10W	6-14	54
284	Melibiose	Disaccharide	CA08F	2-31	17
285	Melibiose	Disaccharide	CDR10	2-37	20
286	Met Enkephalin	Peptide	CHP2MGM	5-12	47
287	Methionine	Amino acid	CK10U	2-1	7
288	Methyl alcohol	Alcohol	CK08EH	2-15	11
289	Methylamine	Amine	SCK01	3-2	22
290	Met-Leu-Tyr	Peptide	CQH3PS	4-26	38
291	Mouse brain sap	Mouse brain sap	ProtEx-DEAE	4-12	33
292	Myoglobin	Protein	CQP30	4-3	29
293	Myoglobin	Protein	ProtEx-DEAE	4-5	31
294	Myoglobin	Protein	CQA31S	4-18	35
295	Myoglobin	Protein	CQA35S	4-18	35
296	Myoglobin	Protein	CQA31S	4-19	35
297	Myoglobin	Protein	CQA35S	4-19	35
298	Myoglobin	Protein	CQK30S	4-20	36
299	Myoglobin	Protein	CQK31S	4-20	36
300	Myoglobin	Protein	PrepEx-DEAE	7-3	60
301	Nitrate	Anion	SCA04	3-12	24
302	Nitrate	Anion	SCA04	3-13	24
303	Nitrate	Anion	SCA04	3-14	25
304	Nitrate	Anion	SCA04	3-15	25
305	Nitrate	Anion	SCA04	3-17	26
306	Nitrite	Anion	SCA04	3-12	24
307	Nitrite	Anion	SCA04	3-13	24
308	D-Norleucine	D-Amino acid	CRS10W	6-3	53
309	L-Norleucine	L-Amino acid	CRS10W	6-3	53
310	D-Norvaline	D-Amino acid	CRS10W	6-3	53
311	L-Norvaline	L-Amino acid	CRS10W	6-3	53
312	Oligosaccharide	Dp1-Dp9	CK04S	2-21,25	15,16
313	Oligosaccharide	Dp1-Dp13	CK04SS	2-22,26	15,16
314	Oligosaccharide	Dp1-Dp15	CK02A	2-23	15
315	Oligosaccharide	Dp1-Dp23	CK02AS	2-24,27	15,16
316	Orotic acid	Uric acid related compound	CHP10M	5-6	44
317	Ovalbumin	Protein	CQP30	4-3	29
318	Ovalbumin	Protein	CQA31S	4-18	35
319	Ovalbumin	Protein	CQA35S	4-18	35
320	Oxalic acid	Carboxylic acid	CK08EH	2-12	11

	Compound	Classification	MCI® GEL Column	Figure	Page
321	Oxine copper	Agricultural chemical	CHP2MG	5-18	50
322	Pancreatin	Digestive enzyme	ProtEx-DEAE	4-17	34
323	Penicillin G	Penicillin antibiotic	CHP55A	7-8	62
324	Penicillin V	Penicillin antibiotic	CHP55A	7-8	62
325	Phenacetin	Ingredients of medicine	CHP2MGM	5-10	46
326	D-Phenylalanine	D-Amino acid	CRS10W	6-2	53
327	L-Phenylalanine	L-Amino acid	CRS10W	6-2	53
328	D-Phenylalanine	D-Amino acid	CRS10W	6-9	53
329	L-Phenylalanine	L-Amino acid	CRS10W	6-9	53
330	D- -Phenylglycine	D-Amino acid	CRS10W	6-13	54
331	L- -Phenylglycine	L-Amino acid	CRS10W	6-13	54
332	Phenobarbital	Anticonvulsant	CHP5C	5-17	49
333	Phosphate	Anion	SCA04	3-12	24
334	Polyethyleneglycol	Water soluble polymer	CQP30	4-2	29
335	Potassium ion	Cation	SCK01	3-1	22
336	Potassium ion	Cation	SCK01	3-3	22
337	Potassium ion	Cation	SCK01	3-4	22
338	Potassium ion	Cation	SCK01	3-5	22
339	Procainamide	Procainamide	CHP2MGM	5-14	48
340	Procaine	Procaine	CHP2MGM	5-14	48
341	Proline	Amino acid	CK10U	2-1	7
342	D-Proline	D-Amino acid	CRS10W	6-2	53
343	L-Proline	L-Amino acid	CRS10W	6-2	53
344	n-Propyl alcohol	Alcohol	CK08EH	2-14	11
345	Propyleneglycol	Polyalcohol	CK08EH	2-15	11
346	Quinine	Cinchona alkaloid	CHP10M	5-5	43
347	Rhamnose	Sugar	CA08F	2-31	17
348	Rhamnose	Sugar	CDR10	2-37	20
349	Ribonuclease A	Protein	ProtEx-SP	4-10	32
350	Ribonuclease A	Protein	CQK30S	4-20	36
351	Ribonuclease A	Protein	CQK31S	4-20	36
352	Ribonuclease A	Protein	CQK31S	4-21	36
353	Ribonuclease A	Protein	CQK30S	4-21	36
354	Ribonuclease A	Protein	CHP2MGM	5-13	47
355	Ribonuclease A	Protein	CHP2MG	5-20	50
356	Ribonuclease A	Protein	PrepEx-SP	7-1	59
357	Ribonuclease A	Protein	CQH3BP	7-5	61
358	Ribonuclease A	Protein	CQH3PP	7-6	61
359	Ribose	Sugar	CK08EC	2-4	9
360	Ribose	Sugar	CA08F	2-31	17
361	Ribose	Sugar	CDR10	2-37	20
362	RNA type III from bakers yeast	RNA	ProtEx-DEAE	4-13	33
363	Rubidium ion	Cation	SCK01	3-1	22
364	Salicin	Phenol glycoside	CK08EC	2-4	9
365	Senoside A	Crude drug	CHP5C	7-10	63
366	Senoside A	Crude drug	CHP20A	7-10	63
367	Senoside A	Crude drug	CHP20Y	7-10	63
368	Senoside A	Crude drug	CHP20Y	7-11	63
369	Senoside B	Crude drug	CHP5C	7-10	63
370	Senoside B	Crude drug	CHP20A	7-10	63
371	Senoside B	Crude drug	CHP20Y	7-10	63
372	Serine	Amino acid	CK10U	2-1	7
373	D-Serine	D-Amino acid	CRS10W	6-4	53
374	L-Serine	L-Amino acid	CRS10W	6-4	53
375	Serotonin	Catecholamine	CHP10M	5-2	42
376	Sodium ion	Cation	SCK01	3-1	22
377	Sodium ion	Cation	SCK01	3-3	22
378	Sodium ion	Cation	SCK01	3-4	22
379	Sodium ion	Cation	SCK01	3-5	22
380	Sorbitol	Sugar alcohol	CK08EC	2-5	9
381	Sorbitol	Sugar alcohol	CK08EH	2-15	11
382	Stachyose	Tetrasaccharide	CK08EC	2-4	9
383	Strontium ion	Cation	SCK01	3-7	23
384	Strontium ion	Cation	SCK01	3-8	23
385	Succinylsulfathiazole	Sulfa drugs	CHP2MGM	5-11	46

Column selection guide

1

Ion exchange columns

2

Packed columns for ion chromatography

3

Bioseparation columns

4

Reversed phase columns

5

Chiral separation columns

6

Chromatography media for preparative use

7

MCI® GEL columns

8

MCI® GEL chromatography media

9

Compounds Index

10

	Compound	Classification	MCI® GEL Column	Figure	Page
386	Sucrose	Disaccharide	CK08EC	2-5	9
387	Sucrose	Disaccharide	CK08EC	2-7	10
388	Sucrose	Disaccharide	CK08EC	2-8	10
389	Sucrose	Disaccharide	CK04S	2-30	16
390	Sulfamerazine	Sulfa drugs	CHP2MGM	5-11	46
391	Sulfanilamide	Sulfa drugs	CHP2MGM	5-11	46
392	Sulfate	Anion	SCA04	3-12	24
393	Sulfate	Anion	SCA04	3-13	24
394	Sulfate	Anion	SCA04	3-14	25
395	Sulfate	Anion	SCA04	3-15	25
396	Sulfate	Anion	SCA04	3-16	26
397	Sufathiazole	Sulfa drugs	CHP2MGM	5-11	46
398	Tartalic acid	Carboxylic acid	CK08EH	2-12	11
399	Tartalic acid	Carboxylic acid	CA08F	2-32	17
400	Tetracycline	Antibiotics	CHP5C	5-16	49
401	Terephthalic acid	Aromatic carboxylic acid	CHP2MG	5-19	50
402	Theobromine	Purine alkaloid	CHP10M	5-4	43
403	Theophylline	Purine alkaloid	CHP10M	5-4	43
404	Theophylline	Purine alkaloid	CHP10M	5-6	44
405	Thiocyanic ion	Anion	SCA04	3-16	26
406	Thiosulfuric ion	Anion	SCA04	3-16	26
407	Thiuram	Agricultural chemical	CHP2MG	5-18	50
408	Threonine	Amino acid	CK10U	2-1	7
409	Thymine	Nucleic base	CDR10	2-35	19
410	-Tocopherol	Vitamin	CHP2MGY	7-13	64
411	-Tocopherol	Vitamin	CHP2MGY	7-13	64
412	-Tocopherol	Vitamin	CHP2MGY	7-13	64
413	-Tocotrienol	Vitamin	CHP2MGY	7-13	64
414	p-Toluic acid	Aromatic carboxylic acid	CHP2MG	5-19	50
415	TPN	Nucleotide	CDR10	2-35	19
416	Transferrin	Protein	CQA35S	4-19	35
417	Transferrin	Protein	CQA31S	4-19	35
418	Transferrin	Protein	CHP2MG	5-20	50
419	Transferrin	Protein	CQH3BP	7-5	61
420	Transferrin	Protein	CQH3PP	7-6	61
421	Trichloroacetic acid	Carboxylic acid	CK08EH	2-17	12
422	Triethyleneglycol	Polyalcohol	CK08EH	2-16	12
423	Trimethylamine	Amine	SCK01	3-2	22
424	Tyrosine	Amino acid	CK10U	2-1	7
425	D-Tyrosine	D-Amino acid	CRS10W	6-2	53
426	L-Tyrosine	L-Amino acid	CRS10W	6-2	53
427	Trypsin Inhibitor	Enzyme	ProtEx-DEAE	4-5	31
428	Trypsin Inhibitor	Enzyme	CQA31S	4-18	35
429	Trypsin Inhibitor	Enzyme	CQA35S	4-18	35
430	Trypsin inhibitor	Enzyme	PrepEx-DEAE	7-3	60
431	Trypsinogen	Enzyme	CQK31S	4-21	36
432	Trypsinogen	Enzyme	CQK30S	4-21	36
433	Tryptophan	Amino acid	CHP10M	5-2	42
434	D-Tryptophan	D-Amino acid	CRS10W	6-10	53
435	L-Tryptophan	L-Amino acid	CRS10W	6-10	53
436	3'-UMP	Nucleotide	CDR10	2-36	19
437	5'-UMP	Nucleotide	CDR10	2-36	19
438	Uracil	Nucleic base	CDR10	2-35	19
439	Uric acid	Uric acid	CHP10M	5-6	44
440	Uridine	Nucleoside	SCK01	3-6	22
441	Urine	Urine	CDR10	2-38	20
442	Ursodeoxycholic acid	Bile acid	CHP10M	5-7	44
443	Valine	Amino acid	CK10U	2-1	7
444	Valine	Amino acid	CK10U	2-2	8
445	D-Valine	D-Amino acid	CRS10W	6-2	53
446	L-Valine	L-Amino acid	CRS10W	6-2	53
447	Vitamin B3	Water soluble vitamin	CHP2MGM	5-15	48
448	Vitamin B6	Water soluble vitamin	CHP2MGM	5-15	48
449	Vitamin B12	Water soluble vitamin	CHP2MGM	5-15	48
450	Vitamin C	Water soluble vitamin	CHP2MGM	5-15	48

	Compound	Classification	MCI® GEL Column	Figure	Page
451	Xanthine	Uric acid related compound	CHP10M	5-6	44
452	Xylitol	Sugar alcohol	CK08EH	2-4	9
453	Xylitol	Sugar alcohol	CK08EC	2-15	11
454	Xylose	Sugar	CA08F	2-31	17
455	Xylose	Sugar	CDR10	2-37	20
456	Zinc ion	Cation	SCK01	3-8	23

Limited warranty

Mitsubishi Chemical Corporation warrants that its pre-packed columns (including separation media products) shall meet published specifications at the time of shipment from Mitsubishi Chemical Corporation. Because of the susceptibility of these products to deterioration, all warranty claims must be made within the stipulated in the listed sales office. All claims shall be deemed waived in the event the purchaser fails to notify the company within the period.

Conditions

- A. The products in this brochure are for laboratory or manufacturing use. They are not intended for drug, medicinal, food additive or household use. Compliance with local and government regulations concerning their use is the responsibility of the purchaser.
- B. Voiding of warranty :
This warranty is null and void if any product has been (1) altered or modified such that its stability or reliability is any way affected ; (2) misused ; or (3) damaged by abuse, negligence or accident. The term "misuse" includes, but is not limited to, use not in compliance with the "Column Handling Instructions".
- C. Limitations and Exclusions :
All recommendations, information and descriptions supplied by Mitsubishi Chemical Corporation with respect to any product in this brochure are believed to be accurate and reliable, but do not constitute warranties. The sole liability of Mitsubishi Chemical Corporation for any breach of warranty is limited to replacement, or at the sole option of Mitsubishi Chemical Corporation a refund of the purchase price.

Changes

All specifications, quantities, designs and prices are subject to change without notice.

General

Neither this publication, nor any products in this brochure shall be construed as recommending the infringement of any patent, nor extending any license, express or implied, nor assuming any liability under any issued or pending patent. The data presented herein have been carefully compiled from our records which we believe to be accurate and reliable. We make, however, no warranties or representations with respect hereto, nor is freedom from any patent to be inferred.

MITSUBISHI
CHEMICAL



please visit

<http://www.diaion.com>

Sales Office

Mitsubishi Chemical Corporation
Environmental Solution Division
Separation Materials Department
14-1, Shiba 4-chome, Minato-ku, Tokyo 108-0014 Japan
Phone: +81-3-6414-3390 Fax: +81-3-6414-3407

1768-5CC0-018
0803MA1000