2.1.3 C18-Phasen

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C18

2.1.3 C18 Phases

2.1 Packed HPLC Columns and Bulk Material

ProntoSIL C18 ace-EPS

The ProntoSIL C18 ace-EPS belongs to the new group of stationary RP-supports with polar embedded groups. The packing is very stable over a wide pH range (pH 1-10). In addition, it offers a maximum of hydrophobicity combined with a maximum of polar selectivity. The silanophilic activity of the support is very low. Ultra strong basic compounds such as amitriptyline can be eluted from the column at neutral pH values with excellent symmetrical peak shapes.

The main application area of these packings is the pharmaceutical industry, where analytes often have basic or acidic groups. For the separation of these compounds these supports exhibit an enhanced polar selectivity. That means: In comparison to a classical bonded C18 column acidic compounds show a higher retention whereas basic compounds show a slight decrease of retention on an embedded polar column. The C18 ace-EPS -bonding type is available in several particle and pore sizes.

- Bulk material available in 10 g, 100 g and 1000 g quantities.
- Columns with ID 100 mm and Length up to 500 mm on request.


RP-Phases with a polar embedded group

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


2.1.3 C18 Phases

2.1 Packed HPLC Columns and Bulk Material

Chromatographic Tests of Chemically bonded ProntoSIL

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not yet specified

2.1.3 C18 Phases

2.1 Packed HPLC Columns and Bulk Material

ProntoPEARL_{sub2} TPP C18 ace-EPS

ProntoPEARL_{sub2} packings are new stationary phases with particle sizes of < 2 µm. The application area is mainly Fast HPLC since high efficiencies can be achieved in a short time. Compared to conventional HPLC columns you can save up to 90% of analysis time. ProntoPEARL_{sub2} packings are available as TPP (Total Porous Phases) or NPP (Non Porous Phases). Non porous packings are mainly used in separation of polymers or biopolymers (Proteins & Peptides). Porous packings are recommended for the separation of small molecules.

Due to the manufacturing process both packings show a narrow particle size distribution. This leads to extremely high efficiencies in packed columns. Therefore we are able to guarantee a minimum plate number of 200.000 plates per meter for each packed ProntoPEARL_{sub2} column.

Using ProntoPEARL_{sub2} columns which are giving lower dilution (band broadening) leads to a higher mass sensitivity and therefore to a better detection limit. (See next page).


Packing Code	Packing	Phase	Particle Size	Form	Pore Size	Surface Area	%C	Endcapping	Manufacturer
.... C18ATP018	ProntoPEARL _{sub2} TPP C18 ace-EPS	C18	1.8 µm	-	60 Å	300 m ² /g	16	no	BISCHOFF

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
2.1.3 C18 Phases

2.1 Packed HPLC Columns and Bulk Material


ProntoSIL C18 H

ProntoSIL C18 H is the standard C18 support of the Alliance for Chromatography. It is applicable in a wide range of RP-chromatography. The packing is fully endcapped and possesses all of the excellent properties a new generation stationary phase can offer. Keeping in line with all ProntoSIL products, this support is based on an ultra pure silica. The C18 H-bonding type is available in several particle and pore sizes. The wide pore supports show excellent properties for the separation of biomolecules such as proteins and peptides.

- Bulk material available in 10 g, 100 g and 1000 g quantities.
- Columns with ID 100 mm and Length up to 500 mm on request.

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


2.1.3 C18 Phases

2.1 Packed HPLC Columns and Bulk Material

ProntoSIL C18 SH

ProntoSIL C18 SH is the stationary phase in the ProntoSIL line with the highest carbon load. It is fully endcapped. Due to the carbon load it shows an excellent shape selectivity and stability even at pH 1. The C18 SH-bonding type is available with a pore size 120 Å and in several particle sizes.

- Bulk material available in 10 g, 100 g and 1000 g quantities.
- Columns with ID 100 mm and Length up to 500 mm on request.

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


L1

2.1.3 C18 Phases

2.1 Packed HPLC Columns and Bulk Material

ProntoSIL C18 AQ

ProntoSIL C18 AQ with its unique bonding technology has been especially developed for the use in aqueous mobile phases with an organic content below 10%. Standard stationary phases with conventional bondings give very poor peak shapes under these chromatographic conditions due to the collapse of C18-brushes in aqueous eluents. ProntoSIL C18 AQ gives excellent peak shapes in these mobile phases resulting in enhanced selectivities. The advantages of the AQ-packings can be demonstrated in applications of polar analytes. The C18 AQ-bonding type is fully endcapped and available in several particle sizes and pore sizes of 120 Å and 200 Å.


Packing Code	Packing	Phase	Particle Size	Form	Pore Size	Surface Area	%C	Endcapping	Manufacturer
.... F184PS030	ProntoSIL 120-3-C18 AQ	C18 AQ	3 µm	spherical	120 Å	300 m ² /g	14	yes	BISCHOFF
.... H184PS030	ProntoSIL 200-3-C18 AQ	C18 AQ	3 µm	spherical	200 Å	200 m ² /g	9	yes	BISCHOFF
.... F184PS050	ProntoSIL 120-5-C18 AQ	C18 AQ	5 µm	spherical	120 Å	300 m ² /g	14	yes	BISCHOFF
.... H184PS050	ProntoSIL 200-5-C18 AQ	C18 AQ	5 µm	spherical	200 Å	200 m ² /g	9	yes	BISCHOFF
.... F184PS100	ProntoSIL 120-10-C18 AQ	C18 AQ	10 µm	spherical	120 Å	350 m ² /g	14	yes	BISCHOFF

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


L1

2.1.3 C18 Phases

2.1 Packed HPLC Columns and Bulk Material

ProntoSIL C18 AQ PLUS

In addition to ProntoSIL C18 AQ, ProntoSIL C18 AQ PLUS has been developed for the use in aqueous mobile phases with an organic content below 10%. In comparison to ProntoSIL C18 AQ, ProntoSIL C18 AQ PLUS shows an enhanced stability at low pH's down to pH 1. Also the packing shows excellent peak shapes in pure aqueous eluents but differs in shape selectivity compared to ProntoSIL C18 AQ. The application field for this support is mainly in combinatorial chemistry where the standard separation conditions are fast gradients from 0 - 100% organic and where the mobile phases include 0.1 %TFA.

Packing Code	Packing	Phase	Particle Size	Form	Pore Size	Surface Area	%C	Endcapping	Manufacturer
.... F183PS050	ProntoSIL 120-5-C18 AQ PLUS	C18 AQ PLUS	5 µm	spherical	120 Å	300 m ² /g	17	yes	BISCHOFF

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L1

2.1.3 C18 Phases

2.1 Packed HPLC Columns and Bulk Material

ProntoSIL EuroBOND C18

ProntoSIL Eurobond C18 Phase is the best value stationary phase of the ProntoSIL product line. The optimized synthesis procedure for this classical bonded C18 phase leads to this best value high quality product. Additionally to that the low price of these columns results from the fact that this stationary phase is offered in two column dimensions (125 x 4.0 mm and 250 x 4.0 mm) only as well as in the pore size 120 Å and in the particle sizes 5 µm and 10 µm. The selectivity of the ProntoSIL Eurobond C18 phase is in between the selectivity offered by the two stationary phases ProntoSIL 120-5-C18 H Phase and ProntoSIL 120-5-C18 SH. The ProntoSIL Eurobond C18 is fully endcapped and can be used in the brought range of RP applications.


Packing Code	Column Dimensions L x ID in mm	Packing	Phase	Particle Size	Form	Pore Area	Surface	%C	Endcapping	Manufacturer
.... F181PS050	125 x 4.0	ProntoSIL 120-5-EuroBOND C18	C18	5 µm	spherical	120 Å	300 m ² /g	17	yes	BISCHOFF
.... F181PS050	250 x 4.0	ProntoSIL 120-5-EuroBOND C18	C18	5 µm	spherical	120 Å	300 m ² /g	17	yes	BISCHOFF
.... F181PS100	125 x 4.0	ProntoSIL 120-10-EuroBOND C18	C18	10 µm	spherical	120 Å	300 m ² /g	17	yes	BISCHOFF
.... F181PS100	250 x 4.0	ProntoSIL 120-10-EuroBOND C18	C18	10 µm	spherical	120 Å	300 m ² /g	17	yes	BISCHOFF

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
2.1.3 C18 Phases

2.1 Packed HPLC Columns and Bulk Material


ProntoSIL HyperSorb ODS

ProntoSIL HyperSorb ODS is the recommended alternative to Hypersil ODS. It is based on a type A silica and has the same bonding density as the original RP support. Due to a slightly lower acidity of the silica the silanophilic activity of the packing is a little bit lower compared to Hypersil ODS. The polarity of the surface is comparable. In the application of these stationary phases the same selectivity but better peak shape for the elution of basic compounds are observed.

- Bulk material available in 10 g, 100 g and 1000 g quantities.
- Columns with ID 100 mm and Length up to 500 mm on request.

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


2.1.3 C18 Phases

ProntoSIL Spheribond ODS 1

ProntoSIL Spheribond ODS 1 is the recommended alternative for Waters Spherisorb ODS 1. It is based on a type A silica and has the same physical chemical properties and bonding density like the original RP support. The selectivity of the two packings is the same. The typical application field of this classical packing with low carbon content and high silanophilic activity is the separation of hydrophobic compounds like fat soluble vitamins.

Packing Code	Packing	Phase	Particle Size	Form	Pore Size	Surface Area	%C	Endcapping	Manufacturer
.... D180SB030	ProntoSIL Spheribond ODS 1	C18	3 µm	spherical	80 Å	220 m ² /g	7	no	BISCHOFF
.... D180SB050	ProntoSIL Spheribond ODS 1	C18	5 µm	spherical	80 Å	220 m ² /g	7	no	BISCHOFF

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
2.1.3 C18 Phases

ProntoSIL Spheribond ODS 2

ProntoSIL Spheribond ODS 2 is the recommended alternative for Waters Spherisorb ODS 2. It is based on a type A silica and has the same physical chemical properties and bonding density like the original RP support. The selectivity of the two packings is the same. The typical application field of this classical packing with moderate silanophilic activity is pharmaceutical analytes. The Waters Spherisorb ODS 2 column is the most popular spherical silica support in pharmaceutical applications.

- Bulk material available in 10 g, 100 g and 1000 g quantities.
- Columns with ID 100 mm and Length up to 500 mm on request.

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



L1

2.1.3 C18 Phases

2.1 Packed HPLC Columns and Bulk Material

Eurosil Bioselect 300 C18

Eurosil Bioselect 300 has been specifically developed for the determination and purification of peptides, proteins and oligonucleotides. The Eurosil Bioselect sorbent is characterized by high loadability and a long lifetime.

Packing Code	Packing	Phase	Particle Size	Form	Pore Size	Surface Area	%C	Endcapping	Manufacturer
.... K181EB050	Eurosil Bioselect 300 C18	C18	5 µm	spherical	300 Å	100 m ² /g	6	yes	KNAUER

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2.1.3 C18 Phases

2.1 Packed HPLC Columns and Bulk Material

Eurospher 100 C18

Eurospher C18 packing material can be universally used in different application areas and has been on the market for more than 15 years. The most often used material is the endcapped Eurospher 100 C18 for the reversed phase mode, followed by the less hydrophobic Eurospher 100 C8. Both can be considered for the determination of water soluble vitamins, steroids, catecholamines, sedatives and in many other application fields in the analytical as well as preparative range. The stationary phase is stable in a pH range of 2 - 8.5. Eurospher C18 is the workhorse of reversed phase applications.

- Bulk material available in 10 g and 100 g quantities.
- Columns with ID 100 mm and Length up to 500 mm on request.

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L1

2.1.3 C18-Phasen

2.1 Packed HPLC Columns and Bulk Material

Kromasil C18

High performance spherical silica for analytical to process scale liquid chromatography. RP Kromasil is manufactured using monofunctional silanes, and all derivatized Kromasil is fully endcapped. This gives high reproducibility and chemical stability.

Packing Code	Packing	Phase	Particle Size	Form	Pore Size	Surface Area	%C	Endcapping	Manufacturer
.... E181KS035	Kromasil 100-3.5-C18	C18	3.5 µm	spherical	110 Å	330 m ² /g	20	yes	EKA CHEMICALS
.... E181KS050	Kromasil 100-5-C18	C18	5 µm	spherical	110 Å	330 m ² /g	20	yes	EKA CHEMICALS
.... E181KS070	Kromasil 100-7-C18	C18	7 µm	spherical	110 Å	330 m ² /g	20	yes	EKA CHEMICALS
.... E181KS100	Kromasil 100-10-C18	C18	10 µm	spherical	110 Å	330 m ² /g	20	yes	EKA CHEMICALS
.... E181KS130	Kromasil 100-13-C18	C18	13 µm	spherical	110 Å	330 m ² /g	20	yes	EKA CHEMICALS
.... E181KS160	Kromasil 100-16-C18	C18	16 µm	spherical	110 Å	330 m ² /g	20	yes	EKA CHEMICALS
.... K181KS050	Kromasil 300-5-C18	C18	5 µm	spherical	300 Å	110 m ² /g	8.1	yes	EKA CHEMICALS
.... K181KS100	Kromasil 300-10-C18	C18	10 µm	spherical	300 Å	110 m ² /g	8.1	yes	EKA CHEMICALS
.... K181KS160	Kromasil 300-16-C18	C18	16 µm	spherical	300 Å	110 m ² /g	8.1	yes	EKA CHEMICALS

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L1

2.1.3 C18 Phases

2.1 Packed HPLC Columns and Bulk Material

LiChrospher RP18 and RP18 endcapped

LiChrospher RP-18 and LiChrospher RP-18 endcapped are reliable and versatile traditionally produced spherical silica support with reversed phase properties.

They are well suited for the chromatography of acidic, neutral and weakly basic compounds, substances found frequently in all analytical fields.

A selectivity and performance guarantee makes sure that these parameters remain constant batch to batch and from year to year.

Both stationary phases differ in polarity.

Superspher RP18 and RP18 endcapped

spherical 4 µm silica particles with octadecyl modification, available with or without endcapping.

Packing Code	Packing	Phase	Particle Size	Form	Pore Size	Surface Area	%C	Endcapping	Manufacturer
.... E181LS050	LiChrospher 100 RP 18e	C18	5 µm	spherical	100 Å	350 m ² /g	21.6	yes	MERCK
.... E181LS100	LiChrospher 100 RP 18e	C18	10 µm	spherical	100 Å	350 m ² /g	21.6	yes	MERCK
.... E189LS050	LiChrospher 100 RP 18	C18	5 µm	spherical	100 Å	350 m ² /g	21	no	MERCK
.... E189LS100	LiChrospher 100 RP 18	C18	10 µm	spherical	100 Å	350 m ² /g	21	no	MERCK
.... E189SS040	Superspher 100 RP 18	C18	4 µm	spherical	100 Å	350 m ² /g	21	no	MERCK
.... E181SS040	Superspher 100 RP 18e	C18	4 µm	spherical	100 Å	350 m ² /g	21.6	yes	MERCK

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- Bulk material available in 10 g and 100 g quantities, Superspher in 10 g only.
- Columns with ID 100 mm and Length up to 500 mm on request.

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L1

2.1.3 C18-Phasen

2.1 Packed HPLC Columns and Bulk Material

Nucleosil C18

Nucleosil C18 phases are suitable for reversed phase chromatography and ion-pair chromatography for the separation of nonpolar to moderately polar compounds such as fatty acids, glycerides, polycyclic aromatics, esters (phthalates), vitamins, steroids, prostaglandins, PTH amino acids, etc.

Packing Code	Packing	Phase	Particle Size	Form	Pore Size	Surface Area	%C	Endcapping	Manufacturer
.... E180NS030	Nucleosil 100-3-C18	C18	3 µm	spherical	100 Å	350 m ² /g	14	yes	MACHEREY-NAGEL
.... E180NS050	Nucleosil 100-5-C18	C18	5 µm	spherical	100 Å	350 m ² /g	14	yes	MACHEREY-NAGEL
.... E180NS070	Nucleosil 100-7-C18	C18	7 µm	spherical	100 Å	350 m ² /g	14	yes	MACHEREY-NAGEL
.... E180NS100	Nucleosil 100-10-C18	C18	10 µm	spherical	100 Å	350 m ² /g	14	yes	MACHEREY-NAGEL
.... F180NS030	Nucleosil 120-3-C18	C18	3 µm	spherical	120 Å	200 m ² /g	11	yes	MACHEREY-NAGEL
.... F180NS050	Nucleosil 120-5-C18	C18	5 µm	spherical	120 Å	200 m ² /g	11	yes	MACHEREY-NAGEL
.... F180NS070	Nucleosil 120-7-C18	C18	7 µm	spherical	120 Å	200 m ² /g	11	yes	MACHEREY-NAGEL
.... F180NS100	Nucleosil 120-10-C18	C18	10 µm	spherical	120 Å	200 m ² /g	11	yes	MACHEREY-NAGEL

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L1

2.1.3 C18 Phases

2.1 Packed HPLC Columns and Bulk Material

Nucleosil C18 AB

C18 phases are suited for reversed phase chromatography and ion-pairing chromatography for separation of nonpolar to moderately polar compounds such as fatty acids, glycerides, polycyclic aromatics, esters (phthalates), vitamins, steroids, prostaglandins, PTH amino acids etc.. In contrast to classical Nucleosil 100-5-C18 packing the Nucleosil 100-5-C18 AB has a high bonding density and therefore provides a shape selectivity.

Packing Code	Packing	Phase	Particle Size	Form	Pore Size	Surface Area	%C	Endcapping	Manufacturer
.... E181NS050	Nucleosil 100-5-C18 AB	C18	5 µm	spherical	100 Å	350 m ² /g	24	yes	MACHEREY-NAGEL

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L1

2.1.3 C18-Phasen

2.1 Packed HPLC Columns and Bulk Material

LiChrosorb RP18

LiChrosorb has been one of the most successful and reliable packing materials used in HPLC for more than 25 years and has been referenced in the literature in the form of several thousand applications. The irregular LiChrosorb RP 18 material is not endcapped. The completely porous irregular particles are classified in the 5, 7 and 10 µm range.

LiChrosorb packing materials offers the complete program of non-polar derivatives (RP8, RP18, RP select B) polar derivatives (Si 60 and Si 100) and derivatives of medium polarity (NH₂, CN and Diol).

Packing Code	Packing	Phase	Particle Size	Form	Pore Size	Surface Area	%C	Endcapping	Manufacturer
.... E689LB050	LiChrosorb RP 18	C18	5 µm	irregular	100 Å	300 m ² /g	16.2	no	MERCK
.... E689LB070	LiChrosorb RP 18	C18	7 µm	irregular	100 Å	300 m ² /g	16.2	no	MERCK
.... E689LB100	LiChrosorb RP 18	C18	10 µm	irregular	100 Å	300 m ² /g	16.2	no	MERCK

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L1

2.1.3 C18 Phases

2.1 Packed HPLC Columns and Bulk Material

Hipak ODS AB

Hipak is the typical second generation packing of Bischoff Chromatography. It is based on a type A silica but has a base deactivation. This packing is our recommended alternative to Hypersil BDS. Both a C18 and a C8 bonding are offered. The packing is available in the particle sizes 3 µm and 5 µm.

Packing Code	Packing	Phase	Particle Size	Form	Pore Size	Surface Area	%C	Endcapping	Manufacturer
.... F181HP030	Hipak ODS AB	C18	3 µm	spherical	130 Å	170 m ² /g	11	yes	BISCHOFF
.... F181HP050	Hipak ODS AB	C18	5 µm	spherical	130 Å	170 m ² /g	11	yes	BISCHOFF

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L1

2.1.3 C18 Phases

2.1 Packed HPLC Columns and Bulk Material

Partisil ODS

Partisil was one of the first commercial available type A silica packings on the HPLC market. The benefit of these packings compared to others was the good batch to batch reproducibility of those stationary phases. It is an irregular shaped support. Several different bonding types are existing. There are three different C18 packings available. Whereas Partisil ODS is a low density monomeric C18 phase with no endcapping, the Partisil ODS(2) is polymeric bonded and Partisil ODS(3) has an additional endcapping combined with a polymeric bonding. The Partisil PAC phase is an alkyl bonded phase with a mixed Amino/Cyano endcapping.

Packing Code	Packing	Phase	Particle Size	Form	Pore Size	Surface Area	%C	Endcapping	Manufacturer
.... D681PA050	Partisil ODS3	C18	5 µm	irregular	85Å	350 m ² /g	10.5	yes	Whatman
.... D681PA100	Partisil ODS3	C18	10 µm	irregular	85Å	350 m ² /g	10.5	yes	Whatman
.... D689PA100	Partisil ODS2	C18	10 µm	irregular	85Å	350 m ² /g	15	no	Whatman
.... D680PA100	Partisil ODS	C18	10 µm	irregular	85Å	350 m ² /g	5	no	Whatman

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2.1.3 C18 Phases

2.1 Packed HPLC Columns and Bulk Material

PolyEncap C18

PolyEncap stationary phases are a family of bonded phases where the silica is polymeric encapsulated by acrylic acid ester polymers that are additionally covalent bonded to the silica surface. This type of stationary phases offers different selectivities compared to classical bonded brush type stationary phases. The PolyEncap phases show a slightly higher pH stability in basic mobile phases. Three different bonding chemistries are offered: These are two RP packings with C4 and C18 and one weak cation exchanger (WCX). The C4 also called PolyEncap A has a comparable selectivity to a brush type C8 packing, whereas the corresponding C18 (PolyEncap B) is similar to a brush type C18 packing. The WCX packing is able to separate mono valent and bivalent cations in a single run.

The two RP bondings are available in four particle sizes (3 µm, 5 µm, 7 µm and 10 µm) and two pore sizes (120 Å and 300 Å) whereas the WCX packing is based on a 5 µm, 120 Å silica.

Packing Code	Packing	Phase	Particle Size	Form	Pore Size	Surface Area	%C	Endcapping	Manufacturer
.... F180PE030	PolyEncap B	C18	3 µm	spherical	120 Å	200 m ² /g	14	no	BISCHOFF
.... F180PE050	PolyEncap B	C18	5 µm	spherical	120 Å	200 m ² /g	14	no	BISCHOFF
.... F180PE070	PolyEncap B	C18	7 µm	spherical	120 Å	200 m ² /g	14	no	BISCHOFF
.... F180PE100	PolyEncap B	C18	10 µm	spherical	120 Å	200 m ² /g	14	no	BISCHOFF
.... K180PE030	PolyEncap B	C18	3 µm	spherical	300 Å	100 m ² /g	7	no	BISCHOFF
.... K180PE050	PolyEncap B	C18	5 µm	spherical	300 Å	100 m ² /g	7	no	BISCHOFF
.... K180PE070	PolyEncap B	C18	7 µm	spherical	300 Å	100 m ² /g	7	no	BISCHOFF
.... K180PE100	PolyEncap B	C18	10 µm	spherical	300 Å	100 m ² /g	7	no	BISCHOFF

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
2.1.3 C18 Phases

Techsphere ODS

The Techsphere range of spherical, silica based materials was introduced because of the problems some manufacturers have producing reproducible media. Techsphere ODS is the material for standard reversed phase separations.

Packing Code	Packing	Phase	Particle Size	Form	Pore Size	Surface Area	%C	Endcapping	Manufacturer
... E180TE030	Techsphere 100 ODS	C18	3 µm	spherical	100 Å	200 m ² /g	11	no	Grace Vydac
... E180TE050	Techsphere 100 ODS	C18	5 µm	spherical	100 Å	200 m ² /g	11	no	Grace Vydac

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L1

2.1.3 C18 Phases

2.1 Packed HPLC Columns and Bulk Material

Zorbax ODS

The classical Zorbax product line is typical for family of spherical first generation (type A) packings. The stationary phases show a high hydrophobicity and the support has an excellent mechanical strength. AfC offers the complete product line of the classical Zorbax supports in all column dimensions.

Packing Code	Packing	Phase	Particle Size	Form	Pore Size	Surface Area	%C	Endcapping	Manufacturer
.... D180ZX050	Zorbax ODS	C18	5 µm	spherical	70 Å	300 m ² /g	20	yes	Agilent Technologies

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Zorbax is a registered trademark of DuPont.

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L1

2.1.3 C18 Phases

2.1 Packed HPLC Columns and Bulk Material

Sphere-Image 80 ODS1 / ODS 2

The technical specifications of the Sphere-Image 80 series of packing materials are nearly equivalent to those of Spherisorb 80, a high-performance packing material which has existed for decades and is very well-known. The Sphere-Image 80 packing material is also available with all of the typical modifications. The cation and anion exchangers Sphere-Image SCX and SAX are also available. The most popular materials in the Sphere-Image series are ODS 1 (not endcapped) and end-capped ODS 2. All of these packing materials exhibit excellent batch-to-batch reproducibility.

Packing Code	Packing	Phase	Particle Size	Form	Pore Size	Surface Area	%C	Endcapping	Manufacturer
.... D181SI030	Sphere-Image ODS2	C18	3 µm	spherical	80 Å	200 m ² /g	12	yes	Phenomenex
.... D181SI050	Sphere-Image ODS2	C18	5 µm	spherical	80 Å	200 m ² /g	12	yes	Phenomenex
.... D181SI100	Sphere-Image ODS2	C18	10 µm	spherical	80 Å	200 m ² /g	12	yes	Phenomenex
.... D189SI030	Sphere-Image ODS1	C18	3 µm	spherical	80 Å	200 m ² /g	7	no	Phenomenex
.... D189SI050	Sphere-Image ODS1	C18	5 µm	spherical	80 Å	200 m ² /g	7	no	Phenomenex
.... D189SI100	Sphere-Image ODS1	C18	10 µm	spherical	80 Å	200 m ² /g	7	no	Phenomenex

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• Bulk material available in 10 g and 100 g quantities.

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not yet specified

2.1.3 C18 Phases

2.1 Packed HPLC Columns and Bulk Material

ProntoPEARL_{sub2} NPP C18

ProntoPEARL_{sub2} packings are new stationary phases with particle sizes of < 2 µm. The application area is mainly Fast HPLC since high efficiencies can be achieved in a short time. Compared to conventional HPLC columns you can save up to 90% of analysis time. ProntoPEARL_{sub2} packings are available as TPP (Total Porous Phases) or NPP (Non Porous Phases). Non porous packings are mainly used in separation of polymers or biopolymers (Proteins & Peptides). Porous packings are recommended for the separation of small molecules.

Due to the manufacturing process both packings show a narrow particle size distribution. This leads to extremely high efficiencies in packed columns. Therefore we are able to guarantee a minimum plate number of 200.000 plates per meter for each packed ProntoPEARL_{sub2} column.

Using ProntoPEARL_{sub2} columns which are giving lower dilution (band broadening) leads to a higher mass sensitivity and therefore to a better detection limit.

Packing Code	Packing	Phase	Particle Size	Form	Pore Size	Surface Area	%C	Endcapping	Manufacturer
.... Z189NP015	ProntoPEARL _{sub2} NPP C18	C18	1.5 µm	-	-	3 m ² /g	0.5	no	BISCHOFF

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2.1.3 C18 Phases

Micra NPS ODS

Micra NPS phases are based on non porous silica. In contrast to conventional porous silica non porous silica has a smaller surface area. Non porous means that the mass transfer of those packings is accelerated. This fact results in higher efficiencies of the corresponding HPLC columns. Due to lower surface area compounds show lower retention on non porous packings compared to porous ones. To achieve the same retention in the chromatography of small molecules about 20% less organic has to be used with RP packings. Therefore the use of non porous packings is in those areas where high organic mobile phases are required. Columns packed with non porous stationary phases are the columns of choice for the chromatography of macromolecules like proteins, polynucleotides and polymers. On porous packings those molecules have to move by diffusion inside the pores. Due to the diffusion coefficient of macromolecules is much lower compared to small molecules the efficiencies of columns packed with porous particles is much lower compared to those that can be achieved on columns packed with non porous stationary phases. There are several bonding types available: two non endcapped C18 bondings called Micra NPS ODS 1 and Micra NPS ODS 2, an endcapped C18 packing Micra NPS ODS 3E and a C30 packing called TAS.

Packing Code *	Packing	Phase	Particle Size	Form	Pore Size	Surface Area	%C	Endcapping	Manufacturer
.... ODS101.5	Micra NPS ODS 1	C18	1.5 µm	spherical	-	3m ² /g	0.5	no	Eprogen
.... ODS201.5	Micra NPS ODS 2	C18	1.5 µm	spherical	-	3m ² /g	0.5	no	Eprogen
.... ODS3E1.5	Micra NPS ODS 3E	C18	1.5 µm	spherical	-	3m ² /g	0.5	yes	Eprogen

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• Bulk material not available.

* The ordering number of all Micra phases is a special packing code.

C18