

Exploiting Selectivity in HPLC and UHPLC With Rational Stationary Phase Design

Thomas J. Waeghe, Ph.D.
twaghe@mac-mod.com

MAC-MOD Analytical, Inc.
800-441-7508



www.mac-mod.com




Outline

- ◆ Chromatographic **selectivity**
- ◆ Stationary phase **design** concepts
- ◆ The unique **ACE[®] C18-AR** and **ACE[®] C18-PFP** phases
- ◆ Introducing the **NEW ACE[®] Excel[™]** UHPLC products
- ◆ Examples
- ◆ Conclusions

Chromatographic Peak Resolution

Efficiency Selectivity Retention

$$R_s = \frac{\sqrt{N}}{4} \cdot \frac{\alpha - 1}{\alpha} \cdot \frac{k}{1+k}$$



The Importance of N , k and α For Resolution

Typical separation:

$N = 10,000$ plates

$k = 3.8 / 4.2$ (4.0 mean)

$\alpha = 1.1$

$$R_s = \frac{\sqrt{N}}{4} \frac{\alpha - 1}{\alpha} \frac{k}{1+k}$$

$$R_s = \frac{1}{4} \sqrt{10,000} \left(\frac{1.1 - 1}{1.1} \right) \left[\frac{4}{1 + 4} \right]$$

$$R_s = 1.8$$

Which looks like



The Importance of N , k and α For Resolution


Double Efficiency (e.g., 5 μm \rightarrow 2.5 μm):

$N = 10,000 \rightarrow 20,000$ plates

$$R_s = \frac{\sqrt{N}}{4} \cdot \frac{\alpha-1}{\alpha} \cdot \frac{k}{1+k}$$

$$R_s = \frac{1}{4} \sqrt{20,000} \left(\frac{1.1-1}{1.1} \right) \left[\frac{4}{1+4} \right]$$

$$R_s = 2.6$$



Opportunity to optimize further, e.g. reduce column length to speed up

The Importance of N , k and α For Resolution


Double Retention Factor (e.g., decrease solvent strength):

$k = 4 \rightarrow 8$

$$R_s = \frac{\sqrt{N}}{4} \cdot \frac{\alpha-1}{\alpha} \cdot \frac{k}{1+k}$$

$$R_s = \frac{1}{4} \sqrt{10,000} \left(\frac{1.1-1}{1.1} \right) \left[\frac{8}{1+8} \right]$$

$$R_s = 2.0$$



Slight improvement in resolution has led to increased analysis time

The Importance of N , k and α For Resolution

Increase Selectivity (e.g., change column):

$\alpha = 1.1 \rightarrow 1.2$

$$R_s = \frac{\sqrt{N}}{4} \cdot \frac{\alpha - 1}{\alpha} \cdot \frac{k}{1 + k}$$

$$R_s = \frac{1}{4} \sqrt{10,000} \left(\frac{1.2 - 1}{1.2} \right) \left[\frac{4}{1 + 4} \right]$$

$$R_s = 3.3$$



$R_s = 1.8$

Can now shorten run time with shorter column and/or faster flow rate



$R_s = 3.3$

~80% Increase

Significant opportunity to speed up for modest change in selectivity

Selectivity: The Key to Chromatographic Peak Resolution


$$R_s = \frac{\sqrt{N}}{4} \cdot \frac{\alpha - 1}{\alpha} \cdot \frac{k}{1+k}$$

Efficiency Selectivity Retention

~40% ~80% ~10%



From the examples, selectivity has the greatest impact on increasing peak resolution



Which Factors^a Affect Selectivity?

**MOST
Influence**



**LEAST
Influence**

Isocratic Separations

- ◆ Column stationary phase
- ◆ Organic modifier
- ◆ pH (ionised analytes only)
- ◆ % Organic modifier
- ◆ Buffer selection
- ◆ Column temperature
- ◆ Buffer concentration

Gradient Separations


- ◆ All parameters for isocratic
- ◆ Gradient steepness
- ◆ k^*
- ◆ Dwell volume
- ◆ Column dimensions

$$k^* = \frac{85 \times t_G \times F}{\Delta\Phi \times V_m \times S}$$


^a Adapted from 'Introduction to Modern Liquid Chromatography', 3rd Edition, Snyder, Kirkland, Dolan, 2010, p.29, Wiley & Sons

Influencing Selectivity – Bonded Phase Effects / Basic Analytes


ACE C18 – Increase Retention




ACE CN – Elution Order




ACE C8 (start point)



ACE C4 – Decrease Retention



ACE Phenyl – Elution Order



Use **ultra high purity silica** for **good chromatography** and **reproducibility**


HPLC End User Surveys^a ...Listening To The Analyst

- ◆ **Column reproducibility** and **column lifetime** are major factors for analysts
 - Have been the **top 2** feedback points since **2007**
 - **Critical** in **pharmaceutical** and other **major industries** for **method transfers / consistency** and **long term** performance

- ◆ **Reversed-phase** is the **dominant** separation mode
 - **C18 & C8 = 60%**; **Phenyl = 16%**; **CN = 9.5%**; **Fluorinated = 5.9%**
 - **92% analysts** use **C18** at some time in their work...they typically meet the above criteria

 - **BUT limited** selectivity

16 Pharmaceutically Relevant Analytes – C18 Columns



C18 phases show 'similar' selectivity...

The Challenge...

- ◆ To **engineer** new phases with **alternative** selectivity but with the robust properties of the **C18** ligand
 - **Reproducible** (column-to-column & batch-to-batch)
 - **Excellent column lifetime**
 - **Superb efficiency** provided by **ultra-inert, ultra-pure** silica particle
 - **Low MS bleed**
 - Usable in **100% aqueous** eluents

- ◆ **Available** for **HPLC & UHPLC** separations

- ◆ **Available** as a '**Phase III Ready**' product family
 - ◆ Globally available, supply chain, reproducible, multiple batches etc



Aromatic Functionality – Engineering New Stationary Phases

- ◆ Phases with aromatic functionality include **phenyl** and pentafluorophenyl (**PFP**) based ligands
- ◆ Advantages
 - Aromatic functionality potentially offer **unique interactions** with analytes (c.f. C18) giving **alternative** selectivity
 - Provides **enhanced retention** of polar compounds
 - Many aromatic functionality-based phases can be used in **100% aqueous** eluents
- ◆ Disadvantages
 - Phenyl / PFP phases **may suffer** phase bleed
 - **Batch-to-batch reproducibility** & **robustness** may be weak


Aromatic Functionality: π – π Interactions

- ◆ A type of electron **donor-acceptor** interaction
- ◆ Originates from π systems in **unsaturated functional groups** on analytes and the stationary phase
- ◆ Types of π - π interaction can be **manipulated** for maximum effect (**orthogonality**) in phase design
 - eg phenyl: **electron rich** ring on the stationary phase also acts as **π -base** and interacts well with **electron poor analytes**
 - eg PFP: **electron poor** ring on the stationary phase also acts as **π -acid** and interacts well with **electron rich analytes**


The Power of π ...Scientific Led Stationary Phase Design



Electron Donating Groups
eg NH_2 , NR_2 , alkyl, OCH_3
 OR , CH_3 , Ar etc




Electron **Rich** Ring




Activity: π -donor (π -base)

Electron Withdrawing Groups
eg NO_2 , halides, NR_3^+ , CO_2H ,
 CN , CO_2R , SO_3H , COH etc




Electron **Deficient** Ring



Activity: π -acceptor (π -acid)

Classic π - π
interaction



How do we exploit these properties for new stationary phases?

C18+Phenyl = ACE® C18-AR

C18+PFP = ACE® C18-PFP

Uniquely Designed Stationary Phases

- ◆ **ACE[®] C18-AR (USP L1)**
 - Ligand has C18 **hydrophobic** element PLUS **phenyl** character
- ◆ **ACE[®] C18-PFP (USP L1)**
 - Ligand has C18 **hydrophobic** element PLUS **PFP** character
- ◆ **Ultra-inert, ultra-pure** silica particle technology as used in all ACE[®] products for **high** peak efficiency
- ◆ Available in **3, 5 & 10 μ m**, (ACE[®]) and **2 μ m** (ACE[®] *Exce/™*)

Multi-mode interaction mechanisms result in **enhanced** chromatographic **selectivity** giving the analyst **new options** for method development

ACE[®] C18-AR: Multi-Mode Separation Mechanisms

- Combining the character of **C18+phenyl** into a single individual phase harnesses **the best** of both phases for **unique** selectivity





Separation mechanism	Typical C18	Typical Phenyl	ACE [®] C18-AR
Hydrophobicity	++++	+ / ++	++++
π - π Interaction	-	+++	+++
Dipole - Dipole	-	+	+
Hydrogen Bonding	-	++	++
Shape Selectivity	++	++	++ / +++

- The predominance of each retention mechanism will be dictated by the analyte's physicochemical properties, its structure and the chromatographic conditions applied

Multi-Mode Interactions Offer the Chromatographer More

ACE® C18-AR Aromatic Selectivity

- ♦ Illustrating **hydrophobicity** and **π -base character / aromatic selectivity** with a **simple example** using substituted aromatics



Log P*

1.2

1.6

1.9

2.7

π -acidity (order)

1


>

2

>

3

-



* Predicted data from ACD Labs software, 30May12

ACE[®] C18-PFP: Multi-Mode Separation Mechanism

- Combining the character of **C18+PFP** into a single individual phase harnesses **the best** of both phases for **unique** selectivity

Separation mechanism	Typical C18	Typical PFP	ACE [®] C18-PFP
Hydrophobicity	++++	+ / ++	++++
π - π Interaction	-	+++	+++
Dipole - Dipole	-	++++	++++
Hydrogen Bonding	-	+++	+++
Shape Selectivity	++	+++	++++

- The predominance of each retention mechanism will be dictated by the analyte's physicochemical properties, its structure and the chromatographic conditions applied


Multi-Mode Interactions Offer the Chromatographer More

ACE[®] C18-PFP Selectivity*

Peak Number:	1	2	3	4	5	6	7	8
	1,2,3-TMB	1,2,4-TMB	1,2-DMB	1,4-DMB	MB	1,3-DMB	1,3,5-TMB	Tol
Log P:	1.7	1.6	1.7	2.1	2.2	2.2	1.6	2.7
π -basicity (order):	1	1	2	2	3	2	1	-

- ◆ Elution / retention **not simply** a function of π -basicity and Log P
- ◆ Retention mechanism for C18-PFP **multi-modal**

ACE® C18-PFP Selectivity





- ◆ **C18 or PFP mechanisms alone not enough** to fully resolve the methoxybenzene **isomers**
- ◆ **ACE C18-PFP mechanism combines hydrophobicity, shape selectivity, dipole-dipole and π - π interactions**
- ◆ **Elution order, retention and selectivity all seen to differ**
- ◆ **Powerful positional isomer and shape selectivity**


1) 1,2,3-trimethoxybenzene, 2) 1,2,4-trimethoxybenzene, 3) 1,2-dimethoxybenzene, 4) 1,4-dimethoxybenzene 5) methoxybenzene, 6) 1,3-dimethoxybenzene, 7) 1,3,5-trimethoxybenzene, 8) toluene (ref) Mobile phase 50:50 v/v MeOH / H₂O; Column= 150 x 4.6 mm id; 1.00 ml/min; 40C; 254 nm

ACE® Phase Comparisons With The Selectivity Descriptor*

Selectivity = 100 x $\sqrt{1 - R^2}$



Selectivity = 100 x $\sqrt{1 - R^2}$
= 100 x $\sqrt{1 - 0.9888}$
= 10.6



* Neue, O’Gara, Méndez “Selectivity in Reversed-Phase Separations: Influence of the Stationary Phase”, J. Chromatogr. A 1127 (2006), 161-174

Ranking ACE® Phase Orthogonality With MeOH and MeCN

- ◆ For the 102 acidic, basic and neutral analytes assessed

MeOH

Column 1	Column 2	Selectivity 'S'
C18	C18-AR	12
C18	C18-PFP	11
C18-AR	C18-PFP	10

MeCN

Column 1	Column 2	Selectivity 'S'
C18	C18-AR	8
C18-AR	C18-PFP	8
C18	C18-PFP	7

MeOH	MeCN	Selectivity Value
C18-PFP	C18	19
C18-AR	C18	18
C18-AR	C18-PFP	18
C18-PFP	C18-AR	18
C18-PFP	C18-PFP	18
C18	C18-AR	17
C18	C18-PFP	17
C18	C18	15
C18-AR	C18-AR	15




Shows value of using the 3 phases in a 2 solvent screen for method development work

What Do I Use These Novel Phases For: ACE C18-PFP?

- ◆ Useful for analytes that contain **electron donating** moieties
eg -NH₂, -NR₂, -OCH₃, -OH, -alkyl, -Ar etc
- ◆ eg nucleotides, nucleosides, nucleobases, halogenated aryl /
aromatics, catecholamines, tetracyclines, beta blockers,
structural isomers, coumarins etc
- ◆ **Excellent** shape and **positional isomer** selectivity

What Do I Use These Novel Phases For: ACE C18-PFP?

- ◆ **C18-PFP: chloroacetophenone halogenated isomers separation**




**Regioisomer
Selectivity**

What Do I Use These Novel Phases For: ACE C18-AR?

- ◆ Useful for analytes that contain **electron withdrawing** moieties eg -NO₂, -halides, -NR₃⁺, -SO₂, -CO₂H, -SO₃H, -CO₂R, -CHO etc
- ◆ eg aromatic compounds, anthocyanins, steroids, analgesics, phenolics, water soluble vitamins, sulphur containing compounds, quinolones, positional isomers etc
- ◆ **Moderate** shape selectivity

What Do I Use These Novel Phases For: ACE C18-AR?




Steroid Mixture


Columns: 4.6 x 150 mm, 3 μ m
A: Water
B: ACN
Temp. 20°C
Flow Rate: 1 mL/min
Gradient: 25 to 46%B, 24 min

- 1) 17 α Estradiol, 2) Prednisolone, 3) Corticosterone, 4) 19-Norethindrone, 5) Estrone, 6) 17 α Ethynylestradiol,
 7) Cortisone 21-acetate, 8) 21-Hydroxyprogesterone, 9) Estriol, 10) 17 β Estradiol, 11) Hydrocortisone,
 12) Cortisone


Steroid Sample




Corticosterone




Cortisone 21-acetate




19-Norethindrone




17 alpha-Ethynylestradiol




17 alpha-Estradiol




17 beta-Estradiol




Prednisolone




Hydrocortisone




Cortisone



Estriol



Estrone




21-Hydroxyprogesterone

Combining Selectivity With The NEW ACE® *Excel*™ Format

- ◆ **NEW high efficiency, ultra-inert** 2µm silica particles suitable for UHPLC at 1000bar (15,000psi)
- ◆ **Nine selectivities** – including the **unique C18-AR** and **C18-PFP**
- ◆ **High reproducibility:** column-to-column and batch-to-batch
- ◆ **Ultra-robust** phases: **NEW** low dispersion column **hardware** and **NEW** High Stability Column (**HSC™**) packing technology
- ◆ Engineered with **lower back pressures (~30%)** compared to other < 2µm phases, due to **2µm** particle size and **frit technology**
- ◆ **Fully scalable** to **ACE® 3µm, 5µm and 10µm** phases
- ◆ **Fully compatible** with all commercial HPLCs and UHPLCs

ACE® *Excel*™ UHPLC Column Robustness


- ◆ **1000 bar for ~2000 gradient runs**
 - ◆ Isocratic efficiency assessments every ~100 runs ← more demanding!



100x2.1mm; MPA 0.1% FA (aq); MPB: 0.1% FA in MeOH; 0.73mL/min; gradient: 20-90%B in 6 mins.

NEW High Stability Column (HSC™) Packing Technology
Significantly Improves UHPLC Column Robustness

ACE[®] Excel[™] UHPLC Columns – Scalability & Reproducibility



UHPLC

HPLC

Prep LC






ACE[®] Excel[™] Has Typically Lower Back Pressure For UHPLC


Specifically engineered for **lower UHPLC** backpressures

ACE Excel 2 C18-AR




P_{max}: 364 bar

Waters Acquity 1.7um BEH C18




P_{max}: 581 bar

Phenomenex Kinetex 1.7um C18



P_{max}: 540 bar

Agilent ZORBAX Eclipse 1.8 um XDB C18



P_{max}: 540 bar

Conditions: A = 5mM formic acid (aq); B = 5mM formic acid in MeOH; tg= 3 to 100%B in 5 min; 0.6 ml/min; 40C; 254nm


1. Paracetamol 2. Hydrochlorothiazide 3. Methylphenylsulphoxide 4. Methylphenylsulphone 5. Aspirin 6. Phenacetin 7. 1,3-dinitrobenzene 8. 1,2,4-Trimethoxybenzene

9. Ethyl benzoate 10. Nimesulide 11. Ibuprofen 12. Indomethacin 13. Mefenamic acid


All trademarks are recognised...comparative separations may not be representative of all applications





Selectivity Changes by Changing Organic Modifier



5) benzoic acid, 6) myrecetin, 7) p-cresol, 8) propranolol, 9) ethylparaben, 10) furosemide, 11) anisole, 12) 1,3,5-trinitrobenzene, 13) toluene, 14) nimesulide, 15) mefenamic acid, 16) 1,2,3-trichlorobenzene









4) phenacetin, 5) p-cresol, 6) acetophenone, 7) quercetin, 8) carvedilol, 9) 1,3-dinitrobenzene, 10) tolmetin 11) 3,4-dichlorobenzoic acid, 12) ethyl benzoate 13) toluene 14) ibuprofen 15) indomethacin, 16) 1,2-dichlorobenzene






Gradient Column Phase Screening: Optimum selectivity + UHPLC efficiency




0–80% ACN/0.1% HCOOH, 10 mM NH₄COO, 12-minutes





Selectivity Changes by Changing Organic Modifier



5) benzoic acid, 6) myrecetin, 7) p-cresol, 8) propranolol, 9) ethylparaben, 10) furosemide, 11) anisole, 12) 1,3,5-trinitrobenzene, 13) toluene, 14) nimesulide, 15) mefenamic acid, 16) 1,2,3-trichlorobenzene



4) phenacetin, 5) p-cresol, 6) acetophenone, 7) quercetin, 8) carvedilol, 9) 1,3-dinitrobenzene, 10) tolmetin 11) 3,4-dichlorobenzoic acid, 12) ethyl benzoate 13) toluene 14) ibuprofen 15) indomethacin, 16) 1,2-dichlorobenzene




Selectivity, Speed & Scaling
Isocratic & Gradient HPLC → UHPLC

ACE[®] Excel[™] C18-PFP Selectivity & Throughput (Isocratic)


Aim: obtain $R_s \geq 1.7$ in shortest possible time for mixture

Waters XBridge 5 μ m C18
 150 x 4.6 mm
 1.00 ml/min
 163 bar



HPLC: 5 μ m C18
13.6 mins
 (L/d_p = 3)

Waters Acquity 1.7 μ m BEH C18
 50 x 2.1 mm
 0.21 ml/min
 246 bar




UPLC: < 2 μ m C18
4.75 mins
 (L/d_p = 2.9)

To maintain R_s and reduce run time, keep L / d_p ratio constant

~ x23 Quicker

ACE Excel 2 μ m C18-PFP
 30 x 2.1 mm
 1.30 ml/min
 492 bar



UHPLC: ACE Excel 2 μ m C18-PFP
< 1 min

Using UHPLC and selectivity, it is possible to dramatically improve resolution allowing shorter columns & increased flow rates


Sample: 1) 1,2-dimethoxybenzene, 2) 1,3-dimethoxybenzene, 3) 1,3,5-trimethoxybenzene, 4) toluene (reference).
 Mobile phase 50:50 MeOH / H₂O; Temperature 40°C; 254 nm

All trademarks are recognised...comparative separations may not be representative of all applications


ACE[®] Excel[™] C18-PFP Selectivity & Throughput (Gradient)

Aim: obtain $R_s \geq 1.7$ in shortest possible time for mixture


ACE 5 μ m C18
 100 x 4.6 mm
 1 ml/min, $t_G = 29$ min
 max pressure: 92 bar
 40 min cycle time



ACE Excel 2 μ m C18
 50 x 2.1 mm
 0.6 ml/min, $t_G = 5$ min
 max pressure: 367 bar
 9 min cycle time




ACE Excel 2 μ m C18-PFP
 30 x 2.1 mm
 2.5 ml/min, $t_G = 0.7$ min
 max pressure: 914 bar
 1 min cycle time



~ x25 Quicker



$R_s \geq 1.7$




1, aspirin; 2, phenacetin; 3, 1,3-dinitrobenzene ; 4, ethyl benzoate; 5, nimesulide; 6, ibuprofen; 7, indomethacin.

Pressure Effects
HPLC ↔ UHPLC

Background

- ◆ Pressure is a **complex** physical parameter that affects **many elements** of a chromatography system
- ◆ Chromatographic **selectivity** and **retention** changes at **elevated pressures** have been investigated and reported^a
- ◆ Observations are **highly dependent** upon the analytes and may be seen with any manufacturer phases operated under UHPLC conditions
- ◆ Changes are typically **not helpful** for HPLC ↔ UHPLC activities

Effect of Pressure on Selectivity and Retention Factor



Initial **2 μ m** and **3 μ m** data are **similar** (A, B)

◆ **Scalability** looks **good**


Agilent 1290, 50 x 2.1 mm (constant flow and restrictor capillary used)

Mobile phase: A=0.1% FA in water: B=0.1% FA in MeOH (51:49 v/v)

Flow Rate: 0.21 ml/min, Temperature: 40 °C

K= Ketoprofen; S= Sulindac; N=Naproxen

Effect of Pressure on Selectivity and Retention Factor



Initial **2 μ m** and **3 μ m** data are **similar** (A, B)

◆ **Scalability** looks **good**

Retention and **selectivity** seen to **change** with pressure (B \rightarrow E)

Agilent 1290, 50 x 2.1 mm (constant flow and restrictor capillary used)

Mobile phase: A=0.1% FA in water: B=0.1% FA in MeOH (51:49 v/v)

Flow Rate: 0.21 ml/min, Temperature: 40 °C

K= Ketoprofen; S= Sulindac; N=Naproxen

Summary: Unwanted Selectivity Changes

- ◆ Pressure **induced k** and **α** changes may be seen for **any manufacturer** phases under UHPLC conditions
- ◆ Changes in selectivity and retention may be **significant** with **ionised analytes** and **large MW** analytes^a, but the impact on neutral molecules is **typically smaller**
- ◆ **Current** discussions / theory focus on changes in analyte **molar volume** as the **principle** cause for **changes** in **k** and **α** observed
- ◆ Successful HPLC ↔ UHPLC **possible**...the analyst just needs to be **vigilant**


Connections : Losses in N and A_s
Peak Dispersion

Background

- ◆ **UHPLC / optimised HPLC** instruments are **very sensitive** to the introduction of **extra column volume**
- ◆ Any time you **install** a column (from **any manufacturer**) it is vital to ensure **good** connections
- ◆ **Aim** for a '**fresh connection**' every time to ensure a **snug fit** between tubing and column and **reduce** the likelihood of an **unwanted** gap and / or tubing slippage
- ◆ **Free movement** of the **ferrule** and **nut** when installing the column gives you a **fresh connection**


Losses in Performance Due to Incorrect Column Fitting

Correctly fitted column



◆ **Correctly fitted columns make the most of your column and system**

Incorrectly fitted column



◆ **Incorrectly connected columns lead to reduced efficiency, reduced peak symmetry, and possibly, leaks**

◆ **Loss of ~23% for N**
Loss of ~11% for A_s

Summary: Column Connections

- ◆ Extra column volume **reduces** peak efficiency and **increases** asymmetry
- ◆ Make a **fresh connection every time** you install **any** column
- ◆ ACE recommends **reusable fittings** for a **fresh connection every time**
- ◆ All ACE® *Exce/™* columns have a **FREE** ‘Making Great UHPLC Connections’ leaflet in **every** box

Making Good UHPLC Column Connections (avoiding extra column volume and unwanted peak dispersion)

High resolution UHPLC systems are extremely sensitive to the introduction of extra column volume (dead volume), which will adversely affect your chromatography and result in unwanted peak dispersion.

Therefore, when installing any UHPLC column, it is vital to ensure that the high pressure inlet tubing is fitted into the column port to the correct depth to avoid the introduction of any extra column volume.

Figure (i) and figure (ii) below illustrate correctly and incorrectly fitted tubing.




Figure (i). Tubing correctly fitted – no gap between tubing and column frit – no extra volume introduced.




Figure (ii). Tubing incorrectly fitted (or tubing that has experienced slippage due to incorrect tightening) – gap between tubing and column frit – extra volume introduced.

Figure (iii) below illustrates the effect that correctly and incorrectly fitted tubing has on a typical UHPLC chromatogram. In this example it can be seen that a 0.5mm gap results in a 23% loss in efficiency and an 11% deterioration in asymmetry.




Figure (iii). Chromatogram illustrating the effect of correct and incorrect positioning of the inlet tubing, ferrule and nut. In this example a 0.5 mm gap results in a 23% loss in efficiency and a 11% deterioration in asymmetry.

Installing UHPLC columns correctly

For all UHPLC columns, correct installation is vital to get the most out of your column.

When you install your UHPLC column (and every time you change a UHPLC column) a ‘fresh connection’ should always be made between the tubing and the UHPLC connector. Avoid UHPLC connectors that have already been swapped onto the tubing, as these do not allow free movement between the tubing and the connector and may lead to a gap between the tubing and the column and the introduction of extra column volume (dead volume) into the system (as described previously).

For a ‘fresh connection’ every time, we recommend the use of ACE® UHPLC reusable column connectors (Part no. EXL-CC (10k) or EXL-CC10 (10µk)). These fittings are suitable for use with 1/16" OD tubing and virtually all manufacturers’ UHPLC systems and all brands of UHPLC columns. They are rated for use up to 25,000psi (1,720 bar) and are also suitable for use up to 100°C. Correctly used, these fittings will provide a ‘fresh connection’ for approximately 10 column installations.

The following simple guidelines will help you install your UHPLC column correctly.

Installation using ACE® UHPLC reusable column connectors (Part no. EXL-CC)

1. Slide the fitting onto the UHPLC tubing (approximately 5mm (1/4") from the end).
2. Insert the assembly into the receiving female port of the column, pushing the tubing in until it bottoms out.
3. Whilst maintaining pressure on the tubing (to ensure that it continues to ‘bottom out’), finger tighten the fitting into the column until snug.
4. Using the ACE® Torque Wrench (Part no. EXL-TW), tighten the fitting to the correct torque (while continuing to maintain pressure on the tubing).

Instructions for the correct use of the ACE® Torque Wrench are provided with the torque wrench.

A starter kit containing 1 x ACE® Torque Wrench and 4 x ACE® UHPLC reusable column connectors is available (Part no. EXL-CCSK).

Installation Using Other UHPLC Connectors

When installing UHPLC columns using other manufacturers’ UHPLC connectors, it is also vital that you ensure that the tubing bottoms out in the column port to avoid dead volume (as previously described). Please refer to the manufacturers’ instructions for further details.



Making Good UHPLC Column Connections


UHPLC column outlet end connections

ACE® UHPLC reusable column connectors (EXL-CC) are suitable for use at both the column inlet and the column outlet, to make good UHPLC connections. Alternatively, at the column outlet end an ACE® reusable PEEK finger-tight column connector (ACE-FT, suitable for use up to 5,000psi, 345 bar) can be used to make a good connection.

To avoid **extra column volume** and unwanted peak dispersion, it is additionally important to make good connections at the outlet end of your UHPLC column. Similar principles as previously described for the inlet end connection should be employed.

ACE® UHPLC Reusable Column Connectors

- Compact, one-piece, easy to use connector
- Pressure rated to 25,000psi (1,720 bar)
- Usable at elevated temperatures (up to 100°C)
- Fits virtually all manufacturers' UHPLC systems
- Fits all brands of UHPLC columns
- Reusable, non-permanent swaging design
- 10 make and break cycles*



Part Numbers

Description	Part No.
UHPLC column inlet / outlet connectors	
ACE® UHPLC reusable column connectors (1/pk) suitable for use up to 25,000psi (1,720 bar)	EXL-CC
ACE® UHPLC reusable column connectors (10/pk) suitable for use up to 25,000psi (1,720 bar)	EXL-CC10
ACE® UHPLC reusable column connectors starter kit (contains 1 x Part no. EXL-TW and 4 x Part no. EXL-CC)	EXL-CCSK
ACE® Torque Wrench	EXL-TW

Description	Part No.
UHPLC column outlet connectors	
ACE® reusable PEEK finger-tight column connector (1/pk) (suitable for use up to 5,000psi, 345 bar)	ACE-FT
ACE® reusable PEEK finger-tight column connector (10/pk) (suitable for use up to 5,000psi, 345 bar)	ACE-FT10

Also downloadable from the ACE website:



www.ace-hplc.com

Overall Summary & Conclusions




- ◆ Understanding the **properties** of building blocks in stationary phase design led to these **unique ACE[®]** products
- ◆ **ACE[®] C18-AR** and **ACE[®] C18-PFP** are powerful tools for method development due to **unique** but **complementary** selectivities
- ◆ These **unique** phases are available for HPLC as the **ACE[®]** range and also for UHPLC as the **NEW ACE[®] Excel[™] 2 μm** format
- ◆ These phases **meet** analyst **demands** of **reproducibility**, **robustness** & **low** phase **bleed** with **excellent** peak efficiency
- ◆ Operating at **high pressures** can deliver **excellent results** but remain vigilant - **selectivity** and **retention** may be affected...and even **column connections** become critical!

Full Information On All ACE Products Available


Unique Selectivities



Free Guides



**MACMOD Analytical = <http://www.mac-mod.com/>
ACT = <http://www.ace-hplc.com>**



Thank You For Your Attention

twaghe@mac-mod.com
info@mac-mod.com

www.mac-mod.com



ACE[®] is a registered trademark of Advanced Chromatography Technologies Ltd. ACE Excel[™] and HSC[™] are trademarks of Advanced Chromatography Technologies Ltd.

UPLC, Xbridge, ACQUITY are trademarks of Waters Corporation; ZORBAX, Eclipse are trademarks of Agilent Technologies Inc.; Kinetex is a trademark of Phenomenex Inc.; GOLD is a trademark of Thermo Fisher Scientific.