


## Magic additives

The majority of all separations performed on the **CHIRAL-AGP** column have been performed using simple mobile phases consisting of a buffer and an uncharged organic modifier (2-propanol, acetonitrile etc.). However, for some compounds, the addition of a charged modifier may induce or improve the enantioselectivity and/or improve the separation efficiency and the resolution. Examples of charged modifiers that have been used are octanoic, hexanoic and heptanoic acid, N,N-dimethyloctylamine (DMOA), tetraethyl- and tetrapropyl-ammonium bromide. Advice on when to use different types of charged modifiers can be found in the Method Development Scheme shipped with each column.

A very illustrative example on the effect that can be obtained with charged modifiers is the addition of **DMOA** to the mobile phase when chromatographing naproxen.



Column: **CHIRAL-AGP** 100x4.0 mm

Mobile phase: 10 mM sodium phosphate buffer pH 7.0


Flow rate: 0.9 ml/min

Detection: UV 225 nm

Note the extreme effect on the second eluted enantiomer of naproxen.

**DMOA** strongly affects the enantioselectivity of all the profens and other acidic compounds, however the most dramatic effect is obtained for naproxen. Another example is ketoprofen, where **DMOA** is essential for the separation:

with 2.5 mM DMOA




Column: **CHIRAL-AGP** 100x4.0 mm

Mobile phase: 10 mM sodium phosphate buffer pH 7.0

Flow rate: 0.9 ml/min Detection: UV 225 nm

**DMOA** can also be used to improve the chromatography of basic compounds on **CHIRAL-AGP**. An example is shown below:



Column: **CHIRAL-AGP** 100x4.0 mm

Mobile phase: 16% acetonitrile and 1 mM DMOA in 10 mM ammonium acetate

Flow rate: 0.8 ml/min

Detection: UV 225 nm

Without the addition of **DMOA** to the mobile phase the separation of clopidogrel is very poor or non-existing, depending on the mobile phase.

