


ANALYSIS OF EPIRUBICIN AND METABOLITES IN HUMAN PLASMA

MICRA NPS[®] is a breakthrough in fast HPLC. *NPS* is ultra-pure, highly uniform non porous silica spheres which provide the LC chromatographer greatly improved mass transfer and lower detection limits. Coupled with enhanced stability and dramatically reduced solvent usage, *NPS* is the ideal column to meet the ever increasing demands placed on today's analytical labs - Improved productivity at a lower cost.

This paper was originally presented at the 5th International Symposium on Pharmaceutical and Biomedical Analysis, in Stockholm, Sweden, September 21-24, 1994 by Professor K.S. Boos and Anne Rudolphi (Institute für Klinische Chemie Ludwig-Maximilians-Universität München).

Injection of biological samples, such as serum or plasma, onto an analytical column packed with porous or non-porous RP-bonded silica increases back pressure and causes a loss of retentive power. This is due to irreversible adsorption and denaturation of matrix proteins on the silica surface. This denaturation by the macromolecules of biological samples can be avoided by a coupled-column HPLC system consisting of a pre-column and the analytical column investigated. The pre-column serves as an on-line sample processing unit that serves two major functions: (a) dynamic extraction and concentration of the analytes of interest by hydrophobic interaction, and (b) eliminates the matrix proteins by size-exclusion. The advantages of using the MICRA *NPS* 1.5 μ ODS-I analytical column are a reduction in analysis time, improved LOQ and LOD, and enhanced efficiency.

Figure 1. Single-column Separation of Epirubicin and Metabolites



20 μ L Aqueous Solution:

Epirubicin (1)	62 ng/mL
Epirubicinol-Aglycon (2)	164 ng/mL
Epirubicinol (3)	56 ng/mL
Epirubicin-Aglycon (4)	89 ng/mL
7-Deoxy-epirubicinol-Aglycon (5)	40 ng/mL

Operating Conditions - Fluorescence Detector; ex: 445nm, em: 560nm

	MICRA NPS ODS-I	LiChrospher	Superspher
Column	33x4.6mm ID	250x4mm ID	250x4mm ID
Aqueous/Organic	85/15	70/30	70/30
Flow Rate	0.7 mL/min	1.0mL/min	1.0mL/min
N-alkyl-ligand	C-18	C-8	C-8
Injection Volume (all)	20 μ L		
Mobile Phase (all)	H ₂ O (0.1% TEA, pH 2.0 w/ TCA), Acetonitrile (85/15, v/v)		

Think small

Think fast

Think **NPS**[®]

Compared to Porous Column Supports, **NPS** Showed:

- 2x Increased Efficiency
- 40% Higher LOD & LOQ
- 50% Reduction in Analysis Time

Figure 2. Structures of Epirubicin and Metabolites




Figure 3. Coupled-column Separation of Spiked Human Plasma


20 μ L Human Plasma Spiked with:

Epirubicin (1)	31 ng/mL
Epirubicinol-Aglycon (2)	82 ng/mL
Epirubicinol (3)	28 ng/mL
Epirubicin-Aglycon (4)	44.5 ng/mL
7-Deoxy-epirubicinol-Aglycon (5)	20 ng/mL

MICRA NPS ODS-I

LiChrospher 60 RP-Select B

Superspher 60 RP-Select B



Eprogen
8205 S. Cass Avenue
Suite #106
Darien, Illinois 60561
USA
Phone: (630) 963-1481
Fax: (630) 963-6432
E-mail: info@eprogen.com
www.eprogen.com

Coupled-column HPLC analysis consisting of a restricted access "Alkyl-Diol" silica (ADS) pre-column followed by the analytical column. Further information regarding this analysis may be obtained by contacting Eprogen, Inc.