


## Separation of 17 Explosives on HALO C18, 2.7 µm Column



### PEAK IDENTITIES:

1. HMX
2. RDX
3. 1,3,5-Trinitrobenzene
4. 1,3-Dinitrobenzene
5. 3,5-Dinitroaniline
6. Nitrobenzene
7. Nitroglycerin
8. Tetryl
9. 2,4,6-Trinitrotoluene
10. 2-Amino-4,6-Dinitrotoluene
11. 4-Amino-2,6-Dinitrotoluene
12. 2,4-Dinitrotoluene
13. 2,6-Dinitrotoluene
14. 2-Nitrotoluene
15. 4-Nitrotoluene
16. 3-Nitrotoluene
17. PETN (pentaerythritol tetranitrate)

### TEST CONDITIONS:

Column: 4.6 x 150 mm, HALO C18, 2.7 µm

Part Number: 92814-702

Mobile Phase: A= water, B= Methanol

Gradient: Time	%B
0.0	25
14.0	35
20.0	62

Flow Rate: 1.5 mL/min.

Pressure: 366 bar to start, maximum: 405 bar

Temperature: 43°C

Detection: UV 220 nm, VWD

Injection Volume: 40 µL

Sample Solvent: 50/50: Water/methanol

Response Time: 0.02 sec.


Data rate: 25 Hz

Flow Cell: 2.5 µL semi-micro


LC System: Shimadzu Prominence UFLC XR

ECV: ~14 µL


### STRUCTURES:




HMX




Nitroglycerin



RDX



Tetryl




Pentaerythritol Tetranitrate

The determination of explosives in the environment is outlined in EPA method 8330B and under the conditions recommended, requires two column phases to determine 17 compounds. However, all 17 explosive compounds can be separated on a HALO C18, 2.7 µm column in less than 20 minutes using a water/methanol gradient.



## HPLC Separation of Explosives: Comparison of HALO<sup>®</sup> to a Competitor SPP Column

238-EX



### PEAK IDENTITIES

1. HMX
2. RDX
3. 1,3,5-Trinitrobenzene
4. 1,3-Dinitrobenzene
5. 3,5-Dinitroaniline
6. Nitrobenzene
7. Nitroglycerin
8. Tetryl
9. 2,4,6-Trinitrotoluene
10. 2-Amino-4,6-dinitrotoluene
11. 4-Amino-2,6-dinitrotoluene
12. 2,4-Dinitrotoluene
13. 2,6-Dinitrotoluene
14. 2-Nitrotoluene
15. 4-Nitrotoluene
16. 3-Nitrotoluene
17. PETN (pentaerythritol tetranitrate)

### TEST CONDITIONS:

**Column:** HALO 90 Å C18, 2.7 μm, 4.6 x 150 mm

**Part Number:** 92814-702

**Competitor Column:** SPP 100 Å C18, 2.6 μm, 4.6 x 150 mm

**Mobile Phase A:** Water

**Mobile Phase B:** Methanol

Gradient:	Time	%B
	0.0	25
	14.0	35
	20.0	62

**Flow Rate:** 1.5 mL/min

**Initial HALO<sup>®</sup> Back Pressure:** 441 bar

**Initial Competitor Back Pressure:** 490 bar

**Temperature:** 43°C

**Detection:** 220 nm

**Injection Volume:** 0.2 μL

**Sample Solvent:** Methanol

**Data Rate:** 100 Hz

**LC System:** Shimadzu Nexera X2


The determination of explosives in the environment is outlined in EPA method 8330B. 17 explosive compounds are separated on a HALO 90 Å C18 column in less than 20 minutes using a water/methanol gradient. These compounds are either used in the manufacture of explosives or propellants. The impurities or degradation of these compounds could be found in water, soil, or sediment samples. Baseline resolution is maintained on the HALO<sup>®</sup> column while there are peak co-elutions with a similar superficially porous particle column.





## UHPLC Separation of Explosives on 2 µm HALO® C18

239-EX



### PEAK IDENTITIES

1. HMX
2. RDX
3. 1,3,5-Trinitrobenzene
4. 1,3-Dinitrobenzene
5. 3,5-Dinitroaniline
6. Nitrobenzene
7. Nitroglycerin
8. Tetryl
9. 2,4,6-Trinitrotoluene
10. 2-Amino-4,6-dinitrotoluene
11. 4-Amino-2,6-dinitrotoluene
12. 2,4-Dinitrotoluene
13. 2,6-Dinitrotoluene
14. 2-Nitrotoluene
15. 4-Nitrotoluene
16. 3-Nitrotoluene
17. PETN (pentaerythritol tetranitrate)

### TEST CONDITIONS:

**Column:** HALO 90 Å C18, 2 µm, 3.0 x 100 mm

**Part Number:** 91813-602

**Mobile Phase A:** Water

**Mobile Phase B:** Methanol

Gradient:	Time	%B
	0.0	25
	6.9	35
	9.9	62

**Flow Rate:** 0.85 mL/min

**Initial Back Pressure:** 571 bar

**Temperature:** 43°C

**Detection:** 220 nm

**Injection Volume:** 0.2 µL

**Sample Solvent:** Methanol

**Data Rate:** 100 Hz

**LC System:** Shimadzu Nexera X2

The determination of explosives in the environment is outlined in EPA method 8330B. 17 explosive compounds are separated on a HALO 90 Å 2 µm C18 column in less than 10 minutes using a water/methanol gradient. These compounds are either used in the manufacture of explosives or propellants. The impurities or degradation of these compounds could be found in water, soil, or sediment samples.

