


Kromasil application guide



Kromasil®

*The way to peak performance
in liquid chromatography*

Kromasil applications – from our lab and the literature

The Kromasil packings and columns have been used over the years for demanding separations all over the world. In this guide we have collected examples of a variety of chromatographic separations, from small synthetic pharmaceuticals, up to peptides and larger molecules.

We hope this guide will be a useful tool when developing new HPLC separation methods in your lab.

Content

Guidelines	3
Applications	6
Amino acids	6
Drugs and metabolites	9
Environmental	23
Food and nutrition	26
Natural products	28
Peptides	31
Vitamins	33
Other	34
Literature references	41
Substance index	44
Kromasil availability	47

Cover figure shows a 3D chromatographic separation, with UV absorption at different wavelengths vs. time.

Copyright

Copyrighted elements of this publication may not be reproduced in any way without the consent of their respective owner.

The Kromasil packings

– some hints for the best performance

The Kromasil family of packing materials is developed to be the perfect choice from analytical to process scale. Kromasil is presently available as bare silica, C4, C8, C18, NH2, CN or as Kromasil Chiral for separation of optical isomers. Pore sizes are 60 Å, 100 Å, and 300 Å, and particle sizes 3.5, 5, 7, 10, 13 and 16 µm. Slurry-packed columns are available from analytical up to 2" inner diameter, all with analytical efficiencies. For larger preparative and industrial scale columns bulk packing is provided.

To learn more about the properties of Kromasil silica please consult our other technical information.

Choice of mobile phase

Normal Phase conditions

Choose mixtures of hexane or heptane, and polar modifiers like alcohols, ethyl acetate, methylene chloride, etc. to adjust retention. The optimum retention factor range is normally $2 \leq k \leq 5$ for a two-component sample, but can be wider for a multi-component sample.

Acidic and basic additives can improve the chromatographic performance. In most cases small amounts of acetic acid or formic acid (0.05 – 0.10%) improve peak shape for acidic or basic solutes. In some cases the combination of acid and an organic amine (e.g. triethylamine) is necessary, for difficult

solutes. The acid should always be in excess relative to the amine, in order to operate at a pH where the silanol groups on the silica are protonated.

Reversed Phase conditions

Choose mixtures of water or buffer, and water miscible solvents like alcohols, acetonitrile, THF, etc. to adjust retention factor k to an optimal range. The pH can be controlled by using a buffer, and in order to minimize the ionization of the silica and the solutes. In order to control peak shape for very basic solutes an additive like TEA (triethylamine) can be added, if necessary. Kromasil is a fully hydroxylated ultra-pure silica, making the surface less acidic, resulting in good peak shape also for basic compounds.

For the C4, C8 and C18 phases, due to the very hydrophobic nature of the surface, it is important to **always keep at least 4 – 5% of organic in the mobile phase**, both when flushing or running the chromatographic separation. The reason is that in the case of a 100% aqueous mobile phase there is a risk that the surface within the porous system in the Kromasil particles is “dewetted”, resulting in a total or partial loss of retention.

This phenomenon of dewetting is more pronounced for high quality, high coverage materials, where the bonding procedure has been very efficient. This will result not only in higher retention times because




Figure 1 | Influence of the mobile phase additives on the separation of mefloquine.

Conditions: Column: 4.6 × 250 mm, Kromasil CHI-DMB, 5 µm
Flow rate: 2 ml/min. Detection: UV 280 nm




Figure 2 | The retention factor, k , vs. organic content. General retention behaviour at low organic content using high density RP phases.

Column length (mm)	Particle size (μm)	Flow rate (ml/min.)	Relative time	Relative R_s	Relative ΔP
250	10	0.5	1	1	1
125	5	1.0	0.25	1	4
87.5	3.5	1.43	0.125	1	8

Table 1 | Relation between optimal flow rate, analysis time and back pressure, ΔP , when going to smaller particles, and shorter columns. The comparisons are for a constant resolution, R_s .

of a higher coverage and hence hydrophobicity, but also in a higher hydrolytic stability, and a longer life-time of the column.

If a 100% aqueous mobile phase has been used accidentally, and the stationary phase has been dewetted, the column can easily be regenerated. Just flush the column with a mobile phase consisting of 40 – 50% or more of organic for 2 – 5 column volumes. After this the column can be equilibrated again with the mobile phase, and the original retention times should be seen.

We also recommend to always filter buffer solutions in order to remove small particulates. It is also preferable to premix aqueous/organic solutions, in order to avoid problems with gas formation in the mobile phase, or a temperature increase or decrease as an effect of endo- or exothermic mixing heats.

The recommended pH range for our RP phases is between pH 1.5 up to 9.5. However, in some applications mobile phases with pH above 11 have been used for continuous chromatography, for several thousands of column volumes.

How to improve speed of separation

There is today a strong driving force towards faster separations, and hence smaller particles and shorter columns. A smaller particle will give a higher efficiency at a higher flow rate; for example will a 5 μm particle give twice the efficiency compared to 10 μm , at twice the flow rate. And if the resolution is to be kept constant one can also reduce column length by 50%. Table 1 gives the relation between the critical parameters when going to smaller particles, and shorter columns.

It can be seen that the combination of smaller particles, shorter columns and higher flow rate will result in much faster analyses. The only drawback is the back pressure, which will increase significantly as particle size goes down.

All in all, one will save a factor 2 in analysis time by going from 5 to 3.5 μm for example, but also experience twice the back pressure.

How to improve resolution

A good resolution in a short period of time is usually a requirement in analytical HPLC. One has essentially three ways to improve the resolution, as can be seen in the equation below:

$$R_s = \frac{1}{4} (\alpha - 1) \cdot \sqrt{N} \cdot \left(\frac{k_1}{1 + k_1} \right)$$

1. **The separation factor, α , can be increased.** This can be done by optimizing stationary and mobile phase, i.e. choosing the best column and mobile phase composition for the specific application
2. **The number of theoretical plates can be increased.** This can be done by:
 - a. Increasing the column length
 - b. Decreasing the particle size
 - c. Optimizing the flow rate. The optimum for small particles is at a higher flow rate than for larger particles (inverse linear relationship)
3. **The retention factor, k , can be increased,** if it is too low. This can be done by adjusting the elution strength of the mobile phase

Scale-up

The analytical separation is very often the start for a scale-up to semi-prep, prep, or large industrial scale chromatography. In the case of Kromasil there is the possibility of seamless scale-up from analytical to semi-prep and prep, and even large diameter Dynamic Axial Compression (DAC) columns. All Kromasil columns independent of diameter are delivered with the same, high efficiency guarantee, and even the large DAC columns can be packed giving the same performance.

We recommend that the method development for the preparative separation is performed using analytical columns, or possibly 10 mm ID if larger volumes of the sample fractions are needed. A small column will make the development work easier, and since the performance is identical in small and large diameter columns, the result in large scale can easily be predicted from the work in small scale. 10 μm particles are recommended, since the efficiency, back pressure, etc. then will be close to the large scale separation. For large scale 10 – 16 μm particles are usually a good choice. However, the performance using a different particle size can easily be predicted.

The optimal running conditions in large industrial scale can also be found by applying a software program, KromaGuide, developed by the Kromasil group. The KromaGuide optimization is part of our technical support to current Kromasil customers, or potential users.




Figure 3 | Scale-up to an 80 cm ID dynamic axial compression (DAC) column, showing analytical performance. Efficiency was 42,000 plates/meter ($h = 2.38$) and asymmetry_{0.1} was 1.19


Conditions: Packing material: 83 kg Kromasil 10 μm Bed length: 25 cm
Eluent: acetonitrile/water (7/3) Flow rate: 20 lit./min.
Sample: uracil and toluene

Applications

Amino acids


Amino acids, PTC derivatives

18 amino acids as phenylthiocarbamyl (PTC) derivatives. (ref. 7)




Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.6 × 200 mm
 Temperature: 38 °C
 Eluent A: 3% ACN in 0.1 M sodium acetate
 Eluent B: ACN:water (80:20; v:v).


1 = PTC-aspartic acid



2 = PTC-glutamic acid




3 = PTC-serine




4 = PTC-glycine




5 = PTC-histidine



6 = PTC-arginine



7 = PTC-threonine




8 = PTC-alanine



9 = PTC-proline




10 = PTC-tyrosine




11 = PTC-valine




12 = PTC-methionine




13 = PTC-cysteine




14 = PTC-isoleucine




15 = PTC-leucine




16 = PTC-norleucine (l. S.)




17 = PTC-phenylalanine



18 = PTC-tryptophan



19 = PTC-lysine




R = PTC derivative group



Gradient: Linear gradient elution. 0 min. 0% B, 13 min. 7% B, 23 min. 23% B, 29 min. 35% B, 35 min. 40% B, 40 min. 100% B, 45 min. 100% B, 47 min. 0% B
 Flow rate: 1 ml/min.
 Detection: UV 254 nm


Amino acids, Fmoc-derivatives

Amino-acid analysis for protein and peptide hydrolysates with precolumn Fmoc (9-fluorenyl methylchloroformate) derivatization. (ref. 30)




Phase: Kromasil 100 Å, 5 µm, C8
 Column: 4 × 250 mm
 Temperature: 45 °C
 Eluent A: sodium acetate buffer (100 mM, pH 4.4):THF:ACN (75:15:10; v:v:v)
 Eluent B: ACN:THF (85:15; v:v).
 Gradient: 0 – 2.5 min. 0% B, 2.5 – 6.6 min. 7% B, 6.6 – 8.3


1 = Fmoc-arginine




2 = Fmoc-serine




3 = Fmoc-aspartic acid




4 = Fmoc-glutamic acid




5 = Fmoc-threonine



6 = Fmoc-glycine




7 = Fmoc-alanine




8 = Fmoc-proline




9 = Fmoc-tyrosine




10 = Fmoc-methionine




11 = Fmoc-OH




12 = Fmoc-NH₂




13 = Fmoc-valine




14 = Fmoc-phenylalanine




15 = Fmoc-tryptophan




16 = Fmoc-cystine




17 = Fmoc-isoleucine




18 = Fmoc-leucine




19 = Fmoc-histidine




20 = Fmoc-cysteine



21 = Fmoc-lysine



R = Fmoc derivative group




min. 14% B, 8.3 – 8.4 min. 21% B, 8.4 – 10 min. 21% B, 10 – 10.1 min. 17% B, 10.1 – 20 min. 19% B, 20 – 29 min. 55% B, 29 – 30 min. 100% B
 Flow rate: 1.5 ml/min.
 Detection: UV 263 nm

Amino acids

Amino acids


Detection of N-acetylaspartate and N-acetylglutamate in cerebral tissue extracts. (ref. 228)



Phase: Kromasil 100 Å, 5 µm, C18
Column: 4.6 × 250 mm
Temperature: 23 °C
Eluent: 2.8 mM tetrabutylammonium hydroxide,
25 mM KH₂PO₄, 1.25% MeOH (pH 7)
Flow rate: 1 ml/min.
Detection: UV 210 nm

Amino acids, N-acetylated


Separation of N-acetylated amino acids. (ref. 348)



Phase: Kromasil 100 Å, 5 µm, C18
Column: 4.6 × 250 mm
Temperature: 23 °C
Eluent: tetrabutylammonium hydroxide 2.8 mM;
KH₂PO₄ 25 mM and 1.25% MeOH, pH 7
Flow rate: 1 ml/min.
Detection: UV 210 nm

Amino acids, benzoylated

Analysis of benzoylated amino acids. (ref. 51a)



Phase: Kromasil 100 Å, 5 µm, C18
Column: 4 × 250 mm
Eluent: acetonitrile-water mixtures
Gradient: 70 – 95% ACN in 30 min.
Flow rate: 1 ml/min.
Detection: UV 274 nm

Amino acids


Aminosalicylic acids

Determination of 5-aminosalicylic acid and 3-aminosalicylic acid. (ref. 279)

1 = 5-aminosalicylic acid



2 = 3-aminosalicylic acid




Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.6 × 200 mm
 Eluent: MeOH:phosphate buffer (35:65; v:v)
 Flow rate: 1 ml/min.
 Detection: UV 254 nm


Carnitines, aminoanthracene derivatives

Determination of L-carnitine, acetyl-L-carnitine and propionyl-L-carnitine in human plasma by HPLC with post-column derivatization with 1-aminoanthracene. (ref. 66)


1 = L-carnitine 1-aminoanthraceneamide




2 = acetyl-L-carnitine 1-aminoanthraceneamide





3 = methansulfonyl-L-carnitine 1-aminoanthraceneamide



4 = propionyl-L-carnitine 1-aminoanthraceneamide



5 = isobutyryl-L-carnitine 1-aminoanthraceneamide



Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.6 × 250 mm
 Eluent: ACN:ammonium acetate (0.1 M, pH 3.5) (30:70; v:v)
 Flow rate: 1.3 ml/min.
 Detection: spectrofluorimetric (λ_{ex} 248 nm, λ_{em} 418 nm)

Boronophenylalanine

Determination of boronophenylalanine in biological samples after precolumn derivatization with o-phthalaldehyde (OPA). (ref. 237)

1 = OPA-aspartic acid



2 = OPA-glutamic acid



3 = OPA-asparagine



4 = OPA-histidine



5 = OPA-serine




6 = OPA-glutamine



7 = OPA-arginine



8 = OPA-citrulline



9 = OPA-glycine



10 = OPA-threonine




11 = OPA-γ-aminobutyric acid (GABA)




12 = OPA-alanine



13 = OPA-tyrosine



14 = OPA-p-boronophenylalanine



15 = OPA-valine



16 = OPA-phenylalanine



17 = OPA-isoleucine



18 = OPA-leucine




19 = OPA-ornithine



20 = OPA-lysine



R = OPA derivative group



Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.6 × 250 mm
 Temperature: 23 °C
 Eluent A: 50 mM CH₃COONa (pH 7.4) : 50 mM NaHPO₄ (pH 7.4) : MeOH : THF (48:48:2:2; v:v:v:v)
 Eluent B: MeOH:water (65:35; v:v).


Gradient: 80% A in 3 min, 80% – 70% A in 12 min, 70% – 50% A in 15 min, 50% – 45% A in 10 min, 45% – 20% A in 10 min, 20% – 15% A in 5 min, 15% – 10% A in 3 min, 10% – 0% A in 2 min, 0% A in 15 min.
 Flow rate: 1.2 ml/min.
 Detection: spectrofluorimetric (λ_{ex} 330 nm, λ_{em} 430 nm)

Drugs and metabolites



Amoxicillin

Measurement of amoxicillin in gastric tissue samples. (ref. 6)

1 = amoxicillin degradation derivative



2 = ampicillin degradation derivative (I. S.)




Phase: Kromasil 100 Å, 5 µm, C18
 Column: 3.2 × 150 mm
 Temperature: 40 °C
 Eluent: MeOH-water (55:45; v:v)
 Flow rate: 0.4 ml/min.
 Detection: fluorescence (λ_{ex} 365 nm, λ_{em} 445 nm)


Amperozide

Separation of amperozide, derivate and metabolite. (ref. 45)



1 = amperozide's N-de-ethyl metabolite (II)



2 = amperozide (I)



3 = amperozide's N-de-ethyl-N-butyl analogue (III)




Phase: Kromasil 100 Å, 5 µm, C18
 Column: 2.1 × 200 mm
 Eluent: MeOH:ammonium phosphate buffer (pH 7.8) (78:22; v:v)
 Flow rate: 0.2 ml/min.
 Detection: UV 265 nm


Anthocyanidins

Separation of cyanidin from 3-O-β-glycosylated anthocyanidins. (ref. 347)


1 = petunidin-3-O-β-glycoside




4 = malvidin-3-O-β-glycoside




7 = cyanidin




2 = delphinidin-3-O-β-glycoside





5 = cyanidin-3-O-β-glycoside



3 = pelargonidin-3-O-β-glycoside



6 = peonidin-3-O-β-glycoside




Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.6 × 250 mm
 Temperature: 23 °C
 Eluent A: HCOOH:water (1:10; v:v)
 Eluent B: HCOOH:water:MeOH (1:9:10; v:v:v)
 Gradient: 0% – 60% A in 5 min., 60% – 45% A in 5 min., 45% – 0% A in 6 min., 0% A in 10 min.

Flow rate: 1.2 ml/min.
 Detection: 520 nm

Drugs and metabolites

Antibacterial drugs


Determination of metronidazole, clotrimazole and chlorhexidine acetate in Shuangzuo effervescent tablets. (ref. 23)



Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.6 × 250 mm
 Eluent: MeOH:buffer (70:30; v:v) (NaAc 24.4 g, HAc 80 ml, (C₄H₉)NBr 4.83 g in 1000 ml water, pH 3.6)
 Flow rate: 1 ml/min.
 Detection: UV 260 nm

Antibacterial drugs, veterinary


Simultaneous determination of sulfaquinoxaline, sulfamethazine and pyrimethamine. (ref. 246)



Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.6 × 150 mm
 Eluent: 40 mM phosphate buffer (pH 3 containing 10 mM ClO₄⁻) : ACN (65:35; v:v)
 Flow rate: 1.5 ml/min.
 Detection: UV 270 nm

Antibacterials, sulfa drugs


Determination of sulfamethoxypyridazine, sulfamethoxazole, sulfadimethoxine and associated compounds. (ref. 267)



Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.6 × 150 mm
 Eluent: 10 mM citrate buffer (pH 3):MeOH
 Gradient: 0 min. 31% MeOH, 4 min. 69% MeOH, 14 min. 69% MeOH, 16 min. 31% MeOH
 Flow rate: 1 ml/min.
 Detection: UV 255 nm

Antibiotics and intermediates

Determination of ceftriaxone, 7-aminocephalosporanic acid (7-ACA) and 7-amino-3-[[[2,5-dihydro-6-hydroxy-2-methyl-5-oxo-1,2,4-triazin-3-yl]-thio]methyl]-cephalosporanic acid (7-ACT). (ref. 129)




Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.6 × 200 mm
 Eluent: ACN:tetrabutyl ammonium bromide:phosphate buffer (pH 7):water (32:0.32:4.4:63.6; v:v:v:v)
 Flow rate: 1 ml/min.
 Detection: UV 270 nm


Drugs and metabolites

Antibiotics and metabolites

Determination of quinupristin and its main metabolites in human plasma. (ref. 143)




Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.6 × 125 mm
 Eluent A: 0.8 ml of 70% perchloric acid (PCA) / litre water
 Eluent B: ACN
 Gradient: 30% B for 11 min., 32% B from 11.1 to 15 min.,
 40% B from 15.6 to 16 min., 38% B from 16.1 to
 34 min., 80% B from 34.1 to 36 min.



Flow rate: 0 – 11 min: 0.5 ml/min., 11 – 36 min: 1 ml/min.
 Detection: fluorescence (λ_{ex} 360 nm and λ_{em} 410 nm)

Anticonvulsants

Determination of theophylline, phenobarbitone, diphenytoin and carbamazepine. (ref. 301b)




1 = theophylline
 2 = phenobarbitone
 3 = diphenytoin
 4 = carbamazepine

Phase: Kromasil 100 Å, 5 µm, C18
 Column: 0.8 × 150 mm
 Eluent: MeOH:water (70:30; v:v)
 Flow rate: 35 µl/min
 Detection: UV 210 nm

Antidepressants

Determination of antidepressant drugs and metabolites. (ref. 49)




1 = amoxapine
 2 = fluvoxamine
 3 = maprotiline
 4 = trimipramine
 5 = mianserine

Phase: Kromasil 100 Å, 5 µm, C18
 Column: 2.1 × 150 mm
 Eluent: ACN:phosphate buffer (40:60; v:v) (pH 6.5)
 Flow rate: 0.35 ml/min.
 Detection: UV 220 nm

Drugs and metabolites

Antidepressants


Analysis of amitriptyline and nortriptyline in plasma. (ref. 58)



Phase: Kromasil 100 Å, 5 µm, C8
 Column: 4 × 250 mm
 Temperature: ambient
 Eluent: ACN:KH₂PO₄ (0.04 M) (40:60; v:v)
 Flow rate: 1 ml/min.
 Detection: UV 240 nm

Antidepressants and metabolites

Simultaneous determination of citalopram, fluoxetine, paroxetine and their metabolites in plasma. (ref. 309)




Phase: Kromasil 100 Å, 3.5 µm, C18
 Column: 0.32 × 300 mm
 Temperature: gradient: 35 °C (3 min.) prior to ramp of 1.3 °C/min. to 100 °C (10 min.)
 Eluent: ACN:NH₄HCOO (45 mM, pH 4) (25:75; v:v)
 Flow rate: 5 µl/min
 Detection: UV 230 nm

Drugs and metabolites


Antifungals

Determination of terbinafine hydrochloride, chlorhexidine and triamcinolone acetonide acetate. (ref. 110)




Anti-HIV

Simultaneous determination of ritonavir and saquinavir. (ref. 126)




Antimicrobials

Determination of metronidazole and nalidixic acid. (ref. 156)



Antimicrobials


Determination of metronidazole and clindamycin. (ref. 268)



Drugs and metabolites

Antipsychotics


Determination of clozapine and loxapine. (ref. 64)



Phase: Kromasil 100 Å, 5 µm, C8
 Column: 4.6 × 150 mm
 Temperature: 31 °C
 Eluent: ACN:water (70:30; v:v) 25 mg ammonium acetate /100 ml mobile phase
 Flow rate: 1.4 ml/min.
 Detection: UV 210 nm

Catecholamines


Determination of catecholamines in pig liver. (ref. 95a)



Phase: Kromasil 100 Å, 5 µm, C8
 Column: 4.6 × 150 mm
 Eluent: 300 ml MeOH + 1.5 ml 1-octanesulfonic acid (200 mg/ml) + 100 ml 1 M NaAc + about 1 litre water (pH 3.8). Volume adjusted to 2 litres with water.
 Flow rate: 0.6 ml/min.
 Detection: electrochemical potential +0.65 V

Catecholamines


Determination of methoxycatecholamines in pig liver. (ref. 95b)



Phase: Kromasil 100 Å, 5 µm, C8
 Column: 4.6 × 150 mm
 Eluent: 300 ml MeOH + 1.5 ml 1-octanesulfonic acid (200 mg/ml) + 100 ml 1 M NaAc + about 1 litre water (pH 3.8). Volume adjusted to 2 litres with water.
 Flow rate: 1.1 ml/min.
 Detection: electrochemical potential +0.8 V

Ciprofloxacin

Determination of ciprofloxacin in pharmaceutical preparations and biological fluids. (ref. 26)




Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.6 × 250 mm
 Temperature: ambient
 Eluent: ACN:MeOH:acetate buffer (pH 3.6; 50 mM) (10:30:60; v:v:v) containing 1 % v/v HAc
 Flow rate: 0.8 ml/min.
 Detection: fluorescence (λ_{ex} 300 nm, λ_{em} 458 nm)

Drugs and metabolites



Clodronate

Simultaneous determination of clodronate and its partial ester derivative. (ref. 97)

1 = clodronate



2 = clodronate monophenylester




Phase: Kromasil 100 Å, 5 µm, C8
 Column: 4.6 × 250 mm
 Eluent: MeOH:ammonium acetate buffer (0.1 M + 0.23 M butylamine, pH 4.6)
 Gradient: linear gradient elution: methanol from 3 to 40 – 60% for between 1.0 and 6.0 min. (not specified)
 Flow rate: 1.2 ml/min.
 Detection: ELS



DRF-2189

Determination of the insulin sensitizing agent DRF-2189 in rat plasma. (ref. 161)

1 = insulin sensitizing agent DRF-2189



2 = troglitazone




Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.6 × 250 mm
 Eluent: 0.05 M NaH₂PO₄:ACN:MeOH (22.5:37.5:40; v:v:v) (pH 5.0)
 Flow rate: 1 ml/min.
 Detection: fluorescence (λ_{ex} 292 nm and λ_{em} 325 nm)


Cytochrome P450 metabolites

Analysis of cytochrome P450 metabolites of arachidonic acid. (ref. 10)


1 = 14,15-DHET




2 = 11,12-DHET



3 = 8,9-DHET



4 = 5,6-DHET




5 = 20-HETE




6 = 14,15-EET



7 = 11,12-EET




8 = 8,9-EET




9 = 5,6-EET



8 = 8,9-EET



9 = 5,6-EET



DHET = dihydroxy-


eicosatrienoic acids

HETE = hydroxy-

eicosatetraenoic acids

EET = epoxy-

eicosatrienoic acids




Phase: Kromasil 100 Å, 5 µm, C18
 Column: 2 × 250 mm
 Eluent: water/ACN with 0.005% HAc
 Gradient: 0 min. 60% ACN, 30 min. 80% ACN, 35 min. 100% ACN 40 min. 100% ACN
 Flow rate: 0.2 ml/min.
 Detection: ESI-MS


Ecstasy analogues

Identification of a homologue derivative of "ecstasy". (ref. 170)

1 = N-methyl-1-(1,3-benzodioxol-5-yl)-2-propanamine (MDMA)





2 = N-ethyl-1-(1,3-benzodioxol-5-yl)-2-propanamine (MDEA)



3 = N-methyl-1-(1,3-benzodioxol-5-yl)-2-butanamine (MBDB)



4 = N-ethyl-1-(1,3-benzodioxol-5-yl)-2-butanamine (EBDB)




Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.6 × 250 mm
 Temperature: ambient
 Eluent: ACN:0.1 M triethylammonium acetate (aq) pH 7.3
 Gradient: 5% to 80% ACN in 25 min.
 Flow rate: 1 ml/min.
 Detection: UV 280 nm

Drugs and metabolites

5-fluorouracil metabolites


Determination of the main metabolites of 5-fluorouracil in plasma. (ref. 116)



Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.6 × 150 mm
 Temperature: 20 °C (ambient)
 Eluent: MeOH:water (3:97; v:v)
 Flow rate: 0.6 ml/min.
 Detection: UV 275 nm

Fotemustine


Stability study of fotemustine in PVC infusion bags. (ref. 124)



Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.6 × 150 mm
 Temperature: ambient
 Eluent: ACN:ammonium acetate buffer (0.05 M, pH 4.5) (30:70; v:v)
 Flow rate: 1 ml/min.
 Detection: UV 230 nm

Ketoprofen


Determination of ketoprofen in vitro in rat skin. (ref. 247)



Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4 × 250 mm
 Temperature: 40 °C
 Eluent: ACN:0.01 M potassium phosphate (pH 1.5) (60:40; v:v)
 Flow rate: 1 ml/min.
 Detection: UV 260 nm

Leukotrienes, cross-reactive

Determination of cross-reactive leukotrienes in biological matrices. (ref. 71 a)




Phase: Kromasil 100 Å, 5 µm, C4
 Column: 2.1 × 100 mm
 Eluent: ACN:K₂HPO₄ 10 mM (pH 7.4) (30:70; v:v)
 Flow rate: 0.2 ml/min.
 Detection: fluorescence (λ_{ex} 544 nm, λ_{em} 572 nm)

Drugs and metabolites

Leukotrienes, cross-reactive


Determination of cross-reactive leukotrienes in biological matrices. (ref. 71b)



Phase: Kromasil 100 Å, 5 µm, C4
 Column: 2.1 × 100 mm
 Eluent: ACN:K₂HPO₄ 10 mM (pH 7.4) (30:70; v:v)
 Flow rate: 0.2 ml/min.
 Detection: fluorescence (λ_{ex} 544 nm, λ_{em} 572 nm)

Megazol


Analysis of megazol in human plasma. (ref. 113)



Phase: Kromasil 100 Å, 10 µm, C8
 Column: 4 × 250 mm
 Temperature: ambient
 Eluent: phosphate buffer (0.068 M, pH 3):MeOH:ACN (65:20:15; v:v:v)
 Flow rate: 0.7 ml/min.
 Detection: UV 360 nm

Meloxicam


Determination of meloxicam in human plasma. (ref. 283)



Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.6 × 150 mm
 Eluent: MeOH:water:ACN:HAc (600:500:50:20; v:v:v:v) + 1.01 g sodium heptanesulfonate
 Flow rate: 1 ml/min.
 Detection: UV 355 nm

Pain relievers

Determination of acetaminophen, ibuprofen and chlorzoxazone. (ref. 154)




Phase: Kromasil 100 Å, 5 µm, C8
 Column: 4.6 × 250 mm
 Temperature: 20 ± 1 °C
 Eluent: ACN:0.2% triethylamine (pH 3.2) (50:50; v:v)
 Flow rate: 1.5 ml/min.
 Detection: UV 215 nm

Drugs and metabolites

PAT-5A


Determination of PAT-5A (5[4-[N-(20pyridyl)-(2s)-pyrrolidine-2-methoxy]phenylmethylene]-thiazolidine-2,4-dione, maleic acid salt), an insulin sensitizing agent, in rat plasma. (ref. 244)



Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.6 × 250 mm
 Eluent: NaH₂PO₄ (0.05 M, pH 4):MeOH (25:75; v:v)
 Flow rate: 1 ml/min.
 Detection: UV 345 nm

Phenolics


Separation of phenolic compounds and corresponding glucuronides. (ref. 103)



Phase: Kromasil 100 Å, 5 µm, C18
 Precolumn: Nucleosil 5µm, C4
 Column: 4.6 × 100 mm (precolumn: 4.6 × 50 mm)
 Temperature: ambient
 Eluent: 30 mM cetyltrimethylammonium bromide in 0.05 M 6-aminohexanoic acid (pH: 5) and 20% ACN (precolumn 7%) (v:v)
 Flow rate: 1 ml/min.
 Detection: UV 300 nm

Phosphonoformate (foscarnet)

Determination of phosphonoformate (foscarnet) in human serum. (ref. 217)




Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.6 × 150 mm
 Eluent: methanol : 40 mM Na₂HPO₄-buffer, pH 7.6 (adjusted with orthophosphoric acid), containing 0.25 mM THAHSO₄ (25:75; v:v)
 Flow rate: 1 ml/min.
 Detection: electrochemical (potential +1.125 V)

Drugs and metabolites

QC test, tricyclic antidepressants


QC test of Kromasil CN. (ref. 342)



Phase: Kromasil 60 Å, 10 µm, CN
 Column: 4.6 × 250 mm
 Temperature: ambient
 Eluent: MeOH:KH₂PO₄ 25 mM pH 6.0 (80:20; v:v)
 Flow rate: 1 ml/min.
 Detection: UV 215 nm

QC test, tricyclic antidepressants


QC test of Kromasil C4. (ref. 349)



Phase: Kromasil 100 Å, 5 µm, C4
 Column: 4.6 × 250 mm
 Temperature: ambient
 Eluent: MeOH:KH₂PO₄ 25 mM pH 6.0 (80:20; v:v)
 Flow rate: 1 ml/min.
 Detection: UV 215 nm

QC test, tricyclic antidepressants


QC test of Kromasil C8. (ref. 350)



Phase: Kromasil 100 Å, 5 µm, C8
 Column: 4.6 × 250 mm
 Temperature: ambient
 Eluent: MeOH:KH₂PO₄ 25 mM pH 6.0 (80:20; v:v)
 Flow rate: 1 ml/min.
 Detection: UV 215 nm

QC test, tricyclic antidepressants

QC test of Kromasil C18. (ref. 351)




Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.6 × 250 mm
 Temperature: ambient
 Eluent: MeOH:KH₂PO₄ 25 mM pH 6.0 (80:20; v:v)
 Flow rate: 1 ml/min.
 Detection: UV 215 nm

Drugs and metabolites


Quinolones

Determination of quinolones in food. (ref. 119)




Salicin

Determination of salicin in extract of willow bark. (ref. 262)




Sildenafil

Determination of sildenafil (Viagra) and its metabolite (UK 103320) with ASTED equipment. (ref. 98)



Sildenafil


Determination of sildenafil citrate (Viagra). (ref. 254)



Drugs and metabolites


Spiramycin

Determination of spiramycin in pig liver. (ref. 94)




Phase: Kromasil 100 Å, 5 µm, C8
 Column: 4.6 × 250 mm
 Temperature: 60 °C
 Eluent: CH₃CN:sodium phosphate buffer (0.05 M pH 2.3)
 (33:67; v:v) + 6 g/l NaClO₄
 Flow rate: 1.1 ml/min.
 Detection: UV 232 nm

1, 2, 3, 4 and 6 =
 substituted base structure
 according to table




α-mycarose



timonacic



5 = spiramycin S




	R1	R2	R3
1. cysteyle neospiramycin I	H	H	timonacic
2. cysteyle spiramycin I	H	α-mycarose	timonacic
3. spiramycin I + cysteyle neospiramycin III	H	α-mycarose	COH
4. cysteyle spiramycin III	COCH ₂ CH ₃	H	timonacic
5. spiramycin S	COCH ₂ CH ₃	α-mycarose	timonacic
6. spiramycin III	COCH ₂ CH ₃	α-mycarose	COH



Steroids

Analysis of clenbuterol hydrochloride and ambroxol hydrochloride. (ref. 331)

1 = clenbuterol



2 = ambroxol




Phase: Kromasil 60 Å, 10 µm, CN
 Column: 4.6 × 250 mm
 Eluent: 1.8 g sodium decanesulphate + 3 g potassium phosphate monobasic + 600 ml water (pH 3.0) + 200 ml ACN + 200 ml MeOH
 Flow rate: 1.5 ml/min.
 Detection: UV 215 nm



Steroids

Analysis of dexamethasone and betamethasone acetate in bovine liver. (ref. 272a)

1 = dexamethasone



2 = betamethasone acetate




Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4 × 150 mm
 Eluent: MeOH:water (80:20; v:v)
 Flow rate: 1 ml/min.
 Detection: UV 240 nm

Drugs and metabolites

Tetracyclines


Determination of tetracyclines as chelates with aluminum(III). (ref. 273)



Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.6 × 250 mm
 Eluent: ACN:DMF:0.05 M citric acid-sodium citrate buffer (pH 2.5) (5:20:75; v:v:v)
 Flow rate: 0.7 ml/min.
 Detection: fluorescence (λ_{ex} 380 nm and λ_{em} 480 nm)

Tramadol


Determination of tramadol and its active metabolite in human plasma. (ref. 130)



Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4 × 250 mm
 Temperature: 30°C ± 3°C
 Eluent: acetonitrile:water (19:81, v:v) cont. 0.06 M NaH₂PO₄ and 0.05 M triethylamine, adjusted to pH 7.90
 Flow rate: 1 ml/min.
 Detection: fluorescence (λ_{ex} 207 nm and λ_{em} 300 nm)

Tranquilizers, veterinary

Analysis of xylazine, haloperidol, acetopromazine, propionylpromazine and chlorpromazine in bovine liver. (ref. 272b)




Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4 × 150 mm
 Eluent: MeOH:water (80:20; v:v)
 Flow rate: 1 ml/min.
 Detection: UV 240 nm

Environmental

Alkyllead


Determination of tetramethyllead and tetraethyllead. (ref. 198)



Phase: Kromasil 100 Å, 5 µm, C18
 Column: 0.32 × 230 mm
 Temperature: start: 50 °C, ramp: 16 °C/min., hold: 100 °C
 Eluent: ACN
 Flow rate: 10 µl/min
 Detection: ICP-MS

Benzimidazole fungicides


Determination of benzimidazole fungicides in fruits. (ref. 112)



Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4 × 150 mm
 Temperature: 55 °C
 Eluent: MeOH-water (50:50; v:v)
 Flow rate: 1 ml/min.
 Detection: UV 285 nm

Explosives


Sensitive determination of RDX, nitroso-RDX metabolites and other munitions in ground water. (ref. 175)



Phase: Kromasil 100 Å, C8
 Column: 2 × 250 mm
 Temperature: 32 °C
 Eluent: isopropanol:water:0.5 M ammonium formate (pH 8 adjusted by ammonium hydroxide) (20:78:2; v:v:v)
 Flow rate: 0.2 ml/min.
 Detection: UV

Herbicides

Analysis of 4-chloro-2-oxo-3(2H)-benzothiazoleacetic ethyl ester and related compounds. (ref. 286)




Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.6 × 200 mm
 Temperature: 30 °C
 Eluent: MeOH:water:HAc (60:40:1; v:v:v)
 Flow rate: 0.7 ml/min.
 Detection: UV 254 nm

Environmental

Nitrophenols


Determination of toxic nitrophenols in the atmosphere. (ref. 183)



Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.4 × 250 mm
 Eluent: A:B (55:45; v:v) A: 0.005 M KH₂PO₄ (pH 4.5 with H₃PO₄):ACN (90:10; v:v) B: 0.005 M KH₂PO₄ (pH 4.5 with H₃PO₄):MeOH (25:75; v:v)
 Flow rate: 1 ml/min.
 Detection: 230 nm

Organonitrogen pesticides


Determination of organonitrogen pesticides in large volumes of surface water. (ref. 132)



Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.6 × 250 mm
 Eluent: Gradient, ACN:water, 15 – 60% ACN for 50 min, 60% for 15 min
 Flow rate: 1 ml/min.
 Detection: APCI-MS

Pesticides and metabolites


Analysis of polar phenolic compounds, pesticides and metabolites in water. (ref. 167)



Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.6 × 250 mm
 Temperature: 65 °C
 Eluent: ACN:water (pH 3 adjusted with sulfuric acid)
 Gradient: From 15 to 25% ACN in 9.3 min., to 50% ACN in 4.3 min., to 100% ACN in 6 min. and then 2 min. isocratic elution at 100% ACN.
 Flow rate: 1 ml/min.
 Detection: UV 280 or 240 nm

Phenylurea herbicides

Determination of phenylurea herbicides in water. (ref. 32)




Phase: Kromasil 100 Å, 5 µm, C18
 Column: 1 × 300 mm
 Eluent: MeOH:water (75:25; v:v) in 0.01 M lithium perchlorate at pH 5.5 (adjusted with 1% phosphoric acid)
 Flow rate: 20 – 40 µl/min
 Detection: UV 254 nm and electrochemical (potential 1,35 V) respectively for the figures

Environmental

Polycyclic aromatic hydrocarbons

Analysis of polycyclic aromatic hydrocarbons. (ref. 184)



Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.6 × 250 mm
 Temperature: 40 °C
 Eluent: CO₂:ACN
 Gradient: 0 min. 100% CO₂, 20 min. 60% CO₂,
 25 min. 60% CO₂


1 = naphthalene



7 = fluoranthene



12 = benzo(k)fluoranthene




2 = acenaphthylene



8 = pyrene




13 = benzo(a)pyrene




3 = acenaphthene



9 = benz(a)anthracene




14 = dibenzo(a,h)anthracene




4 = fluorene



10 = chrysene




15 = indeno(1,2,3-c,d)pyrene




5 = phenanthrene



11 = benzo(b)fluoranthene



16 = benzo(g,h,i)perylene



6 = anthracene



Flow rate: 3 ml/min.
 Detection: UV 210 nm


Quaternary ammonium herbicides

Determination of quaternary ammonium herbicides. (ref. 201)


1 = diquat





2 = paraquat



3 = ethyl viologen



4 = difenzoquat




Phase: Kromasil 100 Å, 5 µm, C8
 Column: 2.1 × 200 mm
 Temperature: 50 °C
 Eluent: Pentafluoropropionic acid in water (15 mM, pH 3.3) : ACN
 Gradient: 0 min. 2% ACN, 5 min. 8.6% ACN, 5.01 min. 40%
 ACN, 13 min. 40% ACN
 Flow rate: 200 µl/min.
 Detection: UV

Food and nutrition

Antioxidants, lipophilic


Determination of lipophilic antioxidants in plasma. (ref. 53)



Phase: Kromasil 100 Å, 5 µm, C8
 Column: 4.6 × 150 mm
 Eluent: Tris (50 mM), HCl (12 mM) and 65% ACN (pH 8.5)
 Flow rate: 1 ml/min.
 Detection: electrochemical, potential +0.70 V

Irganox


Determination of Irganox (an antioxidant). (ref. 20)



Phase: Kromasil 100 Å, 3.5 µm, C18
 Column: 0.25 × 250 mm
 Temperature: gradient: 5 °C/min. from 5 to 40 °C, 2 °C/min. from 40 to 80 °C, 5 °C/min. from 80 to 90 °C
 Eluent: ACN + 10mM TEA + HCOOH
 Flow rate: 5 µl/min.
 Detection: ELS

Irganox


Determination of Irganox. (ref. 208)



Phase: Kromasil 100 Å, 5 µm, C18
 Temperature: from 7 to 90 °C at 3 °C/min.
 Column: 0.32 × 500 mm
 Eluent: ACN
 Flow rate: 5 µl/min
 Detection: UV 280 nm

Sugars

Analysis of sugars in urine. (ref. 82)




Phase: Kromasil 100 Å, 5 µm, NH₂
 Column: 4.6 × 250 mm
 Temperature: ambient
 Eluent: ACN:water (70:30; v:v)
 Flow rate: 1 ml/min.
 Detection: refractive index

Food and nutrition

Sugars


Analysis of sugars. (ref. 315)



Phase: Kromasil 100 Å, 5 µm, NH₂
 Column: 4.6 × 250 mm
 Eluent: ACN:water (75:25; v:v)
 Flow rate: 1 ml/min.
 Detection: RI

Sugars, phosphorylated


Determination of reducing sugars in beef sirloin, with post-column reduction. (ref. 27)




Phase: Kromasil 100 Å, 5 µm, NH₂
 Column: 4 × 250 mm
 Eluent: ACN:water (85:15; v:v) at pH 4.8
 Flow rate: 1.4 ml/min.
 Post column: Post-column reduction at 95 °C with tetrazolium blue (0.7 mM in distilled water and 0.16 M NaOH, 15% EtOH, 0.047M Na-K-tartrate, pH 12.7) before detection.
 Detection: 550 nm

Sugars and polyols, benzoylated

Analysis of benzoylated sugars and polyols. (ref. 51b)




Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4 × 250 mm
 Eluent: Gradient, ACN:water, 0 min. 70% ACN, 30 min. 95% ACN
 Flow rate: 1 ml/min.
 Detection: UV 228 nm



Natural products

Actarit


Determination of actarit and related compounds. (ref. 274)



Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.6 × 250 mm
 Eluent: MeOH:water (70:30; v:v) + 1 % tetrabutylammonium bromide
 Flow rate: 1 ml/min.
 Detection: UV 245 nm

Gaulthersides


Determination of Gaulthersides in Yunnan wintergreen. (ref. 307)



Phase: Kromasil 100 Å, 5 µm, C18
 Column: 3.9 × 250 mm
 Temperature: ambient
 Eluent: MeOH:ACN:water (25:5:70; v:v:v) pH=3.5 (adjusted with H₃PO₄)
 Flow rate: 0.7 ml/min.
 Detection: UV 220 nm

Caffeine and metabolites

Quantitation of caffeine metabolism products. (ref. 271)




Phase: Kromasil 100 Å, 5 µm, C4
 Column: 4 × 250 mm
 Temperature: ambient
 Eluent: acetate buffer (pH 3.5) : MeOH (97:3; v:v)
 Gradient: 0 min. 3% MeOH, 20 min. 20% MeOH
 Flow rate: 1 ml/min.
 Detection: UV 275 nm

Natural products


Ginkgolides

Determination of ginkgolides. (ref. 277)


1 = bilobalide



2 – 4 = ginkgolide



	R ₁	R ₂	R ₃
2. ginkgolide C	OH	OH	OH
3. ginkgolide A	H	H	OH
4. ginkgolide B	OH	H	OH




Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.6 × 250 mm
 Eluent: water:MeOH (77:33; v:v)
 Flow rate: 1 ml/min.
 Detection: refractive index

Nicotine

Clinical assay of nicotine and its metabolite, cotinine, in body fluids. (ref. 306)



1 = cotinine



2 = nicotine




3 = scopolamine




Phase: Kromasil 100 Å, 5 µm, C8
 Column: 4 × 250 mm
 Temperature: ambient, 22 °C
 Eluent: ammonium acetate (0.05 M):CH₃OH (60:40; v:v)
 Flow rate: 1.4 ml/min.
 Detection: UV 262 nm

Sphingoids


Analysis of sphinganine and sphingosine from urine with precolumn o-phthaldialdehyde (OPA) derivatization. (ref. 87)




1 = OPA-sphingosine



2 = OPA-sphinganine



3 = OPA-C20-sphinganine




Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.6 × 250 mm
 Temperature: 45 °C
 Eluent A: 0.07 M K₂HPO₄ in MeOH (1:9; v:v)
 Eluent B: MeOH
 Gradient: 0 min. 0% B, 10 min. 0% B, 30 min. 40% B, 32 min. 100% B, 42 min. 100% B, 44 min. 0% B, 60 min. 0% B

Flow rate: 1.3 ml/min.
 Detection: fluorescence (λ_{ex} 340 nm, λ_{em} 455 nm)

Natural products

TCM, Traditional Chinese Medicine


Determination of gastrodin, p-hydroxybenzyl alcohol, vanillyl alcohol, p-hydroxybenzaldehyde and vanillin from TCM. (ref. 297)



Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.6 × 150 mm
 Temperature: ambient
 Eluents: Eluent A: water, eluent B: MeOH
 Gradient: 0 min 5% B, 9 min 44% B, 12 min 65% B, 15 min 65% B
 Flow rate: 1 ml/min.
 Detection: UV 270 nm

TCM, Traditional Chinese Medicine


Determination of three components in a Chinese doctor-cough syrup. (ref. 210)



Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.6 × 250 mm
 Temperature: 45 °C
 Eluent: MeOH:water:acetic acid (40:60:2; v:v:v) + 5 mM IPR-B₈
 Flow rate: 1 ml/min.
 Detection: UV 245 nm

TCM, Traditional Chinese Medicine


Analysis of caffeine, antipyrine and sodium salicylate in Satongfeng injection. (ref. 215)



Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.6 × 250 mm
 Eluent: 20 mM potassium dihydrogen phosphate: MeOH:glacial acetic acid (55:25:0.4; v:v:v)
 Flow rate: 1 ml/min.
 Detection: UV 242 nm


TCM, Traditional Chinese Medicine

Determination of four components of Ganmaoling capsules. (ref. 258)





Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.6 × 250 mm
 Temperature: 30 °C
 Eluent: ACN:diammonium hydrogen phosphate (pH 3.1, 0.03 M) (12:88; v:v) containing 0.75 – 5 mM sodium sulfonic heptane
 Flow rate: 1 ml/min.
 Detection: UV 214 nm

Analysis of enkephalin peptides, their metabolites and enzyme inhibitors. (ref. 104)




Analysis of five nonapeptides. (ref. 28)



Peptides

Peptides

Separation of 9 peptides. (ref. 316)





Phase: Kromasil 100 Å, 5 µm, C8
 Column: 4.5 × 250 mm
 Temperature: ambient
 Eluent A: 0.1% TFA in 10% CH₃CN and 90% water
 Eluent B: 0.1% TFA in 90% CH₃CN and 10% water
 Gradient: 0 min. 0% B, 8 min. 25% B, 20 min. 75% B

- 1 = oxytocin
Cys-Tyr-Ile-Gln-Asn-Cys-Pro-Leu-Gly-NH2
 $\begin{array}{c} \text{S} \text{-----} \text{S} \end{array}$
- 2 = bradykinin
Arg-Pro-Pro-Gly-Phe-Ser-Pro-Phe-Arg
- 3 = methionine enkephalin
Tyr-Gly-Gly-Phe-Met
- 4 = angiotensin II
Asp-Arg-Val-Tyr-Ile-His-Pro-Phe
- 5 = leucin enkephalin
Tyr-Gly-Gly-Phe-Leu
- 6 = angiotensin I
Asp-Arg-Val-Tyr-Ile-His-Pro-Phe-His-Leu
- 7 = insulin
- 8 = lysosyme
- 9 = melittine
Gly-Ile-Gly-Ala-Val-Leu-Lys-Val-Leu-Thr-Thr-Gly-Leu-Pro-Ala-Leu-Ile-Ser-Trp-Ile-Lys-Arg-Lys-Arg-Gln-Gln-NH2






Flow rate: 2 ml/min.
 Detection: UV 254 nm

Determination of vitamin E in human plasma. (ref. 108)








Phase: Kromasil 100 Å, 5 µm, C1
Column: 4.6 × 100 mm
Temperature: ambient
Eluent: MeOH:ACN:water (50:35:15; v:v:v)
Flow rate: 1.5 ml/min.
Detection: UV 292 nm

Determination of tocopherols and vitamin A in vegetable oils.
(ref. 188)



Phase: Kromasil 100 Å, 5 µm, C18
 Column: 0.2 × 800 mm
 Temperature: 65°C
 Eluent: CO₂ with 8% MeOH
 Pressure: 180 atm
 Detection: electrochemical (potential + 1.80 V versus Quasi-Reference Electrode)

Analysis of soluble vitamins. (ref. 330)




33

Other

Amines

Determination of amines from fish decomposition by dansylchloride derivatisation. (ref. 73)




Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.6 × 250 mm
 Temperature: 25 °C
 Eluent: ACN:water
 Gradient: 0 min 60% ACN, 6 min 75% ACN, 8 min 75% ACN, 13 min 95% ACN, 20 min 95% ACN, 20.01 min 60% ACN

1 = $\text{NH}_3^+ \text{R}$

2 = methylamine
 NHR


3 = tryptamine




4 = 1,3-diaminopropane




5 = putrescine




6 = cadaverine




7 = histamine




8 = tyramine




9 = spermidine



10 = spermine




dansyl derivative group



Flow rate: 1 ml/min.
 Detection: UV 254 nm

Amino alcohols

Separation of derivates of 1-ethylamino-3-phenoxy-propan-2-ol. (ref. 38)



Phase: Kromasil 100 Å, 5 µm, C8
 Column: 0.2 × 900 mm
 Eluent: ACN:ammonium acetate(5 mM) (55:45; v:v)
 Flow rate: 0.95 µl/min.
 Detection: ESI-MS

Aroma extracts in alcoholic beverages

Separation of aroma extracts found in wine and other alcoholic beverages. (ref. 209)


1 = furfural



2 = sotolon




3 = vanillin





4 = guaiacol




5 = 2-phenylethanol



6 = 2-phenylethyl acetat




4 = guaiacol



Phase: Kromasil 100 Å, 5 µm, C18
 Column: 10 × 250 mm
 Eluent: water:ethanol
 Gradient: 0 min. 100% water, 8 min. 80% water, 28 min. 50% water, 40 min. 0% water
 Flow rate: 2 ml/min.
 Detection: UV 220 nm

Aromatics


Separation of mixtures of nitrobenzoic acid and aminobenzoic acid isomers. (ref. 214)



Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.6 × 200 mm
 Temperature: 35 °C
 Eluent: MeOH:water:THF (55:44:1; v:v:v) with β-cyclodextrin at pH 3.0
 Flow rate: 0 – 4 min. 2 ml/min., 4 – 10 min. 2.6 ml/min.
 Detection: UV 254 nm

Aromatics


HPLC analysis of isomers of nitrotoluene and nitrobenzoic acid. (ref. 213)



Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.6 × 200 mm
 Temperature: 35 °C
 Eluent: MeOH:water:THF (55:44:1; v:v:v) with β-cyclodextrin at pH 3.0
 Flow rate: 0 – 4 min. 2 ml/min., 4 – 10 min. 2.6 ml/min.
 Detection: UV 254 nm

Aromatics


Determination of sulfonurea, benzene, toluene, naphtalene, biphenyl, phenanthrene, anthracene. (ref. 301a)



Phase: Kromasil 100 Å, 5 µm, C18
 Column: 0.8 × 150 mm
 Eluent: MeOH:water (80:20; v:v)
 Flow rate: 38 µl/min.
 Detection: UV 254 nm

Aromatics

Separation of benzene and pyridine derivatives. (ref. 40)




Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.6 × 150 mm
 Eluent: ACN:water (56.9:43.1; v:v)
 Flow rate: 1 ml/min.
 Detection: UV 254 nm

Other


ATP degradation products

Determination of ATP degradation products from fish decomposition. (ref. 159)




Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.6 × 250 mm
 Temperature: 25 °C
 Eluents: Eluent A, ACN and eluent B, phosphate buffer (pH 7.00, 60 mM K₂HPO₄ + 40 mM KH₂PO₄)


1 = IMP




2 = ATP




3 = ADP



4 = AMP




5 = hypoxanthine



6 = purine



7 = inosine



Gradient: 0 min. 100% B, 4 min. 98% B, 5 min. 97% B, 8 min. 96% B, 15 min. 96% B, 15.01 min. 100% B
 Flow rate: 1 ml/min.
 Detection: UV 254 nm

Benzene, substituted

Separation of substituted benzene. (ref. 1)

1 = nitrobenzene



2 = fluorobenzene




3 = chlorobenzene



4 = toluene



5 = brombenzene



Phase: Kromasil 100 Å, 5 µm, C8
 Column: 4.6 × 250 mm
 Temperature: 20 °C
 Eluent: ACN:water (60:40; v:v)
 Flow rate: 1 ml/min.
 Detection: UV 210 nm

Chlorinated benzenes

Determination of chlorobenzene and derivatives. (ref. 301c)

1 = chloro-benzene



2 = p-dichloro-benzene



3 = m-dichloro-benzene



4 = 1,2,3-trichloro-benzene




5 = 1,2,4-trichloro-benzene




6 = 1,3,5-trichloro-benzene




7 = 1,2,3,4-tetrachloro-benzene



8 = 1,2,4,5-tetrachloro-benzene




9 = pentachloro-benzene




10 = hexachloro-benzene



11 = 1,2,3,4,5-pentachloro-benzene



12 = 1,2,3,4,6-pentachloro-benzene



13 = 1,2,3,4,5,6-hexachloro-benzene



14 = 1,2,3,4,5,6-hexachloro-benzene




15 = 1,2,3,4,5,6-hexachloro-benzene



Phase: Kromasil 100 Å, 5 µm, C18
 Column: 0.8 × 150 mm
 Eluent: eluent A: ACN, eluent B: water
 Gradient: 0 min. 80% A, 5 min. 80% A, 10 min. 100% A
 Flow rate: 32 µl/min.
 Detection: UV 220 nm

Fatty acids


Analysis of plasma fatty acids as their phenacyl esters. (ref. 193)



Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.6 × 250 mm
 Temperature: ambient
 Eluent: MeOH:water (91:9; v:v)
 Flow rate: 1.15 ml/min.
 Detection: UV 254 nm

Lauric acid


Detection of ester of lauric acid. (ref. 35)



Phase: Kromasil 100 Å, 7 µm, C18
 Column: 4.6 × 150 mm
 Eluent: ACN:MeOH:water (55:10:35; v:v:v)
 0.2% phosphoric acid added
 Flow rate: 1 ml/min.
 Detection: fluorescence (λ_{ex} 357.5 nm and λ_{em} 482 nm)

Organic acids

Separation of formic acid, benzoic acid, lactic acid, acetic acid. (ref. 344)




Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4.6 × 250 mm
 Eluent: KH_2PO_4 -buffer (10 mM, pH 2.5):ACN (95:5; v:v)
 Flow rate: 38 µl/min.
 Detection: UV 254 nm

Other

QC test, neutral compounds


QC test of Kromasil CN. (ref. 341)



Phase: Kromasil 60 Å, 10 µm, CN
 Column: 4.6 × 250 mm
 Eluent: hexane:ethylacetate (90:10; v:v)
 Flow rate: 2 ml/min.
 Detection: UV 254 nm

QC test, neutral compounds


QC test of Kromasil SIL. (ref. 346)



Phase: Kromasil 60 Å, 5 µm, SIL
 Column: 4.6 × 250 mm
 Eluent: hexane:ethylacetate (85:15; v:v)
 Flow rate: 2 ml/min.
 Detection: UV 254 nm


QC test, silanophilic compounds

QC test of Kromasil SIL. (ref. 345)



Phase: Kromasil 60 Å, 5 µm, SIL
 Column: 4.6 × 250 mm
 Eluent: MeCl₂:MeOH (98:2; v:v)
 Flow rate: 2 ml/min.
 Detection: UV 254 nm


QC test, substituted aromatic compounds

QC test of Kromasil NH₂. (ref. 343)


Phase: Kromasil 100 Å, 5 µm, NH₂
 Column: 4.6 × 250 mm
 Eluent: hexane:MeCl₂ (97:3; v:v)
 Flow rate: 1 ml/min.
 Detection: UV 254 nm

Flavonoid glycosides

Analysis of flavonoid glycosides. (ref. 100)




Phase: Kromasil 100 Å, 5 µm, C18
 Column: 3.2 × 250 mm
 Eluent: ACN:water
 Gradient: 0 min. 20% ACN, 10 min. 20% ACN, 18 min. 40% ACN, 28 min. 75% ACN, 30 min. 100% ACN, 37 min. 100% ACN



Flow rate: 0.75 ml/min.
 Detection: UV 280 nm

Organometallic catalysts

Purity testing of organometallic catalysts. (ref. 248)




1 = CpMn(CO)₃
 2 = η⁵-CpMoH(CO)₂PMe₃
 3 = η⁵-CpW(CO)₃Me
 4 = η⁵-CpMoH(CO)₂PPh₃

Phase: Kromasil 100 Å, 5 µm, C18
 Column: 0.32 × 450 mm
 Temperature: 60 °C
 Eluent: carbon dioxide
 Flow rate: 7.2 µl/min.
 Pressure: 100 bar (hold 10 min.) then 10 bar/min. until 180 bar (hold 1 min.), then 10 bar/min. until 300 bar (hold 1 min.), then 10 bar/min. until 400 bar (hold 10 min.)
 Detection: FID

Surfactants

Determination of caproic acid, n-octyl alcohol and n-octyl caproate. (ref. 285)




1 = caproic acid
 2 = n-octyl alcohol
 3 = n-octyl caproate

Phase: Kromasil 100 Å, 5 µm, C18
 Temperature: 30 °C
 Column: 4.6 × 150 mm
 Eluent: MeOH:water (95:5; v:v)
 Flow rate: 1 ml/min.
 Detection: refractive index

Other

Triacylglycerols


Analysis of seven triacylglycerols. (ref. 139)



Phase: Kromasil 100 Å, 5 µm, C18
 Column: 0.7 × 120 mm
 Eluent: (A):ACN, (B):acetone
 Gradient: stepwise: 0 – 5 min. 90% A, 5 – 25 min. 70% A, after 25 min. 40% A.
 Flow rate: 5 – 100 µl/min (not specified)
 Detection: ELS

Triclosan


Determination and stability tests of triclosan in disinfectants. (ref. 8)



Phase: Kromasil 100 Å, 7 µm, C18
 Column: 4.6 × 200 mm
 Eluent: MeOH:ACN:water (40:40:20; v:v:v) containing 0.02 M KH₂PO₄ (pH 2.7)
 Flow rate: 1 ml/min.
 Detection: UV 280 nm

Triglycerides

Analysis of triglyceride profiles in Cretan olive oils. (ref. 96)



Phase: Kromasil 100 Å, 5 µm, C18
 Column: 4 × 250 mm
 Temperature: 40 °C
 Eluent: acetone:ACN (60:40; v:v)
 Flow rate: 0.7 ml/min.
 Detection: refractive index

Literature references

- 1 A. Granfors, Masters Project, Dept. Anal. & Marine Chem., Univ. of Gothenburg (2002)
- 6 J. I. D. Wibawa, D.Fowkes D, P.N. Shaw, D. A. Barrett, J. Chromatogr. B 774(2) (2002) 141-148 ¹
- 7 J. Yang, L.-G. Sun, X.-Z. Bai, H.-T. Zhou, Chin. J. Chromatogr. 20(4) (2002) 369-371 ¹⁰
- 8 J. Li, H.-Y. Zhao, G.-K. Peng, P.-Y. Tian, C.-Q. Li, Chin. J. Chromatogr. 20(4) (2002) 372-374 ¹⁰
- 10 K. Nithipatikom, A. J. Grall, B. B. Holmes, D. R. Harder, J. R. Falck, W. B. Campbell, Anal. Biochem. 298(2) (2001) 327-336 ¹
- 20 P. Molander, A. Holm, E. Lundanes, D.R. Hegna, E. Ommundsen, T. Greibrokk, Analyst 127(7) (2002) 892-897 ⁹
- 23 X. Ouyang, S. Ren, Chin. J. Pharm. Anal. 22(1) (2002) 78-80 ¹¹
- 26 A. Zotou, N. Miltiadou, J. Pharm. Biomed. Anal. 28(3-4) (2002) 559-568 ¹
- 27 M. Aliani, L. J. Farmer, J. Agric. Food Chem. 50(10) (2002) 2760-2766 ¹³
- 28 K. Larsson, W. Hermann, P. Möller, D. Sanchez, J. Chromatogr. A 450(1) (1988) 71-80 ¹
- 30 B. Gustavsson, I. Betnér, J. Chromatogr. A 507 (1990) 67-77 ¹
- 32 R. Boussenadji, P.Dufek, M. Porthault, LC-GC Int. 5(10) (1992) 40-43 ⁷
- 35 M. Chelminska-Bertilsson, S. Allenmark, J. Chromatogr. B 575 (1992) 237-242 ¹
- 38 K.-E. Karlsson, J. Chromatogr. A 647(1) (1993) 31-38 ¹
- 40 H. A. Claessens, E. A. Vermeer, C. A. Cramers, LC-GC Int. 6(11) (1993) 692-693, 696-698, 700 ⁷
- 45 B. Lindegård, H. Björk, J. A. Jönsson, L. Mathiasson, A.-M. Olsson, Anal. Chem. 66(24) (1994) 4490-4497 ⁴
- 49 S. Joron, H. Robert, Biomed. Chromatogr. 8 (1994) 158-164 ²
- 51 J. Oehlke, M.Brudel, I. E. Blasig, J. Chromatogr. B 655(1) (1994) 105-111 ¹
- 53 M. Hedenmo, B.-M. Eriksson, J. Chromatogr. A 661(1-2) (1994) 153-159 ¹
- 58 J. Atta-Politou, K. Tsarpalis, A. Koutselinis, J. Liq. Chromatogr. 17(18) (1994) 3969-3982 ¹²
- 64 B. Jeanny, J.L. Vaillau, P. d'Athis, V. Lahet, Analisis 23(1) (1995) 7-11 ³
- 66 A. Longo, G. Bruno, S.Curti, A. Mancinelli, G. Miotto, J. Chromatogr. B 686(2) (1996) 129-139 ¹
- 71 A. J. Oosterkamp, H. Irth, L. Heintz, G. Marko-Varga, U.R. Tjaden, J. van der Greef, Anal. Chem. 68 (1996) 4101-4106 ⁵
- 73 M. Vallé, P. Malle, S. Bouquelet, J. AOAC Int. 79(5) (1996) 1134-1140 ⁶
- 82 K. Miki, R. Butler, D. Moore, G. Davidson, Clin. Chem. 42(1) (1996) 71-75 ⁸
- 87 M. Castegnaro, L. Garren, I. Gaucher, C.P. Wild, Natural Toxins 4(6) (1996) 284-290 ¹⁴
- 94 P. Mourier, A. Brun, J. Chromatogr. B 704(1-2) (1997) 197-205 ¹
- 95 M. Hay, P. Mormède, J. Chromatogr. B 703(1-2) (1997) 15-23 ¹
- 96 E. Stefanoudaki, F. Kotsifaki, A. Koutsaftakis, Food Chem. 60(3) (1997) 425-432 ¹
- 97 R. Niemi, H. Taipale, M. Ahlmark, J. Vepsäläinen, T. Järvinen, J. Chromatogr. B 701(1) (1997) 97-102 ¹
- 98 J. D. H. Cooper, D.C. Muirhead, J. E. Taylor, P.R. Baker, J. Chromatogr. B 701(1) (1997) 87-95 ¹
- 100 H. L. Constant, K. Slowing, J. G. Graham, J. M. Pezzuto, G. A. Cordell, C. W. W. Beecher, Phytochem. Anal. 8(4) (1997) 176-180 ¹⁵
- 103 B. Lu, M. Koimur, D. Westerlund, Chromatographia 46(1-2) (1997) 72-78 ¹⁶
- 104 V.-P. Ranta, A. Urtti, S. Auriola, J. Chromatogr. A 766(1-2) (1997) 85-97 ¹
- 108 J. D. H. Cooper, R. Thadwal, M. J. Cooper, J. Chromatogr. B 690(1-2) (1997) 355-358 ¹
- 110 X. Ji, K. Yu, W. Zhang, Chin. J. Pharm. Anal. 17(6) (1997) 363-365 ¹¹
- 112 N. Tharsis, J.L. Portillo, F. Broto-Puig, L. Comellas, J. Chromatogr. A 778 (1997) 95-101 ¹
- 113 B. Enanga, C. Labat, H. Boudra, G. Chauvière, M. Keita, B. Bouteille, M. Dumas, G. Houin, J. Chromatogr. B 696(2) (1997) 261-266 ¹
- 116 J. M. Joulia, F. Pinguet, P.Y. Grosse, C. Astre, F.Bressolle, J. Chromatogr. B 692(2) (1997) 427-435 ¹
- 119 M. D. Rose, J. Bygrave, G.W.F. Stubbings, Analyst 123(12) (1998) 2789-2796 ⁹
- 124 T. Dine, F.Khalfi, B.Gressier, M. Luyckx, C.Brunet, L. Ballester, F.Goudaliez, J. Kablan, M.Cazin, J.C. Cazin, J. Pharm. Biomed. Anal. 18(3) (1998) 373-381 ¹
- 126 S. Frappier, D. Breilh, E. Diarte, B. Ba, D. Ducint, J. L. Pellegrin, M.C. Saux, J. Chromatogr. B 714(2) (1998) 384-389 ¹
- 129 W.-X. Zhang, C.-G. Jiang, G.-L. Guo, Chin. J. Pharm. Anal. 18(5) (1998) 314-316 ¹¹
- 130 M. A. Campanero, B. Calahorra, E. Garcia-Quetglás, M.Escolar, J. Honorato, Chromatographia 48(7-8) (1998) 555-560 ¹⁶
- 132 H. Sabik, R. Jeannot, J. Chromatogr. A 818(2) (1998) 197-207 ¹
- 139 M. B. O. Andersson, L. G. Blomberg, J. Microcolumn Sep. 10(3) (1998) 249-254 ¹⁷
- 143 A. Le Liboux, O. Pasquier, G. Montay, J. Chromatogr. B 708(1-2) (1998) 161-168 ¹
- 154 S. Ravisankar, M. Vasudevan, M. Gandhimathi, B. Suresh, Talanta 46(6) (1998) 1577-1581 ¹
- 156 S. Ravisankar, M. Vasudevan, M. J. Nanjan, M. Gandhimathi, B. Suresh, Indian Drugs 35(6) (1998) 359-363 ¹⁹
- 159 M. Vallé, P. Malle, S. Bouquelet, J. AOAC Int. 81(3) (1998) 571-575 ⁶
- 161 H. K. Jajoo, N.V.S. Rao Mamidi, K. Kasiram, A. S. Prakash, V.V.S. Swaroop Kumar, P. Bheema Rao, V. Bhushan, S. Subramaniam, J. Chromatogr. B 707(1-2) (1998) 241-246 ¹
- 167 N. Masqué, M. Galià, R. M. Marcé, F. Borrull, J. High Resolut. Chromatogr. 22(10) (1999) 547-552 ²⁰
- 170 Ph. Baudot, A. Vicherat, M.-L. Viriot, M.-C. Carré, Analisis 27(6) (1999) 523-532 ³
- 175 D. A. Cassada, S.J. Monson, D.D. Snow, R.F. Spalding, J. Chromatogr. A 844(1-2) (1999) 87-95 ¹
- 183 R. Belloli, B. Barletta, E. Bolzacchini, S. Meinardi, M. Orlandi, B. Rindone, J. Chromatogr. A 846(1-2) (1999) 277-281 ¹
- 184 E. Lesellier, Analisis 27(3) (1999) 241-249 ⁵
- 188 F. J. Señoráns, K. E. Markides, L. Nyholm, J. Microcolumn Sep. 11(5) (1999) 385-391 ¹⁷
- 193 I. Zaitseva, M. Ajmal, E. Cersosimo, J. Chromatogr. B 727(1-2) (1999) 15-22 ¹
- 198 R. Trones, A. Tangen, W. Lund, T. Greibrokk, J. Chromatogr. A 835(1-2) (1999) 105-112 ¹

- 201 R. Castro, E. Moyano, M. T. Galceran, J. Chromatogr. A 830(1) (1999) 145-154 ¹
- 208 P. Molander, K. Haugland, D. R. Hegna, E. Ommundsen, E. Lundanes, T. Greibrokk, J. Chromatogr. A 864(1) (1999) 103-109 ¹
- 209 V. Ferreira, P. Hernandez-Orte, A. Escudero, R. Lopez, J. Cacho, J. Chromatogr. A 864(1) (1999) 77-88 ¹
- 210 X.-F. Wei, Y.-Q. Du, Acad. J. of Guangdong Coll. of Pharm. 15(3) (1999) 164-166 ²¹
- 213 Y.-C. Liang, Z.-Z. Yi, Z.-H. Cai, G.-Y. Zhu, Chin. J. Chromatogr. 17(4) (1999) 397-398 ¹⁰
- 214 Y.-C. Liang, Z.-Z. Yi, Z.-H. Cai, G.-Y. Zhu, Anal. Lab. 18(3) (1999) 55-57 ²²
- 215 B. Liao, C. Xin, J. Chin. Pharm. Univ. 30(2) (1999) 121-123 ²³
- 217 B. Ba. Boubakar, A.-G. Corniot, D. Ducint, D. Breilh, J. Grellet, M.-C. Saux, J. Chromatogr. B 724 (1999) 127-136 ¹
- 228 B. Tavazzi, R. Vagnozzi, D. Di Perro, A. M. Amorini, G. Fazzina, S. Signoretti, A. Marmarou, I. Caruso, G. Lazzarino, Anal. Biochem. 277(1) (2000) 104-108 ¹
- 237 D. Di Pierro, G. Lazzarino, F. S. Pastore, B. Tavazzi, F. Del Bolgia, A. M. Amorini, G. Fazzina, R. Giuffrè, Anal. Biochem. 284 (2000) 301-306 ¹
- 244 N. V. S. Rao Mamidi, V. V. S. Swaroop Kumar, K. Katneni, M. Rao Chaluvasi, S. Shreeram, S. Gangadhar, B. Nataraj, P. Hari Kishore, V. Bhushan, S. Subramaniam, J. Pharm. Biomed. Anal. 22(2) (2000) 251-255 ¹
- 246 J. J. Berzas Nevado, G. Castañeda Peñalvo, F. J. Guzmán Bernardo, J. Chromatogr. A 870(1-2) (2000) 169-177 ¹
- 247 E. G. de Jalón, M. Josa, M. A. Campanero, S. Santoyo, P. Ygartua, J. Chromatogr. A 870(1-2) (2000) 143-149 ¹
- 248 I. Bruheim, I. L. Skuland, E. Lundanes, T. Greibrokk, J. Chromatogr. A 868(2) (2000) 261-268 ¹
- 254 M. Lee, D. I. Min, Ther. Drug Monit. 23(1) (2001) 21-26 ²⁴
- 258 X. Ran, D. Hu, J. Wang, Chin. J. Pharm. Anal. 21(5) (2001) 369-370 ¹¹
- 262 L.-S. Li, W.-D. Huang, Q. He, S. Ye, Chin. J. Chromatogr. 19(5) (2001) 446-448 ¹⁰
- 267 J. J. Berzas Nevado, G. Castañeda Peñalvo, F. J. Guzmán Bernardo, Analytica Chimica Acta 442 (2001) 241-248 ¹
- 268 H.-Q. Zhu, J.-L. Ming, Chin. J. Pharm. Anal. 21(3) (2001) 198-200 ¹¹
- 271 K. A. Georga, V. F. Samanidou, I. N. Papadoyannis, J. Chromatogr. B 759(2) (2001) 209-218 ¹
- 272 M. C. Quintana, M. H. Blanco, J. Lacal, L. Hernandez, J. Liq. Chromatogr. Relat. Technol. 24(5) (2001) 735-745 ²⁵
- 273 X.-S. Zhang, Y.-C. Liu, L. Wang, C.-S. Lin, J. China Univ. Sci. Tech. 31(2) (2001) 213-217
- 274 F. Fang, G.-J. Xu, B. Lu, J. Liq. Chromatogr. Relat. Technol. 24(7) (2001) 1021-1027 ²⁵
- 277 Y.-M. Guan, H.-L. You, J. Wang, X. Lin, Chin. J. Pharm. Anal. 21(2) (2001) 107-109 ¹¹
- 279 S. Zhang, J. Yi, Z. Tan, H. Zhou, Zhongguo Yiyuan Yaoxue Zazhi 21(3) (2001) 139-141 ²⁶
- 283 Y. Li, G.-F. Wang, Y.-F. Wu, L.-H. Long, Chin. J. Pharm. Anal. 21(1) (2001) 33-36 ¹¹
- 285 G.-W. Zhou, X.-R. Huang, Y.-Z. Li, G.-Z. Li, W. Hu, S.-F. Song, Fine chemicals 18(1) (2001) 29-30, 33 ²⁷
- 286 T. Shao, Y.-D. Long, T.-B. Huang, Chin. J. Anal. Chem. 29(1) (2001) 74-76 ²⁸
- 297 C.-L. Liu, M.-C. Liu, P.-L. Zhu, Chromatographia 55(5-6) (2001) 317-320 ¹⁶
- 301 Y.-H. Tang, D.-Q. Zhu, Y.-F. Guan, Chin. J. Chromatogr. 19(4) (2001) 289-292 ¹⁰
- 306 I. N. Papadoyannis, V. F. Samanidou, P. G. Stefanidou, J. Liq. Chromatogr. Relat. Technol. 25(13-15) (2001) 2315-2335 ²⁵
- 307 Y.-J. Zhao, Z.-T. Han, W.-Z. Wang, X.-J. Ma, J.-H. Zheng, Chin. J. Anal. Chem. 30(9) (2001) 1109-1111 ²⁸
- 309 P. Molander, A. Thomassen, L. Kristoffersen, T. Greibrokk, E. Lundanes, J. Chromatogr. B 766 (2001) 77-87 ¹
- 315 EKA, Kromasil Application Note
- 316 EKA, Kromasil Application Note
- 330 EKA, Kromasil Application Note
- 331 EKA, Kromasil Application Note
- 341 EKA, Kromasil Application Note
- 342 EKA, Kromasil Application Note
- 343 EKA, Kromasil Application Note
- 344 EKA, Kromasil Application Note
- 345 EKA, Kromasil Application Note
- 346 EKA, Kromasil Application Note
- 347 G. Lazzarino, A.M. Amorini, B. Tavazzi, G. Galvano, Reprinted by permission of Professor G. Lazzarino, Dept. Chem. Sci., Lab. Biochem., Univ. Catania, Catania, Italy
- 348 G. Lazzarino, A.M. Amorini, G. Fazzina, B. Tavazzi, Reprinted by permission of Professor G. Lazzarino, Dept. Chem. Sci., Lab. Biochem., Univ. Catania, Catania, Italy
- 349 EKA, Kromasil Application Note
- 350 EKA, Kromasil Application Note
- 351 EKA, Kromasil Application Note

Copyright footnotes:

1. Reprints from Analytical Biochemistry, Journal of Chromatography (A & B), Journal of Pharmaceutical and Biomedical Analysis, Analytica Chimica Acta, Food Chemistry & Talanta with permission from Elsevier
2. Reprint from Biomedical Chromatography reproduced with permission of John Wiley & Sons Limited
3. Reprints from Analysis with permission from WILEY-VCH
4. Reprint from Analytical Chemistry with permission from the American Chemical Society and Prof Jan Åke Jönsson, Analytical Chemistry, Lund University, Sweden.
5. Reprint from Analytical Chemistry with permission from the American Chemical Society and Prof H. Irth, Analytical Chemistry, Leiden University, The Netherlands
6. Reprints from Journal of AOAC International with permission from AOAC International
7. Reprinted from LC-GC International by permission of Advanstar Communications UK Ltd
8. Reprints from Clinical Chemistry (Washington, D.C.) with permission from the American Association for Clinical Chemistry
9. Reproduced from Analyst by permission of The Royal Society of Chemistry
10. Reprints from Chinese Journal of Chromatography by permission of Sepu
11. Reprints from Chinese Journal of Pharmaceutical Analysis by permission of Yowu Fenxi Zazhi
12. Reprinted by permission of Journal of Liquid Chromatography
13. Reprinted by permission of Journal of Agricultural and Food Chemistry
14. Reprinted by permission of Natural Toxins
15. Reprinted by permission of Phytochemical Analysis
16. Reprinted by permission of Chromatographia
17. Reprinted by permission of Journal of Microcolumn Separation
19. Reprinted by permission of Indian Drugs
20. Reprinted by permission of Journal of High Resolution Chromatography
21. Reprinted from Academic Journal of Guangdong College of Pharmacy by permission of Guangdong Yaoxueyuan Xuebao
22. Reprinted from Analytical Laboratory by permission of Fenxi Shiyanshi
23. Reprinted from Journal of the Chinese Pharmaceutical University by permission of Zhongguo Yaoke Daxue Xuebao
24. Reprinted by permission of Therapeutic Drug Monitoring
25. Reprinted by permission of Journal of Liquid Chromatography & Related Technologies
26. Reprinted by permission of Zhongguo Yiyuan Yaoxue Zazhi
27. Reprinted by permission of Jingxi Huagong
28. Reprinted by permission of Fenxi Huaxue

Substance index

substance	page		
1-[Ethyl-(3-methanesulfonyl-propyl)-amino]-3-phenoxy-propan-2-ol	34	ascorbic acid	33
1-{[Ethyl-(2-hydroxy-3-(propane-1-sulfonyl)-propyl)-amino]-3-phenoxy-propan-2-ol}	34	asparagine, OPA-	8
1-{[Ethyl-(3-(2-hydroxy-propane-1-sulfonyl)-propyl)-amino]-3-phenoxy-propan-2-ol}	34	aspartate	7
3-{3-[Ethyl-(2-hydroxy-3-phenoxy-propyl)-amino]-propane-1-sulfinyl}-propan-1-ol	34	aspartate, N-acetyl-	7
3-{3-[Ethyl-(2-hydroxy-3-phenoxy-propyl)-amino]-propane-1-sulfonyl}-propan-1-ol	34	aspartic acid, Fmoc-	6
ACA, 7-	10	aspartic acid, OPA-	8
acenaphthene	25	aspartic acid, PTC-	6
acenaphthylene	25	ATP	36
acetaminophen	17	atrazine	24
acetic acid	37	benz(a)anthracene	25
acetoacetanilide	38	benzene	35
acetophenone	38	benzo(a)pyrene	25
acetopromazine	22	benzo(b)fluoranthene	25
acetylaminophenylacetic acid, 4-	28	benzo(g,h,i)perylene	25
acetyl-L-carnitine 1-amino-anthraceneamide	8	benzo(k)fluoranthene	25
ACT, 7-	10	benzoic acid	37
actarit	28	benzothiazoleacetic acid, 4-chloro-2-oxo-3(2H)-	23
adenosine	7	benzothiazoleacetic ethyl ester, 4-chloro-2-oxo-3(2H)-	23
ADP	36	benzylmandelate	38
alanine	7	bestatin	31
alanine, Fmoc-	6	betamethasone acetate	21
alanine, PTC-	6	bilobalide	29
ambroxol	21	biochanin-A	39
ametryn	24	biphenyl	35
aminobenzoic acid, m-	35	boronophenylalanine, OPA-p-	8
aminobenzoic acid, o-	35	bradykinin	32
aminobenzoic acid, p-	35	brombenzene	36
aminobuturic acid, OPA- γ -	8	bromhexine	10
aminocephalosporanic acid, 7-	10	butanamine (EBDB), N-ethyl-1-(1,3-benzodioxol-5-yl)-2-	15
aminosalicylic acid, 3-	8	butanamine (MBDB), N-methyl-1-(1,3-benzodioxol-5-yl)-2-	15
aminosalicylic acid, 5-	8	butylbenzene	38
amitriptyline	12	cadaverine	34
amitriptyline	19	caffeine	28
amoxapine	11	caffeine	30
amoxicillin degradation derivative	9	caffeine	38
AMP	36	caproic acid	39
amperozide (I)	9	carbamazepine	11
amperozide's N-de-ethyl metabolite (II)	9	carbendazim	23
amperozide's N-de-ethyl-N-butyl analogue (III)	9	carnitine 1-aminoanthraceneamide,L-	8
ampicillin degradation derivative	9	ceftriaxone	10
angiotensin I	32	cephalosporanic acid, 7-amino-3-[[[(2,5-dihydro-6-hydroxy-2-methyl-5-oxo-1,2,4-trizin-3-yl)-thio]methyl]-	10
angiotensin II	32	chlorhexidine	10
anthracene	25	chlorhexidine	13
anthracene	35	chlorobenzene	36
anthranilic acid	14	chlorobenzothiazole, 2-amino-4-	23
antipyrine	30	chlorpheniramine maleate	30
apigenin	39	chlorpromazine	22
arachidonic acid (phenacyl ester derivative)	37	chlortetracycline	22
arginine, Fmoc-	6	chlorzoxazone	17
arginine, OPA-	8	chrysene	25
arginine, PTC-	6	chrysin	39
		ciprofloxacin	14
		citalopram	12
		citalopram N-oxide	12
		citrulline, OPA-	8
		clenbuterol	21
		clindamycin	13
		clodronate	15
		clomipramine	12
		clomipramine, 8-OH-	12
		clomipramine, 8-OH-desmethyl-	12
		clomipramine, desmethyl-	12
		clorodolate monophenylester	15
		clotrimazole	10
		clozapine	14
		codeine	30
		cotinine	29
		cyanazine	24
		cyanidin	9
		cyanidin-3-O- β -glycoside	9
		cyclohexylalanine	7
		cystein-conjugate (RPR 100391)	11
		cysteine, Fmoc-	6
		cysteine, PTC-	6
		cysteyl neospiramycin I	21
		cysteyl neospiramycin III	21
		cysteyl spiramycin I	21
		cysteyl spiramycin III	21
		cystine, Fmoc-	6
		daidzein	39
		D-Ala-Gly-Phe-Met-NH ₂	31
		DEA	24
		DEA (desethylatrazine)	24
		delphinidin-3-O- β -glycoside	9
		demethyltramadol, o-	22
		desethylatrazine	24
		desisopropylatrazine	24
		desmethylcitalopram	12
		dexamethasone	21
		DHET, 11,12-	15
		DHET, 12,15-	15
		DHET, 14,15-	15
		DHET, 5,6-	15
		DHET, 8,9-	15
		DIA	24
		DIA (desisopropylatrazine)	24
		dibenzo(a,h)anthracene	25
		dichlorobenzene, m-	36
		dichlorobenzene, p-	36
		dichlorobenzothiazole, 2,4-	23
		didesmethylcitalopram	12
		diethylaniline, N,N-	38
		difenzoquat	25
		dihydroxybenzylamine hydrobromide, 3-4-	14
		dimethylamino-3-propyl thiomethylene-5-virginiamycin S (internal standard)	11
		dimethyluric acid, 1,3-	28
		dimethyluric acid, 1,7-	28
		dinitrophenol, 2,3-	24
		dinitrophenol, 2,4-	24
		dinitrophenol, 2,5-	24
		dinitrophenol, 2,6-	24
		dinitrophenol, 3,4-	24
		diphentoin	11
		diquat	25
		diuron	24
		DNB	23
		DNT	23
		DNT, 2A-	23
		dopamine	14
		doxepine	12

EET, 11,12-	15	hypoxanthine	7	methyl tricarbonyl cyclopentadienyl	
EET, 14,15-	15	hypoxanthine	36	tungsten	39
EET, 5,6-	15	ibuprofen	16	methylamine	34
EET, 8,9-	15	ibuprofen	17	methylbenzene	35
elaïdic acid (phenacyl ester derivative)	37	imipramine	19	methylbenzoate	38
enoxacin	20	IMP	36	methyl-D-glucose,3-O-	27
epinephrine	14	indeno(1,2,3-c,d)pyrene	25	methyldopamine,3-O-	14
ethyl viologen	25	inosine	7	methyl-p-hydroxybenzoate	35
ethylcinnamate	38	inosine	36	methyluric acid, 1-	28
fluconazole	22	insulin	32	methyluric acid, 7-	28
fluoranthene	25	insulin sensitizing agent DRF-2189	15	methylxanthine, 1-	28
fluorene	25	iodouracil, 5-	16	methylxanthine, 3-	28
fluorenyl methylchloroformate, 9-	6	Irganox 1076	26	methylxanthine, 7-	28
fluoro-2'-deoxyuridine, 5-	16	isobutyryl-L-carnitine		metolachlor	24
fluorobenzene	36	1-aminoanthraceneamide	8	metoxuron	24
fluorouridine, 5-	16	isoleucine	7	metribuzin	24
fluoxtetine	12	isoleucine, Fmoc-	6	metronidazole	10
fluvoxamine	11	isoleucine, OPA-	8	metronidazole	13
Fmoc-NH ₂	6	isoleucine, PTC-	6	mianserine	11
Fmoc-OH	6	ketoprofen	16	minocycline	22
formic acid	37	ketoprofen	17	MXN	23
formononetin	39	lactic acid	37	myristic acid (phenacyl ester derivative)	37
foscarnet	18	lactose	27	N-acetylaspargate (NAA)	7
fotemustine	16	lactose, α -	27	N-acetylglutamate (NAG)	7
fructose	27	lactulose	26	nalidixic acid	13
furazolidone	13	lauric acid, fluorescent ester of	37	naphtalene	25
furfural	34	leucin enkephalin	32	naphtalene	35
GABA, OPA-alanine, OPA-	8	leucine	7	naringenin	39
galactose	27	leucine, Fmoc-	6	naringin	39
gastrodin	30	leucine, OPA-	8	NBA, 4-	23
gaulthereside D1	28	leucine, PTC-	6	neburon	24
gaulthereside D2	28	leukotriene LTC ₄	16	NH ₃ ⁺ R	34
gaulthereside D3	28	leukotriene LTD ₄	16	nicotinamide	33
genistein	39	leukotriene LTE ₄	16	nicotine	29
ginkgolide A	29	leukotriene LTE ₄	17	nitrobenzene	36
ginkgolide B	29	leukotriene LTE ₄ -sulfone	17	nitrobenzene	38
ginkgolide C	29	leukotriene LTE ₄ -sulfoxide	17	nitrobenzoic acid, o-	35
glucose	26	leukotriene N-acetyl LTE ₄	17	nitrobenzoic acid, m-	35
glucose	27	linoleic acid (phenacyl ester derivative)	37	nitrobenzoic acid, o-	35
glucosepentabenzoate	38	linolenic acid (phenacyl ester derivative)	37	nitrobenzoic acid, p-	35
glutamate	7	linuron	24	nitrophenol, 2-	24
glutamate, N-acetyl-	7	lomefloxacin	20	nitrophenol, 3-	24
glutamic acid	7	loxapine	14	nitrophenol, 4-	18
glutamic acid, Fmoc-	6	luteolin	39	nitrophenol, 4-	24
glutamic acid, OPA-	8	lysine	7	nitrophenyl- β -D-glucuronide, 4-	18
glutamic acid, PTC-	6	lysine, Fmoc-	6	nitrophenylacetic acid, 4-	28
glutamine, OPA-	8	lysine, PTC-	6	nitrophenylacetoneitrile, 4-	28
glutathione-conjugate (RP 69012)	11	lysine, OPA-	8	nitrotoluene, m-	35
glycine	7	lysozyme	32	nitrotoluene, o-	35
glycine, Fmoc-	6	maltose	27	nitrotoluene, p-	35
glycine, OPA-	8	malvidin-3-O- β -glycoside	9	norepinephrine	14
glycine, PTC-	6	mannitol	26	norfluoxtetine	12
guaiacol	34	mannose	27	norleucine, PTC-	6
H290-51	26	maprotyline	11	normetanephrine	14
H290-39	26	megazol	17	nortriptyline	12
haloperidol	22	melittin	32	nortriptyline	19
heptylpyridine	35	meloxicam	17	nortriptyline, E-10-OH-	12
hesperetin	39	metanephrine	14	nortriptyline, Z-10-OH-	12
HETE, 20-	15	methacycline	22	octyl alcohol, n-	39
hexabenzoyl mannitol	27	methansulfonyl-L-carnitine		octyl caproate, n-	39
hexachlorobenzene	36	1-aminoanthraceneamide	8	ofloxacin	20
hexylpyridine	35	methionine	7	oleic acid (phenacyl ester derivative)	37
histamine	34	methionine enkephalin	32	ornitine, OPA-	8
histidine, Fmoc-	6	methionine enkephalin (ME)	31	oxamyl	24
histidine, N-acetyl-	7	methionine enkephalinamide (MEA)	31	oxytetracycline	22
histidine, OPA-	8	methionine, Fmoc-	6	oxytocin	31
histidine, PTC-	6	methionine, PTC-	6	oxytocin	32
HMX	23	methomyl	24	palmitic acid (phenacyl ester derivative)	37
hydroxybenzyl alcohol, p-	30	methoxybenzylamine hydrochloride,		palmitoleic acid (phenacyl ester derivative)	37
hydroxyl-4-chlorobenzothiazole, 2-	23	4-hydroxy-3-	14	paracetamol	28
hydroxylbenzaldehyde, p-	30			paracetamol	30

paraquat	25	sildenafil (viagra)	20	tyramine	34
paraxanthine	28	simazine	24	Tyr-D-Ala	31
paroxetine	12	sodium phenobarbital	16	Tyr-D-Ala-Gly	31
PAT-5A	18	sodium salicylate	30	Tyr-D-Ala-Gly-Phe	31
pelargonidin-3-O- β -glycoside	9	sotolon	34	tyrosine	7
pentabenzoyl mannitol	27	spermidine	34	tyrosine, Fmoc-	6
pentabenzoyl-D-glucose	27	spermine	34	tyrosine, OPA-	8
pentachlorobenzene	36	sphinganine, OPA-	29	tyrosine, PTC-	6
peonidin-3-O- β -glycoside	9	sphingosine, OPA-	29	UK 103320 (viagra metabolite)	20
petunidin-3-O- β -glycoside	9	spiramycin I	21	UK114542-27	20
Phe	31	spiramycin III	21	uric acid	7
Phe-Met	31	spiramycin S	21	valine	7
phenanthrene	25	stearic acid (phenacyl ester derivative)	37	valine, Fmoc-	6
phenanthrene	35	sucrose	27	valine, OPA-	8
phenobarbitone	11	sulfadimethoxine	10	valine, PTC-	6
phenol	24	sulfamethazine	10	vanillin	30
phenylalanine	7	sulfamethoxazole	10	vanillin	34
phenylalanine, Fmoc-	6	sulfamethoxypyridazine	10	vanillyl alcohol	30
phenylalanine, OPA-	8	sulfaquinoxaline	10	vasopressin, 8-Arg-	31
phenylalanine, PTC-	6	sulfoarea	35	vasopressin, d-4-Asn-8-D-	31
phenylephrine	30	tangeretin	39	vasopressin, d-8-D-Arg-	31
phenylethanol, 2-	34	terbinafine	13	vasotocin, 8-Arg-	31
phenylethyl acetat, 2-	34	tetrabenzoyl-D-glucose	27	viagra	20
phenylpropanolamine	19	tetrachlorobenzene, 1,2,3,4-	36	vitamin A	33
phenylpropanolamine hydrochloride	30	tetrachlorobenzene, 1,2,4,5-	36	vitamin B ₁	33
phosphonoformate	18	tetraethyllead	23	vitamin B ₂	33
piperidine, (-)-trans-4-(4-fluorophenyl)-	12	tetramethyllead	23	vitamin B ₃	33
3-(3-hydroxy-4-methoxyphenoxyethyl)-	12	theobromine	28	vitamin B ₆	33
piperidine, (-)-trans-4-(4-fluorophenyl)-	12	theophylline	11	vitamin C	33
3-(4-hydroxy-3-methoxyphenoxyethyl)-	12	theophylline	28	vitamin E	33
proline, Fmoc-	6	thiabendazole	23	xanthine	7
proline, PTC-	6	thiamine chloride	33	xylazine	22
prometryn	24	thiazolidine-2,4-dione,		xylose	27
propanamine (MDEA), N-ethyl-		5[4-[N-(20pyridyl)-(2s)-pyrrolidine-			
1-(1,3-benzodioxol-5-yl)-2-	15	2-methoxy]phenylmethylene],			
propanamine (MDMA), N-methyl-		maleic acid salt	18		
1-(1,3-benzodioxol-5-yl)-2-	15	thiazolidinedione	18		
propane, 1,3-diamino-	34	threonine, Fmoc-	6		
propazine	24	threonine, OPA-	8		
propionyl-L-carnitine 1-amino-		threonine, PTC-	6		
anthraceneamide	8	tinidazol	17		
propionylpromazine	22	TNB	23		
propyl-p-hydroxybenzoate	35	TNT	23		
protriptyline	12	TNX	23		
pseudoephedrine	30	tocopherol, δ -	33		
purine	36	tocopherol, α -	33		
putrescine	34	tocopherol, γ -	33		
pyrene	25	toluene	19		
pyridoxine	33	toluene	35		
pyridylacetonitrile, 3-	38	toluene	36		
pyrimethamine	10	toluene	38		
quercetin	39	tramadol	22		
quinoxaline	38	triacylglycerols	40		
quinupristin (RP 57669)	11	triamcinolone acetonide acetate	13		
RDX	23	tricarboxyl cyclopentadienyl manganese	39		
resorcinol	24	trichlorobenzene, 1,2,3-	36		
retinol	33	trichlorobenzene, 1,2,4-	36		
rhamnose	27	trichlorobenzene, 1,3,5-	36		
rhamnose, L-	26	triclosan	40		
riboflavine	33	triglycerides from olive oil	40		
ribose	27	trimethoprim	10		
ritonavir	13	trimethylphosphine hydride dicarbonyl			
rutin	39	cyclopentadienyl molybdenum	39		
salicin	20	trimipramine	11		
saquinavir	13	triphenylphosphine hydride dicarbonyl			
SCH 39370 (neutral metallo-		cyclopentadienyl molybdenum	39		
endopeptidase inhibitor)	31	troglitazone	15		
scopolamine	29	tryptamine	34		
serine, Fmoc-	6	tryptophan, PTC-	6		
serine, OPA-	8	tryptophane, Fmoc-	6		
serine, PTC-	6	Tyr	31		

Availability of Kromasil

Kromasil 60 Å bulk packings

Phases	Particle sizes, μm					
	3.5	5	7	10	13	16
SIL	□	■	■	■	■	■
CN	□	■	□	■	□	■

Kromasil 100 Å bulk packings

Phases	Particle sizes, μm					
	3.5	5	7	10	13	16
SIL	■	■	■	■	■	■
C4	■	■	■	■	■	■
C8	■	■	■	■	■	■
C18	■	■	■	■	■	■
NH2	■	■	■	■	■	■
Chiral DMB	□	■	□	■	□	■
Chiral TBB	□	■	□	■	□	■

Kromasil 300 Å bulk packings

Phases	Particle sizes, μm					
	3.5	5	7	10	13	16
SIL	□	■	□	■	□	■
C4	□	■	□	■	□	■
C8	□	■	□	■	□	■
C18	□	■	□	■	□	■

■ = available as standard product □ = please inquire!

Kromasil HPLC columns

Kromasil high pressure slurry-packed columns are available in dimensions from 2.1 mm up to 50.8 mm (2") inner diameter, all columns packed with analytical performance. For detailed information on availability please consult our column brochure, or contact us directly.

The moment you adopt our Kromasil High Performance Concept, you join thousands of chromatographers who share a common goal: to achieve better separations when analyzing or isolating pharmaceuticals or other substances.

Not only will you benefit from our patented silica technology, but you gain a strong partner with a reliable track record in the field of silica products. For the past 60 years, Eka Chemicals has pioneered new types of silica. Our long experience in the field of silica chemistry is the secret behind the development of Kromasil, and the success of our Separation Products Group.

Kromasil is available in bulk, or in high-pressure slurry-packed columns. The development, production and marketing of Kromasil are ISO 9001 certified.

Eka Chemicals is a global company with 3,000 people in 30 countries. It is a business unit within Akzo Nobel, one of the world's largest chemical groups, with more than 67,000 employees in 80 countries.

Eka Chemicals AB

Separation Products,

SE-445 80 Bohus, Sweden.

Telephone + 46 31 58 70 00

Fax + 46 31 58 77 27

For NAFTA countries:

Eka Chemicals

7 Livingstone Avenue, Dobbs Ferry,

NY 10522-3401, U S A.

Telephone + 1 914 674 5019

Fax + 1 914 693 3654

kromasil@eka.com

www.kromasil.com

Kromasil®

