

SUGAR COLUMNS


43



SUCREBEAD I

Sucrosebead I is targeted exclusively to carbohydrate analysis. It is based on styrene/divinylbenzene polymer support. Its mono-dispersion character was given by the unique two-step swelling method used for the polymerization, and is advantageous in chromatographic separation.

- Excellent durability
- Efficient peaks at low pressure




Mono-dispersion particles

Sucrosebead I, in combination with the pulsed amperometric detector (PAD), provides a high-sensitivity carbohydrate analysis.


Monosaccharides

Column	: SUCREBEAD I 2.0 mm i.d. x 250 mm
Mobile Phase	: 200mmol/L NaOH
Flow Rate	: 100 μ L/min
Temperature	: 35°C
Detection	: PAD
Sample	: 1. Fucose 2. Mannose 3. Glucose 4. Galactose 5. Ribose



Sugar alcohols

Column	: SUCREBEAD I 2.0 mm i.d. x 250 mm
Mobile Phase	: 200mmol/L NaOH
Flow Rate	: 100 μ L/min
Temperature	: 35°C
Detection	: PAD
Sample	: 1. myo-Inositol 2. meso-Erythritol 3. Xylitol 4. Sorbitol 5. Mannitol 6. Maltitol



SUCREBEAD II

Sucrosebead II is developed to analyze carbohydrates by using anion-exchanging polymer as a stationary phase. Sucrosebead II enables operation under high pH range and high selectivity with carbohydrates.

Optimum for analyzing oligosaccharides and polysaccharides

