

## SEPARATION OF FREE AMINO ACIDS AND PRIMARY AMINES USING DAICEL CROWN ETHER COLUMNS: CROWNPAK CR-I(+)<sup>+</sup> AND CROWNPAK CR-I(-)<sup>-</sup>

### APPLICATION NOTE

#### INTRODUCTION

Daicel Corporation recently introduced a new generation of CROWNPAK<sup>®</sup> chiral selectors that can be used for the separation of free amino acids and primary amines: CROWNPAK CR-I(+)<sup>+</sup> and CR-I(-)<sup>-</sup>. These crown ether selectors are immobilized on 5- $\mu$ m silica support. Immobilization allows use of organic solvents in a wider range for both reversed-phase and normal-phase chromatography modes, thus, enhancing enantioselective resolution of chiral compounds in a shorter analysis time.

The CROWNPAK CR-I(+)<sup>+</sup> and CR-I(-)<sup>-</sup> chiral selectors are complementary to our CHIRALPAK<sup>®</sup> ZWIX chiral selectors. The complementarity provides a total solution for enantio-recognition of a wide variety of amino acids. For example, CR-I selectors afford challenging separations of asparagine, glutamine and serine.


One important feature of both CHIRALPAK ZWIX and CROWNPAK CR-I chiral selectors is the ability to control the elution order. Typically, use of ZWIX(+)<sup>+</sup> and CR-I(+)<sup>+</sup> columns would lead to the reversal of the elution order for free amino acids eluted from ZWIX(-)<sup>-</sup> and CR-I(-)<sup>-</sup> columns.

#### EXPERIMENTAL AND DISCUSSION

A CROWNPAK CR-I(+)<sup>+</sup> column, 3.0 mm i.d. x 150 mm, packed with 5- $\mu$ m particles was used to develop the separation of DL-serine. The mobile phase was a mixture of perchloric acid and acetonitrile.

The CR-I(+)<sup>+</sup> and CR-I(-)<sup>-</sup> selectors are packed in columns of 3.0 mm i.d. and 150 mm long.

#### SEPARATION OF DL-SERINE ON CROWNPAK CR-I(+)<sup>+</sup>



#### CHROMATOGRAPHIC CONDITIONS

- Column Size:** Daicel CROWNPAK CR-I(+)  
3.0 mm i.d. x 150 mm long, 5- $\mu$ m
- Mobile Phase:** HClO<sub>4</sub> a.q.(pH1.0) / ACN=85/15(v/v)
- Flow Rate:** 0.1 ml/min
- UV Detection:** 200 nm
- Column Temperature:** 25 °C

*Note: Recent scientific studies have demonstrated that the brain of Alzheimer's disease patients contain unusually high levels of D-serine. The potential association of the D-serine level with cognitive decline in the patients may lead to the development of a novel and effective biomarker for early detection of the disease.*



#### CHIRAL TECHNOLOGIES, INC.

800 North Five Points Road  
West Chester, PA 19380 USA  
Tel: 610-594-2100  
Fax: 610-594-2325  
www.chiraltech.com  
email: chiral@chiraltech.com

#### CHIRAL TECHNOLOGIES EUROPE

Parc d'Innovation, Bd Gonther d'Andernach  
67404 Illkirch Cedex, France  
Tel: +33 (0) 388 79 52 00  
Fax: +33 (0) 388 66 71 66  
www.chiraltech.com  
e-mail: cte@chiral.fr

