

Fast Online Analysis of Free Urinary Catecholamines Using MICRA NPS[®] and Restricted Access Affinity Precolumns

MICRA NPS[®] is a breakthrough in fast HPLC. *NPS* is ultra-pure, highly uniform non porous silica spheres which provide the LC chromatographer greatly improved mass transfer and lower detection limits. Coupled with enhanced stability and dramatically reduced solvent usage, *NPS* is the ideal HPLC column to meet the ever increasing demands placed on today's analytical labs - Improved productivity at a lower cost.

This paper was originally presented at HPLC '95, in Innsbruck, Austria entitled; "Coupled-Column LC-Analysis of Catecholamines in Urine using a MICRA NPS Analytical Column" by A. Rudolphi, K.-S. Boos, K. D. Bischoff, and T.J. Barder.


A fast coupled column analysis for routine online determination of free urinary catecholamines has been developed using a MICRA NPS column and a Restricted Access Affinity precolumn^{1,2}. Selective extraction of the catecholamines from urine was performed using alkaline pH and boronic acid-modified copolymers supports. Analytical separation and quantitation was performed with a 4.6 x 53 mm MICRA 1.5 μ m NPS ODS-I column using ion-pair RP HPLC with no organic modifier in the mobile phase.

Figure 1. Separation on NPS 1.5 μ m

LOD's well below 5 ng/mL urine

LC Conditions²:

| | |
|--------------------------|---|
| Mobile phase: | 10 mM NaH ₂ PO ₄ + 0.1 mM 1-dodecanesulfonic acid, pH = 2.5 with H ₃ PO ₄ |
| Flow rate: | 0.5 mL/min |
| Detection: | FD; ex = 275 nm, em = 330 nm, time const. = 3 s, 12 μ L flow cell |
| Injection volume: | 50 μ L |



Coupled-column analysis of norepinephrine (1), epinephrine (2) and dopamine (3)

- (A) 50 μ L aqueous standard solution (1, 2, 3 = 1.99, 1.09, 7.03 nmol/mL)
- (B) 50 μ L Human urine from healthy person (1,2,3 = 1.77,0.73, 14.34 nmol/mL)
- (C) 50 μ L Human urine from pheochromocytoma patient (1,2,3 = 2.16, 2.73, 6.71 nmol/mL)
- (D) 50 μ L blank urine after oxidative destruction of endogenous catecholamines (as described in paper)

METHOD DEVELOPMENT

Think small

Think fast

Think **NPS**[®]

Table 1. Linearity & Precision

| Analyte | LINEARITY | | PRECISION | |
|---------|------------------|---------------------|-----------------------------|-----------------------------|
| | Range pmol/mL | corr. coeff. [r] | CV within run (n=10) [%] | CV within day (n=10) [%] |
| NE | 48.9 - 9784 | 0.9994 | 0.97 | 1.68 |
| E | 54.6 -10916 | 0.9983 | 2.25 | 1.90 |
| DA | 105.5 - 8701 | 0.9997 | 3.74 | 2.51 |

Table 2. LOQ & LOD

| Analyte | LOQ | LOD | Mean Recoveries, % |
|---------|----------------|----------------|--------------------|
| | pmol on column | pmol on column | nmol/mL [%RSD] |
| NE | 5.57 | 2.43 | 98.18 [2.8] |
| E | 1.75 | 0.76 | 102.00 [1.8] |
| DA | 36.81 | 15.61 | 101.12 [1.8] |

NE = norepinephrine, E = epinephrine, DA = dopamine

MICRA NPS Performance

Combined with the precolumn, the MICRA NPS column showed:

1. No degradation in performance over the life of the study (5 months);
2. 50% reduction in run time compared to a typical porous column assay³;
3. >70% reduction in solvent consumption³; and
4. Quantitative recovery of catecholamines.

Eprogen
8205 S. Cass Avenue
Suite #106
Darien, Illinois 60561
USA
Phone: (630) 963-1481
Fax: (630) 963-6432
E-mail: info@eprogen.com
www.eprogen.com

References:

- 1 Boos, K.-S., et. al., *J. Chromatogr.* (1988) **456**, 83-104
- 2 Rudolphi, A., et. al., *Chromatographia* (1995) **41**, 645-650
- 3 K.-S. Boos, et. al., *Chromatographia* (1987) **24**, 363-370

#AP17

Rev. 040309
©2009 Eprogen