

# HALO<sup>®</sup>

## C30


ISOMERS HAVE MET THEIR MATCH



# HALO<sup>®</sup> C30

## INTRODUCING THE NEW HALO<sup>®</sup> C30!

Built on proven Fused-Core<sup>®</sup> particle technology, the HALO<sup>®</sup> C30 is designed to deliver fast separations ideal for lipids and isomers compared to your C18.



## FEATURES OF HALO<sup>®</sup> C30


- High shape selectivity for hydrophobic, long-chain structurally related isomers
- Highly reproducible bonded phase coverage resulting in fast, highly efficient, rugged separations
- 100 % Aqueous Compatibility

### Best Applications:

- Fat/Water Soluble Vitamins
- Carotenoids
- Lipids
- Steroids

## QUALITY YOU CAN COUNT ON

The HALO<sup>®</sup> C30 exhibits reliable stability providing you assurance in your separations – injection to injection.




### TEST CONDITIONS

Column: HALO 160 Å C30, 2.7µm, 2.1 x 150mm  
Isocratic: 50-50 ACN/H<sub>2</sub>O  
Flow Rate: 1.1 mL/min.  
Back Pressure: 602 bar  
Temperature: 60°C  
Injection Volume: 1 µl test mix containing uracil, phenol, 1-chloro-4-nitrobenzene, and naphthalene  
Instrument: Shimadzu Nexera  
Detection: UV 254 nm, PDA

## ISOMERS HAVE MET THEIR MATCH

Tocopherols, a form of vitamin E, known for their antioxidant properties contain many isomers which are easily baseline separated on the HALO® C30 compared to a C18 due to the shape selectivity property of the phase.




### PEAK IDENTITIES

1. δ-tocopherol
2. γ-tocopherol
3. β-tocopherol
4. α-tocopherol

### TEST CONDITIONS

Columns: HALO 160 Å C30 and 90 Å C18, 2.7 μm, 4.6 x 150mm  
 Mobile Phase A: Water  
 Mobile Phase B: Methanol  
 Isocratic: 95% B  
 Flow Rate: 1.5 mL/min  
 Temperature: 10°C  
 Injection Volume: 1.5 μL  
 Instrument: Agilent 1200 SL  
 Detection: UV 290 nm, PDA


### STRUCTURE



Tocopherol	R1	R2
Alpha (α)	CH <sub>3</sub>	CH <sub>3</sub>
Beta (β)	CH <sub>3</sub>	H
Gamma (γ)	H	CH <sub>3</sub>
Delta (δ)	H	H

## CAROTENOID SEPARATION

These seven carotenoids and related isomers from a commercially available vitamin formulation are easily resolved under simple mobile phase conditions.



### TEST CONDITIONS

Column: HALO 160 Å C30, 2.7 μm, 3.0 x 150 mm  
 Mobile Phase A: Methanol  
 Mobile Phase B: Ethanol  
 Gradient: 100% A with gradient to 40% B at 20 min.  
 Flow Rate: 0.65 mL/min  
 Temperature: 38°C  
 Data acquisition rate: 2.5 Hz  
 Injection volume: 0.60 μL  
 Instrument: Agilent 1100  
 Detection: UV 471 nm, PDA


### PEAK IDENTITIES

1. Lutein
2. cis-carotenoid 1
3. cis-carotenoid 2
4. α-carotene
5. β-carotene
6. cis-lycopene
7. Lycopene

Data courtesy of Nature's Sunshine Products.

## STEROID SEPARATIONS USING HALO® C30

Glucocorticoids are a powerful form of steroid hormone both naturally produced and prescribed as a means to control a number of conditions caused by inflammation. In the example below six glucocorticoids are separated within three minutes and demonstrate excellent peak shape.



### TEST CONDITIONS


Column: HALO 160 Å C30, 2.7 μm, 4.6 x 150 mm  
 Mobile Phase A: Water  
 Mobile Phase B: 50/50 Acetonitrile/ Methanol  
 Isocratic: 50% B  
 Flow Rate: 1.5 mL/min  
 Back Pressure: 309 bar  
 Temperature: 60°C  
 Injection Volume: 0.5 μL  
 Instrument: Shimadzu Nexera X2  
 Detection: UV 220 nm, PDA

### PEAK IDENTITIES

1. Prednisone
2. Cortisone
3. Prednisolone
4. Hydrocortisone
5. Dexamethasone
6. Corticosterone

## ENHANCED SELECTIVITY AND RESOLUTION FOR FAT SOLUBLE VITAMINS

For a panel of fat soluble vitamins, the HALO® C30 demonstrates superior resolution and selectivity in comparison with a commonly employed RP-Amide phase.



### PEAK IDENTITIES


1. Retinyl acetate (A)
2. Delta tocopherol (E)
3. Ergocalciferol (D2)
4. Cholecalciferol (D3)
5. Alpha tocopherol (E)
6. DL-alpha-tocopherol acetate (E)
7. 2,3-trans-phyloquinone (K)
8. Retinyl palmitate (A)

### TEST CONDITIONS

Columns: HALO C30, RP-Amide, 2.7  $\mu\text{m}$ , 4.6 x 150 mm  
 Isocratic: 100% Methanol  
 Flow Rate: 1.5 mL/min.  
 Temperature: 30 °C  
 Injection: 2  $\mu\text{l}$   
 Instrument: Nexera 062  
 Detection: UV 280 nm, PDA

## COMPETITIVE ADVANTAGE

A HALO® C30 column shows increased retention and baseline resolution of vitamin K1 trans and cis isomers compared to a SPP C30 column. Since the cis isomer of K1 is biologically inactive, it is important to know how much of each isomer is present in vitamin enriched products.



### PEAK IDENTITIES

1. Menadione (K3)
2. Menaquinone 4 (K2)
3. 2,3-trans-phyloquinone (K1)
4. cis-phyloquinone (K1)

### TEST CONDITIONS

Columns: HALO 160 Å C30, 2.7  $\mu\text{m}$ , 4.6 x 150 mm; SPP C30, 2.6  $\mu\text{m}$ , 4.6 x 150 mm  
 Mobile Phase A: Water  
 Mobile Phase B: Methanol  
 Isocratic: 95% B  
 Flow Rate: 1.5 mL/min  
 Temperature: 25 °C  
 Injection Volume: 1  $\mu\text{l}$   
 Instrument: Shimadzu Nexera  
 Detection: UV 280 nm, PDA



# SPECIFICATIONS

Ligand: Triacetyldimethyl  
 Particle Size: 2.7  $\mu\text{m}$   
 Pore Size: 160  $\text{\AA}$

USP Designation: L62  
 Carbon Load: 4.5 %  
 Surface Area: 90  $\text{m}^2/\text{g}$

Endcapped: Yes  
 Low pH Limit /Max T: 2/60°C  
 High pH Limit/Max T: 9/40°C

## PART NUMBERS


ANALYTICAL COLUMNS	
Dimensions: ID x Length (in mm)	Part Number
2.1 x 20	92112-230
2.1 x 30	92112-330
2.1 x 50	92112-430
2.1 x 75	92112-530
2.1 x 100	92112-630
2.1 x 150	92112-730
2.1 x 250	92112-930
3.0 x 20	92113-230
3.0 x 30	92113-330
3.0 x 50	92113-430
3.0 x 75	92113-530
3.0 x 100	92113-630
3.0 x 150	92113-730
3.0 x 250	92113-930
4.6 x 20	92114-230
4.6 x 30	92114-330
4.6 x 50	92114-430
4.6 x 75	92114-530
4.6 x 100	92114-630
4.6 x 150	92114-730
4.6 x 250	92114-930
10.0 x 50	92110-402
10.0 x 75	92110-502
10.0 x 100	92110-602
10.0 x 150	92110-702

CAPILLARY COLUMNS	
Dimensions: ID x Length (in mm)	Part Number
0.075 x 50	91219-402
0.075 x 100	91219-602
0.075 x 150	91219-702
0.1 x 50	91218-402
0.1 x 100	91218-602
0.1 x 150	91218-702
0.2 x 50	91217-402
0.2 x 100	91217-602
0.2 x 150	91217-702
0.3 x 50	91216-402
0.3 x 100	91216-602
0.3 x 150	91216-702
0.5 x 50	91215-402
0.5 x 100	91215-602
0.5 x 150	91215-702
1.0 x 30	92111-302
1.0 x 50	92111-402
1.0 x 75	92111-502
1.0 x 100	92111-602
1.0 x 150	92111-702

GUARD COLUMNS	
Guard columns, 3-pack	
Dimensions: ID x Length (in mm)	Part Number
2.1 x 5	92112-130
3.0 x 5	92113-130
4.6 x 5	92114-130
Guard Column Holder	94900-001



# HALO<sup>®</sup>



advancedmaterialstechnology

FIND YOUR DISTRIBUTOR AT

[advanced-materials-tech.com](http://advanced-materials-tech.com)