

#### Abstract

There is considerable interest in HPLC columns for performing very fast separations either for increasing sample throughput or for use in multidimensional (2-D) separations. One approach for achieving fast separations is to use very small particles for the column packing. However, columns of sub-2-µm particles must be operated at very high back pressures for expected efficiency. In some cases special HPLC apparatus is required for optimum performance with these very small particles. Unique silica-based particles now have been developed for packed HPLC columns that allow ultrafast separations at operating pressures available with most existing instruments. These new "fused core" particles are 2.7- $\mu$ m overall with a 1.7- $\mu$ m solid core and a 0.5- $\mu$ m porous outer shell of ~9 nm pores. This configuration provides a fortunate compromise of characteristics that result in very rapid separations with modest operating pressures. Columns of these new particles exhibit unusual efficiency, with reduced plate heights of 1.5 for small molecules under optimum conditions (> 200,000 plates/ meter). This high level of efficiency is believed to be a result of the very narrow particles size distribution (5%, 1-sigma from the average) and higher particle density because of the solid core. Relative to totally porous sub-2-µm particles, columns of the new particles generate a competitive number of theoretical plates in the same time frame, but with greatly reduced back pressures. Compared to columns of sub-2-µm particles, the short diffusion paths of the thin porous crust and reduced back pressure of the 2.7-µm particles allow mobile phase velocities at the plate height minimum for highest efficiency per unit time while staying within the 400 bar pressure limit of most HPLC instruments. Different from sub-2-µm particles, the 2.7-µm particles can be packed into stable columns of expected efficiency. These 2.7-µm particles are more easily packed efficiently and require inlet frits or retainers that have larger pores that are less likely to be plugged with use. Allowed longer columns with the new particles also exhibit reduced problems with extra-column band broadening effects. Surface areas of ~150 m2/g for the new fused core particles are suitable for separating small molecules with good sample loading characteristics. The new particles exhibit the strength and stability required for operation at pressures of at least 600 bar.

#### Study Purpose

Develop unique high-purity silica particles for stable high-efficiency packed columns that are capable of very fast separations with modest back pressures.






### Specifications for Fused-core Particles

- 2.7 µm overall diameter Highly purified Type
- 1.7 µm solid core
- 0.5-µm-thick porous outer shell

B silica


• 9 nm pores in porous shell 600 bar




# New Fused-Core<sup>TM</sup> Particles for Very Fast HPLC Separations


Joseph J. Kirkland, Timothy J. Langlois & Joseph J. DeStefano Advanced Materials Technology, 3521 Silverside Road, Suite 1-K, Quillen Bldg. Wilmington, DE 19810


• Surface area ~150 m2/g Pressures to at least













Ascentis C8 (3.1µm) ; HALO C8 and HALO C18 (2.7µm) Columns: 50 x 4.6 mm;Temperature: 30°C; Agilent 1100 Mobile phase: 35% ACN/65% potassium phosphate buffer, pH 3.0 Data fitted to Knox equation

> Acentis C8 Halo C18 Halo C8 3 4 Mobile Phase Velocity, mm/sec

Separation of Explosives **Column:** 2.1  $\times$  50 mm Halo C18 Mobile Phase: 27% Methanol / 73% Water Temperature: 40°C Flowrate: 0.85 ml/min Detector: UV @ 254 nm Pressure: 310 bar Instrument: Agilent 1100 8.4-amino-2,6-dinitrotoluene 9. 2-amino-4,6-dinitrotoluene **10.** 2,6-dinitrotoluene II.2,4-dinitrotoluene 12.2-nitrotoluene **13.** 4-nitrotoluene 14.3-nitrotoluene

12 13 14 0 0.5 1 1.5 2 2.5 3 3.5 Time, minutes





## Conclusions Regarding Columns of 2.7-µm Fused-core Particles

- Produce unusual column efficiency; reduced plate heights *h* of ~1.5 for small molecules
- Efficiency and separation speed rivals that of sub-2-µm columns, but with about one-half the pressure drop
- Thin porous outer shell results in superior mass transfer kinetics and better efficiency at high mobile phase velocities, especially for larger molecules
- 2-µm porosity frits allow the convenience and ruggedness of columns with 5 µm particles
- Sample loading competitive with conventional totally porous particles
- Strong particles allow pressures of at least 600 bar