


## > Selection Guide according to the "Neutral" Tanaka Test



N°	Manufacturer - Brand	N°	Manufacturer - Brand	N°	Manufacturer - Brand	N°	Manufacturer - Brand
1	SAF - Supelco ABZ+	11	Thermo Scientific - Hypersil Gold	21	Phenomenex - Gemini	31	Phenomenex - Synergy Fusion RP
2	Waters - Xbridge Shield	12	Waters - Symetry C18	22	Phenomenex - Synergi MaxRP 4µm	32	Interchim - Uptisphere MM1
3	Waters - Xterra C18	13	Phenomenex - Prodigy ODS 3	23	Grace - Altima HP C18	33	Macherey Nagel - Nucleosil HD
4	Interchim - Uptisphere TF	14	Interchim - Uptisphere Strategy C18-2	24	Interchim - Uptisphere ODB	34	Interchim - Uptisphere HDO
5	Interchim - Uptisphere Strategy C12	15	Phenomenex - Luna C18(2)	25	Macherey Nagel - Nucleodur C18 Pyramid	35	YMC - Pack ODS-AQ
6	Waters - Xbridge C18	16	Waters - Xterra MS	26	Merck - Superspher E	36	Phenomenex - Synergy Hydro RP
7	Macherey Nagel - Nucleodur C18 Gravity	17	Thermo Scientific - Hypersil Gold AQ	27	GL Sciences - InertSill ODS 3	37	Agilent - Zorbax RX C18
8	Agilent - Zorbax Extend C18	18	Interchim - Uptisphere Strategy C18-3	28	YMC - PRO C18	38	Merck - Purospher RP18
9	Eka Nobel - Kromasil C18	19	Interchim - Uptisphere PLP	29	Interchim - Uptisphere Strategy RP	39	Merck - Superspher 100A RP18
10	Agilent - Zorbax Bonus RP	20	Agilent - Zorbax Eclipse XDB	30	Interchim - Uptisphere HSC	40	Interchim - Uptisphere NEC

Tanaka Test at 30°C

## > Selection Guide for Biomolecules

Uptisphere	Small molecules & Bio-Drugs	Weak hydrophobic peptides	Very hydrophobic peptides	Insulin	Metabolites isolation	Weak hydrophobic proteins >15 kDa	Moderately hydrophobic proteins >15 kDa	Very hydrophobic proteins >15 kDa
X-Serie 100 Å	OD2							
X-Serie 200 Å	C18	C4	C8	C18AQ				
300 Å					WOD	WOD	WC4	WC4



211 bis avenue J.F. Kennedy  
BP 1140 - 03100 Montluçon  
France  
Tél +33 4 70 03 88 55  
Fax + 33 4 70 03 82 60  
[www.interchim.com](http://www.interchim.com)


+33 4 70 03 73 09 [interchrom@interchim.com](mailto:interchrom@interchim.com) [www.interchim.com](http://www.interchim.com)

USA  
Interchim® Inc  
1536 West 25th Street, Suite 452  
San Pedro CA 90732  
Tel. 800 560 8262  
Fax 310 802 3878  
[web.interchiminc.com](http://web.interchiminc.com)

UK  
Cheshire Sciences UK Ltd  
67 Bridge Street Rows.  
Chester - CH1 1NW  
Tel. 01244 620486  
Fax 01244 620378  
[germany@cheshiresciences.com](mailto:germany@cheshiresciences.com)

Benelux  
Belgium  
Netherlands  
[benelux@interchim.com](mailto:benelux@interchim.com)

Germany & Austria  
Germany  
Austria  
[germany@interchim.com](mailto:germany@interchim.com)



# UHPLC, HPLC & prep LC Columns

Uptisphere®  
Uptisphere® Strategy™  
Uptisphere® 300 Å  
Uptisphere® X-Serie



ANALYTICAL SCIENCES ITM-2013/06/21-C/P



Hotline +33 4 70 03 73 09 • [interchrom@interchim.com](mailto:interchrom@interchim.com) • [www.interchim.com](http://www.interchim.com)

## > UHPLC, HPLC & prep LC Columns

### Uptisphere®

HPLC & prep LC columns for the identification, separation & purification of small moléculles & pharma compounds.

### Uptisphere® Strategy™

UHPLC, Analytical & prep LC columns with high loadability for the identification, separation & purification of small moléculles & pharma compounds.

### Uptisphere® X-serie™


HPLC & prep LC columns for the identification, analysis & purification of small molecules & bio-médicaments at high & low pH.

### Uptisphere® 300 Å

HPLC & prep LC columns for identification, analysis & purification of Proteins, Peptides & Polypeptides.

Description	UHPLC		HPLC		
	1.7µm	2.2µm	3µm	5µm	10µm
<b>Small Molecules</b>					
Uptisphere® Strategy™ C18-3	⊕	⊕	⊕	⊕	⊕
Uptisphere® Strategy™ C18-2	⊕	⊕	⊕	⊕	⊕
Uptisphere® Strategy™ RP	⊕	⊕	⊕	⊕	⊕
Uptisphere® Strategy™ PLP	⊕	⊕	⊕	⊕	⊕
Uptisphere® Strategy™ C12	⊕	⊕	⊕	⊕	⊕
Uptisphere® Strategy™ RPX	⊕	⊕	⊕	⊕	⊕
Uptisphere® Strategy™ C8-2	⊕	⊕	⊕	⊕	⊕
Uptisphere® Strategy™ HILIC	⊕	⊕	⊕	⊕	⊕
Uptisphere® Strategy™ Si	⊕	⊕	⊕	⊕	⊕
Uptisphere® Strategy™ NH2	⊕	⊕	⊕	⊕	⊕
Uptisphere® HSC	⊕	⊕	⊕	⊕	⊕
Uptisphere® ODB	⊕	⊕	⊕	⊕	⊕
Uptisphere® HDO	⊕	⊕	⊕	⊕	⊕
Uptisphere® NEC	⊕	⊕	⊕	⊕	⊕
Uptisphere® TF	⊕	⊕	⊕	⊕	⊕
Uptisphere® PLP	⊕	⊕	⊕	⊕	⊕
Uptisphere® C8	⊕	⊕	⊕	⊕	⊕
Uptisphere® PH	⊕	⊕	⊕	⊕	⊕
Uptisphere® C4	⊕	⊕	⊕	⊕	⊕
Uptisphere® CN	⊕	⊕	⊕	⊕	⊕
Uptisphere® Diol	⊕	⊕	⊕	⊕	⊕
Uptisphere® Si	⊕	⊕	⊕	⊕	⊕
Uptisphere® MM1	⊕	⊕	⊕	⊕	⊕
Uptisphere® NH2	⊕	⊕	⊕	⊕	⊕
Uptisphere® SCX	⊕	⊕	⊕	⊕	⊕
Uptisphere® SAX	⊕	⊕	⊕	⊕	⊕
<b>Biochromatography</b>					
Uptisphere® X-serie OD2	⊕	⊕	⊕	⊕	⊕
Uptisphere® X-serie C18	⊕	⊕	⊕	⊕	⊕
Uptisphere® X-serie AQ	⊕	⊕	⊕	⊕	⊕
Uptisphere® X-serie C8	⊕	⊕	⊕	⊕	⊕
Uptisphere® X-serie C4	⊕	⊕	⊕	⊕	⊕
Uptisphere® 300 Å WOD	⊕	⊕	⊕	⊕	⊕
Uptisphere® 300 Å WC4	⊕	⊕	⊕	⊕	⊕
Uptisphere® 300 Å WD4	⊕	⊕	⊕	⊕	⊕
Uptisphere® 300 Å WT4	⊕	⊕	⊕	⊕	⊕

## > Selection Guide for Small Molecules



## > Stationary Phases & Features

**Strategy™ C18-3**

100 Å - 425 m<sup>2</sup>/g  
3. 5, 10 & 15 µm  
Bonding: C18 %C: 22  
End-capping: multi-step  
pH stability: 1.0 to 12.0  
The high bonding density of C18-3 facilitates a stronger separation of non polar compounds. Multi step bonding technology guarantees a fully end-capped phase, stable under basic pH conditions. C18-3 is the excellent phase for the integral separation of basic drugs up to pH: 12.

**Uptisphere® C18-ODB**

120 Å - 320 m<sup>2</sup>/g  
3. 5 & 10 µm  
Bonding: C18 %C: 18  
End-capping: one-step  
pH stability: 1.5 to 7.0  
Serves a broad-spectrum of analytical requirements for separating non polar compounds.

**Strategy™ C18-2**

100 Å - 425 m<sup>2</sup>/g  
1.7, 2.2, 3 & 5 µm  
Bonding: C18 %C: 19  
End-capping: multi-step  
pH stability: 1.0 to 10.0  
This utility phase serves many pharmaceutical applications. Its 425 m<sup>2</sup>/g surface area providing excellent loading capacity capabilities.

**Uptisphere® C18-TF**

120 Å - 320 m<sup>2</sup>/g  
3, 5 & 10 µm  
Bonding: C18 %C: 14  
End-capping: one-step  
pH stability: 1.5 to 8.0  
Alternative selectivity for challenging separations & for - but not limited to - aromatic, polyphenol, PAHs etc.

**Strategy™ RP**

100 Å - 425 m<sup>2</sup>/g  
3, 5, 10 & 15 µm  
Bonding: C18 %C: 16  
End-capping: one-step hydrophilic  
pH stability: 1.5 to 8.0  
Suitable for mid & non polar compound separation. RP shows excellent mechanical stability under 100% aqueous mobile phase condition.

**Strategy™ C12**

100 Å - 425 m<sup>2</sup>/g  
3 & 5 µm  
Bonding: C12 %C: 16  
End-capping: one-step  
pH stability: 1.5 to 8.0  
Less retentive than C18 with greater capacity

**Uptisphere® C8**

120 Å - 320 m<sup>2</sup>/g  
3 & 5 µm  
Bonding: C8 %C: 11  
End-capping: none  
pH stability: 1.5 to 6.5  
NEC strongly retains the polar and mid-polar compounds. It overcomes peak tailing with compounds that contain chains and/or carbon cycles combined with numerous polar

Brand Name	Code	USP Code	Pore size	Surface Area	Bonding	Bonding Type	Carbon content	End-Capping	pH range	Utilization Mode
Uptisphere® Strategy™	C18-3	L1	100 Å	425 m <sup>2</sup> /g	C18 - octadecyl	mono-fonctionnel	22%	Multi step	1 - 12	Reverse
Uptisphere® Strategy™	C18-3	L1	100 Å	425 m <sup>2</sup> /g	C18 - octadecyl	mono-fonctionnel	19%	Multi step	1 - 10	Reverse
Uptisphere® Strategy™	RP	L1	100 Å	425 m <sup>2</sup> /g	C18 - octadecyl	mono-fonctionnel	16%	one step hydrophilic	1.5 - 8	Reverse
Uptisphere® Strategy™	NEC	L1	100 Å	425 m <sup>2</sup> /g	C18 - octadecyl	mono-fonctionnel	18%		1.5 - 7	Reverse
Uptisphere® Strategy™	PLP	L60	100 Å	425 m <sup>2</sup> /g	Alkyl chain w/ polar group embedded	poly-fonctionnel	14%	Multi step	2.5 - 7.5	Reverse
Uptisphere® Strategy™	C12	L1	100 Å	425 m <sup>2</sup> /g	C12 - dodecyl	mono-fonctionnel	16%	One step	1.5 - 8	Reverse
Uptisphere® Strategy™	RPX	L1	100 Å	425 m <sup>2</sup> /g	Proprietary	mono-fonctionnel			1.5 - 7	Reverse
Uptisphere® Strategy™	C8-2	L7	100 Å	425 m <sup>2</sup> /g	C8 - octyl	mono-fonctionnel	14%	One step	1.5 - 7	Reverse
Uptisphere® Strategy™	HILIC	L3	100 Å	425 m <sup>2</sup> /g	Proprietary				1.5 - 7	Hilic
Uptisphere® Strategy™	SI	L3	100 Å	425 m <sup>2</sup> /g	Silice				1.5 - 7	Normal
Uptisphere® Strategy™	NH2	L8	100 Å	425 m <sup>2</sup> /g	NH2 - amino	mono-fonctionnel	4%		1.5 - 7	Reverse / Normal / Weak Ion exchange
Uptisphere®	HSC	L1	n.a.	n.a.	C18 - octadecyl	mono-fonctionnel	20%	Multi step	1.5 - 8	Reverse
Uptisphere®	ODB	L1	120 Å	320 m <sup>2</sup> /g	C18 - octadecyl	mono-fonctionnel	18%	One step	1.5 - 7	Reverse
Uptisphere®	HDO	L1	120 Å	320 m <sup>2</sup> /g	C18 - octadecyl	mono-fonctionnel	17%	one step hydrophilic	1.5 - 7	Reverse
Uptisphere®	NEC	L1	120 Å	320 m <sup>2</sup> /g	C18 - octadecyl	mono-fonctionnel	16%		1.5 - 6.5	Reverse
Uptisphere®	TF	L1	n.a.	n.a.	C18 - octadecyl	poly-fonctionnel	14%	One step	1.5 - 8	Reverse
Uptisphere®	PLP	L60	120 Å	320 m <sup>2</sup> /g	Alkyl chain w/ polar group embedded	poly-fonctionnel	14%	Multi step	2.5 - 7.5	Reverse
Uptisphere®	PAH	L1	n.a.	n.a.	C18 - octadecyl	poly-fonctionnel	12%		1.5 - 7	Reverse
Uptisphere®	C8	L7	120 Å	320 m <sup>2</sup> /g	C8 - octyl	mono-fonctionnel	11%	One step	2 - 7	Reverse
Uptisphere®	C8U	L7	120 Å	320 m <sup>2</sup> /g	C8 - octyl	mono-fonctionnel	7%		2 - 6.5	Reverse
Uptisphere®	MM1	L44	120 Å	320 m <sup>2</sup> /g	mixed mode RP/SCX	mono-fonctionnel			2 - 6.5	Reverse / Ion exchange
Uptisphere®	C4	L26	120 Å	320 m <sup>2</sup> /g	C4 - butyl	mono-fonctionnel	7%	One step	2 - 7	Reverse
Uptisphere®	CN	L10	120 Å	320 m <sup>2</sup> /g	CN - cyano	mono-fonctionnel	8%	One step	2 - 7	Reverse / Normal
Uptisphere®	PH	L11	120 Å	320 m <sup>2</sup> /g	PH - phenyl	mono-fonctionnel	9%	One step	2 - 7	Reverse / Normal
Uptisphere®	DNAP	L1	120 Å	320 m <sup>2</sup> /g	DNAP - dinitroanilido phenyl	mono-fonctionnel			2 - 6.5	Reverse / Normal
Uptisphere®	OH	L20	120 Å	320 m <sup>2</sup> /g	OH - diol	mono-fonctionnel	6%		2 - 6.5	Reverse / Normal
Uptisphere®	SI	L3	120 Å	320 m <sup>2</sup> /g	Silice				1.5 - 6.5	Normal
Uptisphere®	NH2	L8	120 Å	320 m <sup>2</sup> /g	NH2 - amino	mono-fonctionnel	5%		2 - 6.5	Reverse / Normal / Weak Ion exchange
Uptisphere®	SCX	L50	120 Å	320 m <sup>2</sup> /g	SCX - strong cation exchanger	mono-fonctionnel			1 - 7.5	Ion exchange
Uptisphere®	SAX	L14	120 Å	320 m <sup>2</sup> /g	SAX - strong anion exchanger	mono-fonctionnel			1 - 7.5	Ion exchange
Uptisphere® X-Serie	OD2	L1	n.a.	n.a.	C18 - octadecyl	poly-fonctionnel type II		Multi-step	1 - 13	Reverse
Uptisphere® X-Serie	C18	L1	n.a.	n.a.	C18 - octadecyl	poly-fonctionnel type II		Multi-step	1 - 13	Reverse
Uptisphere® X-Serie	C8	L7	n.a.	n.a.	C8 - octyl	poly-fonctionnel type II		Multi-step	1 - 13	Reverse
Uptisphere® X-Serie	C4	L26	n.a.	n.a.	C4 - butyl	poly-fonctionnel type II		Multi-step	1 - 13	Reverse
Uptisphere® 300Å	WOD	L1	300 Å	100 m <sup>2</sup> /g	C18 - octadecyl	mono-fonctionnel	10%	One step	1.5 - 7	Reverse
Uptisphere® 300Å	PXP	L1	300 Å	100 m <sup>2</sup> /g	C18 - octadecyl	poly-fonctionnel type I	8%	One step	1 - 10	Reverse
Uptisphere® 300Å	WC8	L7	300 Å	100 m <sup>2</sup> /g	C8 - octyl	mono-fonctionnel	8%	One step	2 - 7	Reverse
Uptisphere® 300Å	WD8	L7	300 Å	100 m <sup>2</sup> /g	C8 - octyl	poly-fonctionnel type I	8%	One step	1.5 - 8	Reverse
Uptisphere® 300Å	WC4	L26	300 Å	100 m <sup>2</sup> /g	C4 - butyl	mono-fonctionnel	4%	One step	2 - 7	Reverse
Uptisphere® 300Å	WD4	L26	300 Å	100 m <sup>2</sup> /g	C4 - butyl	poly-fonctionnel type I	4%	One step	1.5 - 8	Reverse
Uptisphere® 300Å	WT4	L26	300 Å	100 m <sup>2</sup> /g	C4 - butyl	tri-fonctionnel	3%	One step	1.5 - 8	Reverse

The comparative data used in this brochure may not be representative of all applications

**Uptisphere® MM1**

120 Å - 320 m<sup>2</sup>/g  
5 µm  
Bonding: mixed mode - RP/SCX  
End-capping: none  
pH stability: 2.0 to 6.5  
Ion exchange & hydrophobic chains are bonded onto the surface of the silica providing unique selectivity. Compounds that possess basic functionality are retained by ion exchange functionality. An appropriate pH will elute ionizable molecules. An organic solvent elutes retained compounds that result from hydrophobic bonding.

**Uptisphere® NH2**

120 Å - 320 m<sup>2</sup>/g  
3, 5, 10 & 15 µm  
Bonding: amino %C: 5  
End-capping: none  
pH stability: 2.0 to 6.5  
Can be either weak anion exchangers for strong acids, or polar media that can interact with OH, NH, SH ...

**Uptisphere® SCX**

120 Å - 320 m<sup>2</sup>/g  
5 & 10 µm  
Bonding: Strong Cation Exchanger - SCX  
End-capping: none  
pH stability: 1.0 to 7.5  
Strong cation exchange (SCX) contains sulfonic acid used to extract weak basic compounds which have one or more positive charges.

**Uptisphere® PLP**

120 Å - 320 m<sup>2</sup>/g  
3 & 5 µm  
Bonding: Alkyl chain w/ polar group embedded %C: 14  
End-capping: multi-step  
pH stability: 2.5 to 7.5  
100% compatibility with aqueous mobile phases. Excellent peak symmetry with basic compounds. Very good retention of mid & non polar compounds combined with numerous polar

**Strategy™ PLP**

100 Å - 425 m<sup>2</sup>/g  
2.2 µm  
Bonding: Alkyl chain w/ polar group embedded %C: 19  
End-capping: multi-step  
pH stability: 2.5 to 7.5

**Uptisphere® Si**

120 Å - 320 m<sup>2</sup>/g  
3 & 5 µm  
Bonding: Polar link C<sub>n</sub>H<sub>2n+1</sub>  
pH stability: 1.5 to 6.5  
High Flow rate &

**Strategy™ Si**

100 Å - 425 m<sup>2</sup>/g  
2.2, 3, 5 & 10 µm  
pH stability: 1.5 to 6.5  
Higher loading capacity vs. Uptisphere

## > Quality Control Report

**Chromatography Report**

Part Number : US3C182-1000648  
mfg # 130852 Lot C18(2)100-3-1240

**Column Information**

Mobile Phase : Acetonitrile  
Bonding: C18-2  
Pore Size (µm) : 3  
Length (mm) : 100  
I.D. (mm) : 4.6

**Test Conditions**

Mobile Phase (ml/min) : 1.0  
Temperature (°C) : 22  
UV (nm) : 254  
Sample Volume (µL) : 10

**SHIPPING / STORAGE SOLVENT :** Acetonitrile

**Résultats d'intégration**

#	Peak Name	Rt	Tailing	Plates (USP)	Resolution (USP)
1	Urethane	0.85	1.23	15021.73	29.00
2	Toluene	2.07	1.23		
3	Naphthalene	3.43	1.07	15905.55	5.50

**SOMME**

Accepted by :

**interchim** www.interchim.com

Technical Support Hot Line : 04 70 03 73 09  
03103 Montluçon Cedex - France  
Tel: 33 (0)4 70 03 88 85  
Fax: 33 (0)4 70 03 82 60  
e-mail: interchim@interchim.com

31-33 rue de Neuilly  
92110 Clichy sur Seine - France  
Tel: 33 (0)1 41 32 34 40  
Fax: 33 (0)1 41 25 15 90  
e-mail: interchim-paris@interchim.com

## > High Performance Hardware

