

A new Kromasil product, Kromasil Diol, has been developed. This new product has been developed in order to achieve high ligand coverage, and the derivatization utilizing a tri-functional silane has been optimized in order to achieve the same batch to batch reproducibility as with a mono-functional silane. Characterization has been performed using solid state NMR, elemental analysis and chromatographic evaluation.

Background


For normal phase chromatography bare silica is the major product sold on the market. However, alternatives with different selectivity and retentivity for polar compounds have been of interest for many years. In certain cases the aim is also enhanced solubility, which can be achieved by introducing alternative mobile phases. Very polar impurities can also be washed out more easily after the analysis or purification. On the market today, the most common alternative NP phases are cyano- and diol-modified silica.

Silica


Kromasil Diol is based on Kromasil 60 Å which is known to have high chemical purity and excellent mechanical stability.

Derivatization

In order to obtain a stationary phase with a pronounced diol character, a tri-functional silane was used in the derivatization. Figure 1 demonstrates the difference in chromatographic properties between a diol phase with high loading (Kromasil 60-5-Diol), a competitor with low loading, and Kromasil 60 Å silica, in the separation of six different phenols. The Kromasil Diol material with high loading provides a significantly better separation.




During the development ²⁹Si CP/MAS NMR and elemental analysis have been run in order to characterize the bonded phase, together with a chromatographic evaluation. The surface coverage is high, a typical value is 3.5 μmol/m². The aim was to have a diol phase with "full" polymerization, which means no T¹ signals, in order to obtain a stationary phase with low silanol activity and high reproducibility. Figure 2 illustrates that this goal has been achieved due to the complete absence of the T¹ signal. The same manufacturing process is applied for both analytical and preparative materials. This is a critical factor when scaling up from analytical to process scale.




Comparison Study and Applications


A comparison of Kromasil Diol with other commercial diol phases, using a test mixture of five different phenols, clearly shows that Kromasil Diol is superior regarding selectivity and symmetry, see Figure 3.



In order to investigate if Kromasil Diol is a complement to the standard NP phases, Cyano and bare silica, the selectivity and elution order of a mixture of different phenols and anilines were studied. Figure 4 clearly illustrates that each phase has its own unique selectivity. Interesting studies are also reported in the literature that illustrates this phenomenon.^{1,2}



Kromasil Diol is also an interesting stationary phase for RP chromatography, offering alternative selectivity and wettability under 100% aqueous conditions. Figure 5 illustrates an application under RP conditions.



References

- [1] Waksmundzka-Hajnos, M., *J. Liq. Chrom. & Rel. Tech.*, **2004**, 2247-2267.
- [2] Waksmundzka-Hajnos, M., *J. Chromatogr. A.*, **2001**, 919, 39-50.