

Investigating the Use of Achiral SFC as an Evolving Tool in the Support of Discovery Chemistry at Novartis

Joseph Twomey, Thierry Mann

Novartis Institutes for BioMedical Research, Inc.


Waters SFC Symposium on Innovative Achiral SFC Solutions for Analysis and Purification
April 20, 2010



Outline of Presentation


- Background of SFC at NIBR
 - Cambridge support
- Normal Phase LC support
 - Single and parallel techniques
- Achiral SFC support
 - Why SFC?
 - SFC vs. NPLC
- SFC for OA Rxn Monitoring

SFC at NIBR Global (Chiral Separations)



GDC Separations Lab -Our Role

Chiral and Normal Phase Separations AND GDC PROJECT SUPPORT






Separations Group NIBRI Cambridge

- Supports the discovery chemistry groups in Cambridge
 - 170 Chemists working in 2 sites
 - An imbedded part of the chemistry department
- 3 members split across 2 buildings
- Service based group
 - Take on difficult separations
 - Develop enabling technologies
- Predominately chiral support via SFC and LC
 - Chemistry mandate to use more NPLC internally

Normal Phase Separations Support

- Increasing number of requests for analytical and preparative separations
- SFC separations 2010
 - Analytical separations: 41
 - Preparative separations: 17 (50 - 2g+)
- Shifting from LC to SFC methods in progress
 - All current submissions are being run SFC first, traditionally run NPLC
- SFC provides many advantages
 - Faster runs and turn around times
 - Lower solvent usage
 - Integrated into Chiral screens

Single column NPLC systems



- Open access chemistry support tool
- Silica, Cyano, Diol, 2-Ethylpyridine, Amino, Phenyl columns
- 4.6 x 100mm, 5um
- 1-3 ml/min
 - Hep/EtOH
 - 0.2% DEA
 - 5-60% Gradients
- 5-10 min runs

Parallel NPLC System

- 8 columns run in parallel (1.25min / col. / inj.)
- 4.6 x 100mm, 5µm, 10 min run times
 - 1.25ml/min flow rates
- Silica, Cyano, Diol, 2-EP, Amino, Phenyl, PVA-Sil, other
 - Same systems as Chiral LC screen, automatic switching of columns via contact closures
- Gradients with Hep/(EtOH/IPA), +/- 0.2% DEA
 - After development in parallel final methods established via single column systems
- Prep systems available are also used for chiral work

Parallel NPLC Output


OAA STD

5-60%EtOH, 1.25ml/min, 4.6x100mm, UV220, B2

SFCNPLCTEST_056


13-Apr-2010

Diode Array
Range: 2.037




SFCNPLCTEST_055

Diode Array
Range: 1.772




SFCNPLCTEST_054

Diode Array
Range: 2.236




SFCNPLCTEST_053

Diode Array
Range: 2.37




SFCNPLCTEST_052

Diode Array
Range: 2.043




SFCNPLCTEST_051

Diode Array
Range: 4.114




SFCNPLCTEST_050

Diode Array
Range: 2.227



SFCNPLCTEST_049


Diode Array
Range: 1.9



Why use SFC for Prep of Achiral samples?


- SFC is on average 3 times faster than HPLC
 - *1 SFC Instrument does what 3 HPLC instruments do*
- Solvent costs reduced by 60-70 %
 - 1 Liter CO₂ = ~ 1 \$
 - *1 Liter Acetonitrile = ~ 20-30 \$ (around 50'000 L/year)*
 - *1 Liter Heptane = ~ 20 \$ (around 10'000 L/year)*
- Solvent removal reduced by about 70 % (CO₂ evaporates)
- Less risk of degradation of purified compound
 - SFC does not need acidic additive in the mobile phase
- Reduced organic solvent consumption: Green chromatography
 - *Critical VOC situation can be improved*
 - *Solvent restriction issue, especially in Boston area*
- Safety: CO₂, main component of the mobile phase (60-90%), is non flammable (fire extinguisher) and is much less toxicity

SFC Screening Process



- Premade sample lists
- Gradients 5-55% Modifier
- MeOH and IPA
 - 0.2% DEA or no modifier
- Princeton Chromatography Columns
 - 4.6 x 100 mm
 - Silica, CN, Diol, 2-EP
- 5ml/min, 6min runs
- 4 columns x 4 solvent systems

SFC Screening Output




- Screening allows for rapid selection of:
 - Columns
 - Solvents
- Isocratic methods quickly selected and tested
- Integrated into Chiral workflow and screens
 - Achiral part of Chiral screening process
 - Avoids surprises

Prep SFC at NIBRI



- Thar 80 systems (2)
 - Single channel analytical system with 10 column oven (2)
- BDS CO2 delivery system
- Princeton Chromatography Columns
 - 20 x 150mm, 5um
 - Silica, Cyano, Diol, 2-EP
- Find we can run 5-10% less solvent than analytical runs and get similar separations
 - May result in longer run times but better separations
- Potential to stack injections or run Prep gradients

SFC vs NPLC Prep Separation Comparison

- How does it compare using a real sample?
- 2.8 g of sample submitted
- Small impurity to remove
- Split into 2 batches
- MD was run on both SFC and NPLC simultaneously
- Princeton Cyano column used in SFC
 - 20 x 150mm 5um





SFC Prep Separation of Sample



- 78 inj. of 2ml of MeOH
 - ~16mg/inj.
- Run time 8 min. / inj.
 - 10.4 hrs total run time
- 70g/min flow rate
- 10% IPA 0.2% DEA
- 4.4L Solvent Used
 - ~1100ml of fraction collected
- 100% pure, 1.45g recovery

NPLC Separation of Sample




- Princeton Diol column
- 82 inj of 1ml EtOH
 - ~15mg/inj.
- Run time 12min
 - 16.4 hours of run time
- 25ml/min flow rate
- 24.6L solvent used
 - 4.1L of fractions
- 1.14g recovery, 98%+

Highlights of SFC vs NPLC


- SFC was 32% faster in total run time
- Both had roughly the same loading and number of injections made
 - NPLC was run over 3 days, SFC completed in 1 day plus PM run
- SFC collected 75% less solvent
- **SFC used 18% of the solvent of NPLC**
- Recovery effected by division of sample possibly not being equal
- Both resulted in high purity samples

Could the SFC have been faster?




- If this sample had been stacked how much faster would it have been
- 3 minutes blank time
 - ~30% faster if stacked
- Stacking has the potential to save
 - ~3hrs total run time
 - ~1.2L of solvent
- Stacking sets up the same as in chiral separations

SFC Prep Separation Example 2





- 140mg submitted
- Chemist had made several attempts to work up
- Run through initial screen
- Successfully separated on 13/16 combinations using gradients of 5-55%
- 2-EP column selected for Prep runs
 - MeOH without modifier

SFC Prep Separation Example 2 Cont.



- 10 injections
 - 14mg / inj.
- 65g/min flow rate
 - 10% MeOH
- 80 min total run time
 - 8 min cycle time
- 600ml of solvent used
 - F1 ~40ml collected
 - F2 ~80ml collected
- 100% pure product

SFC Prep Separation Example 3



- 180mg submitted
 - 17 inj of 10-11mg
- Analytical Method
 - Diol 4.6 x 100mm
 - 40% IPA 0.2% DEA
- Prep Method
 - Diol 20 x 150mm
 - 35% IPA 0.2% DEA
- 4 min. run, 70g/min
- Totals
 - Time 68min
 - ~1.6L solvent
 - ~400 ml fraction

OA-SFC as a complimentary approach to Reaction monitoring

- Is there an optimum mobile phase?
- Is there a best column chemistry / manufacturer
- Will the data be as rich as reversed phase?
- Are we more likely to precipitate?
- Can we make the interface as friendly as LC-MS
 - Mass lynx OA software on Thar system
- How do we interpret the data when we see differences

SFC OA conditions

- Column 1: 3mmx 50 mm ES Industries 2-EP 5 micron
- Column 2: 3mmx 50 mm ES Industries Chromegabond NPI 5 micron
- 10mM ammonium formate in MeOH
- 10% to 40% in 1.5 minutes hold for 0.5 minutes.
- Total flow-3 ml/ min, ~1mL of MeOH per injection.

RP OA conditions

- Column 1: 3mmx 33 mm GL Sciences Inertsil C8 3 micron
- Column 2: 3mmx 33 mm GL Sciences Inertsil C18 3 micron
- 5mM ammonium formate and MeOH:ACN
- 5% to 95% in 1.6 minutes hold for 0.4 minutes.
- Total flow-2 ml/ min, ~5mL of solvent per injection.


RP analysis on C8 2mL/min

Method: **INERTSIL333_C8_NPH_2ml** Date: **13-Apr-2010**

Time: **15:25:50**

Notes:

Chromatograms:




SFC-UV Analysis of previous sample


AWM:
Mobile Phase: 10-40% MeOH in CO₂
oatest_011

Sample ID: 13222-029
Instrument: sfc-100

Submitter:
Column: NPI
Diode Array
240
Range: 9.52e-2



oatest_010




RP Analysis on C18 polar gradient

Method: 333_C18_APH_POLAR_2ml Date: 14-Apr-2010

Time: 12:01:09

Notes:

Chromatograms:



Peak Table:


#	Time	BPM+	BPM-	UV	ELSD	Target
1	1.17	282.87	280.96	10.8	0.0	No
2	1.24	324.83	322.96	61.7	60.7	No

SFC-UV Analysis of previous sample

AWM:
Mobile Phase: 10-40% MeOH in CO₂
oatest_003


Sample ID: 12886-053-3
Instrument: sfc-100

Submitter:
Column: NPI
Diode Array
240
Range: 7.674e-1



oatest_002

Diode Array
240
Range: 8.028e-1




SFC for OA Reaction Monitoring Summary

- Potential to add valuable information
- ~20% of solvent used compared to RP
- No single column chemistry is best.
- Instrumentation Robust enough?
- Instrumentation is large
 - Space available for OA instruments is limited
- BPR stability could be an issue
- We will move forward with our partners in A.S.

Thar 5x SFC System Evaluation

- Achiral columns
 - Silica, Diol, Cyano, 2-Ethylpyridine, Diol-HL
- MeOH 0.2% DEA
 - 5-55% Gradients
- OA-LCMS samples
- 6 minutes for 5 columns
 - 1.25 min. / col. / inj.
- System is over 7' tall on bench



Where are we going next?

- Evaluating system suitability for use in OA SFC-MS
 - Possible replacement for some older LC-MS systems
 - Use in core labs for faster analysis times and solvent issues
 - Looking for universal conditions
- Prep SFC-MS
 - Can we set up a simple assisted use systems for chemist use
- Parallel SFC
 - Demo of Thar X5 systems under way
 - Can we make our 8 Column system work via SFC
- CO₂ supply issues, will this equipment need to be centralized to supply it?

Acknowledgements

- Thierry Mann
- Melissa Grondine
- John Reilly
- Jay Larrow
- Scale Up / Separations Team